

ORDER FORM




Q200R with Captain's Seat

EIPW35

Customer		Ship To	
Account #		Name	
Date		Address	
PO#		City	
Buyer		State	Zip
Mark For		Phone	

Client Measurements

a) Top of head to seat	in
b) Top of shoulder to seat	in
c) Underarm to seat	in
d) Elbow to seat	in
e) Top of leg to seat	in
f) Shoulder width	in
g) Chest width	in
i) Hip width	in
j) Actual sitting depth	in
m) Foot length	in
n) Knee to foot	in
o) Seat width	in



Notes

Chair Specifications

HCPC	K0848 / K0849 / K0856 / K0857
Weights starting at: (with Batteries)	260 lbs
Power Seating	22° Tilt
Seat Type	Rehab Seat w/ Fold Down Back and 25° Adjustable Recline
Seat Size Range	Width 16" to 20" Depth 16" to 22"
Seat-to-Floor Height Range	17" to 18"

Mid-Wheel Drive Base Specifications

Drive Wheel Position	Rear-wheel drive
User Weight Capacity	300 lbs.
Base Size	Width 23", Length 42"
Minimum Turning Radius	18.5" (excludes footplates)
Drive Control	R-net PM 90A / VR-2 90A
Driving Range	Up to 15.5 miles
Drive Motors	2-pole, 6 mph / 10 kph
Drive Wheels	12" w/ Flat Free Inserts 14" w/ Pneumatic Inserts 14" w/ Flat Free Inserts
Front Casters	9" w/ Pneumatic Inserts 9" w/ Flat Free Inserts
Charger	8A Dual Mode (Compact), 10A
Battery Size	Group 22
Color	7 color accents
Ground Clearance	3"

BASE			
OPTION	HCPC	DESCRIPTION	RETAIL
Group 3 Standard Duty Rear Wheel Drive Base <i>End Item Number EIPW35 (Captain's Seat)</i>			\$
350PC304	K0849	Quickie Q200 R NCC No Power Option, Captain's Seat and Back Included in Base Pricing : Q200 R Base, VR-2 90amp Electronics for No Power Seat Functions, 4-wheel Suspension System, 12" Drive Wheels with Dual Density Inserts, 9" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 10kph / 6mph, Weight Capacity = 300 lbs., STFH = 17".	\$7,793.00
350PC305	K0857	Quickie Q200 R SCC Single Power Option, Captain's Seat and Back Included in Base Pricing : Q200 R Base, VR-2 90amp Electronics for No Power Seat Functions, 4-wheel Suspension System, 12" Drive Wheels with Dual Density Inserts, 9" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 10kph / 6mph, Weight Capacity = 300 lbs., STFH = 17".	\$7,793.00
350PC306	K0848	Quickie Q200 R NCC No Power Option, Captain's Seat Back and Solid Seat Included in Base Pricing : Q200 R Base, VR-2 90amp Electronics for No Power Seat Functions, 4-wheel Suspension System, 12" Drive Wheels with Dual Density Inserts, 9" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 10kph / 6mph, Weight Capacity = 300 lbs., STFH = 17".	\$7,793.00
350PC307	K0856	Quickie Q200 R SCC Single Power Option, Captain's Seat Back and Solid Seat Included in Base Pricing : Q200 R Base, VR-2 90amp Electronics for No Power Seat Functions, 4-wheel Suspension System, 12" Drive Wheels with Dual Density Inserts, 9" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 10kph / 6mph, Weight Capacity = 300 lbs., STFH = 17".	\$7,793.00
User Weight			\$
350WG23	–	Up to 300 lbs.	NC
Motor Package			\$
350PF107	–	2-pole 6 mph / 10 kph	NC
Drive Wheels			\$
350DW77	–	Drive Wheel, 12", Solid, Knobby, Black	NC
350DW78	–	Drive Wheel, 14", Pneumatic, Knobby, Black	\$279.00
350DW79	–	Drive Wheel, 14", Solid, Knobby, Black	NC
Caster Wheels			\$
350CW121	–	Caster Wheels 9", Pneumatic, Black	\$279.00
350CW122	–	Caster Wheels 9", Solid, Black	NC
Color Accents			\$
350EZ20	–	Apple Green	NC
350EZ22	–	Black	NC
350EZ21	–	Blue	NC
350EZ23	–	Orange	NC
350EZ24	–	Red	NC
350EZ19	–	White	NC
350EZ25	–	Silver	NC

BASE			
OPTION	HCPC	DESCRIPTION	RETAIL
ISO Transit Tie Downs			\$
350T050	–	ISO Transit Package	NC
LED Lights			\$
350LI04	–	LED Lights and Indicators	\$1,341.00
350LI02	–	Omit LED Lights	NC
350LI08	–	LED Fender Lights	\$542.00
Seat Heights <i>(Seat height is measured from front of seat pan to floor with 0° frame angle.)</i>			\$
350HT55	–	Seat Height 17" (N/A with Tilt)	NC
350HT56	–	Seat Height 18"	NC
Batteries and Chargers			\$
Group 22 Battery			
350AG31	E2361	GP22 Flat Top Pair w/Straps	\$587.00
350AG35	E2361	GP22/50Ah Flat Top Pair w/Straps <i>(Required for Tilt)</i>	\$652.00
Chargers			
350BJ30	–	10A Dual Mode (Compact)	\$106.00
350BJ19	–	8A Dual Mode (Compact)	NC

SEATING

OPTION	HCPC	DESCRIPTION	RETAIL
Seating Options			\$
Pre Set Seat Frame Angle (No Power Seating Options, Increases Seat to Floor Height)			\$
350FA10	–	0° (Adds 0" to STFH)	NC
350FA3	–	3° (Adds 1" to STFH)	NC
350FA4	–	6° (Adds 2" to STFH)	NC
Backrest Frame Type			\$
350BF08	–	Captain's Seat Frame (Includes back cushion)	NC
Seat Type			\$
350ST80	–	JAY®	NC
350ST77	–	Captain's Seat Cushion	NC
350ST75	–	Omit	NC
Seat Size			\$
Seat Pan Width			
350R53	–	16"	NC
350R55	–	18"	NC
350R57	–	20"	NC
		Available Depths	
		16" to 18"	
		16" to 20"	
		18" to 22"	
Seat Pan Depth (Front of Seat Pan to Back Structure)			
350RF51	–	16"	NC
350RF45	–	17"	NC
350RF46	–	18"	NC
350RF47	–	19"	NC
350RF48	–	20"	NC
350RF49	–	21"	NC
350RF50	–	22"	NC
JAY® Seat Cushion			\$
350CU59	E2607	JAY® Ion®	\$441.00
350CU61	E2607	JAY® Union®	\$498.00
350CU64	E2622	JAY® Fusion®	\$492.00
350CU104	E2622	JAY® Fusion® with Cryo® Technology	\$603.00

POSITIONING SUPPORTS & ACCESSORIES

OPTION	HCPC	DESCRIPTION	RETAIL
Armrest Structure			\$
350EG37	–	Flip Back <i>(Includes Waterfall Arm Pads) (Height Adjustable 11" to 14")</i>	NC
Swing-Away Legrest (Pair)			\$
350H177	–	SEDEO PRO Legrests, Swing Away, 70° <i>(14" to 21" Leg Length)</i>	NC
350H178	E0990	SEDEO PRO Legrests, Manual Elevating, Swing Away <i>(Includes Padded Calf Pads) (14".5 to 21.5" Leg Length)</i>	\$391.00
350H179	E1010	SEDEO PRO Legrests, Power Elevating, Swing Away <i>(Includes Padded Calf Pads) (15.5" to 22.5" Leg Length)</i>	\$2,403.00
350H180	–	SEDEO PRO Hanger Receivers Only	NC
Swing Away Legrest Calf Supports			
350DI42	–	Calf Support, Depth Adjustable, Padded <i>(6.5" to 8" Wide, Depending On Chair Width, x 9.5" Tall)</i>	\$441.00
Footplates			
350F370	K0040	Angle Adjustable Flip Up Footplates	\$189.00
Accessories			
350DI44	E0951	Heel Loops	\$56.00
350KG4	–	Knee Pads, Integral Foam <i>(2.5" Wide x 4.5" Tall)</i>	\$224.00
SEDEO Lower Extremity Swing Away Support			
350VD01	E1020	Lower Extremity Support, Padded, Swing Away, SEDEO, Right	\$441.00
350VD02	E1020	Lower Extremity Support, Padded, Swing Away, SEDEO, Left	\$441.00
Pelvic Positioning Belts			
350A837	–	Positioning Belt, Auto Style Buckle, 50 mm 1600 mm / 2"	NC
350A839	–	Positioning Belt, Aircraft Buckle, 50 mm x 1300 mm / 2"	NC
350A840	E0978	Positioning Belt, Aircraft Buckle, Padded 50 mm x 1300 mm / 2"	\$84.00
Accessories			\$
Backpack			
350A1	–	Black Utility	\$61.00
Captain's Seat Accessories			
350A608	E2207	Cane / Crutch Holder	\$290.00
350A610	K0108	Walker Holder	\$290.00
Medial Thigh Support			
3508A80	E0957	Medial Thigh Support	\$189.00
3508A81	E1028	Pommel Swing-Away Hardware	\$279.00
Misc. Accessories			
3508A52	–	Ram Mount Phone	\$96.00
3508A53	–	Ram Mount Tablet	\$391.00
3508A54	–	Ram Mount Cup Holder	\$96.00
3508A58	–	Rear View Mirror	\$96.00
3509A10	–	Utility Hooks	\$58.00

JOYSTICK / ELECTRONICS

OPTION	HCPC	DESCRIPTION	RETAIL
VR-2 Joystick			\$
350J127	–	VR-2 Joystick	NC
R-net Joystick (Remote)			\$
350J121	–	Standard LED Joystick (w/o Screen Display)	NC
350J138	–	Color Joystick (2.3" Color LCD Screen Display, CJSM1)	\$279.00
350J118	–	High Resolution Color Joystick, (CJSM2+BT) 3.5" Color LCD Display, Paddle Switches, Programmable Shortcut Keys, ICONs, Stereo Jack, IR and dual (IOS/Android) Bluetooth control	\$1,178.00
350J119	–	Omit Remote Control	NC
Joystick Mounting			\$
350JM34	–	Fixed on Armrest	NC
350JM35	E1028	Parallel Swing-Away	\$330.00
Joystick Mount (Location)			\$
350HM16	–	Right Hand Armrest	NC
350HM17	–	Left Hand Armrest	NC
Display Language Selection R-net			\$
350L1	–	English	NC
350L3	–	French	NC
350L7	–	Spanish	NC
Joystick Handle (Optional) N/A with HD Joystick			\$
350JH13	E2323	Ball Handle	\$141.00
350JH15	E2323	T-Handle	\$141.00
Electronic Programmings			\$
350PR10	–	R-net PC Programmer w/ Dongle	\$587.00
R-net Power Modules			\$
350ER12	–	Power Module, VR-2 90A, Non-Expandable	NC
350ER13	–	Power Module, R-net, 90A, Expandable	\$1,335.00
350EX15	E2377	Expandable 90A Controller - R-net	\$704.00
350EX16	E2313	Harness for 90A Expandable Controller	\$631.00

JOYSTICK / ELECTRONICS

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

Seating Control

\$

Powered Seating Control

350E040	–	Omit Thru-Drive <i>(Requires Switch Control and Standard LED Joystick Requires Ctrl+5) (N/A w/ VR-2 electronics)</i>	NC
350E037	E2310	Thru-drive Control for 1 Actuator	\$2,146.00
350E038	E2311	Thru-drive Control for 2+ Actuators	\$3,884.00

Switch Control *(Available with Omit Thru-Drive only) (Requires use of Omni2 or Joystick Jacks)*

350E033	–	Dual Toggle for 1 Actuator	NC
350E034	–	2.5" Single Button Switch for 1 Actuator <i>(Button switches are not mounted)</i>	NC
350E035	–	Switch-It Mini 1 button for 1 Actuator	NC
350E036	–	Switch-It Mini 2 button for 1 Actuator	NC
350E041	–	Ctrl+5 Switch Box for 2 Actuators or more	NC

Additional Power Seating Control

350QG11	–	Ctrl+5 Switch Box with Free Assignable Functions	\$1,397.00
---------	---	--	------------

Dual Toggle Mounting Location

350TN1	–	Toggle Mount, Right	NC
350TN2	–	Toggle Mount, Left	NC

Ctrl+5 Mounting Location

350EY19	–	Right Armrest	NC
350EY20	–	Left Armrest	NC

Ctrl+5 Switch Type

350EY22	–	Push Buttons	NC
350EY23	K0108	Toggle Switches	\$335.00

R-net Programming

\$

Assignable Buttons

Assignable Button feature allows for customization of joystick buttons and/or Ctrl+5 Switches allowing access to various wheelchair and seating functions. When enabled, chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle.

350PM12	–	Include Assignable Buttons	NC
---------	---	----------------------------	----



Q200 R with Captain's Seat

\$

Base Price	
Total of additional equipment	
Packaging and transport	
Total Price	

Date

Signature

Month Day Year