



# JAY® J3 Form Back Support

E2617

January 2024

Mark For: \_\_\_\_\_  
 Date: \_\_\_\_\_  
 Dealer Acct #: \_\_\_\_\_  
 Dealer: \_\_\_\_\_  
 Dealer Contact: \_\_\_\_\_  
 Dealer Address: \_\_\_\_\_  
 Dealer City: \_\_\_\_\_ ST: \_\_\_\_\_ ZIP: \_\_\_\_\_  
 Dealer Phone: ( ) \_\_\_\_\_ Fax: ( ) \_\_\_\_\_  
 Confirmation Email: \_\_\_\_\_  
 Confirm Via:  Fax  Email

Submitting for:  Quote  Order  
 PO#: \_\_\_\_\_

### ADDITIONAL SHIPPING INFORMATION

Ship To: \_\_\_\_\_  
 Attention: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 Ship To City: \_\_\_\_\_ ST \_\_\_\_\_ ZIP: \_\_\_\_\_  
 Ship To Phone: ( ) \_\_\_\_\_ Fax ( ) \_\_\_\_\_

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

## JAY® J3 Form Back Support - E2617

### Specifications

Shell Contour (Without Foam) 2.2"  
 Height 12" to 24.0"  
 Cover Material 3DX™ Spacer w/Lycra® Stretch  
 JAY J3 Hardware SH, TH, FS  
 HCPCS E2617












### Design



- Utilizes the JAY® J3 PA design and hardware
- Three foam-in-place densities: extra soft, soft and medium
- Outer cover with offset band allows for the expansion of foam
- Zipper location at top of cover to allow for access during installation (molding)



### STEP 1 - JAY J3 Form Posterior (PA) Back, Hardware Selection, and Accessories (pages 2-4)

BACKREST WIDTH	14"	16"	18"	20"
BACKCANE WIDTHS	14"-16"	16"-18"	18"-20"	20"-22"
MAX. USER TRUNK WIDTH	13"	15"	17"	19"

Back Height	Mid Thoracic (MT)						
	12.0" (S)	13.5" (M)	15.0" (T)				
		J3CUSTPAMTSxxSH	<input type="checkbox"/> 14 \$1,087	<input type="checkbox"/> 16 \$1,087	<input type="checkbox"/> 18 \$1,087	<input type="checkbox"/> 20 \$1,087	
		J3CUSTPAMTSxxTH	<input type="checkbox"/> 14 \$1,138	<input type="checkbox"/> 16 \$1,138	<input type="checkbox"/> 18 \$1,138	<input type="checkbox"/> 20 \$1,138	
		J3CUSTPAMTSxxFS	<input type="checkbox"/> 14 \$1,138	<input type="checkbox"/> 16 \$1,138	<input type="checkbox"/> 18 \$1,138	<input type="checkbox"/> 20 \$1,138	
		J3CUSTPAMTxxSH	<input type="checkbox"/> 14 \$1,087	<input type="checkbox"/> 16 \$1,087	<input type="checkbox"/> 18 \$1,087	<input type="checkbox"/> 20 \$1,087	
		J3CUSTPAMTxxTH	<input type="checkbox"/> 14 \$1,138	<input type="checkbox"/> 16 \$1,138	<input type="checkbox"/> 18 \$1,138	<input type="checkbox"/> 20 \$1,138	
		J3CUSTPAMTxxFS	<input type="checkbox"/> 14 \$1,138	<input type="checkbox"/> 16 \$1,138	<input type="checkbox"/> 18 \$1,138	<input type="checkbox"/> 20 \$1,138	
			J3CUSTPAUTSxxSH	<input type="checkbox"/> 14 \$1,087	<input type="checkbox"/> 16 \$1,087	<input type="checkbox"/> 18 \$1,087	<input type="checkbox"/> 20 \$1,087
			J3CUSTPAUTSxxTH	<input type="checkbox"/> 14 \$1,138	<input type="checkbox"/> 16 \$1,138	<input type="checkbox"/> 18 \$1,138	<input type="checkbox"/> 20 \$1,138
			J3CUSTPAUTSxxFS	<input type="checkbox"/> 14 \$1,138	<input type="checkbox"/> 16 \$1,138	<input type="checkbox"/> 18 \$1,138	<input type="checkbox"/> 20 \$1,138
			J3CUSTPAUTMxxSH	<input type="checkbox"/> 14 \$1,087	<input type="checkbox"/> 16 \$1,087	<input type="checkbox"/> 18 \$1,087	<input type="checkbox"/> 20 \$1,087
			J3CUSTPAUTMxxTH	<input type="checkbox"/> 14 \$1,138	<input type="checkbox"/> 16 \$1,138	<input type="checkbox"/> 18 \$1,138	<input type="checkbox"/> 20 \$1,138
			J3CUSTPAUTMxxFS	<input type="checkbox"/> 14 \$1,138	<input type="checkbox"/> 16 \$1,138	<input type="checkbox"/> 18 \$1,138	<input type="checkbox"/> 20 \$1,138
		J3CUSTPAUTTxxSH	<input type="checkbox"/> 14 \$1,087	<input type="checkbox"/> 16 \$1,087	<input type="checkbox"/> 18 \$1,087	<input type="checkbox"/> 20 \$1,087	
		J3CUSTPAUTTxxTH	<input type="checkbox"/> 14 \$1,138	<input type="checkbox"/> 16 \$1,138	<input type="checkbox"/> 18 \$1,138	<input type="checkbox"/> 20 \$1,138	
		J3CUSTPAUTTxxFS	<input type="checkbox"/> 14 \$1,138	<input type="checkbox"/> 16 \$1,138	<input type="checkbox"/> 18 \$1,138	<input type="checkbox"/> 20 \$1,138	
		J3CUSTPASHSxxSH	<input type="checkbox"/> 14 \$1,087	<input type="checkbox"/> 16 \$1,087	<input type="checkbox"/> 18 \$1,087	<input type="checkbox"/> 20 \$1,087	
		J3CUSTPASHSxxTH	<input type="checkbox"/> 14 \$1,138	<input type="checkbox"/> 16 \$1,138	<input type="checkbox"/> 18 \$1,138	<input type="checkbox"/> 20 \$1,138	
		J3CUSTPASHSxxFS	<input type="checkbox"/> 14 \$1,138	<input type="checkbox"/> 16 \$1,138	<input type="checkbox"/> 18 \$1,138	<input type="checkbox"/> 20 \$1,138	
	J3CUSTPASHMxxSH	<input type="checkbox"/> 16 \$1,087	<input type="checkbox"/> 18 \$1,087	<input type="checkbox"/> 20 \$1,087			
	J3CUSTPASHMxxTH	<input type="checkbox"/> 16 \$1,138	<input type="checkbox"/> 18 \$1,138	<input type="checkbox"/> 20 \$1,138			
	J3CUSTPASHMxxFS	<input type="checkbox"/> 16 \$1,138	<input type="checkbox"/> 18 \$1,138	<input type="checkbox"/> 20 \$1,138			
	J3CUSTPASHTxxSH	<input type="checkbox"/> 18 \$1,087	<input type="checkbox"/> 20 \$1,087				
	J3CUSTPASHTxxTH	<input type="checkbox"/> 18 \$1,138	<input type="checkbox"/> 20 \$1,138				
	J3CUSTPASHTxxFS	<input type="checkbox"/> 18 \$1,138	<input type="checkbox"/> 20 \$1,138				



### STEP 2 - Foam-In-Place Selection

Note: Foam-In-Place must be ordered in combination with backrest & hardware (STEP 1), separate line during ordering process. Will ship separately from backrest & hardware, any JAY accessories, or bundled wheelchair. For spare parts please see [www.sunparts.us](http://www.sunparts.us)



✓	Part #	Description	Price
<input type="checkbox"/>	251600	SUNMATE FOAM-IN-PLACE X-SOFT	INCLUDED
<input type="checkbox"/>	251601	SUNMATE FOAM-IN-PLACE SOFT	INCLUDED
<input type="checkbox"/>	251602	SUNMATE FOAM-IN-PLACE MEDIUM	INCLUDED

## JAY® J3 Back Accessories

All JAY J3 Plus Back accessories are optional. All accessories will be factory mounted as much as possible to insure proper shipping (If selected below).

### SHIPPING INSTRUCTIONS

	Option #	HCPCS	Definition	Price		Option #	HCPCS	Definition	Price
<input checked="" type="checkbox"/>	N/A		Ship Back on Chair (Chair to be ordered separately)	NC	<input type="checkbox"/>	N/A		Ship accessories mounted on back	NC

### Headrests

Use this section for selecting headrest. Head supports are not available on Mid Thoracic Back Shells

#### STEP 1

#### Select your HEADREST POST

Please select one box from the below **Headrest Mount & Bracket** choices.

<input type="checkbox"/>	J3HEDX	E1028	Axys Mount & Bracket	\$254.00	<input type="checkbox"/>	J3HEDO	E1028	ONYX Mount & Bracket	\$229.00
<input type="checkbox"/>	J3HEDL	E1028	LINX Mount & Bracket	\$324.00	<input type="checkbox"/>	J3HEDC	E1028	Cobra Mount & Bracket	\$320.00
<input type="checkbox"/>	J3HEDL2	E1028	LINX <sup>2</sup> Mount & Bracket	\$355.00	<input type="checkbox"/>	J3HEDF	E1028	Cobra Flip Back Mount and Bracket	\$385.00



#### STEP 2

#### Select your HEADREST PAD

Please select one box from the below **Headrest Pad** choices.

<input type="checkbox"/>	14	E0955	J3PLSH14 - 14" Plush Pad	\$201.00	<input type="checkbox"/>	AAPS	E0955	J3AAPS - Adjust-A-Plush Standard (Lg)	\$301.00
<input type="checkbox"/>	10	E0955	J3PLSH10 - 10" Plush Pad	\$186.00	<input type="checkbox"/>	AAPM	E0955	J3AAPM - Adjust-A-Plush Medium	\$301.00
<input type="checkbox"/>	8	E0955	J3PLSH8 - 8" Plush Pad	\$186.00	<input type="checkbox"/>	AAPN	E0955	J3AAPN - Adjust-A-Plush Narrow (Sm)	\$301.00
<input type="checkbox"/>	HUPL	E0955	J3HUPL - HeadsUp Lrg: 5.3" Lng Arms	\$358.00	<input type="checkbox"/>	CL	E0955	J3CRDLLG - Lrg Contour Cradle Pad	\$179.00
<input type="checkbox"/>	HUPM	E0955	J3HUPM - HeadsUp Med: 4.1" Sht Arms	\$358.00	<input type="checkbox"/>	CM	E0955	J3CRDLMD - Med Contour Cradle Pad	\$166.00
<input type="checkbox"/>	HUPS	E0955	J3HUPS - HeadsUp Sml: 3.5" Sht Arms	\$358.00	<input type="checkbox"/>	PL	E0955	J3CPLULG - Large C-Plush Pad	\$179.00
					<input type="checkbox"/>	PM	E0955	J3CPLUMD - Medium C-Plush Pad	\$179.00



#### Universal Mounting Plate

<input type="checkbox"/>	J3UHM	K0108	Universal Headrest Mounting Plate Only	\$60.00
--------------------------	-------	-------	--	---------

## Lateral Supports

Use this section for selecting laterals and hardware. Lateral supports are not available on Lower Thoracic Back Shells

### Lateral Thoracic Support (Fixed)

*Includes fixed mounting hardware. Lateral pads use 3DX covers with 1/2" Sunmate foam in the core*

<input type="checkbox"/>	J3LAT4X4FX	E0956	LATERAL 4" Depth X 4" Height FIXED	\$287.00	<input type="checkbox"/>	J3LAT6X4FX	E0956	LATERAL 6" Depth X 4" Height FIXED	\$287.00
<input type="checkbox"/>	J3LAT4X5FX	E0956	LATERAL 4" Depth X 5" Height FIXED	\$287.00	<input type="checkbox"/>	J3LAT6X5FX	E0956	LATERAL 6" Depth X 5" Height FIXED	\$287.00
<input type="checkbox"/>	J3LAT4X6FX	E0956	LATERAL 4" Depth X 6" Height FIXED	\$287.00	<input type="checkbox"/>	J3LAT6X6FX	E0956	LATERAL 6" Depth X 6" Height FIXED	\$287.00

### Swing-Away Lateral Thoracic Support

Includes: Swing-away mounting hardware, left and right lateral are independent. Different types and sizes can be selected for each side.

Left Side				
<input type="checkbox"/>	<b>JLATSA44L</b>		<b>LAT 4X4 W/ S/A MOUNTS - LH</b>	<b>\$446.00</b>
	JLAT4X4	E0956	Lateral 4" Depth X 4" Height	\$144.00
	JSAMNTL	E1028	Lateral Swing Away Mount - Left	\$302.00
<input type="checkbox"/>	<b>JLATSA45L</b>		<b>LAT 4X5 W/ S/A MOUNTS - LH</b>	<b>\$446.00</b>
	JLAT4X5	E0956	LATERAL 4" Depth X 5" Height	\$144.00
	JSAMNTL	E1028	Lateral Swing Away Mount - Left	\$302.00
<input type="checkbox"/>	<b>JLATSA46L</b>		<b>LAT 4X6 W/ S/A MOUNTS - LH</b>	<b>\$446.00</b>
	JLAT4X6	E0956	LATERAL 4" Depth X 6" Height	\$144.00
	JSAMNTL	E1028	Lateral Swing Away Mount - Left	\$302.00
<input type="checkbox"/>	<b>JLATSA64L</b>		<b>LAT 6X4 W/ S/A MOUNTS - LH</b>	<b>\$446.00</b>
	JLAT6X4	E0956	LATERAL 6" Depth X 4" Height	\$144.00
	JSAMNTL	E1028	Lateral Swing Away Mount - Left	\$302.00
<input type="checkbox"/>	<b>JLATSA65L</b>		<b>LAT 6X5 W/ S/A MOUNTS - LH</b>	<b>\$446.00</b>
	JLAT6X5	E0956	LATERAL 6" Depth X 5" Height	\$144.00
	JSAMNTL	E1028	Lateral Swing Away Mount - Left	\$302.00
<input type="checkbox"/>	<b>JLATSA66L</b>		<b>LAT 6X6 W/ S/A MOUNTS - LH</b>	<b>\$446.00</b>
	JLAT6X6	E0956	LATERAL 6" Depth X 6" Height	\$144.00
	JSAMNTL	E1028	Lateral Swing Away Mount - Left	\$302.00

Right Side				
<input type="checkbox"/>	<b>JLATSA44R</b>		<b>LAT 4X4 W/ S/A MOUNTS - RH</b>	<b>\$446.00</b>
	JLAT4X4	E0956	Lateral 4" Depth X 4" Height	\$144.00
	JSAMNTR	E1028	Lateral Swing Away Mount - Right	\$302.00
<input type="checkbox"/>	<b>JLATSA45R</b>		<b>LAT 4X5 W/ S/A MOUNTS - RH</b>	<b>\$446.00</b>
	JLAT4X5	E0956	Lateral 4" Depth X 5" Height	\$144.00
	JSAMNTR	E1028	Lateral Swing Away Mount - Right	\$302.00
<input type="checkbox"/>	<b>JLATSA46R</b>		<b>LAT 4X6 W/ S/A MOUNTS - RH</b>	<b>\$446.00</b>
	JLAT4X6	E0956	Lateral 4" Depth X 6" Height	\$144.00
	JSAMNTR	E1028	Lateral Swing Away Mount - Right	\$302.00
<input type="checkbox"/>	<b>JLATSA64R</b>		<b>LAT 6X4 W/ S/A MOUNTS - RH</b>	<b>\$446.00</b>
	JLAT6X4	E0956	Lateral 6" Depth X 4" Height	\$144.00
	JSAMNTR	E1028	Lateral Swing Away Mount - Right	\$302.00
<input type="checkbox"/>	<b>JLATSA65R</b>		<b>LAT 6X5 W/ S/A MOUNTS - RH</b>	<b>\$446.00</b>
	JLAT6X5	E0956	Lateral 6" Depth X 5" Height	\$144.00
	JSAMNTR	E1028	Lateral Swing Away Mount - Right	\$302.00
<input type="checkbox"/>	<b>JLATSA66R</b>		<b>LAT 6X6 W/ S/A MOUNTS - RH</b>	<b>\$446.00</b>
	JLAT6X6	E0956	Lateral 6" Depth X 6" Height	\$144.00
	JSAMNTR	E1028	Lateral Swing Away Mount - Right	\$302.00



### Rotational Swing-Away Lateral Thoracic Support

Includes: Rotational swing-away mounting hardware, left and right lateral are independent. Different types and sizes can be selected for each side.

Left Side				
<input type="checkbox"/>	<b>JLATRSA44L</b>		<b>LAT 4X4 W/ S/A MOUNTS - LH</b>	<b>\$502.00</b>
	JLAT4X4	E0956	Lateral 4" Depth X 4" Height	\$144.00
	JRSAMNTL	E1028	Lateral Swing Away Mount - Left	\$358.00
<input type="checkbox"/>	<b>JLATRSA45L</b>		<b>LAT 4X5 W/ S/A MOUNTS - LH</b>	<b>\$502.00</b>
	JLAT4X5	E0956	LATERAL 4" Depth X 5" Height	\$144.00
	JRSAMNTL	E1028	Lateral Swing Away Mount - Left	\$358.00
<input type="checkbox"/>	<b>JLATRSA46L</b>		<b>LAT 4X6 W/ S/A MOUNTS - LH</b>	<b>\$502.00</b>
	JLAT4X6	E0956	LATERAL 4" Depth X 6" Height	\$144.00
	JRSAMNTL	E1028	Lateral Swing Away Mount - Left	\$358.00
<input type="checkbox"/>	<b>JLATRSA64L</b>		<b>LAT 6X4 W/ S/A MOUNTS - LH</b>	<b>\$502.00</b>
	JLAT6X4	E0956	LATERAL 6" Depth X 4" Height	\$144.00
	JRSAMNTL	E1028	Lateral Swing Away Mount - Left	\$358.00
<input type="checkbox"/>	<b>JLATRSA65L</b>		<b>LAT 6X5 W/ S/A MOUNTS - LH</b>	<b>\$502.00</b>
	JLAT6X5	E0956	LATERAL 6" Depth X 5" Height	\$144.00
	JRSAMNTL	E1028	Lateral Swing Away Mount - Left	\$358.00
<input type="checkbox"/>	<b>JLATRSA66L</b>		<b>LAT 6X6 W/ S/A MOUNTS - LH</b>	<b>\$502.00</b>
	JLAT6X6	E0956	LATERAL 6" Depth X 6" Height	\$144.00
	JRSAMNTL	E1028	Lateral Swing Away Mount - Left	\$358.00

Right Side				
<input type="checkbox"/>	<b>JLATRSA44R</b>		<b>LAT 4X4 W/ S/A MOUNTS - RH</b>	<b>\$502.00</b>
	JLAT4X4	E0956	Lateral 4" Depth X 4" Height	\$144.00
	JRSAMNTR	E1028	Lateral Swing Away Mount - Right	\$358.00
<input type="checkbox"/>	<b>JLATRSA45R</b>		<b>LAT 4X5 W/ S/A MOUNTS - RH</b>	<b>\$502.00</b>
	JLAT4X5	E0956	Lateral 4" Depth X 5" Height	\$144.00
	JRSAMNTR	E1028	Lateral Swing Away Mount - Right	\$358.00
<input type="checkbox"/>	<b>JLATRSA46R</b>		<b>LAT 4X6 W/ S/A MOUNTS - RH</b>	<b>\$502.00</b>
	JLAT4X6	E0956	Lateral 4" Depth X 6" Height	\$144.00
	JRSAMNTR	E1028	Lateral Swing Away Mount - Right	\$358.00
<input type="checkbox"/>	<b>JLATRSA64R</b>		<b>LAT 6X4 W/ S/A MOUNTS - RH</b>	<b>\$502.00</b>
	JLAT6X4	E0956	Lateral 6" Depth X 4" Height	\$144.00
	JRSAMNTR	E1028	Lateral Swing Away Mount - Right	\$358.00
<input type="checkbox"/>	<b>JLATRSA65R</b>		<b>LAT 6X5 W/ S/A MOUNTS - RH</b>	<b>\$502.00</b>
	JLAT6X5	E0956	Lateral 6" Depth X 5" Height	\$144.00
	JRSAMNTR	E1028	Lateral Swing Away Mount - Right	\$358.00
<input type="checkbox"/>	<b>JLATRSA66R</b>		<b>LAT 6X6 W/ S/A MOUNTS - RH</b>	<b>\$502.00</b>
	JLAT6X6	E0956	Lateral 6" Depth X 6" Height	\$144.00
	JRSAMNTR	E1028	Lateral Swing Away Mount - Right	\$358.00

## Spine-Align Positioning Componets

J3 Back LT models come standard with a Small Lumbar, MT,UT and SH models come standard with a small and large lumbar.

### SPINE-ALIGN MINI FIT KIT 1

*Includes: 3 wedges, 2 lateral shims, and 2 posterior/lateral pelvic positioning pads.*

<input type="checkbox"/>	J3SPNK112	K0108	SPINE-ALIGN MINI FIT KIT 1 (12")	\$115.00
<input type="checkbox"/>	J3SPNK114	K0108	SPINE-ALIGN MINI FIT KIT 1 (14")	\$115.00
<input type="checkbox"/>	J3SPNK116	K0108	SPINE-ALIGN MINI FIT KIT 1 (16")	\$115.00
<input type="checkbox"/>	J3SPNK118	K0108	SPINE-ALIGN MINI FIT KIT 1 (18")	\$115.00
<input type="checkbox"/>	J3SPNK120	K0108	SPINE-ALIGN MINI FIT KIT 1 (20")	\$115.00

### LARGE LUMBAR PADS (Back Width) (W" x H" x D")

<input type="checkbox"/>	J3SPNLL12	K0108	LG LUMBAR (12") (3.5" x 7" x 1.6")	\$50.00
<input type="checkbox"/>	J3SPNLL14	K0108	LG LUMBAR (14") (3.5" x 7" x 1.6")	\$50.00
<input type="checkbox"/>	J3SPNLL16	K0108	LG LUMBAR (16") (5.5" x 7" x 1.6")	\$50.00
<input type="checkbox"/>	J3SPNLL18	K0108	LG LUMBAR (18") (7.5" x 7" x 1.6")	\$50.00
<input type="checkbox"/>	J3SPNLL20	K0108	LG LUMBAR (20") (9.5" x 7" x 1.6")	\$50.00

### WEDGES (Back Width) (W" x H" x D")

<input type="checkbox"/>	J3SPNWE12	K0108	WEDGES (12") (3.5" x 6.5" x 1.6")	\$37.00
<input type="checkbox"/>	J3SPNWE14	K0108	WEDGES (14") (3.5" x 6.5" x 1.6")	\$37.00
<input type="checkbox"/>	J3SPNWE16	K0108	WEDGES (16") (5.5" x 6.5" x 1.6")	\$37.00
<input type="checkbox"/>	J3SPNWE18	K0108	WEDGES (18") (7.5" x 6.5" x 1.6")	\$37.00
<input type="checkbox"/>	J3SPNWE20	K0108	WEDGES (20") (9.5" x 6.5" x 1.6")	\$37.00

### SPINE-ALIGN MINI FIT KIT 2

*Includes: 5 wedges, 2 lat shims, 2 post/lat pelvic, 1 lg lumbar, 1 flat shim, and 1 curved*










<input type="checkbox"/>	J3SPNK212	K0108	SPINE-ALIGN MINI FIT KIT 2 (12")	\$159.00
<input type="checkbox"/>	J3SPNK214	K0108	SPINE-ALIGN MINI FIT KIT 2 (14")	\$159.00
<input type="checkbox"/>	J3SPNK216	K0108	SPINE-ALIGN MINI FIT KIT 2 (16")	\$159.00
<input type="checkbox"/>	J3SPNK218	K0108	SPINE-ALIGN MINI FIT KIT 2 (18")	\$159.00
<input type="checkbox"/>	J3SPNK220	K0108	SPINE-ALIGN MINI FIT KIT 2 (20")	\$159.00

### SMALL LUMBAR PADS (Back Width) (W" x H" x D")

<input type="checkbox"/>	J3SPNSL12	K0108	SML LUMBAR (12") (1.8" x 3.2" x 0.8")	\$48.00
<input type="checkbox"/>	J3SPNSL14	K0108	SML LUMBAR (14") (3.5" x 4.3" x 1")	\$48.00
<input type="checkbox"/>	J3SPNSL16	K0108	SML LUMBAR (16") (5.5" x 4.3" x 1")	\$48.00
<input type="checkbox"/>	J3SPNSL18	K0108	SML LUMBAR (18") (7.5" x 4.3" x 1")	\$48.00
<input type="checkbox"/>	J3SPNSL20	K0108	SML LUMBAR (20") (9.5" x 4.3" x 1")	\$48.00

### OTHER INDIVIDUAL SPINE-ALIGN COMPONENTS

<input type="checkbox"/>	J3SPNCR	K0108	Curved Shim	\$28.00
<input type="checkbox"/>	J3SPNBS	K0108	Small Bean Bag	\$31.00
<input type="checkbox"/>	J3SPNFL	K0108	Flat Shim	\$28.00
<input type="checkbox"/>	J3SPNBL	K0108	Large Bean Bag	\$34.00
<input type="checkbox"/>	J3SPNLS	K0108	Lateral Shim	\$54.00
<input type="checkbox"/>	J3SPNLP	K0108	Posterior/Lateral Pelvic	\$54.00


									
	<b>Sm Lumbar</b>	<b>Lg Lumbar</b>	<b>Wedge</b>	<b>Curved Slim</b>	<b>Flat</b>	<b>Lateral Slim</b>	<b>Sm Bean Bag</b>	<b>Lg Bean Bag</b>	<b>Post/Lat Pelvic</b>
<b>W</b>	1.8" - 9.5"	3.5" - 9.5"	3.5" - 9.5"	10"	10"	2.4"	6"	10"	5.3"
<b>H</b>	3.2" - 4.3"	7"	6.5"	4"	4"	4.3"	4.5"	3.2"	3.6"
<b>D</b>	0.8" - 1"	1.6"	1.6"	0.5"	0.3"	0.5"	0.75"	0.75"	0.75"

## Chest Straps, Harnesses and Brackets


Use this section for selecting harnesses. Harnesses are not available on Lower Thoracic Back Shells

### CHEST STRAPS & HARNESSSES (Includes strap guides & mounting hardware)


<input type="checkbox"/>	J3HRNCOS	E0960	Center Opening w/ Buckle (Sml)	\$174.00	<input type="checkbox"/>	J3HRNTDS	E0960	Contoured-Dynamic Pull Down (Sml)	\$174.00
<input type="checkbox"/>	J3HRNCOM	E0960	Center Opening w/ Buckle (Med)	\$174.00	<input type="checkbox"/>	J3HRNTDM	E0960	Contoured-Dynamic Pull Down (Med)	\$174.00
<input type="checkbox"/>	J3HRNCOL	E0960	Center Opening w/ Buckle (Lrg)	\$174.00	<input type="checkbox"/>	J3HRNTDL	E0960	Contoured-Dynamic Pull Down (Lrg)	\$174.00
<input type="checkbox"/>	J3HRNCOX	E0960	Center Opening w/ Buckle (X-Lrg)	\$174.00	<input type="checkbox"/>	J3HRNTDX	E0960	Contoured-Dynamic Pull Down (X-Lrg)	\$174.00
<input type="checkbox"/>	J3HRNCDS	E0960	Classic-Dynamic Pull Down (Sml)	\$174.00	<input type="checkbox"/>	J3HRNTSS	E0960	Contoured-Structured Pull Down (Sml)	\$174.00
<input type="checkbox"/>	J3HRNCDM	E0960	Classic-Dynamic Pull Down (Med)	\$174.00	<input type="checkbox"/>	J3HRNTSM	E0960	Contoured-Structured Pull Down (Med)	\$174.00
<input type="checkbox"/>	J3HRNCDL	E0960	Classic-Dynamic Pull Down (Lrg)	\$174.00	<input type="checkbox"/>	J3HRNTSL	E0960	Contoured-Structured Pull Down (Lrg)	\$174.00
<input type="checkbox"/>	J3HRNCDX	E0960	Classic-Dynamic Pull Down (X-Lrg)	\$174.00	<input type="checkbox"/>	J3HRNTSX	E0960	Contoured-Structured Pull Down (X-Lrg)	\$174.00
<input type="checkbox"/>	J3HRNCSS	E0960	Classic-Structured Pull Down (Sml)	\$174.00	<input type="checkbox"/>	J3HRNPSS	E0960	Padded Anterior Trunk Straps (Sml)	\$162.00
<input type="checkbox"/>	J3HRNCSSM	E0960	Classic-Structured Pull Down (Med)	\$174.00	<input type="checkbox"/>	J3HRNPSSM	E0960	Padded Anterior Trunk Straps (Med)	\$162.00
<input type="checkbox"/>	J3HRNCSSL	E0960	Classic-Structured Pull Down (Lrg)	\$174.00	<input type="checkbox"/>	J3HRNPSSL	E0960	Padded Anterior Trunk Straps (Lrg)	\$162.00
<input type="checkbox"/>	J3HRNCSSX	E0960	Classic-Structured Pull Down (X-Lrg)	\$174.00	<input type="checkbox"/>	J3HRNPSSX	E0960	Padded Anterior Trunk Straps (X-Lrg)	\$162.00



	A	B	C	D	E	F	G	H
S	15"	4.5"	2"	9"	4"	4"	2"	8.5"
M	19"	5.5"	2"	11"	5"	5"	3"	9.5"
L	22.5"	6.5"	2"	13"	6"	6"	3.5"	10.5"
XL	25.5"	6.5"	2.5"	14"	7"	7"	4.5"	11.5"



	A	B	C	D	E	F	G	H
S	15"	4.5"	2"	9"	1"	4"	2"	8.5"
M	18.5"	5.5"	2"	11"	1"	5"	2.5"	9.5"
L	22.5"	6.5"	2"	13"	1.5"	6"	3.5"	10.5"
XL	25.5"	6.5"	2.5"	14"	1.5"	7"	4.5"	11.5"



	A	B	C
S	12.3"	9"	6.5"
M	14.3"	11"	9.8"
L	15.5"	11.5"	9.8"
XL	17.5"	12.8"	10"

### GUIDE BRACKETS & HARDWARE ONLY

<input type="checkbox"/>	J3GPSG	K0108	Shoulder Strap Guides	\$47.00	<input type="checkbox"/>	J3HHP	K0108	Hardware Pack Only	\$23.00
--------------------------	--------	-------	-----------------------	---------	--------------------------	-------	-------	--------------------	---------

### VANITY FLAP

<input type="checkbox"/>	117096	K0108	12" to 14" Wide	\$51.00	<input type="checkbox"/>	117099	K0108	18" to 20" Wide	\$51.00
<input type="checkbox"/>	117097	K0108	14" to 16" Wide	\$51.00	<input type="checkbox"/>	117100	K0108	20" to 22" Wide	\$51.00
<input type="checkbox"/>	117098	K0108	16" to 18" Wide	\$51.00					

## JAY YOUR WAY MODIFICATIONS

To order sizes offered via JAY Your Way add a "M" to the front of the part number and follow the prompts.

For additional modifications visit [www.sunrisemedical.com](http://www.sunrisemedical.com).

### PAINTED BACK SHELL

<input type="checkbox"/>	SBC	Painted Back Shell	NC
Color:		<input type="checkbox"/> Match Quickie or Zippie Frame Color	

*Note: Backshells can be painted any Quickie Or Zippie Frame color at no charge please specify the desired color in the NOTES box below. For a full listing of frame colors available please visit [www.sunrisemedical.com](http://www.sunrisemedical.com)*

### MILITARY SERVICE PATCHES

<input type="checkbox"/>	VAPH	Purple Heart Service Patch	NC
<input type="checkbox"/>	VAM	Marine Service Patch	NC
<input type="checkbox"/>	VAA	Army Service Patch	NC
<input type="checkbox"/>	VAAF	Air Force Service Patch	NC
<input type="checkbox"/>	VAN	Navy Service Patch	NC
<input type="checkbox"/>	VACG	Coast Guard Service Patch	NC

### EMBROIDERY

<input type="checkbox"/>	EMBROIDERY	K0108	Embroidery (32 Character Max)	\$25.00
--------------------------	------------	-------	-------------------------------	---------

Location:  Front

Font:  Block  Script

Embroidery Text:

Embroidery Color:

### COVER OPTIONS - Contact Surface

<input type="checkbox"/>	DBS	Dartex Back Surface	NC
<input type="checkbox"/>	3DXBS	3DX Back Surface	NC
<input type="checkbox"/>	RDBS	Reverse Dartex Back Surface	NC

### COVER OPTIONS - Rear Center Flap

<input type="checkbox"/>	DCC	Digital Camouflage Cover	NC
<input type="checkbox"/>	CBC	Black Cover	NC
<input type="checkbox"/>	CBCI	Black Incontinence Cover	NC

### SPINAL FLUID PAD

<input type="checkbox"/>	SFPS	Small Spinal Fluid Pad (16" Height)	\$57.00
<input type="checkbox"/>	SFPM	Medium Spinal Fluid Pad (18" Height)	\$57.00
<input type="checkbox"/>	SFPT	Tall Spinal Fluid Pad (20" Height)	\$57.00
<input type="checkbox"/>	SFPXT	Extra Tall Spinal Fluid Pad (23" Height)	\$57.00

**Don't see the option you want? NO PROBLEM! Write your request below and we can evaluate it and provide a quote.**

**Notes:**



Sunrise Medical (US) LLC  
2842 Business Park Ave. · Fresno, CA · USA

Customer Service: 800-333-4000 Fax: 800-300-7502 Email: [orders@sunmed.com](mailto:orders@sunmed.com) [www.sunrisemedical.com](http://www.sunrisemedical.com)

Sunrise Medical (US) LLC © 1/2024 MK-129762 REV. G