

GO THAT XTRA MILE

The QUICKIE® Xtender™ electric assist wheelchair accessory extends the range of mobility for manual wheelchair users by adding power assisted wheels to a manual frame. The Xtender increases travel distance and propulsion efficiency and reduces overall fatigue while maintaining the look, control, and functionality of a manual wheelchair.





Quick-Release Wheels and Easy Storage

The Xtender's quick-release wheels and quick-sling battery pack allow you to easily transition between manual and assisted propulsion.



Easy to Use

- No programming required
- Automatic wheel synchronization
- Downhill speed control and rollback assistance
- Bright LED alerts and battery power level indicators
- Easy to remove battery
- Battery includes microprocessor for increased performance
- Nickel or Lithium battery options
- Axle plate adapter is easy to mount and offers more seat-to-floor height adiustment



Benefits of Motorized Wheels

The Xtender's motorized rear wheels increase the force you apply to the handrims. This increases range of mobility, reduces energy expenditure by as much as 45%, and decreases fatique. Research shows that the Xtender electric assist may help reduce or prevent secondary injuries to the wrists, elbows, and shoulders, including Repetitive Strain Injury.

XTENDER SPECS

HCPCS Code F0986

Available on

Quickie 2, 2HP, Q7, 5R, 7R, GT. & GP Series Zippie X'Cape Through Built-4-Me: Ouickie 2 Lite & ORi

Overall Weight (off chair) w/Lithium Battery: 37 lbs w/Nickel Battery: 35 lbs

Weight Capacity 265 lbs

Battery Options Lithium Nickel

Battery Capacity Lithium: 9 hours Nickel: 7 hours

Time to Full Battery Charge Lithium: 3 to 4.5 hours Nickel: 2.5 to 3 hours

Average Range Lithium: 22 miles Nickel: 14 miles

Available Tires Pneumatic

Available Handrims Stainless Steel Plastic Coated



