



ORDER FORM

# PAL Classroom Seat

The perfect balance of postural support and classroom practicality.

Customer			
Account			
Account #		Date	
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #	
Buyer			
Mark For			

Shipping Address			
Name			
Address			
City		State	
ZIP		Phone	
Email / Fax			

The HCPCS CODES herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative HCPCS CODES, including misc. codes to ensure access for their enrollees. The use of HCPCS CODES does not ensure coverage or payment.

## PAL Classroom Seat

Size 1 Weight Capacity: 66 lbs | Size 2 Weight Capacity: 77 lbs | Size 3 Weight Capacity: 88 lbs | Size 4 Weight Capacity: 110 lbs

<input checked="" type="checkbox"/> OPTION	DESCRIPTION	RETAIL
<b>Seat</b>		<b>\$</b>
Includes Lap Strap		
<input type="checkbox"/> 135-1600-03	Size 1 PAL Classroom Seat - Blue	\$1,098
<input type="checkbox"/> 135-1600-01	Size 1 PAL Classroom Seat - Green	\$1,098
<input type="checkbox"/> 135-2600-03	Size 2 PAL Classroom Seat - Blue	\$1,206
<input type="checkbox"/> 135-2600-01	Size 2 PAL Classroom Seat - Green	\$1,206
<input type="checkbox"/> 135-3600-03	Size 3 PAL Classroom Seat - Blue	\$1,318
<input type="checkbox"/> 135-3600-01	Size 3 PAL Classroom Seat - Green	\$1,318
<input type="checkbox"/> 135-4600-03	Size 4 PAL Classroom Seat - Blue	\$1,535
<input type="checkbox"/> 135-4600-01	Size 4 PAL Classroom Seat - Green	\$1,535
<b>Side Pad</b>		<b>\$</b>
Must match seat size.		
<input type="checkbox"/> 135-1635-03	Size 1 PAL Side Pads 1" Thick - Blue (Pair)	\$132
<input type="checkbox"/> 135-1635-01	Size 1 PAL Side Pads 1" Thick - Green (Pair)	\$132
<input type="checkbox"/> 135-2635-03	Size 2 PAL Side Pads 1" Thick - Blue (Pair)	\$132
<input type="checkbox"/> 135-2635-01	Size 2 PAL Side Pads 1" Thick - Green (Pair)	\$132
<input type="checkbox"/> 135-3635-03	Size 3 and 4 PAL Side Pads 1" Thick - Blue (Pair)	\$132
<input type="checkbox"/> 135-3635-01	Size 3 and 4 PAL Side Pads 1" Thick - Green (Pair)	\$132
<input type="checkbox"/> 135-1654-03	Size 1 PAL Side Pads 2" Thick - Blue (Pair)	\$132
<input type="checkbox"/> 135-1654-01	Size 1 PAL Side Pads 2" Thick - Green (Pair)	\$132
<input type="checkbox"/> 135-2654-03	Size 2 PAL Side Pads 2" Thick - Blue (Pair)	\$132
<input type="checkbox"/> 135-2654-01	Size 2 PAL Side Pads 2" Thick - Green (Pair)	\$132
<input type="checkbox"/> 135-3654-03	Size 3 and 4 PAL Side Pads 2" Thick - Blue (Pair)	\$132
<input type="checkbox"/> 135-3654-01	Size 3 and 4 PAL Side Pads 2" Thick - Green (Pair)	\$132
<b>Ramped Seat Cushion</b>		<b>\$</b>
<input type="checkbox"/> 135-1651-03	Size 1 PAL Ramped Seat Cushion - Blue	\$89
<input type="checkbox"/> 135-1651-01	Size 1 PAL Ramped Seat Cushion - Green	\$89
<input type="checkbox"/> 135-2651-03	Size 2 PAL Ramped Seat Cushion - Blue	\$89
<input type="checkbox"/> 135-2651-01	Size 2 PAL Ramped Seat Cushion - Green	\$89
<input type="checkbox"/> 135-3651-03	Size 3 PAL Ramped Seat Cushion - Blue	\$89
<input type="checkbox"/> 135-3651-01	Size 3 PAL Ramped Seat Cushion - Green	\$89
<input type="checkbox"/> 135-4651-03	Size 4 PAL Ramped Seat Cushion - Blue	\$89
<input type="checkbox"/> 135-4651-01	Size 4 PAL Ramped Seat Cushion - Green	\$89
<b>Tray and Tray Pad</b>		<b>\$</b>
Must match seat size.		
<input type="checkbox"/> 135-1620	Size 1 PAL Tray	\$241
<input type="checkbox"/> 135-2620	Size 2 PAL Tray	\$241
<input type="checkbox"/> 135-3620	Size 3 PAL Tray	\$241
<input type="checkbox"/> 135-4620	Size 4 PAL Tray	\$241
<input type="checkbox"/> 135-1656	Size 1 PAL Tray Pad	\$111
<input type="checkbox"/> 135-2656	Size 2 PAL Tray Pad	\$111
<input type="checkbox"/> 135-3656	Size 3 PAL Tray Pad	\$111
<input type="checkbox"/> 135-4656	Size 4 PAL Tray Pad	\$111

## PAL Classroom Seat

✓ OPTION	DESCRIPTION	RETAIL
<b>Pommel</b>		\$
Size 1 Pommel suitable for Size 1 seat, Size 2 Pommel for Size 2, 3 and 4 seats.		
<input type="checkbox"/> 135-1602-03	Size 1 PAL Pommel - Blue	\$174
<input type="checkbox"/> 135-1602-01	Size 1 PAL Pommel - Green	\$174
<input type="checkbox"/> 135-3602-03	Size 2 PAL Pommel - Blue	\$174
<input type="checkbox"/> 135-3602-01	Size 2 PAL Pommel - Green	\$174
<b>Footrest</b>		\$
Must match seat size		
<input type="checkbox"/> 135-1603	Size 1 PAL Footrest	\$132
<input type="checkbox"/> 135-2603	Size 2 PAL Footrest	\$132
<input type="checkbox"/> 135-3603	Size 3 PAL Footrest	\$132
<input type="checkbox"/> 135-4603	Size 4 PAL Footrest	\$132
<b>Rocker / Stabilizer</b>		\$
Not available with PAL Ski with Casters		
<input type="checkbox"/> 135-1604	Size 1 PAL Rocker/Stabilizer (Pair)	\$132
<input type="checkbox"/> 135-2604	Size 2 PAL Rocker/Stabilizer (Pair)	\$132
<input type="checkbox"/> 135-3604	Size 3 PAL Rocker/Stabilizer (Pair)	\$132
<input type="checkbox"/> 135-4604	Size 4 PAL Rocker/Stabilizer (Pair)	\$132
<b>Caster</b>		\$
Not available with Rocker/Stabilizer		
<input type="checkbox"/> 135-1610	Size 1 PAL Ski with Casters (Pair)	\$261
<input type="checkbox"/> 135-2610	Size 2 PAL Ski with Casters (Pair)	\$261
<input type="checkbox"/> 135-3610	Size 3 PAL Ski with Casters (Pair)	\$261
<input type="checkbox"/> 135-4610	Size 4 PAL Ski with Casters (Pair)	\$261



### Technical Specifications - Imperial

Frame Sizing Specification	Size 1	Size 2	Size 3	Size 4
Age	1 - 3 years	2 - 5 years	4 - 8 years	7 - 12 years
Child Height	27.6" - 37.4"	35.4" - 41.3"	39.4" - 51.2"	47" - 63"
Max User Weight	66 lbs	77 lbs	88 lbs	110 lbs
Seat to Footrest	5.5" - 9.4"	6.5" - 11.2"	9.1" - 13.8"	11.2" - 16"
Seat Depth	7.1" - 9.4"	8.7" - 11.4"	9.8" - 13"	12" - 16"
Backrest Height	10.4"	12.4"	15"	17.7"
Backrest Angle	95° - 105°	95° - 105°	95° - 105°	95° - 105°
Seat Width	5.9" - 9.8"	6.7" - 10.6"	9.1" - 13"	11.1" - 15.2"
Seat Weight	9.6 lbs	13.2 lbs	17.9 lbs	22 lbs
Tray Size	20"W x 18"D	21.3"W x 19.7"D	23.6"W x 21.7"D	28.7"W x 23.6"D

### Technical Specifications - Metric

Frame Sizing Specification	Size 1	Size 2	Size 3	Size 4
Age	1 - 3 years	2 - 5 years	4 - 8 years	7 - 12 years
Child Height	700 - 950 mm	900 - 1050 mm	1000 - 1300 mm	1200 - 1600 mm
Max User Weight	30 kg	35 kgs	40 kgs	50 kgs
Seat to Footrest	40 - 240 mm	165 - 285 mm	230 - 350 mm	280 - 400 mm
Seat Depth	80 - 240 mm	220 - 290 mm	250 - 330 mm	300 - 400 mm
Backrest Height	265 mm	315mm	380 mm	450 mm
Backrest Angle	95° - 105°	95° - 105°	95° - 105°	95° - 105°
Seat Width	50 - 250 mm	170 - 270 mm	230 - 330 mm	280 - 380 mm
Seat Weight	4.35kg	6 kgs	8.15 kgs	10 kgs
Tray Size	508mm W x 457mm D	541mm W x 500mm D	599mm W x 551mm D	729mmW x 599mmD

#### Sunrise Medical (US) LLC

28412 Business Park Avenue  
 Fresno, CA 93727 USA  
 Tel.: 1-800-333-4000  
 Fax: 1-800-300-7502  
 E-Mail: webmaster@sunmed.com



If you are visually impaired, this document can be viewed in PDF format at [www.SunriseMedical.com](http://www.SunriseMedical.com)



MK-129852 Rev. F ©01/2024



800.333.4000

[www.SunriseMedical.com](http://www.SunriseMedical.com)