



ORDER FORM

PAL Classroom Seat

EILK03

The perfect balance of postural support and classroom practicality.

Customer			
Account			
Account #		Date	
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #	
Buyer			
Mark For			

Shipping Address			
Name			
Address			
City		State	
ZIP		Phone	
Email / Fax			

The HCPCS CODES herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative HCPCS CODES, including misc. codes to ensure access for their enrollees. The use of HCPCS CODES does not ensure coverage or payment.

PAL Classroom Seat

Size 1 Weight Capacity: 66 lbs | Size 2 Weight Capacity: 77 lbs | Size 3 Weight Capacity: 88 lbs | Size 4 Weight Capacity: 110 lbs

✓ OPTION	DESCRIPTION	RETAIL
Seat		\$
Includes Lap Strap		
<input type="checkbox"/> 151ST92	Size 1 PAL Classroom Seat - Blue	\$1,098
<input type="checkbox"/> 151ST93	Size 1 PAL Classroom Seat - Green	\$1,098
<input type="checkbox"/> 151ST94	Size 2 PAL Classroom Seat - Blue	\$1,206
<input type="checkbox"/> 151ST95	Size 2 PAL Classroom Seat - Green	\$1,206
<input type="checkbox"/> 151ST96	Size 3 PAL Classroom Seat - Blue	\$1,318
<input type="checkbox"/> 151ST97	Size 3 PAL Classroom Seat - Green	\$1,318
<input type="checkbox"/> 151ST98	Size 4 PAL Classroom Seat - Blue	\$1,535
<input type="checkbox"/> 151ST99	Size 4 PAL Classroom Seat - Green	\$1,535
Side Pad		\$
Must match seat size.		
<input type="checkbox"/> 151ZV34	Size 1 PAL Side Pads 1" Thick - Blue (Pair)	\$132
<input type="checkbox"/> 151ZV35	Size 1 PAL Side Pads 1" Thick - Green (Pair)	\$132
<input type="checkbox"/> 151ZV36	Size 2 PAL Side Pads 1" Thick - Blue (Pair)	\$132
<input type="checkbox"/> 151ZV37	Size 2 PAL Side Pads 1" Thick - Green (Pair)	\$132
<input type="checkbox"/> 151ZV38	Size 3 and 4 PAL Side Pads 1" Thick - Blue (Pair)	\$132
<input type="checkbox"/> 151ZV39	Size 3 and 4 PAL Side Pads 1" Thick - Green (Pair)	\$132
<input type="checkbox"/> 151ZV40	Size 1 PAL Side Pads 2" Thick - Blue (Pair)	\$132
<input type="checkbox"/> 151ZV41	Size 1 PAL Side Pads 2" Thick - Green (Pair)	\$132
<input type="checkbox"/> 151ZV42	Size 2 PAL Side Pads 2" Thick - Blue (Pair)	\$132
<input type="checkbox"/> 151ZV43	Size 2 PAL Side Pads 2" Thick - Green (Pair)	\$132
<input type="checkbox"/> 151ZV44	Size 3 and 4 PAL Side Pads 2" Thick - Blue (Pair)	\$132
<input type="checkbox"/> 151ZV45	Size 3 and 4 PAL Side Pads 2" Thick - Green (Pair)	\$132
Ramped Seat Cushion		\$
<input type="checkbox"/> 151CU115	Size 1 PAL Ramped Seat Cushion - Blue	\$89
<input type="checkbox"/> 151CU116	Size 1 PAL Ramped Seat Cushion - Green	\$89
<input type="checkbox"/> 151CU117	Size 2 PAL Ramped Seat Cushion - Blue	\$89
<input type="checkbox"/> 151CU118	Size 2 PAL Ramped Seat Cushion - Green	\$89
<input type="checkbox"/> 151CU119	Size 3 PAL Ramped Seat Cushion - Blue	\$89
<input type="checkbox"/> 151CU120	Size 3 PAL Ramped Seat Cushion - Green	\$89
<input type="checkbox"/> 151CU121	Size 4 PAL Ramped Seat Cushion - Blue	\$89
<input type="checkbox"/> 151CU122	Size 4 PAL Ramped Seat Cushion - Green	\$89
Tray and Tray Pad		\$
Must match seat size.		
<input type="checkbox"/> 151A222	Size 1 PAL Tray	\$241
<input type="checkbox"/> 151A223	Size 2 PAL Tray	\$241
<input type="checkbox"/> 151A224	Size 3 PAL Tray	\$241
<input type="checkbox"/> 151A225	Size 4 PAL Tray	\$241
<input type="checkbox"/> 151A226	Size 1 PAL Tray Pad	\$111
<input type="checkbox"/> 151A227	Size 2 PAL Tray Pad	\$111
<input type="checkbox"/> 151A228	Size 3 PAL Tray Pad	\$111
<input type="checkbox"/> 151A229	Size 4 PAL Tray Pad	\$111

PAL Classroom Seat

✓ OPTION	DESCRIPTION	RETAIL
Pommel		\$
Size 1 Pommel suitable for Size 1 seat, Size 2 Pommel for Size 2, 3 and 4 seats.		
<input type="checkbox"/> 1518A114	Size 1 PAL Pommel - Blue	\$174
<input type="checkbox"/> 1518A115	Size 1 PAL Pommel - Green	\$174
<input type="checkbox"/> 1518A116	Size 2 PAL Pommel - Blue	\$174
<input type="checkbox"/> 1518A117	Size 2 PAL Pommel - Green	\$174
Footrest		\$
Must match seat size		
<input type="checkbox"/> 151F378	Size 1 PAL Footrest	\$132
<input type="checkbox"/> 151F379	Size 2 PAL Footrest	\$132
<input type="checkbox"/> 151F380	Size 3 PAL Footrest	\$132
<input type="checkbox"/> 151F381	Size 4 PAL Footrest	\$132
Rocker / Stabilizer		\$
Not available with PAL Ski with Casters		
<input type="checkbox"/> 151RX19	Size 1 PAL Rocker/Stabilizer (Pair)	\$132
<input type="checkbox"/> 151RX20	Size 2 PAL Rocker/Stabilizer (Pair)	\$132
<input type="checkbox"/> 151RX21	Size 3 PAL Rocker/Stabilizer (Pair)	\$132
<input type="checkbox"/> 151RX22	Size 4 PAL Rocker/Stabilizer (Pair)	\$132
Caster		\$
Not available with Rocker/Stabilizer		
<input type="checkbox"/> 151SY5	Size 1 PAL Ski with Casters (Pair)	\$261
<input type="checkbox"/> 151SY6	Size 2 PAL Ski with Casters (Pair)	\$261
<input type="checkbox"/> 151SY7	Size 3 PAL Ski with Casters (Pair)	\$261
<input type="checkbox"/> 151SY8	Size 4 PAL Ski with Casters (Pair)	\$261



Technical Specifications - Imperial

Frame Sizing Specification	Size 1	Size 2	Size 3	Size 4
Age	1 - 3 years	2 - 5 years	4 - 8 years	7 - 12 years
Child Height	27.6" - 37.4"	35.4" - 41.3"	39.4" - 51.2"	47" - 63"
Max User Weight	66 lbs	77 lbs	88 lbs	110 lbs
Seat to Footrest	5.5" - 9.4"	6.5" - 11.2"	9.1" - 13.8"	11.2" - 16"
Seat Depth	7.1" - 9.4"	8.7" - 11.4"	9.8" - 13"	12" - 16"
Backrest Height	10.4"	12.4"	15"	17.7"
Backrest Angle	95° - 105°	95° - 105°	95° - 105°	95° - 105°
Seat Width	5.9" - 9.8"	6.7" - 10.6"	9.1" - 13"	11.1" - 15.2"
Seat Weight	9.6 lbs	13.2 lbs	17.9 lbs	22 lbs
Tray Size	20"W x 18"D	21.3"W x 19.7"D	23.6"W x 21.7"D	28.7"W x 23.6"D

Technical Specifications - Metric

Frame Sizing Specification	Size 1	Size 2	Size 3	Size 4
Age	1 - 3 years	2 - 5 years	4 - 8 years	7 - 12 years
Child Height	700 - 950 mm	900 - 1050 mm	1000 - 1300 mm	1200 - 1600 mm
Max User Weight	30 kg	35 kgs	40 kgs	50 kgs
Seat to Footrest	40 - 240 mm	165 - 285 mm	230 - 350 mm	280 - 400 mm
Seat Depth	80 - 240 mm	220 - 290 mm	250 - 330 mm	300 - 400 mm
Backrest Height	265 mm	315mm	380 mm	450 mm
Backrest Angle	95° - 105°	95° - 105°	95° - 105°	95° - 105°
Seat Width	50 - 250 mm	170 - 270 mm	230 - 330 mm	280 - 380 mm
Seat Weight	4.35kg	6 kgs	8.15 kgs	10 kgs
Tray Size	508mm W x 457mm D	541mm W x 500mm D	599mm W x 551mm D	729mmW x 599mmD

Sunrise Medical (US) LLC
 2842 N. Business Park Avenue
 Fresno, CA 93727 USA



If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.com



MK-129852 Rev. G © 06/2024

