



A truly transportable tilt wheelchair, specifically designed for families on the move.

EIZ22



Customer				Ship To			
Company				Name			
Account #		Date		Address			
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #		City		State	
Buyer				ZIP		Phone	
Mark For				Email / Fax			

The HCPCS CODES herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative HCPCS CODES, including misc. codes to ensure access for their enrollees. The use of HCPCS CODES does not ensure coverage or payment.

### INCLUDED COMPONENTS

- ✓ Silver & Black Folding Frame
- ✓ 10°, 20°, 30° Seat Tilt Settings
- ✓ Quick Adjust Recline Backrest
- ✓ Foot Operated Wheel Lock
- ✓ 12" Rear Wheels & 6" Casters (Solid)
- ✓ Quick Release Rear Wheels
- ✓ Adjustable Stroller Handle
- ✓ Seat Sling with Depth Adjustability

✓ OPTION	HCPCS	DESCRIPTION	RETAIL
----------	-------	-------------	--------

### BASE

#### Zippie Sphynx

\$

For E1161, please select either 20° or 30° Seat Tilt.

<input type="checkbox"/> 149WM103	E1232   E1161 (20° or 30°)	Zippie Sphynx - Includes upholstery, Transit optional	\$ 2,816.00
<input type="checkbox"/> 149WM104	E1234   E1161 (20° or 30°)	Zippie Sphynx without Upholstery - Transit optional	\$ 2,550.00
<input type="checkbox"/> 149WM106	E1232   E1161 (20° or 30°)	Zippie Sphynx with Transit Included - With upholstery	\$ 3,151.00
<input type="checkbox"/> 149WM111	E1232   E1161 (20° or 30°)	Zippie Sphynx Package - Includes Transit, Upholstery, Canopy, Underseat Basket, Rain Cover	\$ 3,844.00

### FRAME

#### Transit System

\$

Includes: Transit Brackets, and Headrest Extension when ordering Sphynx upholstery.

<input type="checkbox"/> 149TO44	K0108	Transit System	\$ 367.00
<input type="checkbox"/> 149TO40		OMIT Transit - Only available on Transit Optional Sphynx Base	NC

#### Seat Tilt

\$

The adjustable static seat tilt setting can be changed between 10°, 20°, or 30°.

<input type="checkbox"/> 149TA15		10° Seat Tilt	NC
<input type="checkbox"/> 149TA16		20° Seat Tilt	NC
<input type="checkbox"/> 149TA17		30° Seat Tilt	NC

✓ OPTION	DESCRIPTION	RETAIL	✓ OPTION	DESCRIPTION	RETAIL
----------	-------------	--------	----------	-------------	--------

**FRAME**

**Seat Width \$ \$**

Selection Required. Seat Width will have an adjustable Seat Depth range.

<input type="checkbox"/> 149W1	10" Seat Width - Adj Seat Depth of 10" to 15"	NC	<input type="checkbox"/> 149W7	16" Seat Width - Adj Seat Depth of 13" to 18"	NC
<input type="checkbox"/> 149W3	12" Seat Width - Adj Seat Depth of 10" to 15"	NC	<input type="checkbox"/> 149W9	18" Seat Width - Adj Seat Depth of 14" to 19"	NC
<input type="checkbox"/> 149W5	14" Seat Width - Adj Seat Depth of 11" to 16"	NC			

**Seat Depth \$ \$**

Selection Required. Seat Depth is adjustable within the range stated above for Seat Width.

<input type="checkbox"/> 149SD3	10" Seat Depth - Only ava on 10"   12" SW	NC	<input type="checkbox"/> 149SD8	15" Seat Depth	NC
<input type="checkbox"/> 149SD4	11" Seat Depth - Only ava on 10"   12"   14" SW	NC	<input type="checkbox"/> 149SD9	16" Seat Depth - NA on 10"   12" SW	NC
<input type="checkbox"/> 149SD5	12" Seat Depth - Only ava on 10"   12"   14" SW	NC	<input type="checkbox"/> 149SD10	17" Seat Depth - Only ava on 16"   18" SW	NC
<input type="checkbox"/> 149SD6	13" Seat Depth - NA on 18" SW	NC	<input type="checkbox"/> 149SD11	18" Seat Depth - Only ava on 16"   18" SW	NC
<input type="checkbox"/> 149SD7	14" Seat Depth	NC	<input type="checkbox"/> 149SD12	19" Seat Depth - Only ava on 18" SW	NC

**Front Seat Height Guide**

Tilt Angle	Seat Depth									
	10"	11"	12"	13"	14"	15"	16"	17"	18"	19"
10°	18"	18"	18"	18.5"	18.5"	18.5"	19"	19"	19"	19.5"
20°	17.5"	17.5"	18"	18.5"	18.5"	19"	19.5"	19.5"	20"	20.5"
30°	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	21"

**FOOTREST**

**Seat-To-Footplate Range (Lower Leg Length) \$**

Selection Required. Reference table below for achievable Seat-To-Footplate distances relative to the Seat Width and Depth.

<input type="checkbox"/> 149L064	Extra Short Seat-To-Footplate Range	NC
<input type="checkbox"/> 149L061	Short Seat-To-Footplate Range	NC
<input type="checkbox"/> 149L062	Medium Seat-To-Footplate Range	NC
<input type="checkbox"/> 149L063	Long Seat-To-Footplate Range	NC

**Achievable Seat-To-Footplate Distances**

Seat Width	Seat-to-Footplate Range Option	Seat Depth									
		10"	11"	12"	13"	14"	15"	16"	17"	18"	19"
10"	Extra Short	6" - 12"	6" - 11"	6" - 10.5"	6" - 10"	5" - 10"	5" - 10"	NA	NA	NA	NA
	Short	12" - 15.5"	11" - 15"	10.5" - 15"	10" - 14.5"	10" - 14.5"	10" - 14.5"	NA	NA	NA	NA
12"	Medium	16" - 17.5"	15" - 17.5"	14.5" - 17"	14" - 17.5"	14" - 17.5"	13.5" - 16.5"	NA	NA	NA	NA
	Long	18" - 19"	16.5" - 19"	16" - 19"	16" - 19"	16" - 19"	15.5" - 19"	NA	NA	NA	NA
14"	Extra Short	NA	6" - 11"	6" - 10.5"	6" - 10"	5" - 10"	5" - 10"	5" - 10"	NA	NA	NA
	Short	NA	11" - 15"	10.5" - 15"	10" - 14.5"	10" - 14.5"	10" - 14.5"	10" - 14.5"	NA	NA	NA
	Medium	NA	15" - 17.5"	14.5" - 17"	14" - 17.5"	14" - 17.5"	13.5" - 16.5"	13" - 16"	NA	NA	NA
	Long	NA	16.5" - 19"	16" - 19"	16" - 19"	16" - 19"	15.5" - 19"	15" - 19"	NA	NA	NA
16"	Extra Short	NA	NA	NA	6" - 10"	5" - 10"	5" - 10"	5" - 10"	6" - 10"	7" - 10"	NA
	Short	NA	NA	NA	10" - 14.5"	10" - 14.5"	10" - 14.5"	10" - 14.5"	10" - 14"	10" - 14"	NA
	Medium	NA	NA	NA	14" - 17.5"	14" - 17.5"	13.5" - 16.5"	13" - 16"	13" - 15.5"	13" - 15.5"	NA
	Long	NA	NA	NA	16" - 19"	16" - 19"	15.5" - 19"	15" - 19"	15" - 19"	14.5" - 19"	NA
18"	Extra Short	NA	NA	NA	NA	10" - 14.5"	10" - 14"	10" - 14"	10" - 14"	10" - 13.5"	10" - 11.5"
	Short	NA	NA	NA	NA	12" - 15"	11.5" - 15"	11.5" - 14.5"	14.5"	11.5" - 14"	11.5" - 13.5"
	Medium	NA	NA	NA	NA	14" - 17"	13.5" - 16.5"	13" - 16"	13" - 15.5"	13" - 15.5"	13" - 15.5"
	Long	NA	NA	NA	NA	15.5" - 19"	15.5" - 19"	15" - 19"	15" - 19"	14.5" - 19"	14.5" - 19"

NOTE: All Seat-To-Footplate distances based on 10° Tilt Setting (Seat-To-Floor Angle).

✓ OPTION	DESCRIPTION	RETAIL	✓ OPTION	DESCRIPTION	RETAIL
----------	-------------	--------	----------	-------------	--------

### UPHOLSTERY

#### Upholstery for E1232 Sphynx Base \$

Selection Required for Upholstery-Included Sphynx Base (E1232). Includes: Seat and back cushions with tension adjustable straps and seat sling.

<input type="checkbox"/> 149UC37	Sphynx Upholstery - Black	STD	<input type="checkbox"/> 149UC51	Sphynx Upholstery - Green	NC
<input type="checkbox"/> 149UC49	Sphynx Upholstery - Red	NC	<input type="checkbox"/> 149UC52	Sphynx Upholstery - Pink	NC
<input type="checkbox"/> 149UC50	Sphynx Upholstery - Blue	NC	<input type="checkbox"/> 149UC53	OMIT Upholstery	NC

#### Upholstery for E1234 Sphynx Base \$

Selection Required for Upholstery-Included Sphynx Base (E1234). Includes: Seat and back cushions with tension adjustable straps and seat sling.

<input type="checkbox"/> 149UC54	Sphynx Upholstery - Black	\$ 266.00	<input type="checkbox"/> 149UC57	Sphynx Upholstery - Green	\$ 266.00
<input type="checkbox"/> 149UC55	Sphynx Upholstery - Red	\$ 266.00	<input type="checkbox"/> 149UC58	Sphynx Upholstery - Pink	\$ 266.00
<input type="checkbox"/> 149UC56	Sphynx Upholstery - Blue	\$ 266.00	<input type="checkbox"/> 149UC53	OMIT Upholstery	NC

### SEATING

#### JAY Options \$

Required to attach JAY Order Forms for bundled products.

<input type="checkbox"/> 149JC20	JAY Cushion	JAY Part #:		See OF
<input type="checkbox"/> 149JB02	JAY Backrest	JAY Part #:		See OF

### POSITIONING

#### Towel Bar \$

Default Towel Bar is 2" for more compact folded package. Select 4" Towel Bar to accommodate posterior contoured back supports.

<input type="checkbox"/> 149QT02		2" Towel Bar	STD
<input type="checkbox"/> 149QT01		4" Towel Bar	NC

#### Head Supports \$

<input type="checkbox"/> 149HE75	E0966	Headrest Extension - Required for Transit Option & Sphynx Upholstery	\$ 160.00
<input type="checkbox"/> 149ZH25		Padded Headwings - May interfere with Full Torso Support Vest	\$ 266.00

#### Trunk & Arm Positioning \$

<input type="checkbox"/> 149K17	E0966	H-Harness - Requires 3-Point Pelvic Belt & Headrest Extension. NA with Full Torso Support Vest.	\$ 193.00	
<input type="checkbox"/> 149JY06	E0966	JAY Center Opening Dynamic - Requires Headrest Extension. Select harness size below.	\$ 179.00	
Select harness size:	<input type="checkbox"/> 149JZ01 Small	<input type="checkbox"/> 149JZ02 Medium	<input type="checkbox"/> 149JZ03 Large	<input type="checkbox"/> 149JZ04 X-Large
<input type="checkbox"/> 149K19		Full Torso Support Vest - Requires Headrest Extension. NA with H-Harness or Lateral Supports.	\$ 239.00	
<input type="checkbox"/> 149K18		Adjustable Lateral Support & Scoli Strap - Requires Armrests. NA with Full Torso Support Vest.	\$ 193.00	
<input type="checkbox"/> 149AA1	E0973	Flip-Up Armrests (Angle Adjustable)	\$ 293.00	
<input type="checkbox"/> 149A206	E0950	Upper Extremity Support Tray - Grey (Requires Armrests)	\$ 213.00	
<input type="checkbox"/> 149A207	E0950	Upper Extremity Support Tray - Clear (Requires Armrests)	\$ 213.00	

✓ OPTION	HCPCS	DESCRIPTION	RETAIL
----------	-------	-------------	--------

### POSITIONING

#### Pelvic Positioning \$

<input type="checkbox"/> 149ZZ29	E0978	3-Point Pelvic Belt - Requires H-Harness. NA with Medial Thigh Support.	\$ 133.00
<input type="checkbox"/> 149ZZ11	E0978	2-Point Positioning Belt with 4.5" Pad	\$ 106.00
<input type="checkbox"/> 149ZZ24	E0978	V-Contour 4-Point Belt - Small Pad	\$ 186.00
<input type="checkbox"/> 149ZZ23	E0978	V-Contour 4-Point Belt - Medium Pad	\$ 186.00
<input type="checkbox"/> 149KP05	E0956	Hip Guides - Inset Mount (Pair) - Reduces seat width by 1" per side	\$ 266.00
<input type="checkbox"/> 149KP06	E0956	Hip Guides - In-Line Mount (Pair)	\$ 266.00
<input type="checkbox"/> 149CU30	K0669	Cushion 1" Foam Riser	\$ 128.00
<input type="checkbox"/> 149CU12	E2294	Wedge Cushion Insert	\$ 282.00
<input type="checkbox"/> 149KH11	E0957   E1028	Medial Thigh Support - Includes hardware. NA with H-Harness & 3-Point Pelvic Belt.	\$ 382.00

#### Foot & Leg Positioning \$

<input type="checkbox"/> 149F26	K0040   E0990	Angle Adjustable Aluminum Footplates (Pair) - Required	\$ 266.00
<input type="checkbox"/> 149HL5	K0038	Leg Strap	\$ 46.00
<input type="checkbox"/> 149A54		Extra Leg Strap	\$ 41.00
<input type="checkbox"/> 149F304	E0951   E0952	Foot & Ankle Positioners - Small	\$ 176.00
<input type="checkbox"/> 149F305	E0951   E0952	Foot & Ankle Positioners - Medium	\$ 176.00
<input type="checkbox"/> 149F306	E0951   E0952	Foot & Ankle Positioners - Large	\$ 176.00
<input type="checkbox"/> 149F307	E0951   E0952	Foot & Ankle Positioners - X-Large	\$ 176.00

### ACCESSORIES \$

<input type="checkbox"/> 149LF16	E1029	Battery Utility Tray - Limit 50lbs, includes 1 utility hook	\$ 482.00
<input type="checkbox"/> 149LF17	E1029	Ventilator Tray - Vertical mount	\$ 442.00
<input type="checkbox"/> 149LF6		Utility Hooks & Zippie Bag	\$ 115.00
<input type="checkbox"/> 149HE76	K0108	Headrest Cover - Canopy - Includes Hardware	\$ 319.00
<input type="checkbox"/> 149A468		Underseat Storage Basket - NA with Battery Tray	\$ 165.00
<input type="checkbox"/> 149LF5		Cup Holder	\$ 80.00
<input type="checkbox"/> 149A466		Curb Tipper	\$ 106.00
<input type="checkbox"/> 149A467		Rain Cover - Requires Headrest Cover - Canopy	\$ 208.00
<input type="checkbox"/> 149A1		Zippie Backpack	\$ 86.00
<input type="checkbox"/> 149A469		Tool Kit	\$ 46.00

### MEASUREMENTS - NOT REQUIRED

Name			1. Top of Head to Seat		7. Width of Head	
Age			2. Top of Shoulder to Seat		8. Shoulder Width	
Diagnosis			3. Undarm to Seat		9. Chest Width	
Height			4. Trunk Depth		10. Hip Width	
Weight			5. Actual Sitting Depth		11. Left Lower Leg Length	
Please use inches (") for unit of measure			6. Foot Length		12. Right Lower Leg Length	