



ZIPPIE® Sphynx™

May 2022

EIZ22

Mark For: _____ Date: _____ Dealer Acct #: _____ Dealer: _____ Dealer Contact: _____ Dealer Address: _____ Dealer City: _____ ST: _____ ZIP: _____ Dealer Phone: _____ Fax: _____ Confirmation Email: _____ Confirm Via: <input type="checkbox"/> Fax <input type="checkbox"/> Email	Submitting for: <input type="checkbox"/> Quote <input type="checkbox"/> Order PO#: _____ <div style="text-align: center;">ADDITIONAL SHIPPING INFORMATION</div> Ship To: _____ Attention: _____ Address: _____ Address: _____ Ship To City: _____ ST: _____ ZIP: _____ Ship To Phone: _____ Fax: _____
---	---

*The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.*

ZIPPIE Sphynx

Each ZIPPIE Sphynx unit includes the following components:

- Silver & Black Folding Frame
- 10°, 20°, 30° Seat Tilt Settings
- Quick Adjust Recline Backrest
- Foot Operated Wheel Lock
- 12" Rear Wheels & 6" Casters (Solid)
- Quick Release Rear Wheels
- Adjustable Stroller Handle
- Seat Sling with Depth Adjustability

Sphynx Model

Selection Required. For sizes greater than 15" use E1161.

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	149WM103	ZIPPIE Sphynx (Includes upholstery, Transit optional)	E1232, E1161	\$2,643
<input type="checkbox"/>	149WM104	ZIPPIE Sphynx without Upholstery (Transit optional)	E1234, E1161	\$2,393
<input type="checkbox"/>	149WM106	ZIPPIE Sphynx with Transit Included (With upholstery)	E1232, E1161	\$2,958
<input type="checkbox"/>	149WM111	ZIPPIE Sphynx Package (Includes Transit, Upholstery, Canopy, Underseat Basket, Rain Cover)	E1232, E1161	\$3,608

Transit System

Transit System includes Transit Brackets, and Headrest Extension when ordering Sphynx upholstery

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	149TO44	Transit System	K0108	\$315
<input type="checkbox"/>	149TO40	OMIT transit (only available on transit optional Sphynx models)		N/C

Seat Tilt

Select Seat Tilt setting, default is 10°. The adjustable static seat tilt setting can be changed between 10°, 20°, or 30°

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	149TA15	10° Seat Tilt		STD
<input type="checkbox"/>	149TA16	20° Seat Tilt		N/C
<input type="checkbox"/>	149TA17	30° Seat Tilt		N/C

Seat Width

Selection Required. Seat width will have an adjustable Seat Depth range

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	149W1	10" Seat Width (Adjustable Seat Depth of 10" to 15")		N/C
<input type="checkbox"/>	149W3	12" Seat Width (Adjustable Seat Depth of 10" to 15")		N/C
<input type="checkbox"/>	149W5	14" Seat Width (Adjustable Seat Depth of 11" to 16")		N/C
<input type="checkbox"/>	149W7	16" Seat Width (Adjustable Seat Depth of 13" to 18")		N/C
<input type="checkbox"/>	149W9	18" Seat Width (Adjustable Seat Depth of 14" to 19")		N/C

Seat Depth

Selection Required. Seat depth is adjustable within the range stated above for the Seat Width

Item Code	Description	HCPCS	Retail
<input type="checkbox"/> 149SD3	10" Seat Depth (Only available on 10" and 12" Seat Widths)		N/C
<input type="checkbox"/> 149SD4	11" Seat Depth (Only available on 10", 12" and 14" Seat Widths)		N/C
<input type="checkbox"/> 149SD5	12" Seat Depth (Only available on 10", 12" and 14" Seat Widths)		N/C
<input type="checkbox"/> 149SD6	13" Seat Depth (N/A on 18" Seat Width)		N/C
<input type="checkbox"/> 149SD7	14" Seat Depth		N/C
<input type="checkbox"/> 149SD8	15" Seat Depth		N/C
<input type="checkbox"/> 149SD9	16" Seat Depth (N/A on 10" or 12" Seat Widths)		N/C
<input type="checkbox"/> 149SD10	17" Seat Depth (Only available on 16" and 18" Seat Widths)		N/C
<input type="checkbox"/> 149SD11	18" Seat Depth (Only available on 16" and 18" Seat Widths)		N/C
<input type="checkbox"/> 149SD12	19" Seat Depth (Only available on 18" Seat Width)		N/C

Tilt Angle	Front Seat Height									
	Seat Depth									
	10"	11"	12"	13"	14"	15"	16"	17"	18"	19"
10°	18"	18"	18"	18.5"	18.5"	18.5"	19"	19"	19"	19.5"
20°	17.5"	17.5"	18"	18.5"	18.5"	19"	19.5"	19.5"	20"	20.5"
30°	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	21"

Seat-To-Footplate Range (Lower Leg Length)

Selection Required. Reference table below for achievable Seat-To-Footplate distances relative to the Seat Width & Depth.

Item Code	Description	HCPCS	Retail
<input type="checkbox"/> 149LO64	Extra Short Seat-To-Footplate Range		N/C
<input type="checkbox"/> 149LO61	Short Seat-To-Footplate Range		N/C
<input type="checkbox"/> 149LO62	Medium Seat-To-Footplate Range		N/C
<input type="checkbox"/> 149LO63	Long Seat-To-Footplate Range		N/C

Seat Width	Seat-to-Footplate Range Option	Achievable Seat-to-Footplate Distances*									
		Seat Depth									
		10"	11"	12"	13"	14"	15"	16"	17"	18"	19"
10" 12"	Extra Short	6"-12"	6"-11"	6"-10.5"	6"-10"	5"-10"	5"-10"	N/A	N/A	N/A	N/A
	Short	12"-15.5"	11"-15"	10.5"-15"	10"-14.5"	10"-14.5"	10"-14.5"	N/A	N/A	N/A	N/A
	Medium	16"-17.5"	15"-17.5"	14.5"-17"	14"-17.5"	14"-17.5"	13.5"-16.5"	N/A	N/A	N/A	N/A
	Long	18"-19"	16.5"-19"	16"-19"	16"-19"	16"-19"	15.5"-19"	N/A	N/A	N/A	N/A
14"	Extra Short	N/A	6"-11"	6"-10.5"	6"-10"	5"-10"	5"-10"	5"-10"	N/A	N/A	N/A
	Short	N/A	11"-15"	10.5"-15"	10"-14.5"	10"-14.5"	10"-14.5"	10"-14.5"	N/A	N/A	N/A
	Medium	N/A	15"-17.5"	14.5"-17"	14"-17.5"	14"-17.5"	13.5"-16.5"	13"-16"	N/A	N/A	N/A
	Long	N/A	16.5"-19"	16"-19"	16"-19"	16"-19"	15.5"-19"	15"-19"	N/A	N/A	N/A
16"	Extra Short	N/A	N/A	N/A	6"-10"	5"-10"	5"-10"	5"-10"	6"-10"	7"-10"	N/A
	Short	N/A	N/A	N/A	10"-14.5"	10"-14.5"	10"-14.5"	10"-14.5"	10"-14"	10"-14"	N/A
	Medium	N/A	N/A	N/A	14"-17.5"	14"-17.5"	13.5"-16.5"	13"-16"	13"-15.5"	13"-15.5"	N/A
	Long	N/A	N/A	N/A	16"-19"	16"-19"	15.5"-19"	15"-19"	15"-19"	14.5"-19"	N/A
18"	Extra Short	N/A	N/A	N/A	N/A	10"-14.5"	10"-14"	10"-14"	10"-14"	10"-13.5"	10"-11.5"
	Short	N/A	N/A	N/A	N/A	12"-15"	11.5"-15"	11.5"-14.5"	11.5"-14.5"	11.5"-14"	11.5"-13.5"
	Medium	N/A	N/A	N/A	N/A	14"-17"	13.5"-16.5"	13"-16"	13"-15.5"	13"-15.5"	13"-15.5"
	Long	N/A	N/A	N/A	N/A	15.5"-19"	15.5"-19"	15"-19"	15"-19"	14.5"-19"	14.5"-19"

**NOTE: All seat-to-footplate distances based on 10 deg tilt setting (seat to floor angle).*

Sphynx Upholstery for E1232 Sphynx Model*Selection Required for Upholstery-Included Sphynx Model (E1232). Includes seat and back cushions with tension adjustable straps and seat sling.*

<input type="checkbox"/>	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	149UC37	Sphynx Upholstery - Black		STD
<input type="checkbox"/>	149UC49	Sphynx Upholstery - Red		N/C
<input type="checkbox"/>	149UC50	Sphynx Upholstery - Blue		N/C
<input type="checkbox"/>	149UC51	Sphynx Upholstery - Green		N/C
<input type="checkbox"/>	149UC52	Sphynx Upholstery - Pink		N/C
<input type="checkbox"/>	149UC53	OMIT upholstery		N/C

Sphynx Upholstery for E1234 Sphynx Model*Selection Required for Upholstery-Optional Sphynx Model (E1234). Includes seat and back cushions with tension adjustable straps and seat sling.*

<input type="checkbox"/>	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	149UC54	Sphynx Upholstery - Black		\$250
<input type="checkbox"/>	149UC55	Sphynx Upholstery - Red		\$250
<input type="checkbox"/>	149UC56	Sphynx Upholstery - Blue		\$250
<input type="checkbox"/>	149UC57	Sphynx Upholstery - Green		\$250
<input type="checkbox"/>	149UC58	Sphynx Upholstery - Pink		\$250
<input type="checkbox"/>	149UC53	OMIT upholstery		N/C

JAY Options*Required to attach JAY order forms for bundled products*

<input type="checkbox"/>	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	149JC20	JAY Cushion JAY Part #: <input type="text"/>		See OF
<input type="checkbox"/>	149JB02	JAY Backrest JAY Part #: <input type="text"/>		See OF

Towel Bar*Select the towel bar size, default is 2" for more compact folded package. Select 4" Towel Bar to accommodate posterior contoured back supports.*

<input type="checkbox"/>	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	149QT02	2" Towel Bar		STD
<input type="checkbox"/>	149QT01	4" Towel Bar		N/C

Head Supports

<input type="checkbox"/>	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	149HE75	Headrest Extension <i>(Required with transit option & Sphynx Upholstery)</i>	E0966	\$150
<input type="checkbox"/>	149ZH25	Padded Headwings <i>(May interfere with full torso vest)</i>		\$250

Trunk & Arm Positioning


<input type="checkbox"/>	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	149KI7	H-Harness <i>(Requires 3-Point Pelvic Belt. Requires Headrest Ext. w/ std. uph. N/A with Full Torso Vest)</i>	E0960	\$182
<input type="checkbox"/>	149JY06	JAY Center Opening Dynamic <i>(Requires Headrest Ext. & Select harness size below)</i>	E0960	\$168
<input type="checkbox"/>	149JZ01	Small		N/C
<input type="checkbox"/>	149JZ02	Medium		N/C
<input type="checkbox"/>	149JZ03	Large		N/C
<input type="checkbox"/>	149JZ04	X-Large		N/C
<input type="checkbox"/>	149KI9	Full Torso Support Vest <i>(Requires Headrest Ext., N/A with H-Harness or Lateral Supports)</i>		\$225
<input type="checkbox"/>	149KI8	Adjustable Lateral Support & Scoli Strap <i>(Requires Armrests, N/A with Full Torso Vest)</i>		\$182
<input type="checkbox"/>	149AA1	Flip-Up Armrests (Angle Adjustable)	E0973	\$275
<input type="checkbox"/>	149A206	Upper Extremity Support Tray - Grey <i>(Requires armrests)</i>	E0950	\$200
<input type="checkbox"/>	149A207	Upper Extremity Support Tray - Clear <i>(Requires armrests)</i>	E0950	\$200

Pelvic Positioning				
	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	149ZZ29	3-Point Pelvic Belt <i>(Requires H-Harness, N/A with Medial Thigh Support)</i>	E0978	\$125
<input type="checkbox"/>	149ZZ11	2-Point positioning belt w/ 4.5" pad	E0978	\$100
<input type="checkbox"/>	149ZZ24	V-Contour 4-Point Belt- Small Pad	E0978	\$175
<input type="checkbox"/>	149ZZ23	V-Contour 4-Point Belt- Med Pad	E0978	\$175
<input type="checkbox"/>	149KP05	Hip Guides - Inset Mount (Pair) <i>(Reduces seat width by 1" per side)</i>	E0956	\$250
<input type="checkbox"/>	149KP06	Hip Guides - In-Line Mount (Pair)	E0956	\$250
<input type="checkbox"/>	149CU30	Cushion 1" Foam Riser	K0669	\$120
<input type="checkbox"/>	149CU12	Wedge Cushion Insert	E2294	\$265
<input type="checkbox"/>	149KH11	Medial Thigh Support <i>(Includes hardware, N/A with H-Harness & 3-Point Pelvic Belt)</i>	E0957, E1028	\$358

Foot & Leg Positioning				
	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	149F26	Angle Adjustable Aluminum Footplates (Pair) Required	K0040, E0990	\$250
<input type="checkbox"/>	149HL5	Leg Strap	K0038	\$43
<input type="checkbox"/>	149A54	Extra Leg Strap		\$38
<input type="checkbox"/>	149F304	Foot & Ankle Positioners - Small	E0951, E0952	\$165
<input type="checkbox"/>	149F305	Foot & Ankle Positioners - Medium	E0951, E0952	\$165
<input type="checkbox"/>	149F306	Foot & Ankle Positioners - Large	E0951, E0952	\$165
<input type="checkbox"/>	149F307	Foot & Ankle Positioners - XL	E0951, E0952	\$165

Accessories				
	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	149LF16	Battery Utility Tray <i>(limit 50lbs, includes 1 utility hook)</i>	E1029	\$452
<input type="checkbox"/>	149LF17	Ventilator Tray <i>(Vertical mount)</i>	E1029	\$415
<input type="checkbox"/>	149LF6	Utility Hooks & ZIPPPIE Bag		\$108
<input type="checkbox"/>	149HE76	Headrest Cover - Canopy <i>(Includes hardware)</i>	K0108	\$300
<input type="checkbox"/>	149A468	Underseat Storage Basket <i>(N/A with Battery Tray)</i>		\$155
<input type="checkbox"/>	149LF5	Cup Holder		\$75
<input type="checkbox"/>	149A466	Curb Tipper		\$100
<input type="checkbox"/>	149A467	Rain Cover <i>(Requires Headrest Cover - Canopy)</i>		\$195
<input type="checkbox"/>	149A1	ZIPPPIE Backpack		\$81
<input type="checkbox"/>	149A469	Tool Kit		\$43

Child's Measurements (Not Required)



Child's Name: _____

Age: _____

Diagnosis: _____

Height: _____

Weight: _____

- | | |
|-------------------------------|-------------------------------|
| 1). Top of Head to Seat: | 7). Width of Head: |
| 2). Top of Shoulder to Seat*: | 8). Shoulder Width: |
| 3). Under Arm to Seat: | 9). Chest Width: |
| 4). Trunk Depth: | 10). Hip Width: |
| 5). Actual Sitting Depth*: | 11). Left Lower Leg Length*: |
| 6). Foot Length: | 12). Right Lower Leg Length*: |

*Critical Dimension



Sunrise Medical (US) LLC2842 Business Park Ave. · Fresno, CA 93727 · USA