



Mygo Stander

February 2023

Mark For:

Date: _____

Dealer Acct #: _____

Dealer: _____

Dealer Contact: _____

Dealer Address: _____

Dealer City: _____ ST: _____ ZIP: _____

Dealer Phone: _____ Fax: _____

Confirmation Email: _____

Confirm Via: Fax Email
Submitting for:
 Quote **Order**

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____

Attention: _____

Address: _____

Address: _____

Ship To City: _____ ST: _____ ZIP: _____

Ship To Phone: _____ Fax: _____

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

Mygo Stander
HCPCS: E0641
Size 1 Weight Capacity: 110lbs | Size 2 Weight Capacity: 132lbs
Packages

Package includes Mygo Stander, Stander Covers, Shoulder Pad and Cover, Contoured Headrest Hardware and Cushion - Black, Fixed Chest Laterals, Chest Wraparound Harness, Supine Derotation Hip Belt, Knee Cup Height Extension Brackets, Basic Knee Straps, Tray, and Sandals with Straps. Substitutions are not available for Packages.

	Item Code	Description	HCPCS	Retail
<input type="radio"/>	PMSTB1	Size 1 Mygo Stander Package - Blue <i>(Includes Medium Sandals)</i>	E0641	\$7,611
<input type="radio"/>	PMSTP1	Size 1 Mygo Stander Package - Pink <i>(Includes Medium Sandals)</i>	E0641	\$7,611
<input type="radio"/>	PMSTB2	Size 2 Mygo Stander Package - Blue <i>(Includes Large Sandals)</i>	E0641	\$7,611
<input type="radio"/>	PMSTP2	Size 2 Mygo Stander Package - Pink <i>(Includes Large Sandals)</i>	E0641	\$7,611

Mygo Stander

Selection Required. Includes manual angle adjustment, polyurethane chest pad, hip pad and knee supports, two piece footplate and chassis.

	Item Code	Description	HCPCS	Retail
<input type="radio"/>	143-1600S	Size 1 Mygo Stander		
	143-1600-1	Mygo Stander Frame	E0641	\$1,670
	143-1600-2	Pivot Base		\$1,326
	143-1600-3	Chest Pad Assembly		\$166
	143-1600-4	Hip Assembly		\$166
	143-1600-5	Knee Supports		\$222
	143-1600-6	Positioning Pad Covers		\$393
	143-1600-7	Swivel Casters		\$437
	143-1600-8	Two Piece Footplate		\$386
<input type="radio"/>	143-2600S	Size 2 Mygo Stander		
	143-2600-1	Mygo Stander Frame	E0641	\$1,670
	143-2600-2	Pivot Base		\$1,326
	143-2600-3	Chest Pad		\$166
	143-2600-4	Hip Pad		\$166
	143-2600-5	Knee Supports		\$222
	143-2600-6	Positioning Pad Covers		\$393
	143-2600-7	Swivel Casters		\$437
	143-2600-8	Two Piece Footplate		\$386

Cover

Selection Required. Must match Stander size. Includes chest and hip pad covers and knee support covers.

	Item Code	Description	HCPCS	Retail
<input type="radio"/>	143-1770-08	Size 1 Mygo Stander Covers - Blue		N/C
<input type="radio"/>	143-1770-07	Size 1 Mygo Stander Covers - Orange		N/C
<input type="radio"/>	143-1770-09	Size 1 Mygo Stander Covers - Pink		N/C
<input type="radio"/>	143-2770-08	Size 2 Mygo Stander Covers - Blue		N/C
<input type="radio"/>	143-2770-07	Size 2 Mygo Stander Covers - Orange		N/C
<input type="radio"/>	143-2770-09	Size 2 Mygo Stander Covers - Pink		N/C

Shoulder Support

Cover selection required with Shoulder Pad. Shoulder Support requires Head Support. Shoulder support pad must be selected for supine and head support use.

	Item Code	Description	HCPCS	Retail
<input type="radio"/>	143-1638	Mygo Stander Shoulder Pad		\$381
<input type="radio"/>	143-1768-08	Mygo Stander Shoulder Pad Cover - Blue		N/C
<input type="radio"/>	143-1768-07	Mygo Stander Shoulder Pad Cover - Orange		N/C
<input type="radio"/>	143-1768-09	Mygo Stander Shoulder Pad Cover - Pink		N/C
<input type="radio"/>		No Cover		N/C

Head Support

Cushion selection required with hardware. Head Support requires Shoulder Support.

	Item Code	Description	HCPCS	Retail
<input type="radio"/>	120-859	Contoured Headrest Support Hardware <i>(Requires Cushion)</i>	E1028	\$376
<input type="radio"/>	137-699-05VIN	Contoured Headrest Cushion Infection Control - Black Vinyl	E0955	\$176
<input type="radio"/>	137-800-05	Contoured Headrest Cushion - Black	E0955	\$172
<input type="radio"/>	117-697	Flat Headrest Support Hardware <i>(Requires Cushion)</i>	E1028	\$231
<input type="radio"/>	117-850-05	Flat Headrest Cushion - Black	E0955	\$118
<input type="radio"/>	120-893-05	Flat Headrest Laterals with Covers - Black (Pair) <i>Requires Flat Headrest. Not available with Contoured Headrest.</i>	E0956	\$330

Wraparound Harness

Harnesses may be used as chest or hip harnesses. If used as chest harness, chest laterals are required. If used as hip harness, hip laterals are required.

	Item Code	Description	Quantity 1	Quantity 2	Retail (Ea)	Retail (Pr)
<input type="radio"/>	143-1780-08	Size 1 Chest/Hip Wraparound Harness - Blue	<input type="radio"/>	<input type="radio"/>	\$274	\$548
<input type="radio"/>	143-1780-07	Size 1 Chest/Hip Wraparound Harness - Orange	<input type="radio"/>	<input type="radio"/>	\$274	\$548
<input type="radio"/>	143-1780-09	Size 1 Chest/Hip Wraparound Harness - Pink	<input type="radio"/>	<input type="radio"/>	\$274	\$548
<input type="radio"/>	143-1775-14	Size 1 Chest/Hip Wraparound Harness Infection Control - Black	<input type="radio"/>	<input type="radio"/>	\$274	\$548
<input type="radio"/>	143-2780-08	Size 2 Chest/Hip Wraparound Harness - Blue	<input type="radio"/>	<input type="radio"/>	\$281	\$562
<input type="radio"/>	143-2780-07	Size 2 Chest/Hip Wraparound Harness - Orange	<input type="radio"/>	<input type="radio"/>	\$281	\$562
<input type="radio"/>	143-2780-09	Size 2 Chest/Hip Wraparound Harness - Pink	<input type="radio"/>	<input type="radio"/>	\$281	\$562
<input type="radio"/>	143-2775-14	Size 3 Chest/Hip Wraparound Harness Infection Control - Black	<input type="radio"/>	<input type="radio"/>	\$281	\$562
<input type="radio"/>	143-3780-08	Size 3 Chest/Hip Wraparound Harness - Blue	<input type="radio"/>	<input type="radio"/>	\$290	\$580
<input type="radio"/>	143-3780-07	Size 3 Chest/Hip Wraparound Harness - Orange	<input type="radio"/>	<input type="radio"/>	\$290	\$580
<input type="radio"/>	143-3780-09	Size 3 Chest/Hip Wraparound Harness - Pink	<input type="radio"/>	<input type="radio"/>	\$290	\$580
<input type="radio"/>	143-3775-14	Size 3 Chest/Hip Wraparound Harness Infection Control - Black	<input type="radio"/>	<input type="radio"/>	\$290	\$580

Lateral Support

Chest Laterals required if ordering Wraparound Harness to be used as chest harness. Hip Laterals required if ordering Wraparound Harness to be used as hip harness.

	Item Code	Description	HCPCS	Retail
<input type="radio"/>	143-1619	Size 1 Fixed Chest Laterals (Pair)		\$165
<input type="radio"/>	143-2619	Size 2 Fixed Chest Laterals (Pair)		\$165
<input type="radio"/>	143-1675	Size 1 and 2 Flip-Away Chest Laterals (Pair)		\$530
<input type="radio"/>	143-1663	Hip Laterals (Pair) <i>Not available with Derotation Belt</i>		\$165

Hip Belt*Not available with Hip Laterals or Hip Wraparound Harness*

	Item Code	Description	HCPCS	Retail
○	143-1798-08	Size 1 Prone Derotation Hip Belt with Hardware - Blue		\$311
○	143-1798-07	Size 1 Prone Derotation Hip Belt with Hardware - Orange		\$311
○	143-1798-09	Size 1 Prone Derotation Hip Belt with Hardware - Pink		\$311
○	143-2798-08	Size 2 Prone Derotation Hip Belt with Hardware - Blue		\$352
○	143-2798-07	Size 2 Prone Derotation Hip Belt with Hardware - Orange		\$352
○	143-2798-09	Size 2 Prone Derotation Hip Belt with Hardware - Pink		\$352
○	143-3798-08	Size 3 Prone Derotation Hip Belt with Hardware - Blue		\$387
○	143-3798-07	Size 3 Prone Derotation Hip Belt with Hardware - Orange		\$387
○	143-3798-09	Size 3 Prone Derotation Hip Belt with Hardware - Pink		\$387
○	143-1783-08	Size 1 Supine Derotation Hip Belt with Hardware - Blue		\$311
○	143-1783-07	Size 1 Supine Derotation Hip Belt with Hardware - Orange		\$311
○	143-1783-09	Size 1 Supine Derotation Hip Belt with Hardware - Pink		\$311
○	143-2783-08	Size 2 Supine Derotation Hip Belt with Hardware - Blue		\$352
○	143-2783-07	Size 2 Supine Derotation Hip Belt with Hardware - Orange		\$352
○	143-2783-09	Size 2 Supine Derotation Hip Belt with Hardware - Pink		\$352
○	143-3783-08	Size 3 Supine Derotation Hip Belt with Hardware - Blue		\$387
○	143-3783-07	Size 3 Supine Derotation Hip Belt with Hardware - Orange		\$387
○	143-3783-09	Size 3 Supine Derotation Hip Belt with Hardware - Pink		\$387

Hip Infill

	Item Code	Description	HCPCS	Retail
○	143-1666	Hip Infill Pad		\$340
○	143-1799-08	Size 1 Hip Infill Pad Cover - Blue		\$143
○	143-1799-07	Size 1 Hip Infill Pad Cover - Orange		\$143
○	143-1799-09	Size 1 Hip Infill Pad Cover - Pink		\$143
○	143-2799-08	Size 2 Hip Infill Pad Cover - Blue		\$164
○	143-2799-07	Size 2 Hip Infill Pad Cover - Orange		\$164
○	143-2799-09	Size 2 Hip Infill Pad Cover - Pink		\$164

Knee

	Item Code	Description	HCPCS	Retail
○	143-1690	Knee Cup Height Extension Brackets (Pair) <i>Increases distance from footplate to top of knee by 4"</i>		\$239
○	143-1776-08	Size 1 Basic Knee Straps - Blue (Pair)		\$187
○	143-1776-07	Size 1 Basic Knee Straps - Orange (Pair)		\$187
○	143-1776-09	Size 1 Basic Knee Straps - Pink (Pair)		\$187
○	143-1776-14	Size 1 Basic Knee Straps Infection Control - Black (Pair)		\$187
○	143-2776-08	Size 2 Basic Knee Straps - Blue (Pair)		\$229
○	143-2776-07	Size 2 Basic Knee Straps - Orange (Pair)		\$229
○	143-2776-09	Size 2 Basic Knee Straps - Pink (Pair)		\$229
○	143-2776-14	Size 2 Basic Knee Straps Infection Control - Black (Pair)		\$229
○	143-1784-08	Split Knee Straps - Blue (Pair)		\$344
○	143-1784-07	Split Knee Straps - Orange (Pair)		\$344
○	143-1784-09	Split Knee Straps - Pink (Pair)		\$344

Accessories

	Item Code	Description	HCPCS	Retail	
○	143-1609	Size 1 Mygo Stander Tray		\$449	
○	143-2609	Size 2 Mygo Stander Tray		\$449	
○	143-1633	Size 1 Mygo Stander Tray with Bowl Insert		\$465	
○	143-2633	Size 2 Mygo Stander Tray with Bowl Insert		\$465	
○	143-1649	Grab Rail / Toy Bar for Tray		\$290	
○	152-1600	Small Sandals with Straps (Pair)		\$267	
○	152-2600	Medium Sandals with Straps (Pair)		\$291	
○	152-3600	Large Sandals with Straps (Pair)		\$314	
○	143-1615	Sandal Raiser 1/2" (Each) <i>Requires Sandal Selection</i>	Quantity:	(1 - 6)	\$131

Technical Specifications - Imperial

Specification	Size 1	Size 2
Age (approximate)	4 - 10 years	8 - 14 years
User Height (Min / Max)	42.5" / 57.1"	53.1" / 66.9"
Maximum User Weight	110 lbs	132 lbs
Chest Support Height	29.9" - 42.5"	38.6" - 51.2"
Chest Width	6.3" - 11"	7.9" - 12.6"
Hip Support Height	18.1" - 28.7"	24.4" - 34.4"
Knee Cup Height	11" - 15"	12.25" - 16.25"
Knee Cup Height with Extension Bracket	15" - 19"	16.25" - 20.25"
Hip Width	8.7" - 13"	11" - 15.4"
Footplate Angle: Plantarflexion / Dorsiflexion	10° / 10°	10° / 10°
Maximum Knee Contracture	25°	25°
Tray Height Range	0" - 6"	0" - 6"
Tray Angle Adjustment: Prone / Supine	45° / 45°	45° / 45°
Product Weight	48 lbs	55 lbs
Product Angle Range: Prone / Supine	0° - 80°	0° - 80°
Stander Footprint (length x width)	38.5" x 24.8"	44.5" x 25.6"

Technical Specifications - Metric

Specification	Size 1	Size 2
Age (approximate)	4 - 10 years	8 - 14 years
User Height (Min / Max)	1080 mm / 1450 mm	1350 mm / 1700 mm
Maximum User Weight	50 kgs	60 kgs
Chest Support Height	760 mm - 1080 mm	980 mm - 1300 mm
Chest Width	160 mm - 280 mm	200 mm - 320 mm
Hip Support Height	460 mm - 730 mm	620 mm - 875 mm
Knee Cup Height	280 mm - 380 mm	310mm - 410 mm
Knee Cup Height with Extension Bracket	381 mm - 481 mm	411 mm - 511 mm
Hip Width	220 mm - 330 mm	280 mm - 390 mm
Footplate Angle: Plantarflexion / Dorsiflexion	10° / 10°	10° / 10°
Maximum Knee Contracture	25°	25°
Tray Height Range	0 mm - 150 mm	0 mm - 150 mm
Tray Angle Adjustment: Prone / Supine	45° / 45°	45° / 45°
Product Weight	22 kgs	25 kgs
Product Angle Range: Prone / Supine	0° - 80°	0° - 80°
Stander Footprint (length x width)	978 mm - 630 mm	1130 mm x 650 mm



Sunrise Medical (US) LLC2842 Business Park Ave. • Fresno, CA 93727 • USA

Customer Service: 800-333-4000 Fax: 800-300-7502 Email: orders@sunmed.com www.sunrisemedical.com

© 2/2023 MK-129850 Rev E