



ORDER FORM

# Mygo Stander

Stand out from the crowd in this modular 3-in-1 adjustable standing system with support for hip, knee & ankle contractures

## Customer

Account			
Account #		Date	
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #	
Buyer			
Mark For			

## Shipping Address

Name			
Address			
City		State	
ZIP		Phone	
Email / Fax			

The HCPCS CODES herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative HCPCS CODES, including misc. codes to ensure access for their enrollees. The use of HCPCS CODES does not ensure coverage or payment.

## Mygo Stander

Size 1 Weight Capacity: 110 lbs | Size 2 Weight Capacity: 132 lbs

✓ OPTION	DESCRIPTION	HCPCS	RETAIL
<b>Packages</b>			<b>\$</b>
Package includes Mygo Stander, Stander Covers, Shoulder Pad and Cover, Contoured Headrest Hardware and Cushion - Black, Fixed Chest Laterals, Chest Wraparound Harness, Supine Derotation Hip Belt, Knee Cup Height Extension Brackets, Basic Knee Straps, Tray, Sandals with Straps. Substitutions are not available for Packages.			
<input type="checkbox"/> PMSTB1	Size 1 Mygo Stander Package - Blue [Includes Medium Sandals]	E0641	\$7,789
<input type="checkbox"/> PMSTP1	Size 1 Mygo Stander Package - Pink [Includes Medium Sandals]	E0641	\$7,789
<input type="checkbox"/> PMSTB2	Size 2 Mygo Stander Package - Blue [Includes Large Sandals]	E0641	\$7,905
<input type="checkbox"/> PMSTP2	Size 2 Mygo Stander Package - Pink [Includes Large Sandals]	E0641	\$7,905
<b>Mygo Stander</b>			<b>\$</b>

Selection Required.

Includes manual angle adjustment, polyurethane chest pad, hip pad and knee supports, two piece footplate and chassis.

<input type="checkbox"/> <b>143-1600S</b>	<b>Size 1 Mygo Stander</b>		
143-1600-1	Mygo Stander Frame	E0641	\$1,709
143-1600-2	Pivot Base		\$1,357
143-1600-3	Chest Pad Assembly		\$170
143-1600-4	Hip Assembly		\$170
143-1600-5	Knee Supports		\$227
143-1600-6	Positioning Pad Covers		\$402
143-1600-7	Swivel Casters		\$447
143-1600-8	Two Piece Footplate		\$395
<input type="checkbox"/> <b>143-2600S</b>	<b>Size 2 Mygo Stander</b>		
143-2600-1	Mygo Stander Frame	E0641	\$1,709
143-2600-2	Pivot Base		\$1,357
143-2600-3	Chest Pad		\$170
143-2600-4	Hip Pad		\$170
143-2600-5	Knee Supports		\$227
143-2600-6	Positioning Pad Covers		\$402
143-2600-7	Swivel Casters		\$447
143-2600-8	Two Piece Footplate		\$395

<b>Cover</b>	<b>\$</b>
--------------	-----------

Selection Required. Must match Stander size.

Includes chest and hip pad covers and knee support covers.

<input type="checkbox"/> 143-1770-08	Size 1 Mygo Stander Covers - Blue	N/C
<input type="checkbox"/> 143-1770-07	Size 1 Mygo Stander Covers - Orange	N/C
<input type="checkbox"/> 143-1770-09	Size 1 Mygo Stander Covers - Pink	N/C
<input type="checkbox"/> 143-2770-08	Size 2 Mygo Stander Covers - Blue	N/C
<input type="checkbox"/> 143-2770-07	Size 2 Mygo Stander Covers - Orange	N/C
<input type="checkbox"/> 143-2770-09	Size 2 Mygo Stander Covers - Pink	N/C

<b>Shoulder Support</b>	<b>\$</b>
-------------------------	-----------

Cover selection required with Shoulder Pad. Shoulder Support requires Head Support.

Shoulder support pad must be selected for supine use and head support use.

<input type="checkbox"/> 143-1638	Mygo Stander Shoulder Pad	\$390
<input type="checkbox"/> 143-1768-08	Mygo Stander Shoulder Pad Cover - Blue	N/C
<input type="checkbox"/> 143-1768-07	Mygo Stander Shoulder Pad Cover - Orange	N/C
<input type="checkbox"/> 143-1768-09	Mygo Stander Shoulder Pad Cover - Pink	N/C
<input type="checkbox"/>	No Cover	N/C

## Mygo Stander

✓ OPTION	DESCRIPTION	RETAIL
Head Support		\$

Selection of Headrest Support Hardware [Contoured/Flat] Requires Cushion

Headrest Laterals only available with Flat Headrest

Head Support selection requires Shoulder Support.

<input type="checkbox"/>	120-859	Contoured Headrest Support Hardware	\$385
<input type="checkbox"/>	137-699-05VIN	Contoured Headrest Cushion Infection Control - Black Vinyl	\$205
<input type="checkbox"/>	137-800-05	Contoured Headrest Cushion - Black	\$176
<input type="checkbox"/>	117-697	Flat Headrest Support Hardware	\$236
<input type="checkbox"/>	117-850-05	Flat Headrest Cushion - Black	\$121
<input type="checkbox"/>	120-893-05	Flat Headrest Laterals with Covers - Black [Pair]	\$338

## Wraparound Harness

\$

Harnesses may be used as chest or hip harnesses.

If used as chest harness, chest laterals are required. If used as hip harness, hip laterals are required.

<input type="checkbox"/>	143-1780-08	QTY	Size 1 Chest/Hip Wraparound Harness - Blue	\$280
<input type="checkbox"/>	143-1780-07	QTY	Size 1 Chest/Hip Wraparound Harness - Orange	\$280
<input type="checkbox"/>	143-1780-09	QTY	Size 1 Chest/Hip Wraparound Harness - Pink	\$280
<input type="checkbox"/>	143-1775-14	QTY	Size 1 Chest/Hip Wraparound Harness Infection Control - Black	\$280
<input type="checkbox"/>	143-2780-08	QTY	Size 2 Chest/Hip Wraparound Harness - Blue	\$288
<input type="checkbox"/>	143-2780-07	QTY	Size 2 Chest/Hip Wraparound Harness - Orange	\$288
<input type="checkbox"/>	143-2780-09	QTY	Size 2 Chest/Hip Wraparound Harness - Pink	\$288
<input type="checkbox"/>	143-2775-14	QTY	Size 3 Chest/Hip Wraparound Harness Infection Control - Black	\$288
<input type="checkbox"/>	143-3780-08	QTY	Size 3 Chest/Hip Wraparound Harness - Blue	\$297
<input type="checkbox"/>	143-3780-07	QTY	Size 3 Chest/Hip Wraparound Harness - Orange	\$297
<input type="checkbox"/>	143-3780-09	QTY	Size 3 Chest/Hip Wraparound Harness - Pink	\$297
<input type="checkbox"/>	143-3775-14	QTY	Size 3 Chest/Hip Wraparound Harness Infection Control - Black	\$297

## Lateral Support

\$

Chest Laterals required if ordering Wraparound Harness to be used as chest harness.

Hip Laterals required if ordering Wraparound Harness to be used as hip harness.

<input type="checkbox"/>	143-1619	Size 1 Fixed Chest Laterals [Pair]	\$169
<input type="checkbox"/>	143-2619	Size 2 Fixed Chest Laterals [Pair]	\$169
<input type="checkbox"/>	143-1675	Size 1 and 2 Flip-Away Chest Laterals [Pair]	\$542
<input type="checkbox"/>	143-1663	Hip Laterals [Pair] Not available with Derotation Belt	\$169

## Hip Infill

\$

<input type="checkbox"/>	143-1666	Hip Infill Pad	\$348
<input type="checkbox"/>	143-1799-08	Size 1 Hip Infill Pad Cover - Blue	\$146
<input type="checkbox"/>	143-1799-07	Size 1 Hip Infill Pad Cover - Orange	\$146
<input type="checkbox"/>	143-1799-09	Size 1 Hip Infill Pad Cover - Pink	\$146
<input type="checkbox"/>	143-2799-08	Size 2 Hip Infill Pad Cover - Blue	\$168
<input type="checkbox"/>	143-2799-07	Size 2 Hip Infill Pad Cover - Orange	\$168
<input type="checkbox"/>	143-2799-09	Size 2 Hip Infill Pad Cover - Pink	\$168

## Mygo Stander

✓ OPTION	DESCRIPTION	RETAIL
----------	-------------	--------

### Hip Belt \$

Not available with Hip Laterals or Hip Wraparound Harness

<input type="checkbox"/>	143-1798-08	Size 1 Prone Derotation Hip Belt with Hardware - Blue	\$318
<input type="checkbox"/>	143-1798-07	Size 1 Prone Derotation Hip Belt with Hardware - Orange	\$318
<input type="checkbox"/>	143-1798-09	Size 1 Prone Derotation Hip Belt with Hardware - Pink	\$318
<input type="checkbox"/>	143-2798-08	Size 2 Prone Derotation Hip Belt with Hardware - Blue	\$360
<input type="checkbox"/>	143-2798-07	Size 2 Prone Derotation Hip Belt with Hardware - Orange	\$360
<input type="checkbox"/>	143-2798-09	Size 2 Prone Derotation Hip Belt with Hardware - Pink	\$360
<input type="checkbox"/>	143-3798-08	Size 3 Prone Derotation Hip Belt with Hardware - Blue	\$396
<input type="checkbox"/>	143-3798-07	Size 3 Prone Derotation Hip Belt with Hardware - Orange	\$396
<input type="checkbox"/>	143-3798-09	Size 3 Prone Derotation Hip Belt with Hardware - Pink	\$396
<input type="checkbox"/>	143-1783-08	Size 1 Supine Derotation Hip Belt with Hardware - Blue	\$318
<input type="checkbox"/>	143-1783-07	Size 1 Supine Derotation Hip Belt with Hardware - Orange	\$318
<input type="checkbox"/>	143-1783-09	Size 1 Supine Derotation Hip Belt with Hardware - Pink	\$318
<input type="checkbox"/>	143-2783-08	Size 2 Supine Derotation Hip Belt with Hardware - Blue	\$360
<input type="checkbox"/>	143-2783-07	Size 2 Supine Derotation Hip Belt with Hardware - Orange	\$360
<input type="checkbox"/>	143-2783-09	Size 2 Supine Derotation Hip Belt with Hardware - Pink	\$360
<input type="checkbox"/>	143-3783-08	Size 3 Supine Derotation Hip Belt with Hardware - Blue	\$396
<input type="checkbox"/>	143-3783-07	Size 3 Supine Derotation Hip Belt with Hardware - Orange	\$396
<input type="checkbox"/>	143-3783-09	Size 3 Supine Derotation Hip Belt with Hardware - Pink	\$396

### Knee \$

Knee Cup Height Extension Brackets increase distance from footplate to top of knee by 4"

<input type="checkbox"/>	143-1690	Knee Cup Height Extension Brackets (Pair)	\$245
<input type="checkbox"/>	143-1776-08	Size 1 Basic Knee Straps - Blue (Pair)	\$191
<input type="checkbox"/>	143-1776-07	Size 1 Basic Knee Straps - Orange (Pair)	\$191
<input type="checkbox"/>	143-1776-09	Size 1 Basic Knee Straps - Pink (Pair)	\$191
<input type="checkbox"/>	143-1776-14	Size 1 Basic Knee Straps Infection Control - Black (Pair)	\$191
<input type="checkbox"/>	143-2776-08	Size 2 Basic Knee Straps - Blue (Pair)	\$234
<input type="checkbox"/>	143-2776-07	Size 2 Basic Knee Straps - Orange (Pair)	\$234
<input type="checkbox"/>	143-2776-09	Size 2 Basic Knee Straps - Pink (Pair)	\$234
<input type="checkbox"/>	143-2776-14	Size 2 Basic Knee Straps Infection Control - Black (Pair)	\$234
<input type="checkbox"/>	143-1784-08	Split Knee Straps - Blue (Pair)	\$352
<input type="checkbox"/>	143-1784-07	Split Knee Straps - Orange (Pair)	\$352
<input type="checkbox"/>	143-1784-09	Split Knee Straps - Pink (Pair)	\$352

### Accessories \$

<input type="checkbox"/>	143-1609	Size 1 Mygo Stander Tray	\$460
<input type="checkbox"/>	143-2609	Size 2 Mygo Stander Tray	\$460
<input type="checkbox"/>	143-1633	Size 1 Mygo Stander Tray with Bowl Insert	\$476
<input type="checkbox"/>	143-2633	Size 2 Mygo Stander Tray with Bowl Insert	\$476
<input type="checkbox"/>	143-1649	Grab Rail / Toy Bar for Tray	\$297
<input type="checkbox"/>	152-1600	Small Sandals with Straps (Pair)	\$273
<input type="checkbox"/>	152-2600	Medium Sandals with Straps (Pair)	\$298
<input type="checkbox"/>	152-3600	Large Sandals with Straps (Pair)	\$321
<input type="checkbox"/>	143-1615	QTY Sandal Raiser 1/2" (Each) Requires Sandal Selection.	\$134

## Technical Specifications - Imperial

Frame Sizing Specification	Size 1	Size 2
Age (approximate)	4 - 10 years	8 - 14 years
User Height (Min / Max)	42.5" / 57.1"	53.1" / 66.9"
Maximum User Weight	110 lbs	132 lbs
Chest Support Height	29.9" - 42.5"	38.6" - 51.2"
Chest Width	6.3" - 11"	7.9" - 12.6"
Hip Support Height	18.1" - 28.7"	24.4" - 34.4"
Knee Cup Height	11" - 15"	12.25" - 16.25"
Knee Cup Height with Extension Bracket	15" - 19"	16.25" - 20.25"
Hip Width	8.7" - 13"	11" - 15.4"
Footplate Angle: Plantarflexion / Dorsiflexion	10° / 10°	10° / 10°
Maximum Knee Contracture	25°	25°
Tray Height Range	0" - 6"	0" - 6"
Tray Angle Adjustment: Prone / Supine	45° / 45°	45° / 45°
Product Weight	48 lbs	55 lbs
Product Angle Range: Prone / Supine	0° - 80°	0° - 80°
Stander Footprint (length x width)	38.5" x 24.8"	44.5" x 25.6"

## Technical Specifications - Metric

Frame Sizing Specification	Size 2	Size 3
Age (approximate)	4 - 10 years	8 - 14 years
User Height (Min / Max)	1080 mm / 1450 mm	1350 mm / 1700 mm
Maximum User Weight	50 kgs	60 kgs
Chest Support Height	760 mm - 1080 mm	980 mm - 1300 mm
Chest Width	160 mm - 280 mm	200 mm - 320 mm
Hip Support Height	460 mm - 730 mm	620 mm - 875 mm
Knee Cup Height	280 mm - 380 mm	310mm - 410 mm
Knee Cup Height with Extension Bracket	381 mm - 481 mm	411 mm - 511 mm
Hip Width	220 mm - 330 mm	280 mm - 390 mm
Footplate Angle: Plantarflexion / Dorsiflexion	10° / 10°	10° / 10°
Maximum Knee Contracture	25°	25°
Tray Height Range	0 mm - 150 mm	0 mm - 150 mm
Tray Angle Adjustment: Prone / Supine	45° / 45°	45° / 45°
Product Weight	22 kgs	25 kgs
Product Angle Range: Prone / Supine	0° - 80°	0° - 80°
Stander Footprint (length x width)	978 mm - 630 mm	1130 mm x 650 mm

### Sunrise Medical (US) LLC

2842 Business Park Avenue  
 Fresno, CA 93727 USA  
 Tel.: 1-800-333-4000  
 Fax: 1-800-300-7502  
 E-Mail: webmaster@sunmed.com



If you are visually impaired, this document can be viewed in PDF format at [www.SunriseMedical.com](http://www.SunriseMedical.com)



MK-129850 Rev. F ©01/2024



800.333.4000

[www.SunriseMedical.com](http://www.SunriseMedical.com)