



# Mygo Stander

January 2022

**Mark For:**

Date: \_\_\_\_\_

Dealer Acct #: \_\_\_\_\_

Dealer: \_\_\_\_\_

Dealer Contact: \_\_\_\_\_

Dealer Address: \_\_\_\_\_

Dealer City: \_\_\_\_\_ ST: \_\_\_\_\_ ZIP: \_\_\_\_\_

Dealer Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Confirmation Email: \_\_\_\_\_

Confirm Via:  Fax  Email

Submitting for:  Quote  Order

PO#: \_\_\_\_\_

**ADDITIONAL SHIPPING INFORMATION**

Ship To: \_\_\_\_\_

Attention: \_\_\_\_\_

Address: \_\_\_\_\_

Address: \_\_\_\_\_

Ship To City: \_\_\_\_\_ ST: \_\_\_\_\_ ZIP: \_\_\_\_\_

Ship To Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

**Mygo Stander****HCPCS: E0641**

Size 1 Weight Capacity: 110lbs | Size 2 Weight Capacity: 132lbs

**Packages**

Package includes Mygo Stander, Stander Covers, Shoulder Pad and Cover, Contoured Headrest Hardware and Cushion - Black, Fixed Chest Laterals, Chest Wraparound Harness, Supine Derotation Hip Belt, Knee Cup Height Extension Brackets, Basic Knee Straps, Tray, Sandals with Straps and Product Assembly. Substitutions are not available for Packages.

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	PMSTB1	Size 1 Mygo Stander Package - Blue (Includes Medium Sandals)	E0641	\$7,374
<input type="checkbox"/>	PMSTP1	Size 1 Mygo Stander Package - Pink (Includes Medium Sandals)	E0641	\$7,374
<input type="checkbox"/>	PMSTB2	Size 2 Mygo Stander Package - Blue (Includes Large Sandals)	E0641	\$7,374
<input type="checkbox"/>	PMSTP2	Size 2 Mygo Stander Package - Pink (Includes Large Sandals)	E0641	\$7,374

**Mygo Stander**

Selection Required. Includes manual angle adjustment, polyurethane chest pad, hip pad and knee supports, two piece footplate and chassis.

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	<b>143-1600S</b>	<b>Size 1 Mygo Stander</b>		
	143-1600-1	Mygo Stander Frame	E0641	\$1,590
	143-1600-2	Pivot Base		\$1,263
	143-1600-3	Chest Pad Assembly		\$158
	143-1600-4	Hip Assembly		\$158
	143-1600-5	Knee Supports		\$211
	143-1600-6	Positioning Pad Covers		\$374
	143-1600-7	Swivel Casters		\$416
	143-1600-8	Two Piece Footplate		\$368
<input type="checkbox"/>	<b>143-2600S</b>	<b>Size 2 Mygo Stander</b>		
	143-2600-1	Mygo Stander Frame	E0641	\$1,590
	143-2600-2	Pivot Base		\$1,263
	143-2600-3	Chest Pad		\$158
	143-2600-4	Hip Pad		\$158
	143-2600-5	Knee Supports		\$211
	143-2600-6	Positioning Pad Covers		\$374
	143-2600-7	Swivel Casters		\$416
	143-2600-8	Two Piece Footplate		\$368

**Cover**

*Selection Required. Must match Stander size. Includes chest and hip pad covers and knee support covers.*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	143-1770-08	Size 1 Mygo Stander Covers - Blue		N/C
<input type="checkbox"/>	143-1770-07	Size 1 Mygo Stander Covers - Orange		N/C
<input type="checkbox"/>	143-1770-09	Size 1 Mygo Stander Covers - Pink		N/C
<input type="checkbox"/>	143-2770-08	Size 2 Mygo Stander Covers - Blue		N/C
<input type="checkbox"/>	143-2770-07	Size 2 Mygo Stander Covers - Orange		N/C
<input type="checkbox"/>	143-2770-09	Size 2 Mygo Stander Covers - Pink		N/C

**Shoulder Support**

*Cover selection required with Shoulder Pad. Shoulder Support requires Head Support.*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	143-1638	Mygo Stander Shoulder Pad		\$363
<input type="checkbox"/>	143-1768-08	Mygo Stander Shoulder Pad Cover - Blue		N/C
<input type="checkbox"/>	143-1768-07	Mygo Stander Shoulder Pad Cover - Orange		N/C
<input type="checkbox"/>	143-1768-09	Mygo Stander Shoulder Pad Cover - Pink		N/C
<input type="checkbox"/>		No Cover		N/C

**Head Support**

*Cushion selection required with hardware. Head Support requires Shoulder Support.*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	120-859	Contoured Headrest Support Hardware <i>(Requires Cushion)</i>	E1028	\$358
<input type="checkbox"/>	143-1787	Contoured Headrest Cushion Infection Control - Black	E0955	\$168
<input type="checkbox"/>	137-800-05	Contoured Headrest Cushion - Black	E0955	\$164
<input type="checkbox"/>	117-697	Flat Headrest Support Hardware <i>(Requires Cushion)</i>	E1028	\$220
<input type="checkbox"/>	117-850-05	Flat Headrest Cushion - Black	E0955	\$112
<input type="checkbox"/>	120-893-05	Flat Headrest Laterals with Covers - Black (Pair) <i>Requires Flat Headrest. Not available with Contoured Headrest.</i>	E0956	\$314

**Wraparound Harness**

*Harnesses may be used as chest or hip harnesses. If used as chest harness, chest laterals are required. If used as hip harness, hip laterals are required.*

	Item Code	Description	Quantity 1	Quantity 2	Retail (Ea)	Retail (Pr)
<input type="checkbox"/>	143-1780-08	Size 1 Chest/Hip Wraparound Harness - Blue	<input type="checkbox"/>	<input type="checkbox"/>	\$261	\$522
<input type="checkbox"/>	143-1780-07	Size 1 Chest/Hip Wraparound Harness - Orange	<input type="checkbox"/>	<input type="checkbox"/>	\$261	\$522
<input type="checkbox"/>	143-1780-09	Size 1 Chest/Hip Wraparound Harness - Pink	<input type="checkbox"/>	<input type="checkbox"/>	\$261	\$522
<input type="checkbox"/>	143-1775-14	Size 1 Chest/Hip Wraparound Harness Infection Control - Black	<input type="checkbox"/>	<input type="checkbox"/>	\$261	\$522
<input type="checkbox"/>	143-2780-08	Size 2 Chest/Hip Wraparound Harness - Blue	<input type="checkbox"/>	<input type="checkbox"/>	\$268	\$536
<input type="checkbox"/>	143-2780-07	Size 2 Chest/Hip Wraparound Harness - Orange	<input type="checkbox"/>	<input type="checkbox"/>	\$268	\$536
<input type="checkbox"/>	143-2780-09	Size 2 Chest/Hip Wraparound Harness - Pink	<input type="checkbox"/>	<input type="checkbox"/>	\$268	\$536
<input type="checkbox"/>	143-2775-14	Size 3 Chest/Hip Wraparound Harness Infection Control - Black	<input type="checkbox"/>	<input type="checkbox"/>	\$268	\$536
<input type="checkbox"/>	143-3780-08	Size 3 Chest/Hip Wraparound Harness - Blue	<input type="checkbox"/>	<input type="checkbox"/>	\$276	\$552
<input type="checkbox"/>	143-3780-07	Size 3 Chest/Hip Wraparound Harness - Orange	<input type="checkbox"/>	<input type="checkbox"/>	\$276	\$552
<input type="checkbox"/>	143-3780-09	Size 3 Chest/Hip Wraparound Harness - Pink	<input type="checkbox"/>	<input type="checkbox"/>	\$276	\$552
<input type="checkbox"/>	143-3775-14	Size 3 Chest/Hip Wraparound Harness Infection Control - Black	<input type="checkbox"/>	<input type="checkbox"/>	\$276	\$552

**Lateral Support**

*Chest Laterals required if ordering Wraparound Harness to be used as chest harness. Hip Laterals required if ordering Wraparound Harness to be used as hip harness.*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	143-1619	Size 1 Fixed Chest Laterals (Pair)		\$157
<input type="checkbox"/>	143-2619	Size 2 Fixed Chest Laterals (Pair)		\$157
<input type="checkbox"/>	143-1675	Size 1 and 2 Flip-Away Chest Laterals (Pair)		\$505
<input type="checkbox"/>	143-1663	Hip Laterals (Pair) <i>Not available with Derotation Belt</i>		\$157

**Hip Belt***Not available with Hip Laterals or Hip Wraparound Harness*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	143-1798-08	Size 1 Prone Derotation Hip Belt with Hardware - Blue		\$296
<input type="checkbox"/>	143-1798-07	Size 1 Prone Derotation Hip Belt with Hardware - Orange		\$296
<input type="checkbox"/>	143-1798-09	Size 1 Prone Derotation Hip Belt with Hardware - Pink		\$296
<input type="checkbox"/>	143-2798-08	Size 2 Prone Derotation Hip Belt with Hardware - Blue		\$335
<input type="checkbox"/>	143-2798-07	Size 2 Prone Derotation Hip Belt with Hardware - Orange		\$335
<input type="checkbox"/>	143-2798-09	Size 2 Prone Derotation Hip Belt with Hardware - Pink		\$335
<input type="checkbox"/>	143-3798-08	Size 3 Prone Derotation Hip Belt with Hardware - Blue		\$369
<input type="checkbox"/>	143-3798-07	Size 3 Prone Derotation Hip Belt with Hardware - Orange		\$369
<input type="checkbox"/>	143-3798-09	Size 3 Prone Derotation Hip Belt with Hardware - Pink		\$369
<input type="checkbox"/>	143-1783-08	Size 1 Supine Derotation Hip Belt with Hardware - Blue		\$296
<input type="checkbox"/>	143-1783-07	Size 1 Supine Derotation Hip Belt with Hardware - Orange		\$296
<input type="checkbox"/>	143-1783-09	Size 1 Supine Derotation Hip Belt with Hardware - Pink		\$296
<input type="checkbox"/>	143-2783-08	Size 2 Supine Derotation Hip Belt with Hardware - Blue		\$335
<input type="checkbox"/>	143-2783-07	Size 2 Supine Derotation Hip Belt with Hardware - Orange		\$335
<input type="checkbox"/>	143-2783-09	Size 2 Supine Derotation Hip Belt with Hardware - Pink		\$335
<input type="checkbox"/>	143-3783-08	Size 3 Supine Derotation Hip Belt with Hardware - Blue		\$369
<input type="checkbox"/>	143-3783-07	Size 3 Supine Derotation Hip Belt with Hardware - Orange		\$369
<input type="checkbox"/>	143-3783-09	Size 3 Supine Derotation Hip Belt with Hardware - Pink		\$369

**Hip Infill**

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	143-1666	Hip Infill Pad		\$324
<input type="checkbox"/>	143-1799-08	Size 1 Hip Infill Pad Cover - Blue		\$136
<input type="checkbox"/>	143-1799-07	Size 1 Hip Infill Pad Cover - Orange		\$136
<input type="checkbox"/>	143-1799-09	Size 1 Hip Infill Pad Cover - Pink		\$136
<input type="checkbox"/>	143-2799-08	Size 2 Hip Infill Pad Cover - Blue		\$156
<input type="checkbox"/>	143-2799-07	Size 2 Hip Infill Pad Cover - Orange		\$156
<input type="checkbox"/>	143-2799-09	Size 2 Hip Infill Pad Cover - Pink		\$156

**Knee**

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	143-1690	Knee Cup Height Extension Brackets (Pair) <i>Increases distance from footplate to top of knee by 4"</i>		\$228
<input type="checkbox"/>	143-1776-08	Size 1 Basic Knee Straps - Blue (Pair)		\$178
<input type="checkbox"/>	143-1776-07	Size 1 Basic Knee Straps - Orange (Pair)		\$178
<input type="checkbox"/>	143-1776-09	Size 1 Basic Knee Straps - Pink (Pair)		\$178
<input type="checkbox"/>	143-1776-14	Size 1 Basic Knee Straps Infection Control - Black (Pair)		\$178
<input type="checkbox"/>	143-2776-08	Size 2 Basic Knee Straps - Blue (Pair)		\$218
<input type="checkbox"/>	143-2776-07	Size 2 Basic Knee Straps - Orange (Pair)		\$218
<input type="checkbox"/>	143-2776-09	Size 2 Basic Knee Straps - Pink (Pair)		\$218
<input type="checkbox"/>	143-2776-14	Size 2 Basic Knee Straps Infection Control - Black (Pair)		\$218
<input type="checkbox"/>	143-1784-08	Split Knee Straps - Blue (Pair)		\$328
<input type="checkbox"/>	143-1784-07	Split Knee Straps - Orange (Pair)		\$328
<input type="checkbox"/>	143-1784-09	Split Knee Straps - Pink (Pair)		\$328

Accessories				
-------------	--	--	--	--

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	143-1609	Size 1 Mygo Stander Tray		\$428
<input type="checkbox"/>	143-2609	Size 2 Mygo Stander Tray		\$428
<input type="checkbox"/>	143-1633	Size 1 Mygo Stander Tray with Bowl Insert		\$443
<input type="checkbox"/>	143-2633	Size 2 Mygo Stander Tray with Bowl Insert		\$443
<input type="checkbox"/>	143-1649	Grab Rail / Toy Bar for Tray		\$276
<input type="checkbox"/>	152-1600	Small Sandals with Straps (Pair)		\$254
<input type="checkbox"/>	152-2600	Medium Sandals with Straps (Pair)		\$277
<input type="checkbox"/>	152-3600	Large Sandals with Straps (Pair)		\$299
<input type="checkbox"/>	143-1615	Sandal Raiser 1/2" (Each) <i>Requires Sandal Selection</i>	<b>Quantity:</b> <input style="width: 50px;" type="text"/> (1 - 6) <input style="width: 50px;" type="text"/>	\$125

Assembly				
----------	--	--	--	--

	Item Code	Description		Retail
<input type="checkbox"/>	ASSY-LKY	Product Assembly Fee		\$126

**Technical Specifications - Imperial**

Specification	Size 1	Size 2
Age (approximate)	4 - 10 years	8 - 14 years
User Height (Min / Max)	42.5" / 57.1"	53.1" / 66.9"
Maximum User Weight	110 lbs	132 lbs
Chest Support Height	29.9" - 42.5"	38.6" - 51.2"
Chest Width	6.3" - 11"	7.9" - 12.6"
Hip Support Height	18.1" - 28.7"	24.4" - 34.4"
Knee Cup Height	11" - 15"	12.25" - 16.25"
Knee Cup Height with Extension Bracket	15" - 19"	16.25" - 20.25"
Hip Width	8.7" - 13"	11" - 15.4"
Footplate Angle: Plantarflexion / Dorsiflexion	10° / 10°	10° / 10°
Maximum Knee Contracture	25°	25°
Tray Height Range	0" - 6"	0" - 6"
Tray Angle Adjustment: Prone / Supine	45° / 45°	45° / 45°
Product Weight	48 lbs	55 lbs
Product Angle Range: Prone / Supine	0° - 80°	0° - 80°
Stander Footprint (length x width)	38.5" x 24.8"	44.5" x 25.6"

**Technical Specifications - Metric**

Specification	Size 1	Size 2
Age (approximate)	4 - 10 years	8 - 14 years
User Height (Min / Max)	1080 mm / 1450 mm	1350 mm / 1700 mm
Maximum User Weight	50 kgs	60 kgs
Chest Support Height	760 mm - 1080 mm	980 mm - 1300 mm
Chest Width	160 mm - 280 mm	200 mm - 320 mm
Hip Support Height	460 mm - 730 mm	620 mm - 875 mm
Knee Cup Height	280 mm - 380 mm	310mm - 410 mm
Knee Cup Height with Extension Bracket	381 mm - 481 mm	411 mm - 511 mm
Hip Width	220 mm - 330 mm	280 mm - 390 mm
Footplate Angle: Plantarflexion / Dorsiflexion	10° / 10°	10° / 10°
Maximum Knee Contracture	25°	25°
Tray Height Range	0 mm - 150 mm	0 mm - 150 mm
Tray Angle Adjustment: Prone / Supine	45° / 45°	45° / 45°
Product Weight	22 kgs	25 kgs
Product Angle Range: Prone / Supine	0° - 80°	0° - 80°
Stander Footprint (length x width)	978 mm - 630 mm	1130 mm x 650 mm



**Sunrise Medical (US) LLC2842 Business Park Ave. • Fresno, CA 93727 • USA**