



ORDER FORM

# MyWay+

EILK01

Unique, open frame gait trainer designed to support children in an upright, hands-free position.

Customer			
Account			
Account #		Date	
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #	
Buyer			
Mark For			

Shipping Address			
Name			
Address			
City		State	
ZIP		Phone	
Email / Fax			

The HCPCS CODES herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative HCPCS CODES, including misc. codes to ensure access for their enrollees. The use of HCPCS CODES does not ensure coverage or payment.

## MyWay+

✓ OPTION	DESCRIPTION	RETAIL
----------	-------------	--------

### Packages

Includes Frame, Harness, User Handles, Guide Handles, Headrest, Saddle and Ankle Guides.

Packages for Frame Sizes 1 & 2 include Size 1 Saddle. Packages for Frame Size 3 includes Size 2 Saddle.

<input type="checkbox"/>	122QE02	Size 1 MyWay+ Package with Size 1 Harness	\$6,086
<input type="checkbox"/>	122QE03	Size 1 MyWay+ Package with Size 2 Harness	\$6,091
<input type="checkbox"/>	122QE04	Size 1 MyWay+ Package with Size 3 Harness	\$6,123
<input type="checkbox"/>	122QE05	Size 2 MyWay+ Package with Size 3 Harness	\$6,743
<input type="checkbox"/>	122QE06	Size 2 MyWay+ Package with Size 4 Harness	\$6,782
<input type="checkbox"/>	122QE07	Size 2 MyWay+ Package with Size 5 Harness	\$6,822
<input type="checkbox"/>	122QE08	Size 3 MyWay+ Package with Size 4 Harness	\$6,955
<input type="checkbox"/>	122QE09	Size 3 MyWay+ Package with Size 5 Harness	\$6,995
<input type="checkbox"/>	122QE10	Size 3 MyWay+ Package with Size 6 Harness	\$7,014

### Frame

Selection Required. Includes multi-functional casters.

<input type="checkbox"/>	122FT17	Size 1 MyWay+ Frame	\$4,100
<input type="checkbox"/>	122FT18	Size 2 MyWay+ Frame	\$4,720
<input type="checkbox"/>	122FT19	Size 3 MyWay+ Frame	\$4,890

### Harness

Selection Required.

<input type="checkbox"/>	122NH34	Size 1 Harness - Purple	\$363
<input type="checkbox"/>	122NH35	Size 2 Harness - Green	\$368
<input type="checkbox"/>	122NH36	Size 3 Harness - Blue	\$400
<input type="checkbox"/>	122NH37	Size 4 Harness - Red	\$439
<input type="checkbox"/>	122NH38	Size 5 Harness - Orange	\$479
<input type="checkbox"/>	122NH39	Size 6 Harness - Grey	\$498



Size 1 - Purple	Size 2 - Green	Size 3 - Blue	Size 4 - Red	Size 5 - Orange	Size 6 - Grey
-----------------	----------------	---------------	--------------	-----------------	---------------

**MyWay+**

✓ OPTION	DESCRIPTION	RETAIL	
----------	-------------	--------	--

**User Handles** \$

User Handles MUST match Frame Size. Complex User Handles Ball Grip standard; Cylinder Grip optional.

<input type="checkbox"/> 122HY1	Size 1 User Handles - Pair	\$315
<input type="checkbox"/> 122HY2	Size 2 User Handles - Pair	\$315
<input type="checkbox"/> 122HY3	Size 3 User Handles - Pair	\$315
<input type="checkbox"/> 122HY17	Size 1 Complex User Handles	\$1,300
<input type="checkbox"/> 122HY18	Size 2 Complex User Handles	\$1,350
<input type="checkbox"/> 122HY19	Size 3 Complex User Handles	\$1,400
<input type="checkbox"/> 122HY22	Cylinder Grip - Pair	\$60



**Guide Handles** \$

Use Size 1 Adult Handles for Size 1 Frame. Use Size 2 Adult Handles for Size 2 & 3 Frame.

<input type="checkbox"/> 122HY20	Size 1 Adult Handles - Fits Size 1 Frame	\$380
<input type="checkbox"/> 122HY21	Size 2 Adult Handles - Fits Size 2 & 3 Frame	\$380



**Headrests** \$

<input type="checkbox"/> 122HE49	Headrest - One Size	\$227
<input type="checkbox"/> 122HE30	Complex Headrest - One Size	\$645



**Saddles** \$

Saddle Size 1 recommended for Size 1 & 2 Frames. Saddle size 2 recommended for Size 2 & 3 Frames.

<input type="checkbox"/> 1229A19	Size 1 Saddle	\$355
<input type="checkbox"/> 1229A20	Size 2 Saddle	\$358
<input type="checkbox"/> 1229A21	Size 1 Narrow Saddle	\$485
<input type="checkbox"/> 1229A22	Size 2 Narrow Saddle	\$490



**Foot and Leg Positioning Accessories** \$

Ankle Guides MUST match Frame Size.

<input type="checkbox"/> 122DI60	Size 1 Drop - Fits Size 1 Narrow Saddle ONLY	\$202
<input type="checkbox"/> 122DI61	Size 2 Drop - Fits Size 2 Narrow Saddle ONLY	\$220
<input type="checkbox"/> 122DI62	Size 3 Drop - Fits Size 2 Narrow Saddle ONLY	\$232
<input type="checkbox"/> 122DI63	Size 1 Ankle Guides - Pair	\$346
<input type="checkbox"/> 122DI64	Size 2 Ankle Guides - Pair	\$346
<input type="checkbox"/> 122DI65	Size 3 Ankle Guides - Pair	\$346
<input type="checkbox"/> 122DI66	Hoist Strap - Pair	\$100



### Technical Specifications - Imperial

Frame Sizing Specification	Size 1	Size 2	Size 3
Age (approximate)	1 - 5 years	4 - 11 years	9 - 16 years
User Shoulder Height (to floor)	23 - 34"	34 - 47"	42 - 59"
User Saddle Height (inside leg)	9.5 - 16"	16 - 24"	21.5 - 32"
Maximum User Weight	55 lbs.	110 lbs.	176 lbs.
Frame Weight	26.5 lbs.	14.5 lbs.	21 lbs.
Upper Support Unit Weight	Size 1 Frame and Upper Support Unit do not separate		23 lbs.
Frame Width	22.5"	26.5"	28.5"
Frame Length	28.5"	33.5"	43.5"
Saddle	9" L x 3.25" W		10" L x 3.25" W

  

Harness Sizing Specification	Size 1	Size 2	Size 3	Size 4	Size 5	Size 6
Age (approximate)	1 - 4 years	2 - 6 years	5 - 9 years	8 - 12 years	12 - 15 years	16 - 18 years
Waist at Navel Level	18 - 22"	19.5 - 23.5"	21.5 - 25.5"	23.5 - 27.5"	27.5 - 31.5"	31.5 - 35.5"
Color	Purple	Green	Blue	Red	Orange	Gray
Maximum User Weight	110 lbs.	110 lbs.	110 lbs.	176 lbs.	176 lbs.	176 lbs.

### Technical Specifications - Metric

Frame Sizing Specification	Size 1	Size 2	Size 3
Age (approximate)	1 - 5 years	4 - 11 years	9 - 16 years
User Shoulder Height (to floor)	635 - 920 mm	870 - 1190 mm	1070 - 1490 mm
User Saddle Height (inside leg)	220 - 440 mm	400 - 610 mm	550 - 810 mm
Maximum User Weight	25 kgs	50 kgs	80 kgs
Frame Weight	12 kgs	6.5 kgs	9.5 kgs
Upper Support Unit Weight	Size 1 Frame and Upper Support Unit do not separate		10.5 kgs
Frame Width	570 mm	670 mm	720 mm
Frame Length	720 mm	850 mm	1100 mm
Saddle	225 mm L x 80 mm W		255 mm L x 80 mm W

  

Harness Sizing Specification	Size 1	Size 2	Size 3	Size 4	Size 5	Size 6
Age (approximate)	1 - 4 years	2 - 6 years	5 - 9 years	8 - 12 years	12 - 15 years	16 - 18 years
Waist at Navel Level	460 - 560 mm	500 - 600 mm	550 - 650 mm	600 - 700 mm	700 - 800 mm	800 - 900 mm
Color	Purple	Green	Blue	Red	Orange	Gray
Maximum User Weight	50 kgs	50 kgs	50 kgs	80 kgs	80 kgs	80 kgs

#### Sunrise Medical (US) LLC

2842 Business Park Avenue  
 Fresno, CA 93727 USA  
 Tel.: 1-800-333-4000  
 Fax: 1-800-300-7502  
 E-Mail: webmaster@sunmed.com



If you are visually impaired, this document can be viewed in PDF format at [www.SunriseMedical.com](http://www.SunriseMedical.com)



MK-130103 Rev. B ©09/2023



800.333.4000

[www.SunriseMedical.com](http://www.SunriseMedical.com)

