



**Zippie X'CAPE**  
**March 2020**  
**EIFZ1**

|  |   |
|--|---|
| Date: _____<br>Dealer Acct #: _____<br>Dealer: _____<br>Dealer Contact: _____<br>Dealer Address: _____<br>Dealer City: _____ ST: _____ ZIP: _____<br>Dealer Phone: ( ) _____ Fax: ( ) _____<br>Confirm Via: <input type="checkbox"/> Fax <input type="checkbox"/> Email<br>Mark For: _____ | <b>SUBMITTING FOR:</b> <input type="checkbox"/> Quote <input type="checkbox"/> Order<br>PO#: _____<br>Ship To: _____<br>Attention: _____<br>Address: _____<br>Address: _____<br>Ship To City: _____ ST: _____ ZIP: _____<br>Ship To Phone: ( ) _____ Fax: ( ) _____ |
|--|---|

*The HCPCS CODES herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative HCPCS CODES, including misc. codes to ensure access for their enrollees. The use of HCPCS CODES does not ensure coverage or payment.*

|  |       |            |       |  |       |            |       |
|--|-------|------------|-------|--|-------|------------|-------|
| <input checked="" type="checkbox"/> Option # | HCPCS | Definition | Price | <input checked="" type="checkbox"/> Option # | HCPCS | Definition | Price |
|--|-------|------------|-------|--|-------|------------|-------|

**1. FRAME**

|   |                          |         |   |   |         |                          |         |       |  |         |  |                          |         |       |  |         |                          |         |       |   |         |
|---|--------------------------|---------|---|---|---------|--------------------------|---------|-------|--|---------|--|--------------------------|---------|-------|--|---------|--------------------------|---------|-------|---|---------|
| <p align="center"><b>A. Zippie X'CAPE - TRANSIT INCLUDED</b><br/>165 lbs. Weight Capacity</p> <p align="center"><b>A. FRAME STYLE</b> (Includes True Fit) (N/A w/ Xtender)</p> <table border="1" style="width:100%"> <tr> <td><input type="checkbox"/></td> <td>121WM87</td> <td>E1238</td> <td>X'CAPE - Transit Incl. without seating<br/>(Custom seating opts available)</td> <td>\$2,395</td> </tr> <tr> <td><input type="checkbox"/></td> <td>121WM84</td> <td>E1236</td> <td>X'CAPE - Transit Incl. w/ Seating<br/>(Sling and Back Uph.)</td> <td>\$2,555</td> </tr> </table> <p align="center"><i>Note: HCPCS Codes shown as submitted. Some sizes fit other HCPC codes i.e. 15"-18" widths = K0005</i></p> | <input type="checkbox"/> | 121WM87 | E1238   | X'CAPE - Transit Incl. without seating<br>(Custom seating opts available) | \$2,395 | <input type="checkbox"/> | 121WM84 | E1236 | X'CAPE - Transit Incl. w/ Seating<br>(Sling and Back Uph.) | \$2,555 | <p align="center"><b>B. Zippie X'CAPE - TRANSIT OPTIONAL</b><br/>165 lbs. Weight Capacity</p> <p align="center"><b>B. FRAME STYLE</b> (Includes True Fit)</p> <table border="1" style="width:100%"> <tr> <td><input type="checkbox"/></td> <td>121WM83</td> <td>E1238</td> <td>X'CAPE - Transit Opt. without seating<br/>(Custom seating opts available)</td> <td>\$2,155</td> </tr> <tr> <td><input type="checkbox"/></td> <td>121WM85</td> <td>E1236</td> <td>X'CAPE - Transit Opt. w/ Seating<br/>(Sling and Back Uph.)</td> <td>\$2,325</td> </tr> </table> <p align="center"><i>Note: HCPCS Codes shown as submitted. Some sizes fit other HCPC codes i.e. 15"-18" widths = K0005</i></p> | <input type="checkbox"/> | 121WM83 | E1238 | X'CAPE - Transit Opt. without seating<br>(Custom seating opts available) | \$2,155 | <input type="checkbox"/> | 121WM85 | E1236 | X'CAPE - Transit Opt. w/ Seating<br>(Sling and Back Uph.) | \$2,325 |
| <input type="checkbox"/>  | 121WM87                  | E1238   | X'CAPE - Transit Incl. without seating<br>(Custom seating opts available) | \$2,395   |         |                          |         |       |  |         |  |                          |         |       |  |         |                          |         |       |   |         |
| <input type="checkbox"/>  | 121WM84                  | E1236   | X'CAPE - Transit Incl. w/ Seating<br>(Sling and Back Uph.)                | \$2,555   |         |                          |         |       |  |         |  |                          |         |       |  |         |                          |         |       |   |         |
| <input type="checkbox"/>  | 121WM83                  | E1238   | X'CAPE - Transit Opt. without seating<br>(Custom seating opts available)  | \$2,155   |         |                          |         |       |  |         |  |                          |         |       |  |         |                          |         |       |   |         |
| <input type="checkbox"/>  | 121WM85                  | E1236   | X'CAPE - Transit Opt. w/ Seating<br>(Sling and Back Uph.)                 | \$2,325   |         |                          |         |       |  |         |  |                          |         |       |  |         |                          |         |       |   |         |

| TRANSIT SYSTEM WC-19     |         |       |                               |       |                          |         |              |    |
|--------------------------|---------|-------|-------------------------------|-------|--------------------------|---------|--------------|----|
| <input type="checkbox"/> | 121TO48 | K0108 | Transit Opt. (N/A w/ Xtender) | \$300 | <input type="checkbox"/> | 121TO40 | Omit Transit | NC |

| FRAME TYPE               |         |  |                         |    |                          |         |                              |    |
|--------------------------|---------|--|-------------------------|----|--------------------------|---------|------------------------------|----|
| <input type="checkbox"/> | 121FT64 |  | Zippie-Swing Away Frame | NC | <input type="checkbox"/> | 121FT65 | Zippie-80° Fixed Front Frame | NC |

**2. SIZE/SEATING**

| 2A. FRAME WIDTH                               |        |  |                                  |    |                          |        |                                  |    |
|---|--------|--|----------------------------------|----|--------------------------|--------|----------------------------------|----|
| <i>8", 9" and 10" widths N/A with Xtender</i> |        |  |                                  |    |                          |        |                                  |    |
| <input type="checkbox"/>                      | 121W48 |  | 8" (Frame Width Growth 9"-10")   | NC | <input type="checkbox"/> | 121W53 | 13" (Frame Width Growth 14"-15") | NC |
| <input type="checkbox"/>                      | 121W49 |  | 9" (Frame Width Growth 10"-11")  | NC | <input type="checkbox"/> | 121W54 | 14" (Frame Width Growth 15"-16") | NC |
| <input type="checkbox"/>                      | 121W50 |  | 10" (Frame Width Growth 11"-12") | NC | <input type="checkbox"/> | 121W55 | 15" (Frame Width Growth 16"-17") | NC |
| <input type="checkbox"/>                      | 121W51 |  | 11" (Frame Width Growth 12"-13") | NC | <input type="checkbox"/> | 121W56 | 16" (Frame Width Growth 17"-18") | NC |
| <input type="checkbox"/>                      | 121W52 |  | 12" (Frame Width Growth 13"-14") | NC |                          |        |                                  |    |

**2B/C. SLING DEPTH, FRAME DEPTH AND BACK POST POSITION. See matrix on page 2.**  
*8"-11" Frame Depth with 20", 22" or 24" wheels and Custom Seating must have hublock or Fixed Front frame. When selecting Frame Depth with custom seating, add 1" to accommodate seating. Example: need 12" depth, order 13" depth.*

| 2B. Sling Depth when choosing with Seating Sling must be up to 3" less than or equal to the Frame Depth (C.) |        |  |     |    |                          |         |     |    |
|--|--------|--|-----|----|--------------------------|---------|-----|----|
| <input type="checkbox"/>   | 121GJ1 |  | 8"  | NC | <input type="checkbox"/> | 121GJ7  | 14" | NC |
| <input type="checkbox"/>   | 121GJ2 |  | 9"  | NC | <input type="checkbox"/> | 121GJ8  | 15" | NC |
| <input type="checkbox"/>   | 121GJ3 |  | 10" | NC | <input type="checkbox"/> | 121GJ9  | 16" | NC |
| <input type="checkbox"/>   | 121GJ4 |  | 11" | NC | <input type="checkbox"/> | 121GJ10 | 17" | NC |
| <input type="checkbox"/>   | 121GJ5 |  | 12" | NC | <input type="checkbox"/> | 121GJ11 | 18" | NC |
| <input type="checkbox"/>   | 121GJ6 |  | 13" | NC |                          |         |     |    |

| 2C. Frame Depth and Back Post Position with Growth Note: Must be equal or greater than Sling Depth if chosen with seat Sling |         |  |   |    |                          |         |                                  |    |
|--|---------|--|---|----|--------------------------|---------|----------------------------------|----|
| <input type="checkbox"/>   | 121FD40 |  | 8" (Frame depth growth 9"-11",<br>N/A w/ 24" wheels)  | NC | <input type="checkbox"/> | 121FD46 | 14" (Frame depth growth 15"-17") | NC |
| <input type="checkbox"/>   | 121FD41 |  | 9" (Frame depth growth 10"-12",<br>N/A w/ 24" wheels) | NC | <input type="checkbox"/> | 121FD47 | 15" (Frame depth growth 16"-18") | NC |
| <input type="checkbox"/>   | 121FD42 |  | 10" (Frame depth growth 11"-13")                      | NC | <input type="checkbox"/> | 121FD48 | 16" (Frame depth growth 17"-19") | NC |
| <input type="checkbox"/>   | 121FD43 |  | 11" (Frame depth growth 12"-14")                      | NC | <input type="checkbox"/> | 121FD49 | 17" (Frame depth growth 18"-20") | NC |
| <input type="checkbox"/>   | 121FD44 |  | 12" (Frame depth growth 13"-15")                      | NC | <input type="checkbox"/> | 121FD50 | 18" (Frame depth growth 19"-20") | NC |
| <input type="checkbox"/>   | 121FD45 |  | 13" (Frame depth growth 14"-16")                      | NC |                          |         |                                  |    |

**SEATING**

| Seat Sling               |         |   |     |
|--------------------------|---------|---|-----|
| <input type="checkbox"/> | 121OP01 | Standard Seat Sling   | STD |
| <input type="checkbox"/> | 121OP04 | Omit Seat Sling (available when base is ordered with seating) | NC  |

| Cushions  |         |                       |       |
|---|---------|-----------------------|-------|
| Some limitations may exist on 2" Cushion and smaller widths |         |                       |       |
| <input type="checkbox"/>                                    | 121CU30 | K0669 1" Foam Cushion | \$120 |
| <input type="checkbox"/>                                    | 121CU31 | K0669 2" Foam Cushion | \$150 |

| JAY Cushions (Optional)   |         |  |            |
|---|---------|--|------------|
| To order any other JAY cushions on-chair, please see the JAY price list.<br>Cushion name or part number must be filled out w/ desired width and depth |         |  |            |
| <input type="checkbox"/>  | 121JC20 | JAY Bundled Cushion                        | Price List |
|   |         | Cushion Name or Part#                      |            |
|   |         | Width (in.):                  Depth (in.): |            |

| Sling and 1" Cushion Matrix |     |                                     |    |     |     |     |     |     |     |     |     |
|-----------------------------|-----|-------------------------------------|----|-----|-----|-----|-----|-----|-----|-----|-----|
|                             |     | Width                               |    |     |     |     |     |     |     |     |     |
|                             |     | 8"                                  | 9" | 10" | 11" | 12" | 13" | 14" | 15" | 16" |     |
| Sling Depth                 | 8"  |                                     |    |     |     |     |     |     |     |     | 8"  |
|                             | 9"  |                                     |    |     |     |     |     |     |     |     | 9"  |
|                             | 10" |                                     |    |     |     |     |     |     |     |     | 10" |
|                             | 11" |                                     |    |     |     |     |     |     |     |     | 11" |
|                             | 12" |                                     |    |     |     |     |     |     |     |     | 12" |
|                             | 13" |                                     |    |     |     |     |     |     |     |     | 13" |
|                             | 14" |                                     |    |     |     |     |     |     |     |     | 14" |
|                             | 15" |                                     |    |     |     |     |     |     |     |     | 15" |
|                             | 16" |                                     |    |     |     |     |     |     |     |     | 16" |
|                             | 17" |                                     |    |     |     |     |     |     |     |     | 17" |
| 18"                         |     |                                     |    |     |     |     |     |     |     | 18" |     |
|                             |     | Width                               |    |     |     |     |     |     |     |     |     |
|                             |     | 1" cushion available at these sizes |    |     |     |     |     |     |     |     |     |
|                             |     | 1" cushion N/A at these sizes       |    |     |     |     |     |     |     |     |     |

| Seat Sling to Frame Depth Matrix                 |     |              |    |     |     |     |     |     |     |     |     |     |
|--|-----|--------------|----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| Difference from front of frame to front of sling |     |              |    |     |     |     |     |     |     |     |     |     |
|  |     | Sling Length |    |     |     |     |     |     |     |     |     |     |
|  |     | 8"           | 9" | 10" | 11" | 12" | 13" | 14" | 15" | 16" | 17" | 18" |
| Frame Length                                     | 8"  | 0"           |    |     |     |     |     |     |     |     |     |     |
|  | 9"  | 0"           | 0" |     |     |     |     |     |     |     |     |     |
|  | 10" | 0"           | 0" | 0"  |     |     |     |     |     |     |     |     |
|  | 11" | 1"           | 1" | 1"  | 0"  |     |     |     |     |     |     |     |
|  | 12" |              | 2" | 2"  | 1"  | 0"  |     |     |     |     |     |     |
|  | 13" |              |    | 3"  | 2"  | 1"  | 0"  |     |     |     |     |     |
|  | 14" |              |    |     | 3"  | 2"  | 1"  | 0"  |     |     |     |     |
|  | 15" |              |    |     |     | 3"  | 2"  | 1"  | 0"  |     |     |     |
|  | 16" |              |    |     |     |     | 2"  | 1"  | 0"  | 0"  |     |     |
|  | 17" |              |    |     |     |     |     | 1"  | 0"  | 0"  | 0"  |     |
| 18"  |     |              |    |     |     |     |     | 0"  | 0"  | 0"  | 0"  |     |

| 2" Cushion Matrix |     |                                     |    |     |     |     |     |     |     |     |     |
|-------------------|-----|-------------------------------------|----|-----|-----|-----|-----|-----|-----|-----|-----|
|                   |     | Width                               |    |     |     |     |     |     |     |     |     |
|                   |     | 8"                                  | 9" | 10" | 11" | 12" | 13" | 14" | 15" | 16" |     |
| Sling Depth       | 8"  |                                     |    |     |     |     |     |     |     |     | 8"  |
|                   | 9"  |                                     |    |     |     |     |     |     |     |     | 9"  |
|                   | 10" |                                     |    |     |     |     |     |     |     |     | 10" |
|                   | 11" |                                     |    |     |     |     |     |     |     |     | 11" |
|                   | 12" |                                     |    |     |     |     |     |     |     |     | 12" |
|                   | 13" |                                     |    |     |     |     |     |     |     |     | 13" |
|                   | 14" |                                     |    |     |     |     |     |     |     |     | 14" |
|                   | 15" |                                     |    |     |     |     |     |     |     |     | 15" |
|                   | 16" |                                     |    |     |     |     |     |     |     |     | 16" |
|                   | 17" |                                     |    |     |     |     |     |     |     |     | 17" |
| 18"               |     |                                     |    |     |     |     |     |     |     | 18" |     |
|                   |     | Width                               |    |     |     |     |     |     |     |     |     |
|                   |     | 2" cushion available at these sizes |    |     |     |     |     |     |     |     |     |
|                   |     | 2" cushion N/A at these sizes       |    |     |     |     |     |     |     |     |     |

**3. CASTERS/SEAT HEIGHT**

| Caster Wheel             |         |                                |       |
|--------------------------|---------|--------------------------------|-------|
| <input type="checkbox"/> | 121CW30 | 3" Micro Light Up              | NC    |
| <input type="checkbox"/> | 121CW3  | 4" Low-Profile Polyurethane    | NC    |
| <input type="checkbox"/> | 121CW31 | 4" Micro Light Up              | \$70  |
| <input type="checkbox"/> | 121CW64 | 4" x 1-1/2" Aluminum Soft Roll | \$185 |
| <input type="checkbox"/> | 121CW5  | 5" Low-Profile Polyurethane    | NC    |
| <input type="checkbox"/> | 121CW51 | 5" Micro Light Up              | \$85  |
| <input type="checkbox"/> | 121CW6  | E2219 5" Soft Roll (Grey)      | \$90  |
| <input type="checkbox"/> | 121CW7  | E2219 5"x1-1/2" Semi-Pneumatic | \$85  |
| <input type="checkbox"/> | 121CW66 | 5"x1-1/2" Aluminum Soft Roll   | \$195 |
| <input type="checkbox"/> | 121CW34 | 6" Poly                        | NC    |
| <input type="checkbox"/> | 121CW35 | E2219 6"x1-1/2" Semi-Pneumatic | \$85  |

| Caster Forks             |         |                   |    |
|--------------------------|---------|-------------------|----|
| <input type="checkbox"/> | 121CF65 | Multiple Position | NC |

| Front Seat Heights Matrix  |                                |                                |                                |                                |
|--|--------------------------------|--------------------------------|--------------------------------|--------------------------------|
| Rear Wheel Size  | Front Caster                   |                                |                                |                                |
|  | <input type="checkbox"/> 3"    | <input type="checkbox"/> 4"    | <input type="checkbox"/> 5"    | <input type="checkbox"/> 6"    |
| <input type="checkbox"/> 12"<br><input type="checkbox"/> 16"<br><input type="checkbox"/> 20"<br><input type="checkbox"/> 22" | <input type="checkbox"/> 13.0" | <input type="checkbox"/> 13.0" | <input type="checkbox"/> 14.0" | <input type="checkbox"/> 15.0" |
|  | <input type="checkbox"/> 13.5" | <input type="checkbox"/> 13.5" | <input type="checkbox"/> 14.5" | <input type="checkbox"/> 15.5" |
|  | <input type="checkbox"/> 14.0" | <input type="checkbox"/> 14.0" | <input type="checkbox"/> 15.0" | <input type="checkbox"/> 16.0" |
|  | <input type="checkbox"/> 14.5" | <input type="checkbox"/> 14.5" | <input type="checkbox"/> 15.5" | <input type="checkbox"/> 16.5" |
|  | <input type="checkbox"/> 15.0" | <input type="checkbox"/> 15.0" | <input type="checkbox"/> 16.0" | <input type="checkbox"/> 17.0" |
|  | <input type="checkbox"/> 15.5" | <input type="checkbox"/> 15.5" | <input type="checkbox"/> 16.5" | <input type="checkbox"/> 17.5" |
|  | <input type="checkbox"/> 16.0" | <input type="checkbox"/> 16.0" | <input type="checkbox"/> 17.0" | <input type="checkbox"/> 18.0" |
|  | <input type="checkbox"/> 16.5" | <input type="checkbox"/> 16.5" | <input type="checkbox"/> 17.5" |                                |
|  | <input type="checkbox"/> 17.0" | <input type="checkbox"/> 17.0" |                                |                                |
|  | <input type="checkbox"/> 17.5" | <input type="checkbox"/> 17.5" | N/A                            | N/A                            |
| <input type="checkbox"/> 24"   | N/A                            | N/A                            | <input type="checkbox"/> 14.0" | <input type="checkbox"/> 15.0" |
|  |                                |                                | <input type="checkbox"/> 14.5" | <input type="checkbox"/> 15.5" |
|  | <input type="checkbox"/> 14.0" | <input type="checkbox"/> 14.0" | <input type="checkbox"/> 15.0" | <input type="checkbox"/> 16.0" |
|  | <input type="checkbox"/> 14.5" | <input type="checkbox"/> 14.5" | <input type="checkbox"/> 15.5" | <input type="checkbox"/> 16.5" |
|  | <input type="checkbox"/> 15.0" | <input type="checkbox"/> 15.0" | <input type="checkbox"/> 16.0" | <input type="checkbox"/> 17.0" |
|  | <input type="checkbox"/> 15.5" | <input type="checkbox"/> 15.5" | <input type="checkbox"/> 16.5" | <input type="checkbox"/> 17.5" |
|  | <input type="checkbox"/> 16.5" | <input type="checkbox"/> 16.5" | <input type="checkbox"/> 17.0" | <input type="checkbox"/> 18.0" |
|  | <input type="checkbox"/> 17.0" | <input type="checkbox"/> 17.0" | <input type="checkbox"/> 17.5" |                                |
|  | <input type="checkbox"/> 17.5" | <input type="checkbox"/> 17.5" |                                |                                |
|  | N/A                            | <input type="checkbox"/> 17.5" | N/A                            | N/A                            |

FRONT SEAT HEIGHTS ARE SHOWN IN THE BOXES IN THE FRONT SEAT HEIGHTS MATRIX. DEPENDING ON CG SETTING, SEAT DEPTH, AND FRAME SIZE SOME OF THE SIZES MAY NOT BE AVAILABLE

|  |  |
|--|--|
|  | Standard Frame   |
|  | Hemi Fame (N/A w/ 90° hanger, 80° hanger footplate in lowest position may cause interference w/caster) |
|  | Depends on Caster setting  |

**4.REAR WHEELS/SEAT HEIGHT**

**XTENDER POWER ASSIST OPTION**

22" & 24" wheels only. RSH 14" to 16" for 22", RSH 15.0" to 17" for 24". 0° camber, pneumatic tire. Anti-Tips no charge. Dual drive mode available through B4Me.

**Xtender Power Assist**

|                          |        |       |  |         |
|--------------------------|--------|-------|--|---------|
| <input type="checkbox"/> | 121XT6 | E0986 | 22" Xtender Power Assist JWX-2 NICKEL  | \$8,045 |
| <input type="checkbox"/> | 121XT7 | E0986 | 22" Xtender Power Assist JWX-2 LITHIUM | \$8,845 |
| <input type="checkbox"/> | 121XT4 | E0986 | 24" Xtender Power Assist JWX-2 NICKEL  | \$8,045 |
| <input type="checkbox"/> | 121XT5 | E0986 | 24" Xtender Power Assist JWX-2 LITHIUM | \$8,845 |

Please select one of the following handrim options

|                          |        |  |                 |    |
|--------------------------|--------|--|-----------------|----|
| <input type="checkbox"/> | 121XH5 |  | Plastic Coated  | NC |
| <input type="checkbox"/> | 121XH6 |  | Stainless Steel | NC |

**Xtender - Other Options**

|                          |         |  |                                   |         |
|--------------------------|---------|--|-----------------------------------|---------|
| <input type="checkbox"/> | 121XA10 |  | Spare NiMH Power Pack - NICKEL    | \$1,395 |
| <input type="checkbox"/> | 121XA11 |  | Spare Li-ion Power Pack - LITHIUM | \$1,795 |

**Xtender Wheel Options**

|                          |        |                |                                   |       |
|--------------------------|--------|----------------|-----------------------------------|-------|
| <input type="checkbox"/> | 121XA8 |                | SPARE 22" Lite Spoke Pneu Alum    | NC    |
| <input type="checkbox"/> | 121XA9 | E2211<br>E2213 | SPARE 22" Lite Spoke Pneu Plastic | \$200 |
| <input type="checkbox"/> | 121XA2 |                | SPARE 24" Lite Spoke Pneu Alum    | NC    |
| <input type="checkbox"/> | 121XA4 | E2211<br>E2213 | SPARE 24" Lite Spoke Pneu Plastic | \$200 |

**Rear Wheels**

N/A w/ Xtender & 5 Spoke Mag Wheels. N/A w/ Hub Lock. Std Mag Wheels only available with One Arm Drive. 24" wheels. N/A with 8" and 9" frame depths. 8"-12" Frame Depth with 20-24" Rear Wheels and ELR's is ONLY available with Hub Lock Wheels

|                          |          |  |  |         |
|--------------------------|----------|--|--|---------|
| <input type="checkbox"/> | 121RW48  |  | 12" Mag (Not available with handrims or one arm drive) | \$55    |
| <input type="checkbox"/> | 121RW14  |  | 16" Mag (Not available with handrims or one arm drive) | \$55    |
| <input type="checkbox"/> | 121RW3   |  | 20" Mag Wheel  | NC      |
| <input type="checkbox"/> | 121RW8   |  | 22" Mag Wheel  | NC      |
| <input type="checkbox"/> | 121RW1   |  | 24" Mag Wheel  | NC      |
| <input type="checkbox"/> | 121RW103 |  | 20" Mag Wheel 5 Spoke                                  | NC      |
| <input type="checkbox"/> | 121RW104 |  | 22" Mag Wheel 5 Spoke                                  | NC      |
| <input type="checkbox"/> | 121RW102 |  | 24" Mag Wheel 5 Spoke                                  | NC      |
| <input type="checkbox"/> | 121RW74  |  | 20" Lite Spoke Wheel                                   | \$200   |
| <input type="checkbox"/> | 121RW75  |  | 22" Lite Spoke Wheel                                   | \$200   |
| <input type="checkbox"/> | 121RW76  |  | 24" Lite Spoke Wheel                                   | \$200   |
| <input type="checkbox"/> | 121RW06  |  | 22" Ultra light wt. Spoke Wheel 2.0                    | \$250   |
| <input type="checkbox"/> | 121RW10  |  | 24" Ultra light wt. Spoke Wheel 2.0                    | \$250   |
| <input type="checkbox"/> | 121RW28  |  | 22" Spinergy SPOX Black Spokes                         | \$820   |
| <input type="checkbox"/> | 121RW95  |  | 24" Spinergy SPOX Black Spokes                         | \$820   |
| <input type="checkbox"/> | 121RW32  |  | 22" Spinergy LX Black Spokes                           | \$1,025 |
| <input type="checkbox"/> | 121RW98  |  | 24" Spinergy LX Black Spokes                           | \$1,025 |

**Rear Wheel Tire**

|                          |         |                |                                  |       |
|--------------------------|---------|----------------|----------------------------------|-------|
| <input type="checkbox"/> | 121RT3  |                | Full Poly Tire                   | NC    |
| <input type="checkbox"/> | 121RT6  | E2211<br>E2212 | High Pressure Clincher           | \$95  |
| <input type="checkbox"/> | 121RT27 | E2211<br>E2212 | Knobby V-TRAK by Kenda           | \$95  |
| <input type="checkbox"/> | 121RT38 | E2218          | SHOX Performance Poly            | \$95  |
| <input type="checkbox"/> | 121RT1  |                | Pneumatic                        | NC    |
| <input type="checkbox"/> | 121RT2  | E2211<br>E2213 | Pneumatic Airless Insert         | \$85  |
| <input type="checkbox"/> | 121RT21 | E2211<br>E2212 | Schwalbe Marathon Plus Evolution | \$140 |
| <input type="checkbox"/> | 121RT28 | E2211<br>E2212 | Schwalbe Speed                   | \$140 |
| <input type="checkbox"/> | 121RT4  |                | Low Poly                         | NC    |

| Tire                     | Wheel | 5 Spoke Mag | Std. Mag            | Lite Spoke  | Ultra Light-weight | Spinergy |
|--------------------------|-------|-------------|---------------------|-------------|--------------------|----------|
| Full Poly Tire           |       | 20",22",24" | 20",22",24"         | 20",22",24" |                    |          |
| High Pressure Clincher   |       |             |                     | 22",24"     | 24"                | 22",24"  |
| SHOX Performance Poly    |       | 20",22",24" |                     | 20",22",24" |                    | 22",24"  |
| Knobby V-TRAK            |       |             |                     | 24"         |                    | 24"      |
| Pneumatic                |       | 20",22",24" | 12",16",20",22",24" | 20",22",24" |                    | 22",24"  |
| Pneumatic Airless Insert |       | 20",22",24" | 12",16",20",22",24" | 20",22",24" |                    | 22",24"  |
| Schwalbe Marathon Plus   |       |             |                     | 24"         | 22",24"            | 22",24"  |
| Schwalbe Speed Run       |       |             |                     | 24"         | 24"                | 22",24"  |
| Low Poly                 |       |             | 12"                 | 24"         |                    |          |

22" = 489 mm 24" = 540 mm

**Handrim**

N/A w/ Xtender. Natural Fit N/A w/ 20" and 22" 5 Spoke Mag Wheel

|                          |         |  |  |       |
|--------------------------|---------|--|--|-------|
| <input type="checkbox"/> | 121HR1  |  | Aluminum Anodized  | NC    |
| <input type="checkbox"/> | 121HR4  |  | Plastic Coated (N/A on Ultra and 22" Spinergy LX,SPOX)               | \$110 |
| <input type="checkbox"/> | 121HR22 |  | Natural Fit - Standard Grip (N/A on Ultra Lite or 20" Lite Spoke)    | \$445 |
| <input type="checkbox"/> | 121HR23 |  | Natural Fit - Super Grip (N/A on Ultra Lite or 20" Lite Spoke)       | \$495 |
| <input type="checkbox"/> | 121HR30 |  | Natural Fit LT - Standard Grip (N/A on Ultra Lite or 20" Lite Spoke) | \$445 |
| <input type="checkbox"/> | 121HR31 |  | Natural Fit LT - Super Grip (N/A on Ultra Lite or 20" Lite Spoke)    | \$495 |
| <input type="checkbox"/> | 121HR24 |  | Spinergy Flex Rim (Only on 24" Spinergy LX)                          | \$700 |
| <input type="checkbox"/> | 121HR7  |  | Omit Handrim (Option may not be available on all rear wheels)        | NC    |

**Axle Plate**

|                          |         |  |          |    |
|--------------------------|---------|--|----------|----|
| <input type="checkbox"/> | 121AL33 |  | Standard | NC |
|--------------------------|---------|--|----------|----|

**4. REAR WHEELS/SEAT HEIGHT Cont.**

**One Arm Drive**

NA w/ Xtender. N/A w/ 12"-16" wheels. Only avail. On 12"-16" widths and 20", 22" and 24" mag wheels. N/A with Natural Fit handrims. N/A with SPHA arms

|                          |         |       |  |       |
|--------------------------|---------|-------|--|-------|
| <input type="checkbox"/> | 121AL16 | E0958 | One Arm Drive                                    | \$900 |
| <input type="checkbox"/> | 121OA2  |       | <b>Standard Handrim</b>                          | NC    |
| <input type="checkbox"/> | 121HM1  |       | <input type="checkbox"/> Right Hand Mount        | NC    |
| <input type="checkbox"/> | 121HM2  |       | <input type="checkbox"/> Left Hand Mount         | NC    |
| <input type="checkbox"/> | 121OA1  |       | <b>Small Handrim</b><br>(N/A w/ 20" rear wheels) | NC    |
| <input type="checkbox"/> | 121HM1  |       | <input type="checkbox"/> Right Hand Mount        | NC    |
| <input type="checkbox"/> | 121HM2  |       | <input type="checkbox"/> Left Hand Mount         | NC    |

**Axle**

|                          |         |  |   |       |
|--------------------------|---------|--|---|-------|
| <input type="checkbox"/> | 121AX1  |  | Quick Release (N/A w/ One Arm Drive)                            | NC    |
| <input type="checkbox"/> | 121AX11 |  | Quad Quick Release (Only Available with 12" and 16" Mag Wheels) | \$125 |

**Center of Gravity Preset** (dependent on wheels, casters & depth, From standard frame depth +2" forward and - depth growth up to -3")

|                          |         |  |                                |    |
|--------------------------|---------|--|--------------------------------|----|
| <input type="checkbox"/> | 121CC28 |  | -0.5"                          | NC |
| <input type="checkbox"/> | 121CC26 |  | 0 (Flush w/ backrest)          | NC |
| <input type="checkbox"/> | 121CC16 |  | 0.5" (Only avail 9"-18" deep)  | NC |
| <input type="checkbox"/> | 121CC20 |  | 1" (Only avail 10"-18" deep)   | NC |
| <input type="checkbox"/> | 121CC15 |  | 1.5" (Only avail 10"-18" deep) | NC |
| <input type="checkbox"/> | 121CC17 |  | 2" (Only avail 16"-18" deep)   | NC |
| <input type="checkbox"/> | 121CC18 |  | 2.5" (Only avail 16"-18" deep) | NC |

**Camber**

|                          |        |  |                                      |    |
|--------------------------|--------|--|--------------------------------------|----|
| <input type="checkbox"/> | 121CI1 |  | 0°                                   | NC |
| <input type="checkbox"/> | 121CI2 |  | 3° (N/A w/ Single Post Ht. Adj. Arm) | NC |

**Rear Seat Height (see Allowable Squeeze chart)**

|                          |         |       |    |
|--------------------------|---------|-------|----|
| <input type="checkbox"/> | 121RH32 | 11"   | NC |
| <input type="checkbox"/> | 121RH33 | 11.5" | NC |
| <input type="checkbox"/> | 121RH9  | 12"   | NC |
| <input type="checkbox"/> | 121RH10 | 12.5" | NC |
| <input type="checkbox"/> | 121RH11 | 13"   | NC |
| <input type="checkbox"/> | 121RH12 | 13.5" | NC |
| <input type="checkbox"/> | 121RH13 | 14"   | NC |
| <input type="checkbox"/> | 121RH14 | 14.5" | NC |
| <input type="checkbox"/> | 121RH15 | 15"   | NC |
| <input type="checkbox"/> | 121RH16 | 15.5" | NC |
| <input type="checkbox"/> | 121RH17 | 16"   | NC |
| <input type="checkbox"/> | 121RH18 | 16.5" | NC |
| <input type="checkbox"/> | 121RH19 | 17"   | NC |
| <input type="checkbox"/> | 121RH26 | 17.5" | NC |
| <input type="checkbox"/> | 121RH20 | 18"   | NC |

| Rear Wheel   | Rear Seat Height |       |     |       |     |       |     |       |     |       |     |       |     |       |     |  |
|--------------|------------------|-------|-----|-------|-----|-------|-----|-------|-----|-------|-----|-------|-----|-------|-----|--|
|              | 11"              | 11.5" | 12" | 12.5" | 13" | 13.5" | 14" | 14.5" | 15" | 15.5" | 16" | 16.5" | 17" | 17.5" | 18" |  |
| 12" low poly |                  |       |     |       |     |       |     |       |     |       |     |       |     |       |     |  |
| 12" pneu     |                  |       |     |       |     |       |     |       |     |       |     |       |     |       |     |  |
| 16"          |                  |       |     |       |     |       |     |       |     |       |     |       |     |       |     |  |
| 20"          |                  |       |     |       |     |       |     |       |     |       |     |       |     |       |     |  |
| 22"          |                  |       |     |       |     |       |     |       |     |       |     |       |     |       |     |  |
| 24"          |                  |       |     |       |     |       |     |       |     |       |     |       |     |       |     |  |
| 24" pneu     |                  |       |     |       |     |       |     |       |     |       |     |       |     |       |     |  |

When tires are not listed default to the wheel size. All Heights are subject to vary slightly

**Allowable Squeeze (may vary +/- 1/2")**

| Frame Depth | Max Squeeze |
|-------------|-------------|
| 8"-9"       | 1.0"        |
| 10"-12"     | 1.5"        |
| 13"-15"     | 2.0"        |
| 16"-18"     | 2.5"        |

**5. WHEEL LOCK**

**Hub Lock**

Only available with 12" or 16" mag wheels or Lite Spoke wheels.  
Not available with One Arm Drive.

|                          |         |  |                             |       |
|--------------------------|---------|--|-----------------------------|-------|
| <input type="checkbox"/> | 121LK01 |  | 12" Hub Lock                | \$100 |
| <input type="checkbox"/> | 121LK02 |  | 16" Hub Lock                | \$100 |
| <input type="checkbox"/> | 121LK04 |  | 20"-24" Hub Lock Lite Spoke | \$100 |

**Hub Lock Mount Position**

Foot lock may limit folding width.

|                          |         |  |                                      |    |
|--------------------------|---------|--|--------------------------------------|----|
| <input type="checkbox"/> | 121HH01 |  | Foot lock Mount                      | NC |
| <input type="checkbox"/> | 121HH02 |  | Seat Mount                           | NC |
| <input type="checkbox"/> | 121HH03 |  | Back Mount (back ht. may be limited) | NC |

**Wheel Lock**

High Mount will be mounted in low position on 12 and 16 rear wheels.  
Some may not be available depending on wheel size and position.

|                          |         |  |                 |       |
|--------------------------|---------|--|-----------------|-------|
| <input type="checkbox"/> | 121WL1  |  | High Mount Push | NC    |
| <input type="checkbox"/> | 121WL2  |  | High Mount Pull | NC    |
| <input type="checkbox"/> | 121WL30 |  | Compact         | \$140 |
| <input type="checkbox"/> | 121WL10 |  | Omit            | NC    |

**Extension Handle (Only avail w/ High Mount Wheel Locks)**

|                          |        |       |         |      |
|--------------------------|--------|-------|---------|------|
| <input type="checkbox"/> | 121EH1 | E0961 | 6" Pair | \$75 |
|--------------------------|--------|-------|---------|------|

**Swing Away 8"-12" Frame Depth  
Wheel Lock go/no-go matrix**

| Frame Depth | C.G.  | 20" wheel | 22" wheel | 24" wheel | Frame Depth | C.G.                          | 20" wheel | 22" wheel | 24" wheel |
|-------------|-------|-----------|-----------|-----------|-------------|-------------------------------|-----------|-----------|-----------|
| 8"          | -0.5" |           |           |           | 11          | -0.5"                         |           |           |           |
|             | 0"    |           |           |           |             | 0"                            |           |           |           |
|             | 0.5"  |           |           |           |             | 0.5"                          |           |           |           |
|             | 1"    |           |           |           |             | 1"                            |           |           |           |
| 9"          | -0.5" |           |           |           | 12          | 1.5"                          |           |           |           |
|             | 0"    |           |           |           |             | -0.5"                         |           |           |           |
|             | 0.5"  |           |           |           |             | 0"                            |           |           |           |
|             | 1"    |           |           |           |             | 0.5"                          |           |           |           |
| 10"         | -0.5" |           |           |           | 12          | 1"                            |           |           |           |
|             | 0"    |           |           |           |             | 1.5"                          |           |           |           |
|             | 0.5"  |           |           |           |             | No-go. User must use hublock. |           |           |           |
|             | 1"    |           |           |           |             | Go                            |           |           |           |
|             | 1.5"  |           |           |           |             |                               |           |           |           |

**6. HANGER AND FOOTREST**

**Z-finity Hangers/Footrest System**

*N/A on Fixed Frame.*

*Some interference may exist in certain configurations.*

|                          |         |   |       |
|--------------------------|---------|---|-------|
| <input type="checkbox"/> | 121H162 | Z-finity Hanger   | STD   |
| <input type="checkbox"/> | 121H160 | Elevating Hanger - Ext Mount<br><i>Includes ELR pad. N/A w/ Z-finity. Some limitations may exist.</i> | \$315 |
| <input type="checkbox"/> | 121H144 | Omit Hanger   | NC    |

**Z-finity Footrest Mount (Must Chose see Illustration)**

|                          |         |                                     |    |
|--------------------------|---------|-------------------------------------|----|
| <input type="checkbox"/> | 121DI22 | Direct Hanger Mt/Fixed Frame - High | NC |
| <input type="checkbox"/> | 121DI23 | Direct Hanger Mt/Fixed Frame - Low  | NC |
| <input type="checkbox"/> | 121DI24 | Extension Tube Mount                | NC |

**Z-finity Footrest Pivot Arm (Must Choose, use matrix for leg length)**

|                          |         |             |    |
|--------------------------|---------|-------------|----|
| <input type="checkbox"/> | 121DI20 | Short 2"    | NC |
| <input type="checkbox"/> | 121DI21 | Long 3-1/4" | NC |

**Footplate**

*Some interference may exist in certain configurations.*

|                          |         |   |       |
|--------------------------|---------|---|-------|
| <input type="checkbox"/> | 121F320 | K0040 Multi-Angle Adjustable - Pair<br><i>(N/A w/ 8"-9" widths)</i> | \$185 |
| <input type="checkbox"/> | 121F321 | K0040 Angle Adjustable Platform Footboard                           | \$190 |
| <input type="checkbox"/> | 121F6   | Omit Footplate  | NC    |

**Extension Tube Length (When Choosing Ext. Tube Mount)**

|                          |        |        |    |
|--------------------------|--------|--------|----|
| <input type="checkbox"/> | 121E62 | Hemi   | NC |
| <input type="checkbox"/> | 121E1  | Short  | NC |
| <input type="checkbox"/> | 121E2  | Medium | NC |
| <input type="checkbox"/> | 121E5  | Long   | NC |

**Hanger Mounted Footrest (Footplates Mount to Hanger Tubing)**

*N/A with 8" and 9" width. Footplate mounted in lowest position may cause interference w/ caster.*

|                          |         |                              |    |
|--------------------------|---------|------------------------------|----|
| <input type="checkbox"/> | 121H157 | 80° Swing In/Out Front Mount | NC |
| <input type="checkbox"/> | 121H158 | 90° Swing In/Out Front Mount | NC |

**Extensions (may interfere with some seat heights & casters)**

|                          |        |                  |      |
|--------------------------|--------|------------------|------|
| <input type="checkbox"/> | 121E23 | 2" (not mounted) | \$55 |
| <input type="checkbox"/> | 121E24 | 4" (not mounted) | \$55 |


**Footplate**

|                          |       |  |       |
|--------------------------|-------|--|-------|
| <input type="checkbox"/> | 121F9 | K0040 Kids Angle Adjustable (only available w/ 10"-16" frame widths) | \$185 |
| <input type="checkbox"/> | 121F8 | K0040 Adult Angle Adj. (only available with 10"-16" widths)          | \$185 |

**Swing In/Out Hanger Leg Lengths**

*(Seat sling to footrest measured)*

| 80° Hanger   | Min   | Max   |
|--------------|-------|-------|
| No Ext. Tube | 5.7"  | 9.7"  |
| 2" Ext. Tube | 11.2" | 12.2" |
| 4" Ext. Tube | 11.2" | 14.1" |
| 90° Hanger   | Min   | Max   |
| No Ext. Tube | 6.0"  | 10.0" |
| 2" Ext. Tube | 11.5" | 12.5" |
| 4" Ext. Tube | 11.5" | 14.5" |




**Footrest Options**

|                          |         |   |      |
|--------------------------|---------|---|------|
| <input type="checkbox"/> | 121HL8  | K0038 Leg Strap   | \$40 |
| <input type="checkbox"/> | 121HL11 | E0951 Heel Loops (N/A w/ Kids Angle Adj. Z-finity. N/A on widths below 14") | \$40 |
| <input type="checkbox"/> | 121IG2  | Neoprene Impact Guards (N/A with footrest in highest mounted position)      | \$90 |

**Z-finity Matrix**

**Direct Hanger Mount**

**Extension Tube Mount**



| Direct Hanger Mounted |                          |       |      |
|-----------------------|--------------------------|-------|------|
| Pivot Arm             | Hanger Mounting Position | Min   | Max  |
| Short                 | High Mount               | 0.5"  | 4.5" |
|                       | Low Mount                | 3.5"  | 8.0" |
| Long                  | High Mount               | -1.0" | 5.5" |
|                       | Low Mount                | 2.5"  | 9.0" |

| Extension Tube Mounted |                       |       |       |
|------------------------|-----------------------|-------|-------|
| Pivot Arm              | Extension Tube Length | Min   | Max   |
| Short                  | Hemi Short            | 5.8"  | 14.1" |
|                        | Short                 | 6.5"  | 15.1" |
|                        | Medium                | 8.4"  | 17.0" |
|                        | Long                  | 10.3" | 19.0" |
| Long                   | Hemi Short            | 4.6"  | 15.3" |
|                        | Short                 | 5.2"  | 16.3" |
|                        | Medium                | 7.2"  | 18.3" |
|                        | Long                  | 9.2"  | 20.2" |

**ELR Ext Mt.(MED Ext Tube)**

| Pivot Arm | Min  | Max   |
|-----------|------|-------|
| Short     | 7.0" | 18.0" |
| Long      | 6.0" | 19.0" |

## 7. BACKRESTS

### BACKREST

|                                  |   |    |                                   |       |   |         |
|----------------------------------|---|----|-----------------------------------|-------|---|---------|
| <input type="checkbox"/> 121BT10 | Standard Backrest ( <i>Angle adj./folding</i> ) | NC | <input type="checkbox"/> 121BT101 | E2398 | Dynamic Backrest ( <i>N/A w/ Fixed Back</i> ) | \$1,495 |
|----------------------------------|---|----|-----------------------------------|-------|---|---------|

### BACK TYPE (*Choose Height Adjustable or Fixed Height*)

|                                   |   |    |                                  |              |  |    |
|-----------------------------------|---|----|----------------------------------|--------------|--|----|
| <input type="checkbox"/> 121BT120 | Height Adjustable ( <i>see matrix below</i> ) | NC | <input type="checkbox"/> 121BT53 | Fixed Height |  | NC |
|-----------------------------------|---|----|----------------------------------|--------------|--|----|

#### Height Adjustable (*See Backrest Height Range matrix*)

*Must choose Lower Base Tube, Backrest Style, and Upper Back Tube.*

#### Lower Base Tube (*See Backrest Height Matrix Below*)

|                                  |   |    |
|----------------------------------|---|----|
| <input type="checkbox"/> 121BP58 | Short                                     | NC |
| <input type="checkbox"/> 121BP59 | Medium                                    | NC |
| <input type="checkbox"/> 121BP60 | Tall ( <i>N/A with Dynamic Backrest</i> ) | NC |

#### Backrest Styles (*See Images Below*)

|                                  |   |       |
|----------------------------------|---|-------|
| <input type="checkbox"/> 121PH16 | Straight Non Push Handle  | NC    |
| <input type="checkbox"/> 121PH6  | Integral Push handle  | NC    |
| <input type="checkbox"/> 121PH17 | Folding Push Handle with Ht. Adj. Back (8"-10" widths, stroller handle will face outward) | \$275 |
| <input type="checkbox"/> 121PH15 | Omit Back Tube ( <i>See omit matrix</i> )   | NC    |

#### Upper Back Tube (*See Backrest Height Matrix Below*)

|                                  |        |    |
|----------------------------------|--------|----|
| <input type="checkbox"/> 121B245 | Short  | NC |
| <input type="checkbox"/> 121B5   | Medium | NC |

#### Height Adjustable Backrest

Preset Back height (*Choose for Set Up*)  
*No Upholstery w/ 19"-24"*


|                          |     |                          |     |                          |     |
|--------------------------|-----|--------------------------|-----|--------------------------|-----|
| <input type="checkbox"/> | 9"  | <input type="checkbox"/> | 15" | <input type="checkbox"/> | 20" |
| <input type="checkbox"/> | 11" | <input type="checkbox"/> | 16" | <input type="checkbox"/> | 21" |
| <input type="checkbox"/> | 12" | <input type="checkbox"/> | 17" | <input type="checkbox"/> | 22" |
| <input type="checkbox"/> | 13" | <input type="checkbox"/> | 18" | <input type="checkbox"/> | 23" |
| <input type="checkbox"/> | 14" | <input type="checkbox"/> | 19" | <input type="checkbox"/> | 24" |

OR

#### Fixed Height

*Fixed Height back N/A with upholstery*


|                                  |                                     |       |
|----------------------------------|-------------------------------------|-------|
| <input type="checkbox"/> 121BP9  | Standard (16" Height)               | NC    |
| <input type="checkbox"/> 121BP10 | Medium (20" Height)                 | NC    |
| <input type="checkbox"/> 121BP11 | Tall (23" Height)                   | NC    |
| <input type="checkbox"/> 121BP61 | Standard w/ Folding PH (16" Height) | \$335 |
| <input type="checkbox"/> 121BP62 | Medium w/ Folding PH (20" Height)   | \$335 |
| <input type="checkbox"/> 121BP63 | Tall w/ Folding PH (23" Height)     | \$335 |



### Backrest Height Ranges

#### Upper Back Tube (*Choose One*)

|                                 |          |          |         | Straight non-push handle Short | Straight non-push handle Medium | Integral push handle Short | Integral push handle Medium | Folding Push Handle Short | Folding Push Handle Medium | Omit |
|---------------------------------|----------|----------|---------|--------------------------------|---------------------------------|----------------------------|-----------------------------|---------------------------|----------------------------|------|
|                                 |          | Standard | Dynamic |                                |                                 |                            |                             |                           |                            |      |
| Lower Base Tube<br>(Choose one) | Standard | Short    |         | 11"-15"                        | 15"-19"                         | 11"-15"                    | 15"-19"                     | 12"-16"                   | 16"-20"                    | 9"   |
|                                 |          | Medium   |         | 15"-19"                        | 18"-23"                         | 15"-19"                    | 18"-23"                     | 16"-20"                   | 19"-24"                    | 13"  |
|                                 |          | Tall     |         | 18"-22"                        | 18"-24"                         | 18"-22"                    | 18"-24"                     | 19"-23"                   | 19"-24"                    | 16"  |
|                                 | Dynamic  | Short    |         | 14"-18"                        | 18"-22"                         | 14"-18"                    | 18"-22"                     | 15"-19"                   | 19"-23"                    | 12"  |
|                                 |          | Medium   |         | 18"-22"                        | 21"-24"                         | 18"-22"                    | 21"-24"                     | 19"-23"                   | 22"-24"                    | 16"  |



**7. BACKRESTS Cont.**

**Back Upholstery**

*Uph N/A with Fixed Height back. Std Uph N/A w/Flip-Up, Adj. Locking Flip-Up and Length Adj Locking Flip-up. 3DX N/A on 8"-11" widths. Ballistic N/A on 8"-9" widths.*

|                          |         |       |                                  |       |
|--------------------------|---------|-------|----------------------------------|-------|
| <input type="checkbox"/> | 121BU1  |       | Standard Upholstery              | STD   |
| <input type="checkbox"/> | 121BU71 | E2611 | Backrest Tension Adj. 3DX Vented | \$330 |
| <input type="checkbox"/> | 121BU70 | E2611 | Backrest Tension Adj. Ballistic  | \$315 |
| <input type="checkbox"/> | 121BU66 |       | Omit Back Upholstery             | NC    |

**JAY Back (Optional)**

*To order any other JAY back on-chair, please See the JAY Price Lst. Requires Part #*

|                          |         |                         |               |
|--------------------------|---------|-------------------------|---------------|
| <input type="checkbox"/> | 121JB02 | JAY Bundled Backrest    | Price List    |
|                          |         | Backrest Name or Part # |               |
|                          |         | Width (in.):            | Height (in.): |

**Stabilizer Bar**

*N/A w/ Dynamic back and Fixed Back. Only available w/ Integral push handle.*

|                          |        |                |       |
|--------------------------|--------|----------------|-------|
| <input type="checkbox"/> | 121SZ1 | Stabilizer Bar | \$105 |
|--------------------------|--------|----------------|-------|

**8. ARMRESTS/SIDEGUARDS**

**Armrest and Armpad - Height Adjustable**

*Flip-Up, Adj. Locking Flip-Up and Length Adj. Locking Flip-up N/A with Std Back Upholstery. Low RSH in combination with large wheels may result in arm interference. SPHA arms N/A with One Arm Drive or Xtender.*

|                          |         |       |  |       |
|--------------------------|---------|-------|--|-------|
| <input type="checkbox"/> | 121AA58 | E0973 | T-Post Xtra Tall(with sideguard)   | \$265 |
| <input type="checkbox"/> | 121AA59 | E0973 | T-Post 2X Tall(with sideguard)   | \$265 |
| <input type="checkbox"/> | 121AA3  |       | Flip-Up  | NC    |
| <input type="checkbox"/> | 121AA1  | E0973 | Adjustable Locking Flip-Up   | \$265 |
| <input type="checkbox"/> | 121AA7  | E0973 | Length Adjust. Locking Flip-Up   | \$275 |
| <input type="checkbox"/> | 121AA11 | E0973 | Single-Post Ht Adj. Std.<br><i>(N/A w/ 3° camber or 8"-9" frame depths)</i>    | \$285 |
| <input type="checkbox"/> | 121AA29 | E0973 | Single-Post Height Adj. Low<br><i>(N/A w/ 3° camber or 8"-9" frame depths)</i> | \$285 |
| <input type="checkbox"/> | 121AA14 |       | Omit Armrest   | NC    |

**Side Guard**

*N/A with height adjustable armrests and T-Posts*

|                          |         |       |             |       |
|--------------------------|---------|-------|-------------|-------|
| <input type="checkbox"/> | 121SG33 | K0108 | Plastic Std | \$140 |
| <input type="checkbox"/> | 121SG29 | K0108 | Plastic Low | \$140 |


**Armpad**

*N/A with T-Post and Adj. Locking Flip-Up and Foam only available with Length Adj. Locking Flip-Up.*


|                          |         |  |  |    |
|--------------------------|---------|--|--|----|
| <input type="checkbox"/> | 121AM23 |  | Desk Classic                               | NC |
| <input type="checkbox"/> | 121AM24 |  | Full Classic <i>(N/A w/ 10"-13" Depth)</i> | NC |
| <input type="checkbox"/> | 121AM9  |  | Foam Short <i>(N/A w/ T Post)</i>          | NC |
| <input type="checkbox"/> | 121AM10 |  | Foam Long <i>(N/A w/ T Post)</i>           | NC |

**Single Post Height Adjustable Armrest**

| Pad     | Armrest Size | Transfer Tube Length (D) | Range (H)         | Pad Length (L) |      |
|---------|--------------|--------------------------|-------------------|----------------|------|
|         |              |                          |                   | Desk           | Full |
| Classic | STD          | 8"                       | 9-1/2" to 14-1/2" | 10-1/4"        | 14"  |
|         | LOW          | 5-1/4"                   | 8" to 12"         |                |      |



| Armrest Size | T-POST           |                |
|--------------|------------------|----------------|
|              | Range (H)        | Pad Length (L) |
| XTALL        | 6" to 7"         | 7"             |
| 2X TALL      | 7-1/2" to 8-1/2" |                |





**9. ACCESSORIES**

**Back Pack** (one Backpack Black/Charcoal already included)

|                          |       |  |                              |      |
|--------------------------|-------|--|------------------------------|------|
| <input type="checkbox"/> | 121A1 |  | Black Utility Backpack       | \$75 |
| <input type="checkbox"/> | 121A2 |  | Backpack Black/Charcoal Kids | \$60 |

**Miscellaneous**

|                          |          |       |  |       |
|--------------------------|----------|-------|--|-------|
| <input type="checkbox"/> | 1218A03  |       | Caddy (Only avail w/ Ext. Mt. hangers) | \$50  |
| <input type="checkbox"/> | 1213A01  |       | Tool Kit                               | \$55  |
| <input type="checkbox"/> | 121A2100 | E2207 | Crutch Holder                          | \$195 |

**Positioning Belt**

|                          |         |       |                                |       |
|--------------------------|---------|-------|--------------------------------|-------|
| <input type="checkbox"/> | 121A823 | E0978 | Belt 1" Dual Adj. Sqz Blk Pad  | \$85  |
| <input type="checkbox"/> | 121A824 | E0978 | Belt 2" Aircraft Buckle        | \$110 |
| <input type="checkbox"/> | 121A825 | E0978 | Belt 2" Aircraft Buckle Padded | \$120 |
| <input type="checkbox"/> | 121A826 | E0978 | Belt 1.5" Auto Buckle          | \$70  |
| <input type="checkbox"/> | 121A814 | E0978 | 1" Comfort Fit APB5116         | \$120 |

**Upper Extremity Support Tray**

*N/A with 8"-9" widths*

|                          |        |       |   |       |
|--------------------------|--------|-------|---|-------|
| <input type="checkbox"/> | 121A29 | E0950 | X-Small (Suggested 10"-12" chair width) | \$255 |
| <input type="checkbox"/> | 121A30 | E0950 | Small (Suggested 13"-14" chair width)   | \$255 |
| <input type="checkbox"/> | 121A31 | E0950 | Medium (Suggested 15"-17" chair width)  | \$255 |
| <input type="checkbox"/> | 121A32 | E0950 | Large (Suggested 18"-20" chair width)   | \$255 |

**Canopy (SunShade)**

*Need 2" back post clearance to mount canopy hardware*

|                          |         |       |                                 |       |
|--------------------------|---------|-------|---------------------------------|-------|
| <input type="checkbox"/> | 1218A06 | K0108 | Canopy (Sunshade) 10"-14" Width | \$275 |
| <input type="checkbox"/> | 1218A07 | K0108 | Canopy (Sunshade) 15"-20" Width | \$275 |

**Seat Pouch**

|                          |         |  |       |      |
|--------------------------|---------|--|-------|------|
| <input type="checkbox"/> | 1211A01 |  | Black | \$50 |
|--------------------------|---------|--|-------|------|

**Anti-Tip**

|                          |        |       |   |       |
|--------------------------|--------|-------|---|-------|
| <input type="checkbox"/> | 121AT1 | E0971 | Rear Anti-Tip                             | \$150 |
| <input type="checkbox"/> | 121AO1 | E0971 | Front Anti-Tip<br>(N/A with transit opt.) | \$150 |

**Spoke Guard**

*N/A with 12" and 16" rear wheels or 5 spoke mag or Xtender*

|                          |         |       |              |       |
|--------------------------|---------|-------|--------------|-------|
| <input type="checkbox"/> | 1214A01 | K0065 | Black        | \$115 |
| <input type="checkbox"/> | 1214A02 | K0065 | Clear Zippie | \$115 |

**Designer Spoke Guard**

*Only available with 22" or 24" Spoke rear wheels*

|                          |         |       |                 |       |
|--------------------------|---------|-------|-----------------|-------|
| <input type="checkbox"/> | 1214A10 | K0065 | Basketball      | \$200 |
| <input type="checkbox"/> | 1214A11 | K0065 | Chrome Wheels   | \$200 |
| <input type="checkbox"/> | 1214A12 | K0065 | Graffiti        | \$200 |
| <input type="checkbox"/> | 1214A13 | K0065 | Happy Lions     | \$200 |
| <input type="checkbox"/> | 1214A14 | K0065 | Hot Wheels      | \$200 |
| <input type="checkbox"/> | 1214A15 | K0065 | Mice n Cheese   | \$200 |
| <input type="checkbox"/> | 1214A16 | K0065 | Princess        | \$200 |
| <input type="checkbox"/> | 1214A17 | K0065 | Space Galactica | \$200 |
| <input type="checkbox"/> | 1214A18 | K0065 | Super Hero      | \$200 |
| <input type="checkbox"/> | 1214A19 | K0065 | Zebra           | \$200 |

**10. COLOR** (Standard with frame and cross brace color choices)

**Frame Color (Must choose)**

**Standard Colors**

|                          |         |  |                             |    |
|--------------------------|---------|--|-----------------------------|----|
| <input type="checkbox"/> | 121C140 |  | Aztec Gold                  | NC |
| <input type="checkbox"/> | 121C2   |  | Black (N/A as base w/Zebra) | NC |
| <input type="checkbox"/> | 121C141 |  | Black Cherry                | NC |
| <input type="checkbox"/> | 121C10  |  | Black Opal                  | NC |
| <input type="checkbox"/> | 121C38  |  | Blue Opal                   | NC |
| <input type="checkbox"/> | 121C36  |  | Candy Blue                  | NC |
| <input type="checkbox"/> | 121C14  |  | Candy Purple                | NC |
| <input type="checkbox"/> | 121C35  |  | Candy Red                   | NC |
| <input type="checkbox"/> | 121C142 |  | Electric Blue               | NC |
| <input type="checkbox"/> | 121C143 |  | Evergreen                   | NC |
| <input type="checkbox"/> | 121C126 |  | Glow                        | NC |
| <input type="checkbox"/> | 121C144 |  | Green Apple                 | NC |
| <input type="checkbox"/> | 121C147 |  | Hot Sparkle Pink            | NC |
| <input type="checkbox"/> | 121C148 |  | Mauve Pink                  | NC |
| <input type="checkbox"/> | 121C103 |  | Rootbeer                    | NC |
| <input type="checkbox"/> | 121C114 |  | Sparkle Silver              | NC |
| <input type="checkbox"/> | 121C145 |  | Sunrise Orange              | NC |
| <input type="checkbox"/> | 121C113 |  | Titanium Color Paint        | NC |
| <input type="checkbox"/> | 121C8   |  | Yellow                      | NC |

**Standard Matte Colors**

|                          |         |  |                     |    |
|--------------------------|---------|--|---------------------|----|
| <input type="checkbox"/> | 121C116 |  | Matte Black         | NC |
| <input type="checkbox"/> | 121C151 |  | Matte Black Cherry  | NC |
| <input type="checkbox"/> | 121C152 |  | Matte Purple        | NC |
| <input type="checkbox"/> | 121C150 |  | Matte Electric Blue | NC |
| <input type="checkbox"/> | 121C149 |  | Matte Evergreen     | NC |

**Custom Colors (Must choose cross brace std. color)**

|                          |         |  |  |       |
|--------------------------|---------|--|--|-------|
| <input type="checkbox"/> | 121C43  |  | Mossy Oak Camouflage                                     | \$395 |
| <input type="checkbox"/> | 121C42  |  | Stars and Stripes  | \$395 |
| <input type="checkbox"/> | 121C108 |  | Zebra (Choose base color above:<br>White Std. N/A Black) | \$395 |

**Cross Brace Color**

**Standard Colors**

|                          |         |  |                                    |    |
|--------------------------|---------|--|------------------------------------|----|
| <input type="checkbox"/> | 121C600 |  | Aztec Gold - Cross Brace           | NC |
| <input type="checkbox"/> | 121C601 |  | Black - Cross Brace                | NC |
| <input type="checkbox"/> | 121C602 |  | Black Cherry - Cross Brace         | NC |
| <input type="checkbox"/> | 121C603 |  | Black Opal - Cross Brace           | NC |
| <input type="checkbox"/> | 121C604 |  | Blue Opal - Cross Brace            | NC |
| <input type="checkbox"/> | 121C605 |  | Candy Blue - Cross Brace           | NC |
| <input type="checkbox"/> | 121C606 |  | Candy Purple - Cross Brace         | NC |
| <input type="checkbox"/> | 121C607 |  | Candy Red - Cross Brace            | NC |
| <input type="checkbox"/> | 121C608 |  | Electric Blue - Cross Brace        | NC |
| <input type="checkbox"/> | 121C609 |  | Evergreen - Cross Brace            | NC |
| <input type="checkbox"/> | 121C610 |  | Glow - Cross Brace                 | NC |
| <input type="checkbox"/> | 121C611 |  | Green Apple - Cross Brace          | NC |
| <input type="checkbox"/> | 121C612 |  | Hot Sparkle Pink - Cross Brace     | NC |
| <input type="checkbox"/> | 121C613 |  | Mauve Pink - Cross Brace           | NC |
| <input type="checkbox"/> | 121C614 |  | Rootbeer - Cross Brace             | NC |
| <input type="checkbox"/> | 121C615 |  | Sparkle Silver - Cross Brace       | NC |
| <input type="checkbox"/> | 121C616 |  | Sunrise Orange - Cross Brace       | NC |
| <input type="checkbox"/> | 121C617 |  | Titanium Color Paint - Cross Brace | NC |
| <input type="checkbox"/> | 121C618 |  | Yellow - Cross Brace               | NC |

**Standard Matte Colors**

|                          |         |  |                                   |    |
|--------------------------|---------|--|-----------------------------------|----|
| <input type="checkbox"/> | 121C619 |  | Matte Black - Cross Brace         | NC |
| <input type="checkbox"/> | 121C620 |  | Matte Black Cherry - Cross Brace  | NC |
| <input type="checkbox"/> | 121C621 |  | Matte Purple - Cross Brace        | NC |
| <input type="checkbox"/> | 121C622 |  | Matte Electric Blue - Cross Brace | NC |
| <input type="checkbox"/> | 121C623 |  | Matte Evergreen - Cross Brace     | NC |

**Additional Options Through Built-4 Me**

To inquire about any unique customization for Quickie products, contact: Phone: (800) 333 - 4000, Fax: (800) 333 - 9011

**Commonly Requested Options Include:**

- ✓ Solid Seat
- ✓ Grade Aid
- ✓ Drum Brake



Sunrise Medical (US) LLC  
2842 Business Park Ave. · Fresno, CA 93727 · USA

Customer Service: 800-333-4000 Fax: 800-300-7502 Email: orders@sunmed.com Web: www.sunrisemedical.com

© 3/2020 MK-100020 REV P