



Quickie 2 EZ Sku

July 2026

EIQ2X



The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

Customer				Ship To			
Account				Name			
Account #		Date		Address			
☐ Quote	☐ Order	PO #		City		State	
Buyer				ZIP		Phone	
Mark For				Email / Fax			
☐ Quote has Additional Sunrise Products with different Quote / Order #				Quote / Order #			

✓	Option #	HCPCS	Definition	Price
Base model				
<i>E1236 coding still applies for EIQ2X with sizes 14" or less. 300 lbs weight capacity</i>				
☐	002WM32	K0005	Quickie 2	\$2,467.00
☐	002E2A1		Quickie 2 EZ sku	STD

Frame Type				
☐	002RS1		Rear Side Frame Straight	NC
☐	002RS2		Rear Side Frame Rounded	NC

Heavy Duty Package				
<i>Includes Reinforced Back Upholstery & Stabilizer Bar. 350 lbs weight capacity</i>				
☐	002HV1		Heavy Duty (18"-24" Wide only)	\$583.00

Transit Option				
☐	002TO46	E1022	Transit System WC19	\$374.00

Frame Width						
☐	002W5	14"	NC	☐	002W9	18" NC
☐	002W6	15"	NC	☐	002W10	19" NC
☐	002W7	16"	NC	☐	002W11	20" \$583.00
☐	002W8	17"	NC		(E2201)	

Seat Depth						
☐	002SD7	14"	NC	☐	002SD11	18" NC
☐	002SD8	15"	NC	☐	002SD72	19" \$88.00
☐	002SD9	16"	NC	☐	002SD13	20" \$641.00
☐	002SD10	17"	NC		(E2203)	

Frame Depth						
<i>19"-20" Steel reinforced</i>						
☐	002FD6	15"	NC	☐	002FD9	18" NC
☐	002FD7	16"	NC	☐	002FD10	19" NC
☐	002FD8	17"	NC	☐	002FD36	20" NC

Seat Sling				
☐	002SL1		Standard Seat Sling	STD

JAY seat cushions (On-chair)				
<i>By opting in the below options the cushions will match the chair width/depth.</i>				
☐	002JC15	E2622	JAY X2	\$501.00
☐	002JC17	E2622	JAY Fusion Reduced Profile / Fluid	\$515.00
☐	002JC18	E2607	JAY Union Reduced Profile	\$520.00
☐	002JC10	E2607	JAY ION	\$462.00
☐	002JC06	E2622	JAY LitePro	\$525.00

Front Seat Height						
☐	002SH12	15.0"	NC	☐	002SH18	18.0" NC
☐	002SH13	15.5"	NC	☐	002SH20	18.5" NC
☐	002SH14	16.0"	NC	☐	002SH21	19.0" NC
☐	002SH15	16.5"	NC	☐	002SH23	19.5" NC
☐	002SH16	17.0"	NC	☐	002SH24	20.0" NC
☐	002SH17	17.5"	NC			

Caster Fork				
☐	002CF46		Multi Position Fork-Short	NC
☐	002CF47		Multi Position Fork- Tall	NC
☐	002CF12		5 1/4" Fork	NC
☐	002CF17		6 " Fork	NC

Caster Housing				
☐	002CP2		Standard	STD

✓	Option #	HCPCS	Definition	Price
NEW!	Caster Wheels			
☐	002CW3		4" Low Profile Polyurethane	NC
☐	002CW77		4" x 1.25" Semi Pneumatic	NC
☐	002CW64		4" x 1.5" Aluminum Soft Roll	\$216.00
☐	002CW5		5" Low Profile Polyurethane	NC
☐	002CW78		5" x 1.5" Semi-Pneumatic	NC
☐	002CW66		5" x 1.5" Aluminum Soft Roll	\$227.00
☐	002CW34		6" Polyurethane	NC
☐	002CW33	E2214 E2215	6" Pneumatic	\$88.00
☐	002CW79		6" x 1.5" Semi Pneumatic	NC
☐	002CW68		6" x 1.5" Aluminum Soft Roll	\$239.00
☐	002CW32		8" Polyurethane	NC

Rear Seat Height						
☐	002RH13	14.0"	NC	☐	002RH26	17.5" NC
☐	002RH14	14.5"	NC	☐	002RH20	18.0" NC
☐	002RH15	15.0"	NC	☐	002RH27	18.5" NC
☐	002RH16	15.5"	NC	☐	002RH21	19.0" NC
☐	002RH17	16.0"	NC	☐	002RH37	19.5" NC
☐	002RH18	16.5"	NC	☐	002RH22	20.0" NC
☐	002RH19	17.0"	NC			

Axle Plate				
☐	002AL102		Standard Adjustable (All in-one)	STD

Axle				
☐	002AX1		Quick Release Steel	STD

Rear Wheel Type				
☐	002RW103		20" Mag Wheel 5 Spoke	NC
☐	002RW104		22" Mag Wheel 5 Spoke	NC
☐	002RW102		24" Mag Wheel 5 Spoke	NC
☐	002RW92		24" Lite Spoke Wheel	NC
☐	002RW93		25" Lite Spoke Wheel	NC
☐	002RW95		24" Spinergy SPOX Black	\$985.00
☐	002RW96		25" Spinergy SPOX Black	\$985.00
☐	002RW98		24" Spinergy LX Black	\$1,237.00
☐	002RW99		25" Spinergy LX Black	\$1,237.00
☐	002RW50		24" Mountain	\$292.00

Rear Wheel Tire				
☐	002RT1		Pneumatic	NC
☐	002RT3		Full Poly	NC
☐	002RT24	E2211 E2213	Pneumatic Airless Insert - Infinergy® Lightweight - E-TPU material	\$114.00
☐	002RT38	E2218	SHOX Performance Poly	\$110.00
☐	002RT25	E2218	Primo Orion Poly	\$126.00
☐	002RT21	E2211 E2212	Schwalbe Marathon Plus Evolution	\$164.00
☐	002RT12	E2211 E2212	Mountain Tire	\$110.00

Handrim				
☐	002HR1		Aluminum Anodized	NC
☐	002HR4		Plastic Coated - N/A on every rim	\$129.00
☐	002HR22		Natural Fit - Standard Grip	\$549.00
☐	002HR39		Ellipse 3R	\$500.00
☐	002HR41		Carbolife - Curve L Tetra grip	\$575.00

Center of Gravity Preset						
Setting	☐ 1.0"	☐ 1.5"	☐ 2.0"	☐ 2.5"	☐ 3.0"	☐ 3.5"
Camber set a 2 deg default, 2" CoG default						


Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price							
Wheel Lock				Armrest (cont)											
☐	002WL46		High Mount Push		☐	002AA11	E0973 Single Post Hgt. Adj. Std (8.5"-13.5")	\$333.00							
☐	002WL50		High Mount Pull		☐	002AA29	E0973 Single Post Hgt. Adj. Low (7"-11")	\$333.00							
Extension Handle				☐					002AA65	E0973 Flip Back Height Adj (LITE 8.5"-13.5")	\$367.00				
☐	002EH7	E0961	Extension Handle (Pair)		☐	002AA14		Omit Armrest		NC					
COMING SOON Anti Roll-Back				Armpads											
	new	E2228	Automatic Bracking system - On chair		☐	002AM23		Desk Length - Classic (10")		NC					
PRATIKO Automatic Braking System. Installed by Sunrise Medical. Wheel Lock selection required - If no selection, Push to Lock will be provided. Weight capacity of 250 lbs (113 kg).				☐					002AM24		Full Length - Classic (14")	NC			
Hangers				Side Guard											
☐	002H145		60° Swing In/Out Hanger		☐	002SG33	K0108	Plastic Standard (8.5" High)		\$169.00					
☐	002H139		70° Swing In/Out Hanger		Anti-Tip										
☐	002H141		80° Swing In/Out Hanger		☐	002AT1	E0971	Rear Anti-Tip		\$175.00					
☐	002H160	E0990	Elevating Hanger - Ext Mount		☐	002AT3		Rear Anti-Tip OMIT		NC					
Footplates				Backpack & Pouch											
☐	002F1		Composite (13"-22"w)		☐	002A1		Quickie Backpack Black		\$88.00					
☐	002F271	K0040	Composite Angle Adj-(16"-22"w)		☐	002A2		Quickie Kids Backpack Black		\$70.00					
☐	002F3	K0040	Alum Angle Adj - Ext Mt (12"-22"w)		☐	0021A01		Seat Pouch Black /Charcoal		\$58.00					
Footplate Options				Positioning Belts											
☐	002HL1	E0951	Heel Loop		☐	002A141	E0978	Belt 1.5" Auto Buckle w/Clamps		\$82.00					
Footrest Extension Tubes				☐					002A117	E0978	Belt 2" Aircraft Buckle w/Clamps	\$129.00			
☐	002E1		Short		☐	002A118	E0978	Belt 2" Aircraft Buckle Pad w/Clamp		\$139.00					
☐	002E2		Medium		☐	002A120		Positioning Belt Slide Clamps		\$58.00					
☐	002E5		Long		Standard Frame Colors										
* Select lower leg length setting at least 2" shorter than FSH.				☐					002C214	Black (Gloss)	NC	☐	002C221	Red Pepper	NC
Back Type				☐					002C215	Black (Matte)	NC	☐	002C222	Royal Blue	NC
☐	002BT10		Standard		☐	002C211		Celestial Sapphir	NC	☐	002C223	Salute Green	NC		
☐	002BT20		Angle Adjustable Back (3-5" Pivot)		☐	002C212		Desert Sage	NC	☐	002C224	Scarlet Noir	NC		
NEW! Back Post Options				☐					002C213	Desert Tan	NC	☐	002C225	Silver Ore	NC
☐	002BP1		Push Handle Straight		☐	002C216		Metallic Shadow	NC	☐	002C226	Solar Orange	NC		
☐	002BP24		Push Handle 8°		☐	002C217		Midnight Blue	NC	☐	002C227	Starlight White	NC		
☐	002BP5		Angle Adj. Straight Push Handles		☐	002C218		Mystic Purple	NC	☐	002C228	Storm Gray	NC		
Back Height				☐					002C219	Petal Pink	NC	☐	002C229	Turbo Green	NC
☐	002B5		Medium (12"-15" Std Back height range)		☐	002C220		Pink Power	NC	☐	002C230	Velocity Yellow	NC		
☐	002B9		Tall (16"-19" Std Back height range)		NEW! Legacy - Frame Colors										
Back Height Setting				☐					002C141	Black Cherry	NC	☐	002C14	Candy Purple	NC
☐	002UX14	14"	NC	☐	002UX22	18"	NC	☐	002C10	Black Opal	NC	☐	002C142	Electric Blue	NC
☐	002UX16	15"	NC	☐	002UX24	19"	NC	Personalized Frame Color							
☐	002UX18	16"	NC	☐	002UX26	20"	NC	☐	002C112		Desert Camouflage		\$460.00		
☐	002UX20	17"	NC	☐	002UX28	21"	NC	☐	002C43		Mossy Oak Camouflage		\$460.00		
									☐	002C42		Stars & Stripes		\$460.00	
									☐	002C108		Zebra		\$460.00	
									Anodized Color Package						
									Standard Axle plate, Multi Forks and Aluminum Soft Roll casters only. N/A with Frog Legs.						
									☐	002EB07		Black - Caster only		\$58.00	
									☐	002EB08		Blue		\$187.00	
									☐	002EB13		Gold		\$187.00	
									☐	002EB09		Lime Green		\$187.00	
									☐	002EB10		Purple		\$187.00	
									☐	002EB11		Red		\$187.00	
									☐	002EB12		Ti Grey		\$187.00	
									ADDITIONAL INFORMATION						
									Quickie 2 product page.						
															
NEW! Armrest															
☐	002AA18	E0973	Padded Swing-Away STD (9" to 13")												
☐	002AA176	E0973	Padded Swing-Away with Side guard Standard (9" to 13")												
	Setting	☐ 9"	☐ 10"	☐ 11"	☐ 12"	☐ 13"	☐ 14"								

All 4" Casters						
Fork Size	Stem Size	Bolt Position				
		5	4	3	2	1
Multi Position Short	Std	-	-	16.5	17.0	17.5
	1/2"	-	-	-	17.5	18.0
	1"	-	-	-	18.0	18.5
Multi Position Tall	Std	-	-	18.5	19.0	19.5
	1/2"	-	-	-	19.5	20.0
	1"	-	-	-	20.0	20.5
3"	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
5 1/4"	Std	-	-	17.0	17.5	18.0
	1/2"	-	-	-	18.0	18.5
	1"	-	-	-	18.5	19.0
6"	Std	-	-	17.5	18.0	18.5
	1/2"	-	-	-	18.5	19.0
	1"	-	-	-	19.0	19.5
7"	Std	-	-	18.5	19.0	19.5
	1/2"	-	-	-	19.5	20.0
	1"	-	-	-	20.0	20.5
8" x 2"	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
Frog Legs	Std	-	-	16.5	17.0	17.5
Integral Multi Pos Forks (Low)	Std	-	-	-	16.0	16.5
Integral Multi Pos Forks (Std)	Std	-	-	-	16.5	17.0

All 5" Casters						
Fork Size	Stem Size	Bolt Position				
		5	4	3	2	1
Multi Position Short	Std	-	-	16.5	17.0	17.5
	1/2"	-	-	-	18.0	18.5
	1"	-	-	-	18.5	19.0
Multi Position Tall	Std	-	-	18.5	19.0	19.5
	1/2"	-	-	-	20.0	20.5
	1"	-	-	-	20.5	21.0
3"	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
5 1/4"	Std	-	-	17.5	18.0	18.5
	1/2"	-	-	-	18.5	19.0
	1"	-	-	-	19.0	19.5
6"	Std	-	-	18.0	18.5	19.0
	1/2"	-	-	-	19.0	19.5
	1"	-	-	-	19.5	20.0
7"	Std	-	-	19.0	19.5	20.0
	1/2"	-	-	-	20.0	20.5
	1"	-	-	-	20.5	21.0
8" x 2"	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
Frog Legs	Std	-	-	-	17.5	18.0
Integral Multi Pos Forks (Low)	Std	-	-	-	-	-
Integral Multi Pos Forks (Std)	Std	-	-	17.0	17.5	18.0


All 6" Casters						
Fork Size	Stem Size	Bolt Position				
		5	4	3	2	1
Multi Position Short	Std	-	-	17.5	18.0	18.5
	1/2"	-	-	-	18.5	19.0
	1"	-	-	-	19.0	19.5
Multi Position Tall	Std	18.5	19.0	19.5	20.0	20.5
	1/2"	-	-	-	20.5	21.0
	1"	-	-	-	21.0	21.5
3"	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
5 1/4"	Std	-	-	18.0	18.5	19.0
	1/2"	-	-	-	19.0	19.5
	1"	-	-	-	19.5	20.0
6"	Std	-	18.0	18.5	19.0	19.5
	1/2"	-	-	-	19.5	20.0
	1"	-	-	-	20.0	20.5
7"	Std	-	19.0	19.5	20.0	20.5
	1/2"	-	-	-	20.5	21.0
	1"	-	-	-	21.0	21.5
8" x 2"	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
Frog Legs	Std	-	-	-	-	18.5
Integral Multi Pos Forks (Low)	Std	-	-	-	-	-
Integral Multi Pos Forks (Std)	Std	-	-	-	17.5	18.0

All 8" x 1.5" Casters						
Fork Size	Stem Size	Bolt Position				
		5	4	3	2	1
Multi Position Short	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
Multi Position Tall	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
3"	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
5 1/4"	Std	-	-	-	19.0	19.5
	1/2"	-	-	-	20.0	20.5
	1"	-	-	-	20.5	21.0
6"	Std	-	-	19.5	20.0	20.5
	1/2"	-	-	-	20.5	21.0
	1"	-	-	-	21.0	21.5
7"	Std	-	20.0	20.5	21.0	21.5
	1/2"	-	-	-	21.5	22.0
	1"	-	-	-	22.0	-
8" x 2"	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
Frog Legs	Std	-	-	-	-	-
Integral Multi Pos Forks (Low)	Std	-	-	-	-	-
Integral Multi Pos Forks (Std)	Std	-	-	-	-	-



Rear Side Frame	Min RSH	Max RSH
Curved	13"	21.0"
Straight	13	22.0"
DAB	13	22.0"

Rear Tire	Min RSH	Max RSH
20" Pneumatic, Poly	13"	19.0"
22" Pneumatic, Poly	14.0"	20.0"
22" HP Clincher Blk	13.5"	19.5"
22" Schwalbe	13.5"	19.5"
24" Pneumatic, Poly	15"	21.0"
24" HP Clincher Blk	14.5"	20.5"
24" Knobby	15"	21.0"
24" Schwalbe	14.5"	20.5"
24" SHOX Poly	14.5"	20.5"
24" Mountain	15.0"	21.0"
25" HP Clincher Blk	15.0"	21.0"
25" Schwalbe	15.0"	21.0"
25" SHOX Poly	14.5"	20.5"
26" Pneumatic	16.0"	22.0"
26" HP Clincher Blk	15.5"	21.5"
26" Schwalbe	15.5"	21.5"



Handrims	Wheels				
	5 Spoke Mag	Lite Spoke	Ultra Lightweight	Spinergy	Mountain
Aluminum Anodized	20,22,24	20,22,24,26	22,24,25	22,24,25,26	24
Plastic Coat	20,22,24	22,24	-	24,26	-
Natural Fit - Standard Grip	24	22,24,25,26	-	22,24,25,26	-
Natural Fit - Super Grip	24	22,24,25,26	-	22,24,25,26	-
Q Grip Neoprene Coated	-	24,25	-	24,25	-
Ellipse 3R	24	24,25	-	24,25	-
Gekko	-	24,25	-	24,25	-
Curve L - Tetra Grip	-	24,25	-	24,25	-
Omit Handrim	24	24	24	24	24

		Seat Depth										
		12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"
Frame Depth	13"	S	S									
	14"	S	S	S								
	15"	S	S	S	S							
	16"		S	S	S	S						
	17"			S	S	S	S					
	18"				S	S	S	S				
	19"					S	S	S	S			
	20"						S	S	S	S		
	21"							S	S	S	S	
	22"								S	S	S	S

"S" = Standard Frame Style

"H" = Hemi Frame Style

*Only (6) Frame Depths avail. on Hemi Style Frame



Sunrise Medical (US) LLC
12002 Volunteer Blvd. Mount Juliet, TN 37122