





**FRONT FRAME AND FOOTREST**


Front Frame Options

XY080 Seat taper design

Specify width

mm


Optimise measurements (XY010 Seat Width -100mm)



Footrest Design (Choose One)


XY100 Round design

Castor Size	Max Height	With XB200
3"	115mm	165mm
4"	125mm	175mm
5"	135mm	185mm



XY105 Square design

Castor Size	Max Height	With XB200
3"	135mm	185mm
4"	145mm	195mm
5"	155mm	205mm



XB200 High level footrest (Only with XB010, XB080 & XB150)  
\$135.00



XB200 High level footrest (Only with XB010, XB080 & XB150)  
\$135.00



Footrest Individual Measurements

XY115 Front footrest height


mm

XY120 Rear footrest height

mm (Std = 20 lower than front height)

XY125 Footrest length

mm (Std = 100mm) (160mm and over requires heavy duty upgrade XB150)



Footrest Style and Cover

FXTXB080 Tubular footrest welded w/ Aluminium plate (Supplied with spacer for +20mm adjustment)

FXTXB110 Tubular footrest welded w/ Carbon Fibre plate (Supplied with spacer for +20mm adjustment)  
\$127.00

FXTXB120 Clamp-on aluminium plate footrest (Angle/height adjustable)

Footrest Options

XB165 Rubber grip tape (3M safety walk)  
\$35.00

## CASTOR WHEELS, FORKS & WHEEL LOCKS

### Castor Wheels 100mm (4")

- XC060 Soft roll 4"x1.3"  
*(Black only)*
- XC110 Aluminium centre 4"x1"  
\$122.00 *(Silver hub and black urethane tyre only)*
- XC130 Aluminium centre 4"x1.4"  
\$122.00 *(Silver hub and black urethane tyre only)*

- i**
- Silver hub & black tyre
  - Gold hub & black tyre
  - Purple hub & black tyre
  - Black hub & black tyre
  - Green hub & black tyre
  - Red hub & black tyre
  - Blue hub & black tyre

### 125mm (5")

- XC070 Soft roll 5"x1.3"  
*(Black only)*
- XC120 Aluminium centre 5"x1"  
\$122.00 *(Silver hub and black urethane tyre only)*
- XC100 Pneumatic Aluminium centre 5"x1"  
\$136.00 *(Grey only)*
- XC140 Aluminium centre 5"x1.4"  
\$122.00 *(Silver hub and black urethane tyre only)*

### 150mm (6")

- XC080 Soft roll 6"x1.3"  
*(Black hub / grey tyre only)*
- XC145 Aluminium centre 6"x1.4"  
\$122.00 *(Silver hub and black urethane tyre only)*

### Castor Forks

- XC192 Anti-Flutter castor fork  
*(Std when no option selected)*
- XC210 Froglegs suspension castor fork  
\$443.00

Note: Colour follows accent

- i**  Black  Silver

### Castor Pin Width

- XY330 Castor pin width  
*(XY010 plus 70mm)*



### Wheel Lock Options (Choose one)

- XD020 Underseat Scissor
- XD050 Push to lock
- XD030 Compact scissor  
\$96.00
- XD046 Out Front push to lock  
\$165.00
- XD035 Compact scissor with lightweight fixing  
\$113.00 *(Only with fixed centre of gravity)*
- XD085 Eagle Brakes (Indicate colour)  
\$208.00
- XD086 Eagle Brakes with L/W fixing (Indicate colour)  
\$225.00 *(Only with fixed centre of gravity)*

- XD010 Omit wheel locks

- BLACK
- PURPLE
- RED
- BLUE

Black Std if no option selected

**SIDEGUARDS AND ARMRESTS**

**Sideguard**

**Aluminium Options**

XE080 Lift out aluminium  
\$389.00

FXTXE130 Fold flat aluminium  
\$337.00

**Carbon Fibre Options**

XE100 Lift out carbon fibre  
\$527.00

FXTXE150 Fold flat carbon fibre  
\$546.00

**ABS Plastic Options**


XE120 Lift out ABS plastic  
\$252.00


FXTXE170 Fold flat ABS plastic  
\$337.00


---

**Sideguard Measurements**


**Optimised measurements**


25mm above tyre  
(Follow tyre) 

Individual above tyre  
(Follow tyre) 

Level with tyre  
(Follow tyre) 

**Individual measurements**

Sideguard length 

Sideguard height 

**Sideguard with Fender**

**Aluminium Options**

XE090 Lift out aluminium with fender  
\$422.00

FXTXE140 Fold flat aluminium with fender  
\$370.00


**Carbon Fibre Options**


XE110 Lift out carbon fibre with fender  
\$613.00

FXTXE160 Fold flat carbon fibre with fender  
\$665.00


**Sideguard Measurements**


**Fender height**

10mm above tyre  
(Standard) 

Individual above tyre  
(Min = 10mm) 

**Fender position**

Centre of tyre  
(Standard if no option selected) 


Outside of tyre 

**Armrest Options**

XF010 Lift out armrest with plastic sideguard  
\$434.00  
(Indicate fixed height)

XF020 Lift out armrest with Aluminium sideguard  
\$434.00  
(Indicate fixed height)

Height (Top of seatrail to top of armrest)



**BACKREST**

Backrest Height


XY230 Fixed height backrest  
(Maximum 400mm)

**i** Specify height  mm



XY245 Adjustable height backrest  
\$127.00 (Maximum height 400mm)  
(+/-25mm)

**i** Specify height  mm



Backrest Design

FXTXG060 Folding backrest with Q Lock  
(Black)

Ergonomic Options


XY360 Backrest angle

**i**  -5  -3  0  +3  +5  
Individual  Deg  
(Max -7/+7 Deg)




XG100 Ergonomic backrest  
\$77.00

**i** Vertical height  mm  
**i** Rearward position  mm




XG110 Inward taper  
\$77.00 (inside of backrest posts)

**i** Specify width or per side  mm



XG120 Outward flare  
\$77.00 (inside of backrest posts)

**i** Specify width or per side  mm



Push handle Options

XG070 Screw in push handles  
\$149.00 (N/A with XY245)

XG080 Flip down push handles  
\$229.00 (N/A with XY245)

XG085 Flip down push handles height adjustable  
\$252.00 (N/A with XY245)

XG097 Integrated push handle  
\$251.00 (quickie design)

**CAMBER BAR**

Adjustable Centre of Gravity

XH110 Square back adjustable centre of gravity (25-150mm)

(Max seat width 500mm)

Specify measurement

mm



XH110

XH120 Round back adjustable centre of gravity (60-150mm)

Specify measurement

mm




XH120

Fixed Centre of Gravity

FXTXH090 V-style fixed centre of gravity (25-125mm)  
\$318.00

Specify measurement

mm



FXTXH090

Specify Camber Angle

XY390 Camber angle

0  1  2  3  5

Specify Camber Angle

XY390 Camber angle

0  1  2  3  5

Camber Bar Options

FXTXK227 Standard Axle pin (120 or 130mm)

FXTXK226 100mm Titanium Axle pin \$173.00 (Fixed Axle only)

XK220 Titanium 130mm quick release axle pins (pair) \$109.00 (N/A with XK010, XK070 or FXTXK226)

XK240 Quad release axle nuts (pair) \$149.00 (N/A with FXTXK226)

Anti-Tip Options

XJ010 Swingaway anti tip LEFT side \$136.00 (Only with square back frame. Not with XH040)

XJ020 Swingaway anti tip RIGHT side \$136.00 (Only with square back frame. Not with XH040)

XJ100 Single sport anti tip \$422.00 (centre mounted)

REAR WHEELS	
Rear Wheel Size	
<input type="checkbox"/> XY420 24" Wheels	<input type="checkbox"/> XY430 25" Wheels
Rear Wheel Gap	
<div style="display: flex; align-items: center;"> <div style="margin-right: 20px;"> <input checked="" type="checkbox"/> XY460 Rear wheel gap   15 <input type="checkbox"/> 20 <input type="checkbox"/> 25  <input type="checkbox"/> Individual <input style="width: 100px; border: 1px solid black;" type="text"/> mm             </div> <div style="flex-grow: 1;"> </div> </div>	
Rear Wheel Options	
<input type="checkbox"/> XK010 Fusion 16 spoke <i>(Silver Hub, black spokes)</i>	<input type="checkbox"/> XK320 Spinergy XLX 12 spoke, cross/radial spoked <b>\$759.00</b>
<input type="checkbox"/> XK050 Spinergy LX 12 spoke <b>\$721.00</b>	<input type="checkbox"/> XK060 Spinergy SPOX 18 spoke <b>\$543.00</b>
<input type="checkbox"/> XK150 Spinergy CLX 18 carbon blade spoke <b>\$1,433.00</b> <i>N/A IN SILVER HUB</i>	<input type="checkbox"/> XK070 Spinergy BLXL 18 blade spoke <b>\$759.00</b> <i>N/A IN SILVER HUB</i>
<b>Spinergy Colour Options:</b>	
<input type="checkbox"/> Black spokes	<input type="checkbox"/> White spokes
<input type="checkbox"/> Red spokes	<input type="checkbox"/> Blue spokes
<input type="checkbox"/> Purple spokes	<input type="checkbox"/> Orange spokes
<input type="checkbox"/> Yellow spokes	<input type="checkbox"/> Pink spokes
<input type="checkbox"/> Green spokes	<input type="checkbox"/> XK180 Individual spoke design <b>\$169.00</b> <i>(Specify colour)</i>
<b>Spinergy Anodise Options</b>	
<input type="checkbox"/> XK190 Black hub <i>(std option)</i>	<input type="checkbox"/> XK200 Silver hub <b>\$88.00</b> <i>(n/a with BLXL/CLX)</i>
<input type="checkbox"/> XK202 Purple hub/rims <b>\$263.00</b> <i>(LX, XLX &amp; BLXL only)</i>	<input type="checkbox"/> XK201 Red hub/rims <b>\$263.00</b> <i>(LX, XLX &amp; BLXL only)</i>
<input type="checkbox"/> XK203 Blue hub/rims <b>\$263.00</b> <i>(LX, XLX &amp; BLXL only)</i>	<input type="checkbox"/> XK204 Gold hub/rims <b>\$263.00</b> <i>(LX, XLX &amp; BLXL only)</i>
<input type="checkbox"/> XK205 Individual hub/rims <b>\$304.00</b> Select colour: _____	
<b>Supply Without Wheels</b>	
<input type="checkbox"/> XK040 Supply without Wheels <b>-\$177.00</b>	Wheel Type: <input style="width: 150px;" type="text"/> Tyre Type: <input style="width: 150px;" type="text"/>
Handrim Position / Options	
<input type="checkbox"/> XY480 Narrow set handrim	<input type="checkbox"/> XY485 Wide set handrim
<input type="checkbox"/> XL090 Remove excess mounting tabs <b>\$97.00</b> <i>(recommended for small wheel gaps)</i>	
Handrim Styles	
<input type="checkbox"/> XL010 Aluminium hard anodised 6 tab <i>(Standard if no option selected)</i>	<input type="checkbox"/> XL040 Titanium 6 tab <b>\$231.00</b>
<input type="checkbox"/> XL160 Elipse 6 tab <b>\$324.00</b> <i>(replacement for Surge)</i>	<input type="checkbox"/> XL035 Natural Fit LT <b>\$290.00</b>
<input type="checkbox"/> XL020 Tetra grip 6 tab <b>\$272.00</b>	<input type="checkbox"/> XL150 Carbolife Gekko <input type="checkbox"/> Black <input type="checkbox"/> Silver <b>\$573.00</b>
<input type="checkbox"/> XL100 Para grip 6 tab <b>\$272.00</b>	<input type="checkbox"/> XL170 Carbolife Curve L (Black) <b>\$614.00</b>
Tyre Options	
<input type="checkbox"/> XM010 Kenda <input type="checkbox"/> Grey <input type="checkbox"/> Blue <input type="checkbox"/> Red <i>(Std)</i>	<input type="checkbox"/> XM060 Schwalbe One (24",25" only) <b>\$91.00</b>
<input type="checkbox"/> XM130 Kenda black <b>\$34.00</b>	<input type="checkbox"/> XM050 Kenda K885 mountain bike tyre (24",25" only) <b>\$26.00</b>
<input type="checkbox"/> XM020 Schwalbe Marathon Plus <b>\$66.00</b>	
<input type="checkbox"/> XM040 Shox solid airless <b>\$107.00</b> (Black)	
Inner Tube Valve Options	
<input type="checkbox"/> XY500 Schraeder inner tube valve <b>(Standard if no option selected)</b>	<input type="checkbox"/> XY490 Presta inner tube valve
Wheel Accessories	
<input type="checkbox"/> XK310 Offroad wheel set 24" rim with handrim, MTB tyre and axle pins <b>\$620.00</b> <i>(Replaces 25" wheel with XM010)</i>	<input type="checkbox"/> XK300 Offroad wheel set 22" rim with handrim, MTB tyre and axle pins <b>\$620.00</b> <i>(Replaces 24" wheel with XM010)</i>



**FINISH OPTIONS**

Frame & Component Finish

XN010 Powder coated frame

XN150 Powder coated components (Frame colour if not specified)

Frame Comp.

Gloss

Frame Comp.

Matt

Frame Comp.

Special

- Black gloss 103
- White gloss 100
- Grey gloss 132
- Cream gloss 135
- Red gloss 108
- Yellow gloss 120
- Forest green gloss 143
- Green gloss 107
- Dark blue gloss 151
- Royal blue gloss 154
- Blue gloss 174
- Duckegg blue gloss 141
- Pink gloss 188
- Orange gloss 180
- Purple gloss 157

- Black matt 105
- White matt 101
- Grey matt 131
- Light grey matt 133
- Navy blue matt 152

- Chrome Mirror 300
- Sparkle silver 170
- Metallic black matt 104
- Metallic grey matt 130
- Black sandpaper 201
- Grey sandpaper 202
- White sandpaper 200

XN200 Individual paint colour (State Colour/RAL code)  
\$165.00

Frame Comp

- FXTXN030 Beadblast frame  
\$229.00 (All parts black unless XN160 chosen)
- FXTXN050 Polished frame  
\$631.00 (All parts black unless XN170 chosen)

- XN160 Beadblast components  
\$126.00
- XN170 Polished components  
\$186.00

Accent Colour (Choose One)

**Note:** Accent colour includes forks only

- Black accent
- Silver accent
- Red accent
- Blue accent
- Purple accent

XN250 Individual Accent Colour  
\$165.00

\_\_\_\_\_

**UPHOLSTERY**

Seat Upholstery

XP040 Seat sling strap adjustable

XR050 Carbon fibre seat plate  
\$324.00

Backrest Upholstery

XP050 Airtech backrest strap adjustable

XP060 Backrest upgrade to full leather (Black)  
\$333.00

XP070 Backrest upgrade to leather corners (Black)  
\$115.00

RGK Logo Colour

- White (Std)
- Green
- Blue
- Individual: \_\_\_\_\_
- Red
- Orange
- Purple

Model Logo Colour

- White (Std)
- Green
- Blue
- Individual: \_\_\_\_\_
- Red
- Orange
- Purple

Backrest Accessories

XP110 Backrest with zip pocket  
\$85.00 (N/A with Sub4 backrest)

XP170 Individual backrest embroidery  
\$131.00 PLEASE COMPLETE PAGE 11

**CUSHION**

Cushion foam options

XP080 Cushion foam  
\$80.00 *(select density and size below)*

Soft (CMS040)

Medium (CMS050)

50mm

75mm

100mm

Cushion cover upgrades

XP085 Nylon cushion cover  
*(Black)*

XP090 Air Lite cushion cover  
\$62.00 *(Black)*

XP100 Leather cushion cover  
\$204.00 *(Black)*

Jay Products - PLEASE PROVIDE JAY ORDER FORM

JAY Cushion  
*(Price dependant on size / option)*

Jay Order form attached

JAY Backrest  
*(Price dependant on size / option)*

Jay Order form attached

JAY Whitmyer Headrest  
*(Price dependant on size / option)*

Jay Order form attached

**ACCESSORIES**

XR080 Frontwheel assembly  
\$700.00 *(Only with tubular style footrest)*

XP120 Frame protectors Neoprene zip fasten  
\$69.00 *(Black)*

XP420 Frame protectors Nylon velcro fasten  
\$46.00 *(Black)*

XR010 Seat pouch velcro attach *(Black)*  
\$36.00

XR200 Rucksack with RGK logo (Black)  
\$61.00

XR100 FX carry bag  
\$322.00

XR090 Tool kit

XS060 Calf strap

ADDITIONAL NOTES

INDIVIDUAL BACKREST EMBROIDERY

PLEASE PROVIDE JPEG IMAGE AT TIME OF ORDER AND NOTE POSITION BELOW



SPEC NOTES: