



Mygo Seat

September 2021

Mark For:

Date: _____

Dealer Acct #: _____

Dealer: _____

Dealer Contact: _____

Dealer Address: _____

Dealer City: _____ ST: _____ ZIP: _____

Dealer Phone: _____ Fax: _____

Confirmation Email: _____

Confirm Via: Fax Email

Submitting for: Quote Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____

Attention: _____

Address: _____

Address: _____

Ship To City: _____ ST: _____ ZIP: _____

Ship To Phone: _____ Fax: _____

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

Mygo Seating System

Size 1 Weight Capacity: 110lbs | Size 2 Weight Capacity: 132lbs

Seat Shell

Selection Required

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	117-610S	Size 1 Mygo Seat Shell		
	117-610-1	Seat Base with Depth Adjustment	E2294	\$806
	117-610-2	Seat Back - Fixed	E2293	\$806
	117-610-3	Seat Back Reclining Mechanism	E1014	\$424
	117-610-4	Seat Base Hip Guides / Support	E0956	\$118
	117-610-5	Seat Base Medial Thigh Guides	E0957	\$146
	117-610-6	Seat Base Leg Adductors	K0108	\$118
<input type="checkbox"/>	137-600S	Size 2 Mygo Seat Shell		
	137-600-1	Seat Base with Depth Adjustment	E2294	\$894
	137-600-2	Seat Back - Fixed	E2293	\$894
	137-600-3	Seat Back Reclining Mechanism	E1014	\$424
	137-600-4	Seat Base Hip Guides / Support	E0956	\$118
	137-600-5	Seat Base Medial Thigh Guides	E0957	\$146
	137-600-6	Seat Base Leg Adductors	K0108	\$118

Seat Cover

Selection Required. Must match Seat Shell size.

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	117-1851-08	Size 1 Mygo Seat Cover - Blue		N/C
<input type="checkbox"/>	117-1851-06	Size 1 Mygo Seat Cover - Gray		N/C
<input type="checkbox"/>	117-1851-07	Size 1 Mygo Seat Cover - Orange		N/C
<input type="checkbox"/>	117-1851-09	Size 1 Mygo Seat Cover - Pink		N/C
<input type="checkbox"/>	117-2851-08	Size 2 Mygo Seat Cover - Blue		N/C
<input type="checkbox"/>	117-2851-06	Size 2 Mygo Seat Cover - Gray		N/C
<input type="checkbox"/>	117-2851-07	Size 2 Mygo Seat Cover - Orange		N/C
<input type="checkbox"/>	117-2851-09	Size 2 Mygo Seat Cover - Pink		N/C

Base

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	102-2753	Size 2 Hi Low Chassis (Foot Pedal Operated. Includes Push Handle)	E1399	\$1,990
<input type="checkbox"/>	194-700	Mobility Base Interface Plate (Allows mounting Mygo Seat onto Sunrise Wheelchairs)	K0108	\$465

Pelvic Fixation

Selection of either Pelvic Harness or Pelvic Cradle required.

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	117-801	Pelvic Harness & Brackets Small - Black	E0978	\$144
<input type="checkbox"/>	117-802	Pelvic Harness & Brackets Medium - Black	E0978	\$144
<input type="checkbox"/>	120-789	Pelvic Harness Spacer Pads - Black (Pair)	K0108	\$82
<input type="checkbox"/>	133-728	Pelvic Cradle Small - Black	K0108	\$402
<input type="checkbox"/>	133-729	Pelvic Cradle Medium - Black	K0108	\$446
<input type="checkbox"/>	137-632	Hip Laterals (Pair) <i>Requires Pelvic Cradle</i>	E0956	\$298
<input type="checkbox"/>	117-791	Obliquity Kit	K0108	\$146

Dynamic Backrest Kit

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	137-660	Dynamic Backrest 50N Cylinder	K0108	\$403
<input type="checkbox"/>	137-661	Dynamic Backrest 100N Cylinder	K0108	\$411

Armrest

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	117-1777	Size 1 Mygo Adjustable Armrests (Pair)	E0973	\$294
<input type="checkbox"/>	117-2777	Size 2 Mygo Adjustable Armrests (Pair)	E0973	\$323

Lower Leg Positioning

Footplate selection required with Sandals

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	117-1776	Size 1 Mygo Angled Footplates (Pair)	K0040	\$294
<input type="checkbox"/>	117-2776	Size 2 Mygo Angled Footplates (Pair)	K0040	\$323
<input type="checkbox"/>	152-1600	Small Sandals w/ Straps (Pair)	E0951, E0952	\$241
<input type="checkbox"/>	152-2600	Medium Sandals w/ Straps (Pair)	E0951, E0952	\$263
<input type="checkbox"/>	152-3600	Large Sandals w/ Straps (Pair)	E0951, E0952	\$284
<input type="checkbox"/>	117-787	Sandal Raiser 1" (Pair) <i>Requires Sandal Selection</i>	Quantity: <input type="text"/> (1-3 Pairs)	K0108 \$85
<input type="checkbox"/>	117-765	Knee Pads (Pair)	K0108	\$146

Harness

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	120-799	Small Trunk Harness	E0960	\$261
<input type="checkbox"/>	120-757	Medium Trunk Harness	E0960	\$261
<input type="checkbox"/>	117-761	Large Trunk Harness	E0960	\$261

Lateral Support

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	117-837	Flip-Away Laterals with Cover (Pair)	E0956, E1028	\$512
<input type="checkbox"/>	137-634	Flip-Away Laterals Extended Height with Cover (Pair) <i>Required w/ Size 2 Seat</i>	E0956, E1028	\$629
<input type="checkbox"/>	117-803	Rigid Laterals with Cover (Pair)	E0956	\$293
<input type="checkbox"/>	137-633	Rigid Laterals Extended Height with Cover (Pair)	E0956	\$410
<input type="checkbox"/>	137-637	Flip-Away Complex Lateral Horizontal with Cover Left	E0956, E1028	\$328
<input type="checkbox"/>	137-638	Flip-Away Complex Lateral Horizontal with Cover Right	E0956, E1028	\$328
<input type="checkbox"/>	137-635	Flip-Away Complex Lateral Vertical with Cover Left	E0956, E1028	\$328
<input type="checkbox"/>	137-636	Flip-Away Complex Lateral Vertical with Cover Right	E0956, E1028	\$328
<input type="checkbox"/>	137-808-05	Chest Support - Black <i>(Requires Laterals. Not available with Complex Laterals)</i>	E0960	\$165

Headrest

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	120-859	Contoured Headrest Support Hardware <i>(Requires Cushion)</i>	E1028	\$340
<input type="checkbox"/>	137-800-05	Contoured Headrest Cushion - Black	E0955	\$156
<input type="checkbox"/>	117-697	Flat Headrest Support Hardware <i>(Requires Cushion)</i>	E1028	\$209
<input type="checkbox"/>	117-850-05	Flat Headrest Cushion - Black	E0955	\$106
<input type="checkbox"/>	120-893-05	Flat Headrest Laterals with Covers - Black (Pair) <i>Requires Flat Headrest. Not available with Contoured Headrest.</i>	E0956	\$298
<input type="checkbox"/>	117-745	Removable Head Support Bracket for 1/2" Box Section	E1028	\$65

Tray*Armrest selection required with Tray*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	117-724	Tray - Gray	E0950	\$352
<input type="checkbox"/>	117-270	Padded Tray Insert - Gray	K0108	\$306
<input type="checkbox"/>	117-769	Grab Rail / Toy Bar for Tray	K0108	\$270

Accessories

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	137-621	Shoulder Support Hardware for Mygo Size 1 <i>Shoulder Support Cushion required. Shoulder Support standard on Size 2.</i>	K0108	\$249
<input type="checkbox"/>	137-801-08	Shoulder Support Cushion - Blue	K0108	\$142
<input type="checkbox"/>	137-801-06	Shoulder Support Cushion - Gray	K0108	\$142
<input type="checkbox"/>	137-801-07	Shoulder Support Cushion - Orange	K0108	\$142
<input type="checkbox"/>	137-801-09	Shoulder Support Cushion - Pink	K0108	\$142

Assembly

	Item Code	Description	Retail
<input type="checkbox"/>	ASSY-LKY	Product Assembly Fee	\$120

Universal Specifications		
--------------------------	--	--

Specification	Size 1	Size 2
Age (approximate)	3 - 10 years	8 - 14 years
Backrest Angle: Prone	10°	10°
Backrest Angle: Recline	25°	25°
Footplate: Abduction	8°	8°
Footplate: Adduction	12°	20°
Footplate: Plantar/Dorsiflexion	10° / 10°	10° / 10°

Size 2 Hi Lo Chassis (Excluding Push Bar)		
Tilt in Space: Prone	10°	10°
Tilt in Space: Recline	25°	25°

Technical Specifications - Imperial		
-------------------------------------	--	--

Specification	Size 1	Size 2
User Height	41.3" - 59"	50" - 66"
User Weight	39.6 - 110 lbs	39.6 - 132 lbs
Seat Depth	10.6" - 16.5"	13.8" - 18.5"
Seat Width	8" - 12.8"	8.7" - 13.6"
Backrest Height	14.2" - 18.5"	18.1" - 22.4"
Seat to Sandal	8.5" - 13.8"	12.4" - 18.5"
Chest Width (Distance between laterals)	6.7" - 10.6"	6.7" - 10.6"
Armrest Height	6.3" - 8.3"	8.3" - 10.2"
Knee Width	3.5" - 4.3"	4.7" - 5.5"
Tray Size	21.6" x 18.9"	21.6" x 18.9"
Seat Unit Weight	22 lbs	31.9 lbs

Size 2 Hi Lo Chassis (Excluding Push Bar)		
Top of Seat to Floor	13.4" - 25.8"	13.4" - 25.8"
Base Weight	26.5 lbs	26.5 lbs

Technical Specifications - Metric		
-----------------------------------	--	--

Specification	Size 1	Size 2
User Height	1050 mm - 1500 mm	1270 mm - 1680 mm
User Weight	18 - 50 kg	18 - 60 kg
Seat Depth	270 mm - 420 mm	350 mm - 470 mm
Seat Width	200 mm - 325 mm	220 mm - 345 mm
Backrest Height	360 mm - 470 mm	460 mm - 570 mm
Seat to Sandal	215 mm - 350 mm	315 mm - 470 mm
Chest Width (Distance between laterals)	170 mm - 270 mm	170 mm - 270 mm
Armrest Height	160 mm - 210 mm	210 mm - 260 mm
Knee Width	90 mm - 110 mm	120 mm - 140 mm
Tray Size	550 mm x 480 mm	550 mm x 480 mm
Seat Unit Weight	10 kgs	14.5 kgs

Size 2 Hi Lo Chassis (Excluding Push Bar)		
Top of Seat to Floor	340 mm - 655 mm	340 mm - 655 mm
Base Weight	12 kgs	12 kgs



Sunrise Medical (US) LLC2842 Business Park Ave. • Fresno, CA 93727 • USA