



Quickie GP/GPV

August 2026

EIGP70

Mark For:

Date: _____

Dealer Acct #: _____

Dealer: _____

Dealer Contact: _____

Dealer Address: _____

Dealer City: _____ ST: _____ ZIP: _____

Dealer Phone: () _____ Fax: () _____

Confirmation Email: _____

Confirm Via: ☐ Fax ☐ Email

Submitting for:

☐ Quote

☐ Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____

Attention: _____

Address: _____

Address: _____

Ship To City: _____ ST: _____ ZIP: _____

Ship To Phone: () _____ Fax: () _____

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
---	----------	-------	------------	-------	---	----------	-------	------------	-------

1. FRAME

Base Model

For sizes 14" or less please use HCPCS E1235

☐	024WM23	K0005	Quickie GP-70°	\$2,234.00
☐	024WM23	K0005	Quickie GPV-70° (NA 12"-13" width)	\$2,234.00

Frame Angle

☐	024FA6	60° (N/A 18 ¼ and 19 ¼ frame length)	NC
☐	024FA7	70°	NC
☐	024FA1	75° (N/A with rigid footrest)	NC
☐	024FA8	80° (N/A with rigid footrest)	NC

2. SIZE / SEATING

Seat Width

Seat Width 1/2" Narrower

☐	024W3	12"	NC
☐	024W4	13"	NC
☐	024W5	14"	NC
☐	024W6	15"	NC
☐	024W7	16"	NC
☐	024W8	17"	NC
☐	024W9	18"	NC
☐	024W41	19"	\$147.00
☐	024W11	E2201 20"	\$583.00
☐	024W12	E2201 21"	\$583.00
☐	024W13	E2201 22"	\$583.00

Seat Depth

21 & 22" seat depths only available with 16" or wider chairs

☐	024SD5	12"	NC
☐	024SD6	13"	NC
☐	024SD7	14"	NC
☐	024SD8	15"	NC
☐	024SD9	16"	NC
☐	024SD10	17"	NC
☐	024SD11	18"	NC
☐	024SD72	19"	\$88.00
☐	024SD13	E2203 20"	\$641.00
☐	024SD14	E2203 21"	\$641.00
☐	024SD33	E2204 22"	\$1,051.00

Frame Length

☐	024FL30	18-1/4" (Avail. on 12" Sling Depth)	NC
☐	024FL31	19-1/4" (Avail. on 13" Sling Depth)	NC
☐	024FL32	20-1/4" (Avail. on 14", 15" Sling Depth)	NC
☐	024FL33	21-1/4" (Avail. on 14", 15", 16" Sling Depth)	NC
☐	024FL34	22-1/4" (Avail. on 16", 17" Sling Depth)	NC
☐	024FL35	23-1/4" (Avail. on 17", 18" Sling Depth)	NC
☐	024FL36	24-1/4" (Avail. on 18" Sling Depth)	NC
☐	024FL37	25-1/4" (Avail. on 19" Sling Depth)	NC
☐	024FL38	26-1/4" (Avail. on 20" Sling Depth)	NC
☐	024FL39	27-1/4" (Avail. on 21" Sling Depth)	NC
☐	024FL40	28-1/4" (Avail. on 22" Sling Depth)	NC

Frame Height

See Appendix 11 for height chart. **Short frame N/A w/Xtender**

☐	024FH2	Short (front 9.25"-rear 8.25") N/A with Xtender Power Assist	\$129.00
☐	024FH5	Standard (front 10" - rear 9")	NC
☐	024FH7	Tall (front 11" - rear 10")	\$129.00

Cushion

☐	024CU9	Reinforced Seat Sling	\$36.00
☐	024CU4	Omit Seat Sling	-\$15.00

JAY Bundled Cushions - Optional

To order any JAY cushions on-chair, please see JAY Price List. Cushion name or part number must be filled out with desired width & depth.

☐	024JC20	JAY Bundled Cushion	Price List
		Cushion Name or Part #:	
		Width: _____	Depth: _____

3. CASTERS

Caster Wheels				
☐	024CW1		3" Micro	NC
☐	024CW3		4" Low Profile Polyurethane	NC
☐	024CW31		4" Micro Lighted	\$82.00
☐	024CW6	E2219	5" Soft Roll Semi-Pneumatic	\$105.00
☐	024CW5		5" Low Profile Polyurethane	NC
☐	024CW7	E2219	5" x 1.5" Semi-Pneumatic <i>only avail. For use in bottom fork holes</i>	\$100.00
☐	024CW33	E2214 E2215	6" Pneumatic	\$88.00
☐	024CW34		6" Polyurethane	NC
☐	024CW35	E2219	6" x 1.5" Semi-Pneumatic <i>only avail. For use in bottom fork holes</i>	\$100.00
☐	024CW21	E2214 E2215	8" Pneumatic <i>Only available on 60° frame option</i>	\$100.00

Caster Wheels - Continued				
☐	024CW32		8" Polyurethane <i>Only available on 60° frame option</i>	NC
☐	024CW36	E2214 E2215	8" x 2" Pneumatic <i>Only available on 60° frame option</i>	\$169.00
☐	024CW37	E2214 K0108	8" x 2" Pneumatic w/Airless Insert <i>Only available on 60° frame option</i>	\$263.00

Caster Fork				
☐	024CS5		3/4" Longer Fork Stem Bolt	NC
☐	024CS9		1-1/2" Longer Fork Stem Bolt	NC
☐	024CS3	K0073	Caster Pin Locks	\$227.00

4. REAR WHEELS

Xtender Power Assist				
<i>RSH 15"-19" with 22" wheels, RSH 16"-20" with 24" wheels. Spare batteries must match set up. Anti-Tips included. Dual Drive mode available w/B4ME</i>				
☐	024XT6	E0986	22" Xtender Power Assist JWX-2 NICKEL	\$10,404.00
☐	024XT7	E0986	22" Xtender Power Assist JWX-2 LITHIUM	\$11,390.00
☐	024XT4	E0986	24" Xtender Power Assist JWX-2 NICKEL	\$10,404.00
☐	024XT5	E0986	24" Xtender Power Assist JWX-2 LITHIUM	\$11,390.00
<i>Please select one of the following handrim options</i>				
☐	024XH5		Plastic Coated	NC
☐	024XH6		Stainless Steel	NC

Xtender - Other Options				
☐	024XA10		Spare NiMH Power Pack - NICKEL	\$2,128.00
☐	024XA11		Spare Li-ion Power Pack - LITHIUM	\$2,594.00

Rear Wheel Size				
<i>Spinergy N/A w/ Xtender, spoke guards. Contact Built 4 Me for Spinergy Color Options</i>				
☐	024RW8		22" Mag	NC
☐	024RW9		22" Spoke	\$234.00
☐	024RW1		24" Mag	NC
☐	024RW2		24" Spoke	NC
☐	024RW11		26" Spoke <i>(N/A with Xtender Power Assist)</i>	\$234.00
☐	024RW42		24" Spinergy (540) SPOX Black	\$985.00
☐	024RR25		Omit Rear Wheel and Axle Intended Size for Set Up :	-\$55.00

Rear Wheel Tire				
☐	024RT1		Pneumatic	NC
☐	024RT2	E2211 E2213	Pneumatic Airless Insert <i>N/A with 26" wheels & Xtender</i>	\$100.00
☐	024RT3		Full Poly <i>N/A with 26" wheels</i>	NC
☐	024RT6	E2211 E2212	High-Pressure Clincher	\$110.00

Hand Rim				
☐	024HR1		Aluminum Anodized	NC
☐	024HR4		Plastic Coated <i>-N/A on all rims</i>	\$129.00
☐	024HR7		Omit Hand Rim	NC

Projection Hand Rim				
<i>Projections N/A w/ plastic coated handrims. N/A w/ Spoke guards</i>				
☐	024HR8		8 Oblique (24"/26")	\$379.00

Axle				
☐	024AX1		Quick Release Steel	NC
☐	024AX11		Quad Quick Release	\$147.00

Axle Plate				
☐	024AL2		Standard Adjustable	NC
☐	024AL13		Double Slotted Adjustable	\$70.00
☐	024AL6	E0959	Amputee Adjustable <i>N/A with padded swing-away armrests or stroller handles</i>	\$205.00

5. WHEEL LOCKS

Wheel Lock				
☐	024WL1		High Mount Push	NC
☐	024WL2		High Mount Pull	NC
☐	024WL5		High Mount Scissor <i>N/A w/ extension handles, grade aides</i>	NC
☐	024WL12		Do Not Mount Wheel Locks	NC
☐	024WL10		Omit	NC

Extension Handle				
☐	024EH1	E0961	6" Pair	\$88.00

Grade Aid				
☐	024GA1	E0974	Grade Aid <i>N/A w/ poly, high-pressure clincher tires, or Xtender</i>	\$175.00

6. FOOTREST

Footplates

Heel loop or leg strap not included, See Footplate Options Section

1. GP FOOTREST OPTIONS

☐	024F56	K0040	Rigid N/A w/ Ext. Tubes; N/A with 75° or 80° frame angle	NC
☐	024F50	K0040	Platform N/A w/ Ext. Tubes	\$110.00
☐	024F1	K0040	Composite Footplate	NC
☐	024F3	K0040	Adult Angle Adj.	\$216.00
☐	024F48	K0040	Platform Flip Up	\$221.00

2. GPV FOOTREST OPTIONS

☐	024F56	K0040	Rigid N/A w/ Ext. Tubes; N/A with 75° or 80° frame angle	NC
☐	024F50	K0040	Platform N/A w/ Ext. Tubes	\$110.00
☐	024F48	K0040	Platform Flip Up	\$221.00

Footrest Extension Tubes

☐	024E1		Short (14" to 16-1/2")	NC
☐	024E2		Medium (16-1/2" to 19")	NC
☐	024E5		Long (19" to 21-1/2")	NC
☐	024E16		High Mount Only avail. with adult angle adjustable footplates	\$63.00

Footplate Options (Ext Tube & Direct Mount)

☐	024HL1	E0951	Heel Loop	\$47.00
☐	024HL5	K0038	Leg Strap	\$47.00
☐	024HL4		Omit Heel Loop	NC
☐	024HL6		Omit Leg Strap	NC
☐	024F7		Omit Footplate & Ext Tube	-\$15.00
☐	024IG2		Impact Guard Neoprene	\$105.00

7. BACKREST

Back Type

☐	024BT3		Non-Folding Backrest (Fixed) N/A with push handles or adj. upholstery	
☐	024BB1		8" Back Height	NC
☐	024BB2		9" Back Height	NC
☐	024BB3		10" Back Height	NC
☐	024BB4		11" Back Height	NC
☐	024BB5		12" Back Height	NC
☐	024BB6		13" Back Height	NC
☐	024BB7		14" Back Height	NC
☐	024BB8		15" Back Height	NC
☐	024BB9		16" Back Height	NC
☐	024BT1		Folding Backrest (Adjustable)	
☐	024BA1		11" - 14" Back Height	\$75.00
☐	024BA2		12" - 16" Back Height	\$75.00
☐	024BA3		14" - 18" Back Height	\$75.00
☐	024BA4		16" - 20" Back Height	\$75.00

Back Upholstery

☐	024BU1		Standard	NC
☐	024BU8		Adjustable	\$205.00
☐	024BU16		Reinforced (N/A w/ adj. upholstery)	\$36.00
☐	024BU14		Omit Back Upholstery	-\$15.00

Backrest Options

☐	024BT33		Folding w/ Integral Push Handles	NC
☐	024A85		Bolt-On Push Handle (Folding only)	\$134.00
☐	024SR2		Regular Stroller Handle N/A with amputee axle plate	\$158.00
☐	024SR3		Tall Stroller Handle N/A with amputee axle plate	\$158.00

JAY Back (Optional)

☐	024JB01		J3 Back (See J3 Order Form)	Price List
--------------------------	---------	--	-----------------------------	------------

8. ARMREST / SIDEGUARDS

Single Post

☐	024AA9	E0973	Single Post Ht. Adjust. Low Must be removed to fold back down	\$333.00
☐	024AM1		Desk Length - Classic (10")	NC
☐	024AM2		Full Length - Classic (14")	NC
☐	024AA8	E0973	Single Post Ht. Adjust. Std Must be removed to fold back down	\$333.00
☐	024AM1		Desk Length - Classic (10")	NC
☐	024AM2		Full Length - Classic (14")	NC

Padded Swing Away Armrest

All Padded Swing-Away N/A with amputee axle plate or Xtender Power Assist

☐	024AA19	E0973	Padded S-A Medium	\$263.00
☐	024AA16	E0973	Padded S-A Adult	\$263.00

Side Guard

☐	024SG3	K0108	Plastic Regular	\$169.00
☐	024SG4	K0108	Plastic Kids	\$169.00

9. KEY OPTIONS AND ACCESSORIES

Rear Anti-Tip

☐	024AT1	E0971	Rear Anti-Tip	\$175.00
--------------------------	--------	-------	---------------	----------

Backpack

☐	024A1		Quickie Backpack Black	\$88.00
☐	024A2		Quickie Kids Backpack Black	\$70.00

Seat Pouch

☐	024A8		Quickie Seat Pouch Black	\$58.00
--------------------------	-------	--	--------------------------	---------

Positioning Belts

☐	024A142	E0978	Belt 1.5" Auto Buckle w/Clamps	\$82.00
☐	024A12	E0978	Belt 2" Aircraft Buckle w/Clamps	\$129.00
☐	024A13	E0978	Belt 2" Aircraft Buckle Pad w/Clamp	\$139.00

Spoke Guards

N/A w/ Xtender Power Assist; Only available w/ Mag and Spoke wheels

☐	024A23	K0065	Black Spoke Guards	\$134.00
☐	024A24	K0065	Clear Spoke Guards	\$134.00

Crutch Holder

☐	024A53	E2207	Left Side	\$227.00
☐	024A52	E2207	Right Side	\$227.00

Extra Leg Strap

☐	024A55	K0038	Extra Leg Strap Black	\$36.00
--------------------------	--------	-------	-----------------------	---------

Miscellaneous

☐	024A83		Tool Kit	\$63.00
☐	024A999		Quickie Water Bottle & Holder	\$63.00

Tray Table

☐	024A31	E0950	Tray Table	\$298.00
--------------------------	--------	-------	------------	----------

10. COLOR

NEW! Standard Frame Colors				
☐	024C214		Black (Gloss)	NC
☐	024C215		Black (Matte)	NC
☐	024C211		Celestial Sapphire	NC
☐	024C212		Desert Sage	NC
☐	024C213		Desert Tan	NC
☐	024C216		Metallic Shadow	NC
☐	024C217		Midnight Blue	NC
☐	024C218		Mystic Purple	NC
☐	024C219		Petal Pink	NC
☐	024C220		Pink Power	NC
☐	024C221		Red Pepper	NC
☐	024C222		Royal Blue	NC
☐	024C223		Salute Green	NC
☐	024C224		Scarlet Noir	NC
☐	024C225		Silver Ore	NC
☐	024C226		Solar Orange	NC
☐	024C227		Starlight White	NC
☐	024C228		Storm Gray	NC
☐	024C229		Turbo Green	NC
☐	024C230		Velocity Yellow	NC

Personalized Frame Color				
☐	024C112		Desert Camouflage	\$460.00
☐	024C43		Mossy Oak Camouflage	\$460.00
☐	024C42		Stars & Stripes	\$460.00
☐	024C108		Zebra (Choose above base color) Std with White unless above base paint is chosen N/A w/ Black, Matte Black)	\$460.00

11. APPENDIX

GP Series Rear STF Chart

*All measurements given are an estimate and have been rounded to the nearest .5"

SHORT Frame (9 1/4" Front / 8 1/4" Rear)			Rear Axle Plate Hole Position					
			Standard			Inverted		
			Top	Middle	Bottom	Top	Middle	Bottom
Wheel	Rim	Tire						
22"	Mag / Spoke	Full Poly	14.5	15	15.5	16	16.5	17
22"	Mag / Spoke	Pneum / Insert	14.5	15	16	16	17	17.5
24"	Mag / Spoke	Full Poly	15	15.5	16.5	16.5	17	17.5
24"	Spoke / Spinerger	Clincher	15	16	16.5	16.5	17	18
24"	Mag / Spoke / Spinerger	Pneum / Insert	15.5	16	17	17	17.5	18
26"	Spoke / Spinerger	Clincher	16	17	17.5	17.5	18	19
26"	Spoke / Spinerger	Pneumatic	16.5	17	17.5	18	18.5	19

TALL Frame (11" Front / 10" Rear)			Rear Axle Plate Hole Position					
			Standard			Inverted		
			Top	Middle	Bottom	Top	Middle	Bottom
Wheel	Rim	Tire						
22"	Mag / Spoke	Full Poly	16	16.5	17.5	17.5	18	19
22"	Mag / Spoke	Pneum / Insert	16.5	17	17.5	18	18.5	19
24"	Mag / Spoke	Full Poly	17	17.5	18	18.5	18.5	19.5
24"	Spoke / Spinerger	Clincher	17	17.5	18	18.5	19	19.5
24"	Mag / Spoke / Spinerger	Pneum / Insert	17	17.5	18.5	19	19	20
26"	Spoke	Clincher	18	18.5	19	19.5	20	20.5
26"	Spoke	Pneumatic	18	19	19.5	19.5	20	21

STANDARD Frame (10" Front / 9" Rear)			Rear Axle Plate Hole Position					
			Standard			Inverted		
			Top	Middle	Bottom	Top	Middle	Bottom
Wheel	Rim	Tire						
22"	Mag / Spoke	Full Poly	15	15.5	16.5	16.5	17.5	18
22"	Mag / Spoke	Pneum / Insert	15.5	16	17	17	17.5	18
24"	Mag / Spoke	Full Poly	16	16.5	17	17	18	18.5
24"	Spoke / Spinerger	Clincher	16	16.5	17	17.5	18	18.5
24"	Mag / Spoke / Spinerger	Pneum / Insert	16	16.5	17.5	17.5	18	18.5
26"	Spoke / Spinerger	Clincher	17	17.5	18	18.5	19	19.5
26"	Spoke / Spinerger	Pneumatic	17	18	18.5	18.5	19.5	20

Maximum Dump 2 3/4"

*** Custom Frame Height Through Configurator - REAR 6 1/2 - 13" (Max frame squeeze = 3 1/2") Calculate Accordingly

GP Series Front STF Chart

*All measurements given are an estimate and have been rounded to the nearest .5"


SHORT Frame (9 1/4" Front / 8 1/4" Rear)									
Fork	Spacer	Fork Hole	3"	4" Poly	5" Poly	6" Poly	6" Pneu	8" Poly	8" Pneu
			4" Lighted	5" Soft Roll	5" x 1.5	6" x 1.5			8x2
4 1/4	STD	Bottom	15.50	16.00	16.50	17.00	17.00	x	x
4 1/4	3/4	Bottom	16.50	17.00	17.00	17.50	18.00	x	x
4 1/4	1 1/2	Bottom	17.00	17.50	18.00	18.50	18.50	x	x
5 1/4	STD	Top	x	x	16.50	17.00	17.00	x	x
5 1/4	STD	Bottom	16.50	17.00	17.5	18.00	18.00	19.00	20.5
5 1/4	3/4	Top	x	x	17.50	18.00	18.00	x	x
5 1/4	3/4	Bottom	17.50	18.00	18.50	18.50	19.00	19.50	20.00
5 1/4	1 1/2	Top	x	x	18.00	18.50	18.50	x	x
5 1/4	1 1/2	Bottom	18.00	18.50	19.00	19.50	19.50	20.50	20.50
6	STD	Top	x	x	17.50	18.00	18.00	18.50	19.00
6	STD	Bottom	17.50	18.00	18.50	19.00	19.00	19.50	20.00
6	3/4	Top	x	x	18.00	18.50	18.50	19.50	19.50
6	3/4	Bottom	18.00	18.50	19.00	19.50	19.50	20.50	20.50
6	1 1/2	Top	x	x	19.00	19.50	19.50	20.00	20.50
6	1 1/2	Bottom	19.00	19.50	20.00	20.50	20.50	21.00	21.50
7	STD	Top	x	x	18.50	19.00	19.00	19.50	20.00
7	STD	Bottom	18.50	19.00	19.50	20.00	20.00	20.50	21.00
7	3/4	Top	x	x	19.00	19.50	19.50	20.50	20.50
7	3/4	Bottom	19.00	19.50	20.00	20.50	20.50	21.50	21.50
7	1 1/2	Top	x	x	20.00	20.50	20.50	21.00	21.50
7	1 1/2	Bottom	20.00	20.50	21.00	21.00	21.50	22.00	22.50
8x2	STD	Top	x	x	x	x	x	x	19.00
8x2	STD	Bottom	x	x	x	x	x	x	20.00
8x2	3/4	Top	x	x	x	x	x	x	19.50
8x2	3/4	Bottom	x	x	x	x	x	x	20.50
8x2	1 1/2	Top	x	x	x	x	x	x	20.50
8x2	1 1/2	Bottom	x	x	x	x	x	x	21.00

STANDARD Frame (10" Front / 9" Rear) ADD 3/4" to each measurement above

TALL Frame (11" Front / 10" Rear) ADD 1 3/4" to each measurement above

*** Custom Frame Height through Configurator - Front 9 - 13" - Calculate Accordingly (Minimum Front STF = 15 1/2" / Max = 26")

Frame Dimensions



Sunrise Medical (US) LLC
12002 Volunteer Blvd. Mount Juliet, TN 37122