

ENHANCE YOUR GAME!

Elite The future is here. A quantum leap in design, the Elite X, and Elite CX are the perfect combination of new technology and materials.

LIGHTNESS: 1 2 3 4 5

STIFFNESS: **1 2 3 4 5**

OPTIONS: **1 2 3 4**

Key Features

- A10 heat treated frame
- Tri-Frame design
- Low ride caster forks
- Titanium fasteners
- Ceramic bearings
- Welded sideguards with carbon fiber infill
- Straight carbon fiber axle (Elite X & Elite CX)
- Carbon fiber seating system (Elite CX)



Testimonials



"For over 30 years I've been playing basketball, I've been blessed to have played in some of the best chairs ever made. My latest, the RGK Elite, may be the best one yet."

Patrick Anderson (Canadian National Team and regarded as one of the best wheelchair basketball players of all time).



"The Elite X is so lightweight, and coming with that extra ability to adapt the toe in toe out post build with the adjustable carbon fiber camber bar, has enhanced my game. High quality equipment allows me to focus solely on the sport."

Bill Latham (Australian National Team).



