

Empulse



R20 shown on a QUICKIE® IRIS®

E | R20

POWER FOR THE CAREGIVERS

Empulse R20

E | R20

The only TRUE power assist device designed to support the caregivers.



LIGHTWEIGHT

At just under 10 lbs. the R20 is one of the lightest power assists on the market.

PERFORMANCE

No compromise on performance. The R20 achieves a 9 mile range, climbs up to 10° slopes, and handles a max user weight of 375 lbs.

RELIABILITY

A solid, long lasting tire, lithium-ion battery, and a 3 year motor warranty promise quality and reliability.

Compatible with most wheelchair frames.
(Frame tube diameter: 0.8" to 1.3")



R20 shown on a QUICKIE® 2



Lift Hook

Lift hook supports the wheel off the ground when not in use or during transport.



Convenient Carry Bag

Compact and easy to carry or store.



Simple Attachment and Removal

Attaches easily to brackets mounted on the frame.



Convenient Folding

The only power assist on the market that can be folded while still connected to the wheelchair.



Simple Control

Thumb throttle and three push-button speed settings.

R20 Specifications	
Weight Capacity	375 lbs
Range*	up to 9 miles
Max Slope	10°
Weight (including battery)	9.9 lbs
Frame Compatibility	Tube diameter 0.8" to 1.3"
Adjustable Max Speed	1.2 mph 2.3 mph 3.1 mph
Electronics	Charger 100 - 240V, AC 42V & 2.0A
Battery	Lithium Ion 5.8-6 Ah / 36 V (209 - 216 Wh)
Charging Time	Approx. 4.5 to 5 hours
Motor	250 W / 30 Nm
Wheelchair Seat Height Range	15.5" to 20.5"
Wheelchair Seat Width Range	13" to 23"

See order form as specifications may vary based upon configuration.

* Depending on user weight and topography



R20 Web Page



R20 Order Form



EMPOWER YOUR LIFE



Sunrise Medical (US) LLC, Fresno, CA 93727
MK-130001 Rev.B ©10.2022

800.333.4000 www.SunriseMedical.com/Empulse

