



JAY® GS CUSHION

E2622/ E2624

February 2023

Mark For: _____
 Date: _____
 Dealer Acct #: _____
 Dealer: _____
 Dealer Contact: _____
 Dealer Address: _____
 Dealer City: _____ ST: _____ ZIP: _____
 Dealer Phone: () _____ Fax: () _____
 Confirmation Email: _____
 Confirm Via: Fax Email

Submitting for: Quote Order
 PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____
 Attention: _____
 Address: _____
 Address: _____
 Ship To City: _____ ST: _____ ZIP: _____
 Ship To Phone: () _____ Fax: () _____

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

JAY GS Cushion

JAY GS CUSHION*

✓	Option #	HCPCS	Definition	Actual Width	Actual Depth	Price
<input type="checkbox"/>	JGS10	E2622	JAY GS Cushion 10" X 10" - 13"	10.2"	13.2"	\$627.00
<input type="checkbox"/>	JGS12	E2622	JAY GS Cushion 12" X 12" - 15"	12.1"	15.0"	\$627.00
<input type="checkbox"/>	JGS14	E2622	JAY GS Cushion 14" X 14" - 17"	14.1"	17.0"	\$627.00
<input type="checkbox"/>	JGS16	E2622	JAY GS Cushion 16" X 16" - 19"	15.8"	18.9"	\$627.00
<input type="checkbox"/>	JGS18	E2622	JAY GS Cushion 18" X 17" - 20"	17.7"	19.9"	\$627.00

JAY GS 2" POSITIONING** CUSHION*

<input type="checkbox"/>	JGS10P2	E2624	JAY GS 2" Pos Cushion 10" X 10" - 13"	10.2"	13.2"	\$672.00
<input type="checkbox"/>	JGS12P2	E2624	JAY GS 2" Pos Cushion 12" X 12" - 15"	12.1"	15.0"	\$672.00
<input type="checkbox"/>	JGS14P2	E2624	JAY GS 2" Pos Cushion 14" X 14" - 17"	14.1"	17.0"	\$672.00
<input type="checkbox"/>	JGS16P2	E2624	JAY GS 2" Pos Cushion 16" X 16" - 19"	15.8"	18.9"	\$672.00
<input type="checkbox"/>	JGS18P2	E2624	JAY GS 2" Pos Cushion 18" X 17" - 20"	17.7"	19.9"	\$672.00

JAY GS 3" POSITIONING** CUSHION*

<input type="checkbox"/>	N/A		Not Available 10" Wide w/ 3" Pos Pieces			
<input type="checkbox"/>	N/A		Not Available 12" Wide w/ 3" Pos Pieces			
<input type="checkbox"/>	JGS14P3	E2624	JAY GS 3" Pos Cushion 14" X 14" - 17"	14.1"	17.0"	\$672.00
<input type="checkbox"/>	JGS16P3	E2624	JAY GS 3" Pos Cushion 16" X 16" - 19"	15.8"	18.9"	\$672.00
<input type="checkbox"/>	JGS18P3	E2624	JAY GS 3" Pos Cushion 18" X 17" - 20"	17.7"	19.9"	\$672.00

JAY GS CUSHION* with Extra Cover

<input type="checkbox"/>	JGS10-2	E2622	JAY GS Cushion 10" X 10" - 13"	10.2"	13.2"	\$742.00
<input type="checkbox"/>	JGS12-2	E2622	JAY GS Cushion 12" X 12" - 15"	12.1"	15.0"	\$742.00
<input type="checkbox"/>	JGS14-2	E2622	JAY GS Cushion 14" X 14" - 17"	14.1"	17.0"	\$742.00
<input type="checkbox"/>	JGS16-2	E2622	JAY GS Cushion 16" X 16" - 19"	15.8"	18.9"	\$742.00
<input type="checkbox"/>	JGS18-2	E2622	JAY GS Cushion 18" X 17" - 20"	17.7"	19.9"	\$742.00

JAY GS 2" POSITIONING** CUSHION* with Extra Cover

<input type="checkbox"/>	JGS10P2-2	E2624	JAY GS 2" Pos Cushion 10" X 10" - 13"	10.2"	13.2"	\$788.00
<input type="checkbox"/>	JGS12P2-2	E2624	JAY GS 2" Pos Cushion 12" X 12" - 15"	12.1"	15.0"	\$788.00
<input type="checkbox"/>	JGS14P2-2	E2624	JAY GS 2" Pos Cushion 14" X 14" - 17"	14.1"	17.0"	\$788.00
<input type="checkbox"/>	JGS16P2-2	E2624	JAY GS 2" Pos Cushion 16" X 16" - 19"	15.8"	18.9"	\$788.00
<input type="checkbox"/>	JGS18P2-2	E2624	JAY GS 2" Pos Cushion 18" X 17" - 20"	17.7"	19.9"	\$788.00

JAY GS 3" POSITIONING** CUSHION* with Extra Cover

<input type="checkbox"/>	JGS14P3-2	E2624	JAY GS 3" Pos Cushion 14" X 14" - 17"	14.1"	17.0"	\$788.00
<input type="checkbox"/>	JGS16P3-2	E2624	JAY GS 3" Pos Cushion 16" X 16" - 19"	15.8"	18.9"	\$788.00
<input type="checkbox"/>	JGS18P3-2	E2624	JAY GS 3" Pos Cushion 18" X 17" - 20"	17.7"	19.9"	\$788.00

**All Positioning Cushions Include: Lateral Pelvic(Pr), Lateral Thigh (PR), and Medial Thigh (Ea) Supports.


JAY GS Size Matrix

		WIDTH											
		10	11	12	13	14	15	16	17	18	19	20	
DEPTH	10	JYW	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	
	11	JYW	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	
	12	JYW		JYW	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	
	13	JYW	JYW	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	
	14	N/A	JYW	JYW		JYW	N/A	N/A	N/A	N/A	N/A	N/A	
	15	N/A	JYW	JYW		JYW	N/A	N/A	N/A	N/A	N/A	N/A	
	16	N/A	N/A	N/A	JYW	JYW		JYW	N/A	N/A	N/A	N/A	
	17	N/A	N/A	N/A	JYW	JYW		JYW		JYW		JYW	N/A
	18	N/A	N/A	N/A	N/A	N/A	JYW	JYW		JYW		JYW	N/A
	19	N/A	N/A	N/A	N/A	N/A	JYW		JYW		JYW	N/A	N/A
20	N/A	N/A	N/A	N/A	N/A	N/A	N/A	JYW		JYW	N/A	N/A	

<input type="checkbox"/>	=	Standard Offering
JYW	=	Offered Via Jay Your Way
N/A	=	Not available

JAY GS Measuring Guide

Measurements made from completed cushion with covers.



Width

X = Width

Note: Measurement made from rear bottom edge to bottom edge of the cushion.
Ex: JGS12 X=12"

Depth

Y = Depth

Note: Measurement made from front edge of rear cutout and front lower of cushion.
Ex: JGS12 Y = 15"

Bottom View

*All JAY GS Cushions ship standard with a 3 layer cover. The cover includes a top layer of 4 way stretch fabric, a middle layer of our Ionic™ impregnated fabric, and a bottom layer of 3DX micro climatic fabric.

JAY YOUR WAY MODIFICATIONS

Note: To order sizes offered via JAY Your Way add a "M" to the front of the part number and follow the prompts.

Fluid Volume Change										
<input type="checkbox"/>	OF			Fluid Volume Overfill						NC
<input type="checkbox"/>	<input type="checkbox"/> 5%	<input type="checkbox"/> 10%	<input type="checkbox"/> 15%	<input type="checkbox"/> 20%	<input type="checkbox"/> 25%	<input type="checkbox"/> 30%	<input type="checkbox"/> 35%	<input type="checkbox"/> 40%		
<input type="checkbox"/>	UF			Fluid Volume Underfill (reduction in well only)						NC
<input type="checkbox"/>	<input type="checkbox"/> -5%	<input type="checkbox"/> -10%	<input type="checkbox"/> -20%	<input type="checkbox"/> -40%	<input type="checkbox"/> -50%					
Obliquity Fluid Volume Change										
<input type="checkbox"/>	OOFR			Obliquity Overfill - Right Side						NC
<input type="checkbox"/>	<input type="checkbox"/> 5%	<input type="checkbox"/> 10%	<input type="checkbox"/> 15%	<input type="checkbox"/> 20%	<input type="checkbox"/> 25%	<input type="checkbox"/> 30%	<input type="checkbox"/> 35%	<input type="checkbox"/> 40%		
<input type="checkbox"/>	OOFL			Obliquity Overfill - Left Side						NC
<input type="checkbox"/>	<input type="checkbox"/> 5%	<input type="checkbox"/> 10%	<input type="checkbox"/> 15%	<input type="checkbox"/> 20%	<input type="checkbox"/> 25%	<input type="checkbox"/> 30%	<input type="checkbox"/> 35%	<input type="checkbox"/> 40%		
<input type="checkbox"/>	OUFR			Obliquity Underfill - Right Side						NC
<input type="checkbox"/>	<input type="checkbox"/> -5%	<input type="checkbox"/> -10%	<input type="checkbox"/> -20%	<input type="checkbox"/> -40%	<input type="checkbox"/> -50%					
<input type="checkbox"/>	OUFL			Obliquity Underfill - Left Side						NC
<input type="checkbox"/>	<input type="checkbox"/> -5%	<input type="checkbox"/> -10%	<input type="checkbox"/> -20%	<input type="checkbox"/> -40%	<input type="checkbox"/> -50%					

Size Change			
<input type="checkbox"/>	IW	Increase Width	N/A
<input type="checkbox"/>	1"		
<input type="checkbox"/>	RW	Reduce Width	NC
<input type="checkbox"/>	-1"		
Cover Options			
<input type="checkbox"/>	CVR	Cover Options	NC
<input type="checkbox"/>	Rvrs Dartex	<input type="checkbox"/> Air Exchange	<input type="checkbox"/> Incontinent <input type="checkbox"/> 3DX Only

Don't see the option you want? NO PROBLEM! Write your request below and we can evaluate it and provide a quote.

Notes:

Please visit www.sunparts.us for spare parts information.



Sunrise Medical (US) LLC
2842 Business Park Ave. · Fresno, CA 93727 · USA

Customer Service: 800-333-4000 Fax: 800-300-7502 Email: orders@sunmed.com www.sunrisemedical.com
Sunrise Medical (US) LLC. © 2/2023, MK-100295 REV. E