

Submitting for: Quote Order

Account: _____
 Date: _____
 PO#: _____
 Buyer: _____
 Marked For: _____

Ship To: _____
 Name: _____
 Address: _____
 City,State,Zip: _____
 Phone: _____

The HCPCS CODES herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative HCPCS CODES, including misc. codes to ensure access for their enrollees. The use of HCPCS CODES does not ensure coverage or payment.

a	Option #	HCPCS	Definition	Price	a	Option #	HCPCS	Definition	Price
---	----------	-------	------------	-------	---	----------	-------	------------	-------

1. QUICK SHIP MIGHTY LITE

* The Quick-Ship Mighty Lite includes the "H" Harness & Padded Covers (E0960) and a Three Points Positioning Belt (E0978)

Navy Blue Upholstery

<input type="checkbox"/>	201NB93	E1236	12" Seat Width	\$2,155
<input type="checkbox"/>	201NB94	E1236	14" Seat Width	\$2,155
<input type="checkbox"/>	201NB95	E1236	16" Seat Width	\$2,155
<input type="checkbox"/>	201NB96	E1236	18" Seat Width	\$2,155

Raspberry Upholstery

<input type="checkbox"/>	201NB97	E1236	12" Seat Width	\$2,155
<input type="checkbox"/>	201NB98	E1236	14" Seat Width	\$2,155
<input type="checkbox"/>	201NB99	E1236	16" Seat Width	\$2,155
<input type="checkbox"/>	201NB100	E1236	18" Seat Width	\$2,155

Quick Ship Options

<input type="checkbox"/>	201TO44	K0108	Transit Option	\$200
<input type="checkbox"/>	201HE75	E0966	Headrest Extension	\$100

2. CONFIGURED MIGHTY LITE

1. FRAME

Frame Color

<input type="checkbox"/>	201C2	Black	std
--------------------------	-------	-------	------------

Transit Option

<input type="checkbox"/>	201TO44	K0108	Transit Option	\$200
--------------------------	---------	-------	----------------	-------

2. SIZE/SEATING

Seat Width

<input type="checkbox"/>	201S49	E1236	12" Seat Width	\$1,955
<input type="checkbox"/>	201S50	E1236	14" Seat Width	\$1,955
<input type="checkbox"/>	201S51	E1236	16" Seat Width	\$1,955
<input type="checkbox"/>	201S52	E1236	18" Seat Width	\$1,955

Seat Angle

<input type="checkbox"/>	201UA05	15°	std
<input type="checkbox"/>	201UA06	22°	NC

Seat to Back Angle

<input type="checkbox"/>	201UB23	85°	NC
<input type="checkbox"/>	201UB10	90°	std
<input type="checkbox"/>	201UB24	95°	NC

3. WHEELS/WHEEL LOCKS

Wheels

<input type="checkbox"/>	201RT37	10" Semi-pneumatic (Pr.)	std
--------------------------	---------	--------------------------	------------

Wheel Lock

<input type="checkbox"/>	201WL44	Toggle Wheel Lock (Pr.)	std
--------------------------	---------	-------------------------	------------

Casters

<input type="checkbox"/>	201CW111	8" Semi-pneumatic (Pr.)	std
--------------------------	----------	-------------------------	------------

4. POSITIONING

Pelvic Positioning

<input type="checkbox"/>	201A807	E0978	Three Points Positioning Belt (3rd PT may interfere w/ foam wedge)	\$50
<input type="checkbox"/>	201A808		Foam Wedge (will reduce clearance between tray and legs)	\$127

Foot & Leg Positioning

<input type="checkbox"/>	201F302		Ht. Adj. S/A Footrest (Pr.)	std
<input type="checkbox"/>	201F303	K0040	Angle Adj. Footplates(Pr.)	\$150

Foot & Ankle Positioning Accessories

<input type="checkbox"/>	201F304		Small (Pr.)	\$180
<input type="checkbox"/>	201F305		Medium (Pr.)	\$180
<input type="checkbox"/>	201F306		Large (Pr.)	\$180
<input type="checkbox"/>	201F307		Xtra Large (Pr.) (not avail. w/ 12" or 14" seat widths)	\$180

Medial Thigh

<input type="checkbox"/>	201KH11	E0957	Medial Thigh Support (not avail w/ foam wedge)	\$150
--------------------------	---------	-------	--	-------

Substitute Navy or Raspberry pack-cloth on Above

<input type="checkbox"/>	Navy	NC
<input type="checkbox"/>	Raspberry	NC

Trunk Positioning

<input type="checkbox"/>	201K17	E0960	H Harness & Padded Covers (Cannot order Full Torso Support Vest with this item)	\$150
<input type="checkbox"/>	201K18	E0956	Adj. Lat.Support & Scoli Strap (Cannot order Full Torso Support Vest with this item)	\$162

Substitute Navy or Raspberry pack-cloth on Above

<input type="checkbox"/>	Navy	NC
<input type="checkbox"/>	Raspberry	NC

Trunk Positioning Accessories

<input type="checkbox"/>	201K19		Full Torso Support Vest (cannot order "H" Harness or Adj. Lat. Support & Scoli Strap with this item)	\$250
--------------------------	--------	--	--	-------

Substitute Navy or Raspberry pack-cloth on Above

<input type="checkbox"/>	Navy	NC
<input type="checkbox"/>	Raspberry	NC

Head Support

<input type="checkbox"/>	201ZH23		Padded Headwings (may interfere with Full Torso Vest)	\$191
<input type="checkbox"/>	201HE75	E0966	Headrest Extension	\$100

Substitute Navy or Rasp. pack-cloth for Head Supports

<input type="checkbox"/>	Black	NC
<input type="checkbox"/>	Navy	NC
<input type="checkbox"/>	Raspberry	NC

5. KEY OPTIONS

Accessories

<input type="checkbox"/>	201LF10		Closure Strap	std
<input type="checkbox"/>	201HE76		Headrest Cover (Canopy)	\$190
<input type="checkbox"/>	201LF9	E0950	Upper Extremity Support Tray (Foam Wedge reduces clearance between tray and legs)	\$230
<input type="checkbox"/>	201LF6		Utility Bag	\$65

5. COLOR

Color Upholstery

<input type="checkbox"/>	201UC3		Navy	NC
<input type="checkbox"/>	201UC16		Raspberry	NC

