

A man in a wheelchair is smiling and holding two basketballs high in the air. He is wearing a grey t-shirt with a logo and purple shorts. The background is a gymnasium with bleachers. A large, stylized number '2' is overlaid on the right side of the image.

ONE SIZE DOES
NOT FIT ALL.
**BUT ONE CHAIR
DOES.**

QUICKIE 2

PERSONAL FIT IS ESSENTIAL

The QUICKIE® 2 is the most sought-after folding ultra lightweight manual wheelchair. Featuring a modular design that can change with you, the QUICKIE 2 offers the broadest range of possible configurations in the industry: over 100,000 different complex rehab options and counting. Designed to meet all of your unique needs without compromising fit or performance, only one wheelchair can live up to the name QUICKIE.



STRONG AND LIGHT

The revolutionary engineering of the QUICKIE 2 offers one of the strongest and most durable frames on the market without sacrificing weight savings.

We make the QUICKIE 2 out of the same material used to make rockets and airplanes, leveraging the strength of 7000 series aerospace aluminum for increased durability and reliability.

The design of the QUICKIE 2 crossbrace is something you won't find anywhere else. We take that 7000 series aerospace aluminum and extrude it into diamond-shaped tubing welded in a four-bar pattern. This optimized geometry creates a stronger, more rigid crossbrace at a lighter weight, ensuring long-lasting durability you can rely on.



MULTIPLE WEIGHT CAPACITY OPTIONS

Fit every body safely and comfortably with the QUICKIE 2 Lite (265-lb user weight capacity) or QUICKIE 2 Heavy Duty (350-lb user weight capacity).





LARGE ASSORTMENT OF SEAT WIDTHS AND DEPTHS

Offering 11 seat depths, 8 frame depths, and 13 frame widths, the QUICKIE 2 is designed to fit you both today and tomorrow.

UNIQUE CASTER HOUSING, FORK AND WHEEL SIZES

Achieve a more nimble ride with added stability by leveraging the QUICKIE 2's unique range of caster sizes and options to better navigate your environment.

QUICKIE 2 LITE OPTION

The QUICKIE 2's minimalist counterpart, the QUICKIE 2 Lite, features an extremely lightweight, modular frame with a rounded rear for unique style and either fixed or swing-away front hangers.

PATENTED Q-FIT TECHNOLOGY

Our Q-Fit Technology guarantees a tight fit that can be maintained for the life of your chair.



JAY[®] SUPERIOR CLINICAL SEATING

JAY cushions and backs come in a wide variety of sizes featuring varying levels of stability and pressure management. Prefer air, fluid, or cooling technology in your cushion or back? QUICKIE 2 has you covered with the industry's largest in-house selection of seating options. From JAY Your Way to made to order, we believe that comfort is not optional, it's a must.



MOST COLOR OPTIONS

With more than 24 frame colors, 4 patterned frame wraps, and 7 anodized parts colors, you can design a chair that's truly unique, just like you.



IMPORTANT NOTE: The color options shown are for reference only. They are not intended as accurate representations. Refer to the QUICKIE Color Brochure or a QUICKIE Color Ring for representations of QUICKIE colors.

Zebra Kolorfusion™ pattern is available on any base color. Shown here on Candy Blue.

QUICKIE 2 Specifications

HCPCS Code	K0005
Overall Weight¹	25.4 lbs
Transport Weight¹	13.6 lbs
Weight Capacity	Standard: 300 lbs Heavy Duty: 350 lbs Lite: 265 lbs
Seat Width	12" to 24"
Seat Depth	12" to 22"
Back Angle Adjustment	Fore (-) 3° to Aft (+) 15°
Front Seat Height	13.5" to 22"
Rear Seat Height	13" to 22"
Hanger Angle	60°, 70°, 80°, 90°
Transit Approved	Yes ²

1. Standard Weight: 16" x 16" standard configuration, with wheels, hangers, composite footplates, and padded swing-away hanger. Transport Weight: 16" x 16" standard configuration, without rear wheels or foot supports.

2. Conforms with ANSI/RESNA WC-19 standards. Please refer to Quickie® & Zippie® WC-19 Transit Securement Supplement (MK-100072) for additional information. Supplement is available for download at www.SunriseMedical.com



Sunrise Medical (US) LLC, Fresno, CA 93727
MK-100079 Rev. E ©01.2021

800.333.4000 www.SunriseMedical.com

