



JAY® Fusion™ Core Adjustable Cushion

E2622 / E2623

January 2022

Mark For: _____

Date: _____

Dealer Acct #: _____

Dealer: _____

Dealer Contact: _____

Dealer Address: _____

Dealer City: _____

ST: _____

ZIP: _____

Dealer Phone: () _____

Fax: () _____

Confirmation Email: _____

Confirm Via: Fax Email

Submitting for: Quote Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____

Attention: _____

Address: _____

Address: _____

Ship To City: _____

ST: _____

ZIP: _____

Ship To Phone: () _____

Fax: () _____

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

JAY FUSION CORE CUSHION

JAY Fusion Core Cushion w/ JAY FLOW Fluid Insert

✓	Option #	HCPCS	Definition	Price
<input type="checkbox"/>	JFUSIONXXYYC	E2622	JAY Fusion Core Cushion 14-21 wide	\$458
<i>Note: XX = Width and YY = Depth (Select appropriate assembly # on page 2)</i> Includes: Reverse Dartex fabric outer cover, foam base, and a JAY Flow fluid insert with a 300 lb weight capacity.				
<input type="checkbox"/>	JFUSIONXXYYC	E2623	JAY Fusion Core Cushion 22-24 wide	\$563
<i>Note: XX = Width and YY = Depth (Select appropriate assembly # on page 2)</i> Includes: Reverse Dartex fabric outer cover, foam base, and a JAY Flow fluid insert with a 500 lb weight capacity.				


JAY Fusion Core Size Matrix

		WIDTH											
		14	15	16	17	18	19	20	21	22	23	24	
DEPTH	14						JYW	JYW	JYW	JYW	N/A	N/A	
	15							JYW	JYW	JYW	JYW		
	16							JYW	JYW	JYW	JYW		
	17								JYW	JYW			
	18									JYW	JYW		
	19	JYW											
	20	JYW											
	21	JYW											
	22	N/A	JYW										
	23	N/A	N/A	JYW									
24	N/A	N/A	N/A	JYW									

<input type="checkbox"/>	=	Standard Offering
JYW	=	Offered Via Jay Your Way
N/A	=	Not available

JAY Fusion Core Measuring Guide

Measurements made from completed cushion with covers.



Bottom View

Width

X = Width

Note: Measurement made from outside edge to outside edge of cushion.

Ex: JFUSION1618C

X = 16"

Depth

Y = Depth

Note: Measurement made from front edge of rear cutout and front lower edge of cushion.

Ex: JFUSION1618C

Y = 18"

Note: To order sizes offered via JAY Your Way add a "M" to the front of the part number and follow the prompts. For additional modifications see Page 2 or visit www.sunrisemedical.com to view even more!

Please visit www.sunparts.us for spare parts information.

JAY Fusion Core Cushion w/ Fluid Insert

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
<input type="checkbox"/>	JFUSION1414C	E2622	JAY Fusion Core 14" X 14" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION1918C	E2622	JAY Fusion Core 19" X 18" w/ Fluid	\$458
<input type="checkbox"/>	JFUSION1415C	E2622	JAY Fusion Core 14" X 15" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION1919C	E2622	JAY Fusion Core 19" X 19" w/ Fluid	\$458
<input type="checkbox"/>	JFUSION1416C	E2622	JAY Fusion Core 14" X 16" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION1920C	E2622	JAY Fusion Core 19" X 20" w/ Fluid	\$458
<input type="checkbox"/>	JFUSION1417C	E2622	JAY Fusion Core 14" X 17" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION1921C	E2622	JAY Fusion Core 19" X 21" w/ Fluid	\$458
<input type="checkbox"/>	JFUSION1418C	E2622	JAY Fusion Core 14" X 18" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION1922C	E2622	JAY Fusion Core 19" X 22" w/ Fluid	\$458
<input type="checkbox"/>	JFUSION1514C	E2622	JAY Fusion Core 15" X 14" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION1923C	E2622	JAY Fusion Core 19" X 23" w/ Fluid	\$458
<input type="checkbox"/>	JFUSION1515C	E2622	JAY Fusion Core 15" X 15" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION1924C	E2622	JAY Fusion Core 19" X 24" w/ Fluid	\$458
<input type="checkbox"/>	JFUSION1516C	E2622	JAY Fusion Core 15" X 16" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2015C	E2622	JAY Fusion Core 20" X 15" w/ Fluid	\$458
<input type="checkbox"/>	JFUSION1517C	E2622	JAY Fusion Core 15" X 17" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2016C	E2622	JAY Fusion Core 20" X 16" w/ Fluid	\$458
<input type="checkbox"/>	JFUSION1518C	E2622	JAY Fusion Core 15" X 18" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2017C	E2622	JAY Fusion Core 20" X 17" w/ Fluid	\$458
<input type="checkbox"/>	JFUSION1519C	E2622	JAY Fusion Core 15" X 19" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2018C	E2622	JAY Fusion Core 20" X 18" w/ Fluid	\$458
<input type="checkbox"/>	JFUSION1520C	E2622	JAY Fusion Core 15" X 20" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2019C	E2622	JAY Fusion Core 20" X 19" w/ Fluid	\$458
<input type="checkbox"/>	JFUSION1521C	E2622	JAY Fusion Core 15" X 21" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2020C	E2622	JAY Fusion Core 20" X 20" w/ Fluid	\$458
<input type="checkbox"/>	JFUSION1614C	E2622	JAY Fusion Core 16" X 14" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2021C	E2622	JAY Fusion Core 20" X 21" w/ Fluid	\$458
<input type="checkbox"/>	JFUSION1615C	E2622	JAY Fusion Core 16" X 15" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2022C	E2622	JAY Fusion Core 20" X 22" w/ Fluid	\$458
<input type="checkbox"/>	JFUSION1616C	E2622	JAY Fusion Core 16" X 16" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2023C	E2622	JAY Fusion Core 20" X 23" w/ Fluid	\$458
<input type="checkbox"/>	JFUSION1617C	E2622	JAY Fusion Core 16" X 17" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2024C	E2622	JAY Fusion Core 20" X 24" w/ Fluid	\$458
<input type="checkbox"/>	JFUSION1618C	E2622	JAY Fusion Core 16" X 18" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2117C	E2622	JAY Fusion Core 21" X 17" w/ Fluid	\$458
<input type="checkbox"/>	JFUSION1619C	E2622	JAY Fusion Core 16" X 19" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2118C	E2622	JAY Fusion Core 21" X 18" w/ Fluid	\$458
<input type="checkbox"/>	JFUSION1620C	E2622	JAY Fusion Core 16" X 20" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2119C	E2622	JAY Fusion Core 21" X 19" w/ Fluid	\$458
<input type="checkbox"/>	JFUSION1621C	E2622	JAY Fusion Core 16" X 21" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2120C	E2622	JAY Fusion Core 21" X 20" w/ Fluid	\$458
<input type="checkbox"/>	JFUSION1622C	E2622	JAY Fusion Core 16" X 22" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2121C	E2622	JAY Fusion Core 21" X 21" w/ Fluid	\$458
<input type="checkbox"/>	JFUSION1714C	E2622	JAY Fusion Core 17" X 14" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2122C	E2622	JAY Fusion Core 21" X 22" w/ Fluid	\$458
<input type="checkbox"/>	JFUSION1715C	E2622	JAY Fusion Core 17" X 15" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2123C	E2622	JAY Fusion Core 21" X 23" w/ Fluid	\$458
<input type="checkbox"/>	JFUSION1716C	E2622	JAY Fusion Core 17" X 16" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2124C	E2622	JAY Fusion Core 21" X 24" w/ Fluid	\$458
<input type="checkbox"/>	JFUSION1717C	E2622	JAY Fusion Core 17" X 17" w/ Fluid	\$458	22" & Greater in Width				
<input type="checkbox"/>	JFUSION1718C	E2622	JAY Fusion Core 17" X 18" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2217C	E2623	JAY Fusion Core 22" X 17" w/ Fluid	\$563
<input type="checkbox"/>	JFUSION1719C	E2622	JAY Fusion Core 17" X 19" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2218C	E2623	JAY Fusion Core 22" X 18" w/ Fluid	\$563
<input type="checkbox"/>	JFUSION1720C	E2622	JAY Fusion Core 17" X 20" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2219C	E2623	JAY Fusion Core 22" X 19" w/ Fluid	\$563
<input type="checkbox"/>	JFUSION1721C	E2622	JAY Fusion Core 17" X 21" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2220C	E2623	JAY Fusion Core 22" X 20" w/ Fluid	\$563
<input type="checkbox"/>	JFUSION1722C	E2622	JAY Fusion Core 17" X 22" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2221C	E2623	JAY Fusion Core 22" X 21" w/ Fluid	\$563
<input type="checkbox"/>	JFUSION1723C	E2622	JAY Fusion Core 17" X 23" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2222C	E2623	JAY Fusion Core 22" X 22" w/ Fluid	\$563
<input type="checkbox"/>	JFUSION1814C	E2622	JAY Fusion Core 18" X 14" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2223C	E2623	JAY Fusion Core 22" X 23" w/ Fluid	\$563
<input type="checkbox"/>	JFUSION1815C	E2622	JAY Fusion Core 18" X 15" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2224C	E2623	JAY Fusion Core 22" X 24" w/ Fluid	\$563
<input type="checkbox"/>	JFUSION1816C	E2622	JAY Fusion Core 18" X 16" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2319C	E2623	JAY Fusion Core 23" X 19" w/ Fluid	\$563
<input type="checkbox"/>	JFUSION1817C	E2622	JAY Fusion Core 18" X 17" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2320C	E2623	JAY Fusion Core 23" X 20" w/ Fluid	\$563
<input type="checkbox"/>	JFUSION1818C	E2622	JAY Fusion Core 18" X 18" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2321C	E2623	JAY Fusion Core 23" X 21" w/ Fluid	\$563
<input type="checkbox"/>	JFUSION1819C	E2622	JAY Fusion Core 18" X 19" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2322C	E2623	JAY Fusion Core 23" X 22" w/ Fluid	\$563
<input type="checkbox"/>	JFUSION1820C	E2622	JAY Fusion Core 18" X 20" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2323C	E2623	JAY Fusion Core 23" X 23" w/ Fluid	\$563
<input type="checkbox"/>	JFUSION1821C	E2622	JAY Fusion Core 18" X 21" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2324C	E2623	JAY Fusion Core 23" X 24" w/ Fluid	\$563
<input type="checkbox"/>	JFUSION1822C	E2622	JAY Fusion Core 18" X 22" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2419C	E2623	JAY Fusion Core 24" X 19" w/ Fluid	\$563
<input type="checkbox"/>	JFUSION1823C	E2622	JAY Fusion Core 18" X 23" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2420C	E2623	JAY Fusion Core 24" X 20" w/ Fluid	\$563
<input type="checkbox"/>	JFUSION1824C	E2622	JAY Fusion Core 18" X 24" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2421C	E2623	JAY Fusion Core 24" X 21" w/ Fluid	\$563
<input type="checkbox"/>	JFUSION1915C	E2622	JAY Fusion Core 19" X 15" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2422C	E2623	JAY Fusion Core 24" X 22" w/ Fluid	\$563
<input type="checkbox"/>	JFUSION1916C	E2622	JAY Fusion Core 19" X 16" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2423C	E2623	JAY Fusion Core 24" X 23" w/ Fluid	\$563
<input type="checkbox"/>	JFUSION1917C	E2622	JAY Fusion Core 19" X 17" w/ Fluid	\$458	<input type="checkbox"/>	JFUSION2424C	E2623	JAY Fusion Core 24" X 24" w/ Fluid	\$563

JAY Your Way - Modifications

Fluid Volume Change					
<input type="checkbox"/>	OF	Fluid Volume Overfill			NC
<input type="checkbox"/>	5%	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	10%	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	15%	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	20%	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	25%	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	30%	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	35%	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	40%	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	UF	Fluid Volume Underfill			NC
<input type="checkbox"/>	-5%	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	-10%	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	-15%	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	-20%	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Obliquity Fluid Volume Change					
<input type="checkbox"/>	OOFR	Obliquity Overfill - Right Side			NC
<input type="checkbox"/>	5%	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	10%	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	15%	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	20%	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	25%	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	30%	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	35%	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	40%	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	OOFL	Obliquity Overfill - Left Side			NC
<input type="checkbox"/>	5%	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	10%	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	15%	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	20%	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	25%	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	30%	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	35%	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	40%	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	OUFR	Obliquity Underfill - Right Side			NC
<input type="checkbox"/>	-5%	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	-10%	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	-15%	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	-20%	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	OUFL	Obliquity Underfill - Left Side			NC
<input type="checkbox"/>	-5%	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	-10%	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	-15%	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	-20%	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Length Change					
<input type="checkbox"/>	IL	Increase Length			NC
<input type="checkbox"/>	+1"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	+2"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	+3"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	+4"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	RL	Reduce Length			NC
<input type="checkbox"/>	-1"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	-2"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	IW	Increase Width			NC
<input type="checkbox"/>	+1"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	+2"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	+3"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	+4"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	RW	Reduce Width			NC
<input type="checkbox"/>	-1"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	-2"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Reduced Profile					
<input type="checkbox"/>	RP	Reduced Profile Option			NC

Notes: Removes 5/8" of the top layer soft foam from the cushion for a reduced profile.

Don't see the option you want? NO PROBLEM! Write your request below and we can evaluate it and provide a quote.

Notes:

