



Freestander

August 2023

Mark For: _____ Date: _____ Dealer Acct #: _____ Dealer: _____ Dealer Contact: _____ Dealer Address: _____ Dealer City: _____ ST: _____ ZIP: _____ Dealer Phone: _____ Fax: _____ Confirmation Email: _____ Confirm Via: <input type="checkbox"/> Fax <input type="checkbox"/> Email	Submitting for: <input type="checkbox"/> Quote <input type="checkbox"/> Order PO#: _____ <div style="text-align: center;">ADDITIONAL SHIPPING INFORMATION</div> Ship To: _____ Attention: _____ Address: _____ Address: _____ Ship To City: _____ ST: _____ ZIP: _____ Ship To Phone: _____ Fax: _____
---	---

*The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.*

Freestander	HCPCS: E0638
--------------------	---------------------

Freestander

Selection required. Includes Chest and Hip Pad, Flexible Lateral, PU Knee Supports, Footplate, Tray Cover and Activity Bowl.

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	LFS/1-MW	Size 1 Freestander - Meadow		\$2,753
<input type="checkbox"/>	LFS/2-MB	Size 2 Freestander - Mondeo Blue		\$2,983
<input type="checkbox"/>	LFS/3-MB	Size 3 Freestander - Mondeo Blue		\$3,374

Additional Positioning

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	LFS/1/B-MW	Size 1 Hip Belt - Meadow		\$185
<input type="checkbox"/>	LFS/2/B-MB	Size 2 Hip Belt - Mondeo Blue		\$208
<input type="checkbox"/>	LFS/3/B-MB	Size 3 Hip Belt - Mondeo Blue		\$226
<input type="checkbox"/>	LFS/1/L-MW	Size 1 Solid Pelvic Plate - Meadow <i>Pelvic Band selection required</i>		\$323
<input type="checkbox"/>	LFS/2/L-MB	Size 2 Solid Pelvic Plate - Mondeo Blue <i>Pelvic Band selection required</i>		\$365
<input type="checkbox"/>	LFS/3/L-MB	Size 3 Solid Pelvic Plate - Mondeo Blue <i>Pelvic Band selection required</i>		\$408
<input type="checkbox"/>	LPS/1/B-MW	Size 1 Pelvic Band - Meadow		\$190
<input type="checkbox"/>	LPS/2/B-MB	Size 2 Pelvic Band - Mondeo Blue		\$208
<input type="checkbox"/>	LPS/3/B-MB	Size 3 Pelvic Band - Mondeo Blue		\$226

Additional Knee Supports

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	LPS/1/K-MW	Size 1 Upholstered Knee Cups - Meadow (Pair)		\$129
<input type="checkbox"/>	LPS/2/K-MB	Size 2 Upholstered Knee Cups - Mondeo Blue (Pair)		\$129
<input type="checkbox"/>	LPS/3/K-MB	Size 3 Upholstered Knee Cups - Mondeo Blue (Pair)		\$129
<input type="checkbox"/>	LPS/4/K-MB	Size 4 Upholstered Knee Cups - Mondeo Blue (Pair)		\$129

Additional Foot Positioning

Sandal selection required with Sandal Raiser

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	152-1600	Small Sandals with Straps (Pair)		\$267
<input type="checkbox"/>	152-2600	Medium Sandals with Straps (Pair)		\$291
<input type="checkbox"/>	152-3600	Large Sandals with Straps (Pair)		\$314
<input type="checkbox"/>	117-787	Sandal Raiser 1" (Each) Requires Sandal Selection	Quantity:	(1-6) \$93

Technical Specifications - Imperial

Specification	Size 1	Size 2	Size 3
Age (approximate)	1 1/2 - 5 years	4 - 10 years	9 - 18 years
User Body Height (Min / Max)	31.5" / 45.3"	43.3" / 57"	55.1" / 72.8"
Maximum User Weight	132 lbs	176 lbs	198 lbs
Chest Pad Height (Min / Max)	20.7" / 33.5"	26.8" / 43.3"	38.4" / 56.1"
Width Between Chest Supports (Min / Max)	6" / 10"	8" / 12"	10" / 14"
Activity Tray	16" x 21"	18" x 23"	20" x 25"

Technical Specifications - Metric

Specification	Size 1	Size 2	Size 3
Age (approximate)	1 1/2 - 5 years	4 - 10 years	9 - 18 years
User Body Height (Min / Max)	800 mm / 1151 mm	1100 mm / 1448 mm	1400 mm / 1849 mm
Maximum User Weight	60 kgs	80 kgs	90 kgs
Chest Pad Height (Min / Max)	470 mm / 680 mm	680 mm / 975 mm	975 mm / 1350 mm
Width Between Chest Supports (Min / Max)	152 mm / 254 mm	203 mm / 305 mm	254 mm / 356 mm
Activity Tray	400 mm x 505 mm	450 mm x 575 mm	500 mm x 625 mm



Sunrise Medical (US) LLC2842 Business Park Ave. • Fresno, CA 93727 • USA