



ORDER FORM

# Freestander

Upright accessible standing frame  
from Early Intervention to adult.

Customer			
Account			
Account #		Date	
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #	
Buyer			
Mark For			

Shipping Address			
Name			
Address			
City		State	
ZIP		Phone	
Email / Fax			

The HCPCS CODES herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative HCPCS CODES, including misc. codes to ensure access for their enrollees. The use of HCPCS CODES does not ensure coverage or payment.

## Freestander

Size 1 Weight Capacity: 132 lbs | Size 2 Weight Capacity: 176 lbs | Size 3 Weight Capacity: 198 lbs

✓ OPTION	DESCRIPTION	HCPCS	RETAIL
----------	-------------	-------	--------

Freestander			\$
-------------	--	--	----

Selection Required. Includes Chest and Hip Pad, Flexible Lateral, PU Knee Supports, Footplate, Tray Cover and Activity Bowl.

<input type="checkbox"/>	LFS/1-MW	Size 1 Freestander - Meadow	E0638	\$2,818
<input type="checkbox"/>	LFS/2-MB	Size 2 Freestander - Mondeo Blue	E0638	\$3,053
<input type="checkbox"/>	LFS/3-MB	Size 3 Freestander - Mondeo Blue	E0638	\$3,453

Additional Positioning			\$
------------------------	--	--	----

If Solid Pelvic Plate selected, pelvic band selection required.

<input type="checkbox"/>	LFS/1/B-MW	Size 1 Hip Belt - Meadow		\$189
<input type="checkbox"/>	LFS/2/B-MB	Size 2 Hip Belt - Mondeo Blue		\$213
<input type="checkbox"/>	LFS/3/B-MB	Size 3 Hip Belt - Mondeo Blue		\$231
<input type="checkbox"/>	LFS/1/L-MW	Size 1 Solid Pelvic Plate - Meadow		\$331
<input type="checkbox"/>	LFS/2/L-MB	Size 2 Solid Pelvic Plate - Mondeo Blue		\$374
<input type="checkbox"/>	LFS/3/L-MB	Size 3 Solid Pelvic Plate - Mondeo Blue		\$418
<input type="checkbox"/>	LPS/1/B-MW	Size 1 Pelvic Band - Meadow		\$194
<input type="checkbox"/>	LPS/2/B-MB	Size 2 Pelvic Band - Mondeo Blue		\$213
<input type="checkbox"/>	LPS/3/B-MB	Size 3 Pelvic Band - Mondeo Blue		\$231

Additional Knee Supports			\$
--------------------------	--	--	----

<input type="checkbox"/>	LPS/1/K-MW	Size 1 Upholstered Knee Cups - Meadow (Pair)		\$132
<input type="checkbox"/>	LPS/2/K-MB	Size 2 Upholstered Knee Cups - Mondeo Blue (Pair)		\$132
<input type="checkbox"/>	LPS/3/K-MB	Size 3 Upholstered Knee Cups - Mondeo Blue (Pair)		\$132
<input type="checkbox"/>	LPS/4/K-MB	Size 4 Upholstered Knee Cups - Mondeo Blue (Pair)		\$132

Additional Foot Positioning			\$
-----------------------------	--	--	----

Sandal selection required with Sandal Raiser.

<input type="checkbox"/>	152-1600	Small Sandals with Straps (Pair)		\$273
<input type="checkbox"/>	152-2600	Medium Sandals with Straps (Pair)		\$298
<input type="checkbox"/>	152-3600	Large Sandals with Straps (Pair)		\$321
<input type="checkbox"/>	117-787	Sandal Raiser 1" (Each) Requires Sandal Selection.	QTY	\$95

**Technical Specifications - Imperial**

Frame Sizing Specification	Size 1	Size 2	Size 3
Age (approximate)	1 1/2 - 5 years	4 - 10 years	9 - 18 years
User Body Height (Min / Max)	31.5" / 45.3"	43.3" / 57"	55.1" / 72.8"
Maximum User Weight	132 lbs	176 lbs	198 lbs
Chest Pad Height (Min / Max)	20.7" / 33.5"	26.8" / 43.3"	38.4" / 56.1"
Width Between Chest Supports (Min / Max)	6" / 10"	8" / 12"	10" / 14"
Activity Tray	16" x 21"	18" x 23"	20" x 25"

**Technical Specifications - Metric**

Frame Sizing Specification	Size 1	Size 2	Size 3
Age (approximate)	1 1/2 - 5 years	4 - 10 years	9 - 18 years
User Body Height (Min / Max)	800 mm / 1151 mm	1100 mm / 1448 mm	1400 mm / 1849 mm
Maximum User Weight	60 kgs	80 kgs	90 kgs
Chest Pad Height (Min / Max)	470 mm / 680 mm	680 mm / 975 mm	975 mm / 1350 mm
Width Between Chest Supports (Min / Max)	152 mm / 254 mm	203 mm / 305 mm	254 mm / 356 mm
Activity Tray	400 mm x 505 mm	450 mm x 575 mm	500 mm x 625 mm

**Sunrise Medical (US) LLC**  
 2842 Business Park Avenue  
 Fresno, CA 93727 USA  
 Tel.: 1-800-333-4000  
 Fax: 1-800-300-7502  
 E-Mail: webmaster@sunmed.com



**If you are visually impaired, this document can be viewed in PDF format at [www.SunriseMedical.com](http://www.SunriseMedical.com)**

