



Squiggles Seat

February 2022

Mark For:

Date: _____

Dealer Acct #: _____

Dealer: _____

Dealer Contact: _____

Dealer Address: _____

Dealer City: _____ ST: _____ ZIP: _____

Dealer Phone: _____ Fax: _____

Confirmation Email: _____

Confirm Via: Fax EmailSubmitting for: Quote Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____

Attention: _____

Address: _____

Address: _____

Ship To City: _____ ST: _____ ZIP: _____

Ship To Phone: _____ Fax: _____

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

Squiggles Seating System

Weight Capacity: 48.4 lbs

Seat Shell

Selection Required

Item Code	Description	HCPC	Retail
<input type="checkbox"/> 120-610S	Squiggles Seat Shell		
120-610-1	Seat Base with Depth Adjustment	E2294	\$1,021
120-610-2	Seat Back - Fixed	E2293	\$1,021
120-610-3	Seat Back Reclining Mechanism	E1014	\$446
120-610-4	Seat Base Hip Guides / Support	E0956	\$124

Seat Cover

Selection Required. Includes seat base and backrest cushion covers.

Item Code	Description	HCPC	Retail
<input type="checkbox"/> 120-710-03	Squiggles Seat Covers - Blue		N/C
<input type="checkbox"/> 120-710-01	Squiggles Seat Covers - Green		N/C
<input type="checkbox"/> 120-710-02	Squiggles Seat Covers - Orange		N/C
<input type="checkbox"/> 120-710-04	Squiggles Seat Covers - Pink		N/C

Base

Item Code	Description	HCPC	Retail
<input type="checkbox"/> 120-800	Size 1 Hi Low Chassis (Foot Pedal Operated.)	E1399	\$1,828
122-710	Hi Low Chassis Push Handle		\$266
<input type="checkbox"/> 102-2753	Size 2 Hi Low Chassis (Foot Pedal Operated. Includes Push Handle)	E1399	\$2,094
120-818	Adapter to Mount Squiggles Seat on Size 2 Hi Low Chassis		\$204
<input type="checkbox"/> 194-600	Mobility Base Interface Plate (Allows mounting Squiggles Seat onto Sunrise Wheelchairs)	K0108	\$480

Foot Positioning

Sandals and Footplate Raisers require Squiggles Seat Adjustable Footrest

Item Code	Description	HCPC	Retail
<input type="checkbox"/> 120-639	Squiggles Seat Adjustable Footrest	K0040	\$226
<input type="checkbox"/> 152-1600	Small Sandals w/ Straps (Pair)	E0951, E0952	\$254
<input type="checkbox"/> 152-2600	Medium Sandals w/ Straps (Pair)	E0951, E0952	\$277
<input type="checkbox"/> 120-816	Footplate Raiser 1" (Each)	E0951, E0952	\$86

Quantity: _____ (1 to 3)

Harness

	Item Code	Description	HCPC	Retail
<input type="checkbox"/>	120-844-03	Pelvic Harness Small - Blue	E0978	\$153
<input type="checkbox"/>	120-844-01	Pelvic Harness Small - Green	E0978	\$153
<input type="checkbox"/>	120-844-02	Pelvic Harness Small - Orange	E0978	\$153
<input type="checkbox"/>	120-844-04	Pelvic Harness Small - Pink	E0978	\$153
<input type="checkbox"/>	117-848-03	Pelvic Harness Medium - Blue	E0978	\$153
<input type="checkbox"/>	117-848-01	Pelvic Harness Medium - Green	E0978	\$153
<input type="checkbox"/>	117-848-02	Pelvic Harness Medium - Orange	E0978	\$153
<input type="checkbox"/>	117-848-04	Pelvic Harness Medium - Pink	E0978	\$153
<input type="checkbox"/>	120-789	Pelvic Harness Spacer Pads - Black (Pair)	K0108	\$86
<input type="checkbox"/>	120-799	Small Trunk Harness	E0960	\$275
<input type="checkbox"/>	120-757	Medium Trunk Harness	E0960	\$275

Lateral Support

Cover selection required with Lateral support hardware.

	Item Code	Description	HCPC	Retail
<input type="checkbox"/>	120-855	Flip-Away Lateral Support Hardware (Pair)		\$321
<input type="checkbox"/>	120-838-03	Lateral Support Covers - Blue (Pair)	E0956	\$147
<input type="checkbox"/>	120-838-01	Lateral Support Covers - Green (Pair)	E0956	\$147
<input type="checkbox"/>	120-838-02	Lateral Support Covers - Orange (Pair)	E0956	\$147
<input type="checkbox"/>	120-838-04	Lateral Support Covers - Pink (Pair)	E0956	\$147
<input type="checkbox"/>	120-656	Rigid Lateral Support Hardware (Pair)	E1028	\$195
<input type="checkbox"/>	120-838-03	Lateral Support Covers - Blue (Pair)	E0956	\$147
<input type="checkbox"/>	120-838-01	Lateral Support Covers - Green (Pair)	E0956	\$147
<input type="checkbox"/>	120-838-02	Lateral Support Covers - Orange (Pair)	E0956	\$147
<input type="checkbox"/>	120-838-04	Lateral Support Covers - Pink (Pair)	E0956	\$147
<input type="checkbox"/>	120-954	Laterals w/ Protection Pads - Black (Pair)	E0956	\$412
<input type="checkbox"/>	137-808-03	Chest Support - Blue <i>(Requires Laterals)</i>	E0960	\$180
<input type="checkbox"/>	137-808-01	Chest Support - Green <i>(Requires Laterals)</i>	E0960	\$180
<input type="checkbox"/>	137-808-02	Chest Support - Orange <i>(Requires Laterals)</i>	E0960	\$180
<input type="checkbox"/>	137-808-04	Chest Support - Pink <i>(Requires Laterals)</i>	E0960	\$180

Tray

	Item Code	Description	HCPC	Retail
<input type="checkbox"/>	120-672	Tray - Gray	E0950	\$382
<input type="checkbox"/>	120-688	Extra High Tray +2" Adjustment	E0950	\$483
<input type="checkbox"/>	120-270	Padded Tray Insert - Gray	K0108	\$316
<input type="checkbox"/>	117-769	Grab Rail / Toy Bar for Tray	K0108	\$284

Pommel

	Item Code	Description	HCPC	Retail
<input type="checkbox"/>	120-641	Pommel	E0957	\$108

Headrest

Cushion selection required with Headrest Support hardware

	Item Code	Description	HCPC	Retail
<input type="checkbox"/>	120-859	Contoured Headrest Support Hardware	E1028	\$358
<input type="checkbox"/>	137-800-03	Contoured Headrest Cushion - Blue	E0955	\$164
<input type="checkbox"/>	137-800-01	Contoured Headrest Cushion - Green	E0955	\$164
<input type="checkbox"/>	137-800-02	Contoured Headrest Cushion - Orange	E0955	\$164
<input type="checkbox"/>	137-800-04	Contoured Headrest Cushion - Pink	E0955	\$164
<input type="checkbox"/>	117-697	Flat Headrest Support Hardware	E1028	\$220
<input type="checkbox"/>	117-850-03	Flat Headrest Cushion - Blue	E0955	\$112
<input type="checkbox"/>	117-850-01	Flat Headrest Cushion - Green	E0955	\$112
<input type="checkbox"/>	117-850-02	Flat Headrest Cushion - Orange	E0955	\$112
<input type="checkbox"/>	117-850-04	Flat Headrest Cushion - Pink	E0955	\$112
<input type="checkbox"/>	118-667	Lateral Bracket Assembly Hardware (Pair) <i>Requires Flat Headrest Support Hardware and Cushion</i>	E1028	\$202
<input type="checkbox"/>	137-807-03	Flat Headrest Lateral Cushions - Blue (Pair)	E0956	\$183
<input type="checkbox"/>	137-807-01	Flat Headrest Lateral Cushions - Green (Pair)	E0956	\$183
<input type="checkbox"/>	137-807-02	Flat Headrest Lateral Cushions - Orange (Pair)	E0956	\$183
<input type="checkbox"/>	137-807-04	Flat Headrest Lateral Cushions - Pink (Pair)	E0956	\$183
<input type="checkbox"/>	117-745	Removable Head Support Bracket for 1/2" Box Section	E1028	\$68

Accessories

Cover selection required with Shoulder Support Hardware

	Item Code	Description	HCPC	Retail
<input type="checkbox"/>	120-779-03	1" Back Spacer Pad and Flip-Away Lateral Supports - Blue	E0956, E1297	\$774
<input type="checkbox"/>	120-779-01	1" Back Spacer Pad and Flip-Away Lateral Supports - Green	E0956, E1297	\$774
<input type="checkbox"/>	120-779-02	1" Back Spacer Pad and Flip-Away Lateral Supports - Orange	E0956, E1297	\$774
<input type="checkbox"/>	120-779-04	1" Back Spacer Pad and Flip-Away Lateral Supports - Pink	E0956, E1297	\$774
<input type="checkbox"/>	120-702	Shoulder Support Hardware	K0108	\$262
<input type="checkbox"/>	120-828-03	Shoulder Support Cover - Blue	K0108	\$149
<input type="checkbox"/>	120-828-01	Shoulder Support Cover - Green	K0108	\$149
<input type="checkbox"/>	120-828-02	Shoulder Support Cover - Orange	K0108	\$149
<input type="checkbox"/>	120-828-04	Shoulder Support Cover - Pink	K0108	\$149

Assembly

	Item Code	Description	Retail
<input type="checkbox"/>	ASSY-LKY	Product Assembly Fee	\$126

Technical Specifications

Specification	Imperial	Metric
Age (Approximate)	1 - 5 years	1 - 5 years
Stature	29.5" - 43.3"	750 mm - 1100 mm
Maximum User Weight	48.4 lbs	22 kg
Seat Depth	7.5" - 11.4"	190 mm - 290 mm
Seat Width	6.3" - 10.2"	160 mm - 260 mm
Hip Width	6.3" - 8.7"	160 mm - 220 mm
Backrest Height	12.6" - 15"	320 mm - 380 mm
Backrest Height with Shoulder Support	16.5" - 18.9"	420 mm - 480 mm
Minimum Chest Width	6.3"	160 mm
Maximum Knee Width	3.5"	90 mm
Minimum Distance Between Pommel and Side of Seat	7.5"	190 mm
Minimum Seat to Footplate	4.8"	122.5 mm
Seat to Floor	11.4" - 26.4"	290 mm - 670 mm
Backrest Angle (Prone / Recline)	10° / 25°	10° / 25°
Footplate Angle (Plantar / Dorsiflexion)	10° / 10°	10° / 10°
Tray Size	18.1" x 22"	460 mm x 560 mm
Seat Weight	13.2 lbs	6 kgs
Tray Height	4.5" - 6.7"	115 mm x 170 mm
Squiggles Seat Hi-Lo Chassis Foot Operated (Excluding Push Bar)	Imperial	Metric
Weight	19.8 lbs	9 kgs
Length	25.2"	640 mm
Width	21.3"	540 mm
Minimum Dimensions for Storage	12.4"	315 mm



Sunrise Medical (US) LLC 2842 Business Park Ave. • Fresno, CA 93727 • USA