



ORDER FORM

# Squiggles Seat

# EILK03

Ultimate combination of postural seating and everyday comfort for young children while still maintaining a fun, colorful design.

Customer			
Account			
Account #		Date	
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #	
Buyer			
Mark For			

Shipping Address			
Name			
Address			
City		State	
ZIP		Phone	
Email / Fax			

The HCPCS CODES herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative HCPCS CODES, including misc. codes to ensure access for their enrollees. The use of HCPCS CODES does not ensure coverage or payment.

## Squiggles Seat

Weight Capacity: 48.4 lbs

✓ OPTION	DESCRIPTION	HCPCS	RETAIL
<b>Seat Shell</b>			<b>\$</b>
Selection Required.			
<input type="checkbox"/>	151FT72 Squiggles Seat Shell		
	151BS28 Seat Base with Depth Adjustment	E2294	\$1,097
	151BF38 Seat Back - Fixed	E2293	\$1,097
	151BF39 Seat Back Reclining Mechanism	E1014	\$479
	151GU5 Seat Base Hip Guides / Support	E0956	\$133
<b>Seat Cover</b>			<b>\$</b>
Selection Required. Includes seat base and backrest cushion covers.			
<input type="checkbox"/>	151YJ14 Squiggles Seat Covers - Blue		N/C
<input type="checkbox"/>	151YJ15 Squiggles Seat Covers - Green		N/C
<input type="checkbox"/>	151YJ17 Squiggles Seat Covers - Pink		N/C
<b>Base</b>			<b>\$</b>
<input type="checkbox"/>	151NB116 Size 1 Hi Low Chassis (Foot Pedal Operated.)	E1399	\$1,964
	151PH35 Hi Low Chassis Push Handle		\$286
<input type="checkbox"/>	151NB114 Size 2 Hi Low Chassis (Foot Pedal Operated. Includes Push Handle)	E1399	\$2,251
	151JI04 Adapter to Mount Squiggles Seat on Size 2 Hi Low Chassis		\$219
<input type="checkbox"/>	151NB117 Mobility Base Interface Plate Allows mounting Squiggles Seat onto Sunrise Wheelchairs	K0108	\$516
<b>Foot Positioning</b>			<b>\$</b>
Sandals and Footplate Raisers require Squiggles Seat Adjustable Footrest			
<input type="checkbox"/>	151DI67 Squiggles Seat Adjustable Footrest	K0040	\$243
<input type="checkbox"/>	1518A107 Small Sandals w/ Straps (Pair)	E0951, E0952	\$273
<input type="checkbox"/>	1518A108 Medium Sandals w/ Straps (Pair)	E0951, E0952	\$298
<input type="checkbox"/>	1519A24 <input type="checkbox"/> QTY Footplate Raiser 1" (Each)	E0951, E0952	\$92
<b>Harness</b>			<b>\$</b>
<input type="checkbox"/>	151NH109 Pelvic Harness Small - Blue	E0978	\$165
<input type="checkbox"/>	151NH110 Pelvic Harness Small - Green	E0978	\$165
<input type="checkbox"/>	151NH111 Pelvic Harness Small - Orange	E0978	\$165
<input type="checkbox"/>	151NH112 Pelvic Harness Small - Pink	E0978	\$165
<input type="checkbox"/>	151NH113 Pelvic Harness Medium - Blue	E0978	\$165
<input type="checkbox"/>	151NH114 Pelvic Harness Medium - Green	E0978	\$165
<input type="checkbox"/>	151NH115 Pelvic Harness Medium - Pink	E0978	\$165
<input type="checkbox"/>	151NH108 Pelvic Harness Spacer Pads - Black (Pair)	K0108	\$92
<input type="checkbox"/>	151NH105 Butterfly Harness Small	E0960	\$406
<input type="checkbox"/>	151NH106 Butterfly Harness Medium	E0960	\$406

## Squiggles Seat

✓ OPTION	DESCRIPTION	HCPCS	RETAIL
Lateral Support			\$
Cover selection required with Lateral Support Hardware			
<input type="checkbox"/> 151AW26	<b>Flip-Away Lateral Support Hardware (Pair)</b>		\$345
<input type="checkbox"/>	151IU19 Lateral Support Covers - Blue (Pair)	E0956	\$158
<input type="checkbox"/>	151IU20 Lateral Support Covers - Green (Pair)	E0956	\$158
<input type="checkbox"/>	151IU21 Lateral Support Covers - Orange (Pair)	E0956	\$158
<input type="checkbox"/>	151IU22 Lateral Support Covers - Pink (Pair)	E0956	\$158
<input type="checkbox"/> 151AW27	<b>Rigid Lateral Support Hardware (Pair)</b>	E1028	\$210
<input type="checkbox"/>	151IU19 Lateral Support Covers - Blue (Pair)	E0956	\$158
<input type="checkbox"/>	151IU20 Lateral Support Covers - Green (Pair)	E0956	\$158
<input type="checkbox"/>	151IU21 Lateral Support Covers - Orange (Pair)	E0956	\$158
<input type="checkbox"/>	151IU22 Lateral Support Covers - Pink (Pair)	E0956	\$158
<input type="checkbox"/> 151KE54	<b>Laterals w/ Protraction Pads - Black (Pair)</b>	E0956	\$443
<input type="checkbox"/>	151ZY25 Chest Support - Black (Requires Laterals)	E0960	\$193
Tray			\$
<input type="checkbox"/> 151A230	Tray - Gray	E0950	\$410
<input type="checkbox"/> 151A231	Extra High Tray +2" Adjustment	E0950	\$519
<input type="checkbox"/> 151A216	Padded Tray Insert - Black	K0108	\$386
<input type="checkbox"/> 1518A106	Grab Rail / Toy Bar for Tray	K0108	\$305
Pommel			\$
<input type="checkbox"/> 1518A118	Pommel	E0957	\$116

## Squiggles Seat

✓ OPTION	DESCRIPTION	HCPCS	RETAIL	
Headrest			\$	
Cushion selection required with hardware.				
Lateral Bracket Assembly Requires Flat Headrest Support Hardware and Cushion.				
<input type="checkbox"/>	<b>151AW24</b>	<b>Contoured Headrest Support Hardware</b>	E1028	\$385
<input type="checkbox"/>	151Y040	Contoured Headrest Cushion - Blue	E0955	\$176
<input type="checkbox"/>	151Y042	Contoured Headrest Cushion - Green	E0955	\$176
<input type="checkbox"/>	151Y057	Contoured Headrest Cushion - Orange	E0955	\$176
<input type="checkbox"/>	151Y043	Contoured Headrest Cushion - Pink	E0955	\$176
<input type="checkbox"/>	<b>151AW22</b>	<b>Flat Headrest Support Hardware</b>	E1028	\$236
<input type="checkbox"/>	151Y048	Flat Headrest Cushion - Blue	E0955	\$121
<input type="checkbox"/>	151Y050	Flat Headrest Cushion - Green	E0955	\$121
<input type="checkbox"/>	151Y051	Flat Headrest Cushion - Pink	E0955	\$121
<input type="checkbox"/>	<b>151AW23</b>	<b>Lateral Bracket Assembly Hardware (Pair)</b>	E1028	\$217
		Requires Flat Headrest Support Hardware and Cushion		
<input type="checkbox"/>	151Y048	Flat Headrest Lateral Cushions - Blue (Pair)	E0956	\$121
<input type="checkbox"/>	151Y050	Flat Headrest Lateral Cushions - Green (Pair)	E0956	\$121
<input type="checkbox"/>	151Y051	Flat Headrest Lateral Cushions - Pink (Pair)	E0956	\$121
<input type="checkbox"/>	151JG16	Removable Head Support Bracket for 1/2" Box Section	E1028	\$73
Accessories			\$	
<input type="checkbox"/>	151KE55	1" Back Spacer Pad and Flip-Away Lateral Supports - Blue	E0956, E1297	\$832
<input type="checkbox"/>	151KE56	1" Back Spacer Pad and Flip-Away Lateral Supports - Green	E0956, E1297	\$832
<input type="checkbox"/>	151KE57	1" Back Spacer Pad and Flip-Away Lateral Supports - Pink	E0956, E1297	\$832
<input type="checkbox"/>	<b>151AW28</b>	<b>Shoulder Support Hardware</b>	K0108	\$281
<input type="checkbox"/>	151SS13	Shoulder Support Cover - Blue	K0108	\$160
<input type="checkbox"/>	151SS14	Shoulder Support Cover - Green	K0108	\$160
<input type="checkbox"/>	151SS15	Shoulder Support Cover - Pink	K0108	\$160

## Technical Specifications

Frame Sizing Specification	Imperial	Metric
Age (Approximate)	1 - 5 years	1 - 5 years
Stature	29.5" - 43.3"	750 mm - 1100 mm
Maximum User Weight	48.4 lbs	22 kg
Seat Depth	7.5" - 11.4"	190 mm - 290 mm
Seat Width	6.3" - 10.2"	160 mm - 260 mm
Hip Width	6.3" - 8.7"	160 mm - 220 mm
Backrest Height	12.6" - 15"	320 mm - 380 mm
Backrest Height with Shoulder Support	16.5" - 18.9"	420 mm - 480 mm
Minimum Chest Width	6.3"	160 mm
Maximum Knee Width	3.5"	90 mm
Minimum Distance Between Pommel and Side of Seat	7.5"	190 mm
Minimum Seat to Footplate	4.8"	122.5 mm
Seat to Floor	11.4" - 26.4"	290 mm - 670 mm
Backrest Angle (Prone / Recline)	10° / 25°	10° / 25°
Footplate Angle (Plantar / Dorsiflexion)	10° / 10°	10° / 10°
Tray Size	18.1" x 22"	460 mm x 560 mm
Seat Weight	13.2 lbs	6 kgs
Tray Height	4.5" - 6.7"	115 mm x 170 mm
Squiggles Seat Hi-Lo Chassis Foot Operated (Excluding Push Bar)	Imperial	Metric
Weight	19.8 lbs	9 kgs
Length	25.2"	640 mm
Width	21.3"	540 mm
Minimum Dimensions for Storage	12.4"	315 mm

Sunrise Medical (US) LLC  
2842 N. Business Park Avenue  
Fresno, CA 93727 USA



If you are visually impaired, this document can be viewed in PDF format at [www.SunriseMedical.com](http://www.SunriseMedical.com)



MK-129856 Rev. 1 © 06/2024

