



ORDER FORM

Squiggles Seat

EILK03

Ultimate combination of postural seating and everyday comfort for young children while still maintaining a fun, colorful design.

Customer

Account			
Account #		Date	
☐ Quote	☐ Order	PO #	
Buyer			
Mark For			

Shipping Address

Name			
Address			
City		State	
ZIP		Phone	
Email / Fax			

The HCPCS CODES herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative HCPCS CODES, including misc. codes to ensure access for their enrollees. The use of HCPCS CODES does not ensure coverage or payment.

Squiggles Seat

Weight Capacity: 48.4 lbs

✓	OPTION	DESCRIPTION	HCPCS	RETAIL
Seat Shell				\$
Selection Required.				
☐	151FT72	Squiggles Seat Shell		
	151BS28	Seat Base with Depth Adjustment	E2294	\$1,119
	151BF38	Seat Back - Fixed	E2293	\$1,119
	151BF39	Seat Back Reclining Mechanism	E1014	\$489
	151GU5	Seat Base Hip Guides / Support	E0956	\$136
Seat Cover				\$
Selection Required. Includes seat base and backrest cushion covers.				
☐	151YJ14	Squiggles Seat Covers - Blue		N/C
☐	151YJ15	Squiggles Seat Covers - Green		N/C
☐	151YJ17	Squiggles Seat Covers - Pink		N/C
Base				\$
☐	151NB116	Size 1 Hi Low Chassis (Foot Pedal Operated.)	E1399	\$2,003
	151PH35	Hi Low Chassis Push Handle		\$292
☐	151NB114	Size 2 Hi Low Chassis (Foot Pedal Operated. Includes Push Handle)	E1399	\$2,296
	151JI04	Adapter to Mount Squiggles Seat on Size 2 Hi Low Chassis		\$223
☐	151NB117	Mobility Base Interface Plate Allows mounting Squiggles Seat onto Sunrise Wheelchairs	K0108	\$526
Foot Positioning				\$
Sandals and Footplate Raisers require Squiggles Seat Adjustable Footrest				
☐	151DI67	Squiggles Seat Adjustable Footrest	K0040	\$248
☐	1518A107	Small Sandals w/ Straps (Pair)	E0951, E0952	\$278
☐	1518A108	Medium Sandals w/ Straps (Pair)	E0951, E0952	\$304
☐	1519A24	QTY Footplate Raiser 1" (Each)	E0951, E0952	\$94
Harness				\$
☐	151NH109	Pelvic Harness Small - Blue	E0978	\$168
☐	151NH110	Pelvic Harness Small - Green	E0978	\$168
☐	151NH112	Pelvic Harness Small - Pink	E0978	\$168
☐	151NH113	Pelvic Harness Medium - Blue	E0978	\$168
☐	151NH114	Pelvic Harness Medium - Green	E0978	\$168
☐	151NH115	Pelvic Harness Medium - Pink	E0978	\$168
☐	151NH108	Pelvic Harness Spacer Pads - Black (Pair)	K0108	\$94
☐	151NH105	Butterfly Harness Small	E0960	\$414
☐	151NH106	Butterfly Harness Medium	E0960	\$414

Squiggles Seat

✓ OPTION	DESCRIPTION	HCPCS	RETAIL
----------	-------------	-------	--------

Lateral Support			\$
-----------------	--	--	----

Cover selection required with Lateral Support Hardware

☐ 151AW26	Flip-Away Lateral Support Hardware (Pair)		\$352
☐ 151IU19	Lateral Support Covers - Blue (Pair)	E0956	\$161
☐ 151IU20	Lateral Support Covers - Green (Pair)	E0956	\$161
☐ 151IU22	Lateral Support Covers - Pink (Pair)	E0956	\$161
☐ 151AW27	Rigid Lateral Support Hardware (Pair)	E1028	\$214
☐ 151IU19	Lateral Support Covers - Blue (Pair)	E0956	\$161
☐ 151IU20	Lateral Support Covers - Green (Pair)	E0956	\$161
☐ 151IU22	Lateral Support Covers - Pink (Pair)	E0956	\$161
☐ 151KE54	Laterals w/ Protraction Pads - Black (Pair)	E0956	\$452
☐ 151ZY25	Chest Support - Black (Requires Laterals)	E0960	\$197

Tray			\$
------	--	--	----

☐ 151A230	Tray - Gray	E0950	\$418
☐ 151A231	Extra High Tray +2" Adjustment	E0950	\$529
☐ 151A216	Padded Tray Insert - Black	K0108	\$394
☐ 1518A106	Grab Rail / Toy Bar for Tray	K0108	\$311

Pommel			\$
--------	--	--	----

☐ 1518A118	Pommel	E0957	\$118
-----------------------------------	--------	-------	-------

Headrest			\$
----------	--	--	----

Cushion selection required with hardware.

Lateral Bracket Assembly Requires Flat Headrest Support Hardware and Cushion.

☐ 151AW24	Contoured Headrest Support Hardware	E1028	\$393
☐ 151Y040	Contoured Headrest Cushion - Blue	E0955	\$180
☐ 151Y042	Contoured Headrest Cushion - Green	E0955	\$180
☐ 151Y043	Contoured Headrest Cushion - Pink	E0955	\$180
☐ 151AW22	Flat Headrest Support Hardware	E1028	\$241
☐ 151Y044	Flat Headrest Cushion - Blue	E0955	\$121
☐ 151Y046	Flat Headrest Cushion - Green	E0955	\$121
☐ 151Y047	Flat Headrest Cushion - Pink	E0955	\$121
☐ 151AW23	Lateral Bracket Assembly Hardware (Pair) Requires Flat Headrest Support Hardware and Cushion	E1028	\$221
☐ 151Y048	Flat Headrest Lateral Cushions - Blue (Pair)	E0956	\$123
☐ 151Y050	Flat Headrest Lateral Cushions - Green (Pair)	E0956	\$123
☐ 151Y051	Flat Headrest Lateral Cushions - Pink (Pair)	E0956	\$123
☐ 151JG16	Removable Head Support Bracket for 1/2" Box Section	E1028	\$74

Accessories			\$
-------------	--	--	----

☐ 151KE55	1" Back Spacer Pad and Flip-Away Lateral Supports - Blue	E0956, E1297	\$849
☐ 151KE56	1" Back Spacer Pad and Flip-Away Lateral Supports - Green	E0956, E1297	\$849
☐ 151KE57	1" Back Spacer Pad and Flip-Away Lateral Supports - Pink	E0956, E1297	\$849
☐ 151AW28	Shoulder Support Hardware	K0108	\$287
☐ 151SS13	Shoulder Support Cover - Blue	K0108	\$163
☐ 151SS14	Shoulder Support Cover - Green	K0108	\$163
☐ 151SS15	Shoulder Support Cover - Pink	K0108	\$163

Technical Specifications		
Frame Sizing Specification	Imperial	Metric
Age (Approximate)	1 - 5 years	1 - 5 years
Stature	29.5" - 43.3"	750 mm - 1100 mm
Maximum User Weight	48.4 lbs	22 kg
Seat Depth	7.5" - 11.4"	190 mm - 290 mm
Seat Width	6.3" - 10.2"	160 mm - 260 mm
Hip Width	6.3" - 8.7"	160 mm - 220 mm
Backrest Height	12.6" - 15"	320 mm - 380 mm
Backrest Height with Shoulder Support	16.5" - 18.9"	420 mm - 480 mm
Minimum Chest Width	6.3"	160 mm
Maximum Knee Width	3.5"	90 mm
Minimum Distance Between Pommel and Side of Seat	7.5"	190 mm
Minimum Seat to Footplate	4.8"	122.5 mm
Seat to Floor	11.4" - 26.4"	290 mm - 670 mm
Backrest Angle (Prone / Recline)	10° / 25°	10° / 25°
Footplate Angle (Plantar / Dorsiflexion)	10° / 10°	10° / 10°
Tray Size	18.1" x 22"	460 mm x 560 mm
Seat Weight	13.2 lbs	6 kgs
Tray Height	4.5" - 6.7"	115 mm x 170 mm
Squiggles Seat Hi-Lo Chassis Foot Operated (Excluding Push Bar)	Imperial	Metric
Weight	19.8 lbs	9 kgs
Length	25.2"	640 mm
Width	21.3"	540 mm
Minimum Dimensions for Storage	12.4"	315 mm

Sunrise Medical (US) LLC
12002 Volunteer Blvd
Mount Juliet, TN 37122



If you are visually impaired, this document
can be viewed in PDF format at
www.SunriseMedical.com



MK-129856 Rev. L © 09/2025

