



# Kids Quickie 2 SE

## January 2018

### EIZ10

Date: _____	<b>SUBMITTING FOR:</b> <input type="checkbox"/> Quote <input type="checkbox"/> Order
Dealer Acct #: _____	
Dealer: _____	PO#: _____
Dealer Contact: _____	Ship To: _____
Dealer Address: _____	Attention: _____
Dealer City: _____ ST: _____ ZIP: _____	Address: _____
Dealer Phone: (____) _____ Fax: (____) _____	Address: _____
Confirm Via: <input type="checkbox"/> Fax <input type="checkbox"/> Email	Ship To City: _____ ST: _____ ZIP: _____
Mark For: _____	Ship To Phone: (____) _____ Fax: (____) _____

*The HCPCS CODES provided are based on code verification through PDAC or our interpretation of the code definitions.  
The use of HCPCS CODES does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.*

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
---	----------	-------	------------	-------	---	----------	-------	------------	-------

#### KIDS QUICKIE 2 SE

*K0005 coding still applies for EIZ10 with certain sizes - Standard Frame - 250 lbs. Maximum User Weight - Uses Rear Side Frame Depth Adj. & Depth Adj. Back*

<input type="checkbox"/>	142WM07	E1236	Kids Quickie 2 SE w/ Seating	\$2,395
--------------------------	---------	-------	------------------------------	---------

*Note: HCPCS Codes shown as submitted. Some sizes fit other HCPCS codes i.e. 15"-16" widths = K0005*

<input type="checkbox"/>	142WM08	E1238	Kids Quickie 2 SE w/o Seating	\$2,525
--------------------------	---------	-------	-------------------------------	---------

#### Select - Frame Type

<input type="checkbox"/>	142RS1	Rear Side Straight	NC
--------------------------	--------	--------------------	----

#### Transit System

*Incl. Stabilizer Bar, must order Push Handle. N/A with padded swing-away armrests or stroller handles.*

<input type="checkbox"/>	142TO1	K0108	Transit System	\$315.00
--------------------------	--------	-------	----------------	----------

#### Frame Width

<input type="checkbox"/>	142W2		11" (By Upholstery)	NC
<input type="checkbox"/>	142W3		12"	NC
<input type="checkbox"/>	142W4		13"	NC
<input type="checkbox"/>	142W5		14"	NC
<input type="checkbox"/>	142W37		15"	\$250
<input type="checkbox"/>	142W40		16"	\$250
<input type="checkbox"/>	142W38		17"	\$250
<input type="checkbox"/>	142W39		18"	\$250

#### Select - Frame Type

<input type="checkbox"/>	142FT54	Cross Brace with Seat Rails	NC
<input type="checkbox"/>	142FT55	Cross Brace w/o Seat Rails	NC

#### Frame Length

<input type="checkbox"/>	142FL1	Kids (12"-14")	STD
<input type="checkbox"/>	142FL9	Regular (13"-16")	STD
<input type="checkbox"/>	142FL18	Long (15"-18")	STD

#### Seat Rail Depth (Sling Not Included)

<input type="checkbox"/>	142SN05		10"	NC
<input type="checkbox"/>	142SN06		11"	NC
<input type="checkbox"/>	142SN07		12"	NC
<input type="checkbox"/>	142SN08		13"	NC
<input type="checkbox"/>	142SN09		14"	NC
<input type="checkbox"/>	142SN10		15"	\$250
<input type="checkbox"/>	142SN11		16"	\$250
<input type="checkbox"/>	142SN12		17"	\$250
<input type="checkbox"/>	142SN13		18"	\$250

#### FRAME COLOR

#### Standard Colors

<input type="checkbox"/>	142C140	Aztec Gold	NC
<input type="checkbox"/>	142C2	Black (N/A as base w/ Zebra)	NC
<input type="checkbox"/>	142C141	Black Cherry	NC
<input type="checkbox"/>	142C10	Black Opal	NC
<input type="checkbox"/>	142C38	Blue Opal	NC
<input type="checkbox"/>	142C36	Candy Blue	NC
<input type="checkbox"/>	142C14	Candy Purple	NC
<input type="checkbox"/>	142C35	Candy Red	NC
<input type="checkbox"/>	142C142	Electric Blue	NC
<input type="checkbox"/>	142C143	Evergreen	NC
<input type="checkbox"/>	142C126	Glow	NC
<input type="checkbox"/>	142C144	Green Apple	NC
<input type="checkbox"/>	142C147	Hot Sparkle Pink	NC
<input type="checkbox"/>	142C148	Mauve Pink	NC
<input type="checkbox"/>	142C103	Rootbeer	NC
<input type="checkbox"/>	142C114	Sparkle Silver	NC

#### Standard Colors - Cont.

<input type="checkbox"/>	142C145	Sunrise Orange	NC
<input type="checkbox"/>	142C113	Titanium Color Paint	NC
<input type="checkbox"/>	142C8	Yellow	NC

#### Standard Matte Colors

<input type="checkbox"/>	142C116	Matte Black	NC
<input type="checkbox"/>	142C151	Matte Black Cherry	NC
<input type="checkbox"/>	142C152	Matte Purple	NC
<input type="checkbox"/>	142C150	Matte Electric Blue	NC
<input type="checkbox"/>	142C149	Matte Evergreen	NC

#### Custom Colors

<input type="checkbox"/>	142C43	Mossy Oak Camouflage	\$295
<input type="checkbox"/>	142C42	Stars and Stripes	\$295
<input type="checkbox"/>	142C108	Zebra (Choose base color above: White Std, N/A Black)	\$295

**SEATING**

**JAY Cushions (Optional)**

To order any JAY cushions on-chair, please see the JAY price list. Cushion name or part number must be filled out with desired width and depth

<input type="checkbox"/>	142JC20	JAY Bundled Cushion	Price List
		Cushion Name or Part #	
		Width (in.):	Depth (in.):

**Seat**

<input type="checkbox"/>	142SL1	Sling	STD
<input type="checkbox"/>	142SL16	Omit Sling, Include Rails	NC

**Cushion**

<input type="checkbox"/>	142CU1	K0669	2" Cushion	\$120
--------------------------	--------	-------	------------	-------

**JAY Backs (Optional)**

Ordering JAY backs on-chair, please see the JAY price list. Requires part #

<input type="checkbox"/>	142JB02	JAY Bundled Backrest	Price List
		Backrest Name or Part #	
		Width (in.):	Height (in.):

**FRAME DEPTH**

**Frame Depth - With Seat Rails**

<input type="checkbox"/>	142BL14	1" offset from seat rail	STD
<input type="checkbox"/>	142BL15	2" offset from seat rail	NC

**Frame Depth - Without Seat Rails**

<input type="checkbox"/>	142FD1	10"	NC
<input type="checkbox"/>	142FD2	11"	NC
<input type="checkbox"/>	142FD3	12"	NC
<input type="checkbox"/>	142FD4	13"	NC
<input type="checkbox"/>	142FD5	14"	NC
<input type="checkbox"/>	142FD6	15"	NC
<input type="checkbox"/>	142FD7	16"	NC
<input type="checkbox"/>	142FD8	17"	NC
<input type="checkbox"/>	142FD9	18"	NC
<input type="checkbox"/>	142FD10	19"	NC
<input type="checkbox"/>	142FD11	20"	NC

**BACKREST**

**Back Height**

Standard With Push Handles

<input type="checkbox"/>	142B1	Low (8.5" - 12")	NC
<input type="checkbox"/>	142B5	Medium (12" - 15.5")	NC
<input type="checkbox"/>	142B9	Tall (15.5" - 19")	NC

**Stroller Handles**

<input type="checkbox"/>	142BP27	Depth Adjustable Standard	NC
<input type="checkbox"/>	142BP28	Depth Adjustable Med	NC
<input type="checkbox"/>	142BP29	Depth Adjustable Tall	NC
<input type="checkbox"/>	142BP33	Swing Away Adjustable	\$100

**Upholstery**

<input type="checkbox"/>	142BU1	Standard	STD
<input type="checkbox"/>	142BU8	Adjustable Upholstery	\$150
<input type="checkbox"/>	142BP2	Omit Upholstery	NC

**Extension**

<input type="checkbox"/>		Removable	\$100
--------------------------	--	-----------	-------

**Backrest Options**

<input type="checkbox"/>	142BP1	Push Handle	NC
<input type="checkbox"/>	142BP3	Push Handle 8° (Med or Tall Only)	\$25
<input type="checkbox"/>	142BP2	Omit Push Handle - Straight N/A w/ transit	NC

**ARMRESTS**

**Armrest**

<input type="checkbox"/>	142AA8	E0973	Height Adj. Adult	\$285.00
<input type="checkbox"/>	142AM1		Desk Length (10")	NC
<input type="checkbox"/>	142AM2		Full Length (14")	NC
<input type="checkbox"/>	142AA9	E0973	Height Adj. Kids	\$285.00
<input type="checkbox"/>	142AM1		Desk Length (10")	NC
<input type="checkbox"/>	142AM2		Full Length (14")	NC
<input type="checkbox"/>	142AA14		Omit Armrest	NC

**Armrest**

<input type="checkbox"/>	142AA16	E0973	Padded Swing Away Adult N/A w/ amputee axle plate or transit	NC
<input type="checkbox"/>	142AA19	E0973	Padded Swing Away Medium N/A w/ amputee axle plate or transit	NC
<input type="checkbox"/>	142AA17	E0973	Padded Swing Away Kids N/A w/ amputee axle plate or transit	STD
<input type="checkbox"/>	142AA6	E0973	Flip Back Height Adjustable (N/A w/ Kids Frame)	\$255.00
<input type="checkbox"/>	142AM1		Desk Length (10")	NC
<input type="checkbox"/>	142AM2		Full Length (14")	NC

**FOOTRESTS**

Choose 1 or 2

1. Hanger Type (Mount With Extension Tube)				2. Hanger Type (Direct Mount)				
<b>Hanger</b>				<b>Hanger</b>				
<input type="checkbox"/>	142H1		60° Swing Away Hanger	NC	<input type="checkbox"/>	142H7	70° Swing Away Hangers Front Mount (Only avail. w/ Kids or Adult Angle Adj. or 90° footboard footplates)	NC
<input type="checkbox"/>	142H2		70° Swing Away Hanger	STD	<input type="checkbox"/>	142H8	80° Swing Away Hangers Front Mount (Only avail. w/ Kids or Adult Angle Adj. or 90° footboard footplates)	NC
<input type="checkbox"/>	142H14	E0990	Elevating Hanger (Only avail. w/ alum. Footplates or angle adj. Available on 14"-18" frame widths. N/A w/ locking angle adj. footplates or short ext. tubes.)	\$295	<input type="checkbox"/>	142H6	90° Swing Away Hangers (Only avail. w/ Kids or Adult Angle Adj. or 90° footboard footplates)	NC
<b>Extension Tube (Required for Option 1)</b>				<b>Extensions (Optional)</b>				
<input type="checkbox"/>	142E1		Short (see matrix for length)	NC	<input type="checkbox"/>	142E23	K0040 2" Extensions (N/A w/90° hangers & elevating leg rests)	\$55
<input type="checkbox"/>	142E2		Medium (see matrix for length)	NC	<input type="checkbox"/>	142E24	K0040 4" Extensions (N/A w/90° hangers & elevating leg rests)	\$55
<input type="checkbox"/>	142E5		Long (see matrix for length)	NC	<b>Footplate Type</b>			
<input type="checkbox"/>	142E12		Short ELR (ELR only)	NC	<input type="checkbox"/>	142F9	K0040 Kids Angle Adjustable Footplate (10" - 16" W)	\$185
<input type="checkbox"/>	142E15		Long ELR (ELR only)	NC	<input type="checkbox"/>	142F8	K0040 Adult Angle Adjustable Footplate (Available on 12" - 16" W)	\$185
<input type="checkbox"/>	142E16		High Mount (N/A elevating leg rests. Only avail. w/ adult or kids angle adj. footplates.)	\$55	<b>Footplate</b>			
<input type="checkbox"/>	142F1	E0951	Composite Footplate (12"-16"w)	STD	<input type="checkbox"/>	142F4	K0040 Kids Angle Adj. Footplate	\$185
<input type="checkbox"/>	142F4	K0040	Kids Angle Adj. Footplate	\$185	<input type="checkbox"/>	142f3	K0040 Adult Angle Adj. Footplate (12"-16"w)	\$185
<input type="checkbox"/>	142f3	K0040	Adult Angle Adj. Footplate (12"-16"w)	\$185	<input type="checkbox"/>	142F48	K0040 Platform Flip Up Footplate (N/A with short ext. tubes. Only avail. with 60° or 70°)	\$190
<input type="checkbox"/>	142F48	K0040	Platform Flip Up Footplate (N/A with short ext. tubes. Only avail. with 60° or 70°)	\$190	<input type="checkbox"/>	142F61	K0040 Locking Angle Adjustable. (N/A Below 15" W - Requires Heel Loops)	\$190
<input type="checkbox"/>	142F61	K0040	Locking Angle Adjustable. (N/A Below 15" W - Requires Heel Loops)	\$190	<input type="checkbox"/>	142F5	Aluminum Footplate Available w/ ELRs	NC
<input type="checkbox"/>	142F5		Aluminum Footplate Available w/ ELRs	NC	<b>Footrest Options</b>			
<input type="checkbox"/>	142HL9	E0951	Heel Loops	\$50	<input type="checkbox"/>	142HL8	K0038 Leg Straps	\$35
<input type="checkbox"/>	142HL8	K0038	Leg Straps	\$35	<input type="checkbox"/>	142A54	K0038 Extra Leg Strap	\$30
<input type="checkbox"/>	142A54	K0038	Extra Leg Strap	\$30	<b>Impact Guard</b>			
<input type="checkbox"/>	142IG2		Neoprene (Only avail. with 60° or 70°)	\$85				

CHOOSE 1 OR 2

**CASTERS**

Caster Wheels				Caster Wheels - Continued				
<input type="checkbox"/>	142CW3		4" Polyurethane	NC	<input type="checkbox"/>	142CW39	8" x 1.5" Semi-Pneumatic (only avail. for use in bottom fork holes)	\$85
<input type="checkbox"/>	142CW31		4" Micro-Lighted	\$70	<input type="checkbox"/>	142CW36	E2214 E2215 8" x 2" Pneumatic (N/A w/ 90° elevating leg rests. N/A w/kids frame length. Requires selection of Caster Tire Tube/Insert)	\$145
<input type="checkbox"/>	142CW5		5" Low Profile Polyurethane	NC	<input type="checkbox"/>	142CW12	Omit Caster Wheels	NC
<input type="checkbox"/>	142CW7	E2219	5" x 1.5" Semi-Pneumatic (Only avail. for use in bottom fork holes)	\$85	<b>Caster Tire Tube/Insert</b>			
<input type="checkbox"/>	142CW6	E2219	5" Soft Roll Semi-Pneumatic	\$90	Requires a selection of 6", 8" or 8" x 2" Pneumatic Tire.			
<input type="checkbox"/>	142CW33	E2214 E2215	6" Pneumatic (Requires selection of Caster Tire Tube/Insert)	\$75	<input type="checkbox"/>	142CZ2	K0064 Airless Insert (Only avail. With 8"x2" Pneumatic)	\$65
<input type="checkbox"/>	142CW34		6" Polyurethane	STD	<input type="checkbox"/>	142CZ1	Caster Tube (all sizes)	\$20
<input type="checkbox"/>	142CW35	E2219	6" x 1.5" Semi-Pneumatic (only avail. for use in bottom fork holes)	\$85	<b>Caster Options</b>			
<input type="checkbox"/>	142CW21	E2214 E2215	8" Pneumatic (N/A w/ elevating leg rests. N/A w/kids frame length. Requires selection of Caster Tire Tube/Insert)	\$85	<input type="checkbox"/>	142CS5	3/4" Longer Fork Stem Bolt	NC
<input type="checkbox"/>	142CW32		8" Polyurethane (N/A w/ elevating leg rests or w/kids frame length.)	NC	<input type="checkbox"/>	142CS10	1-1/2" Longer Fork Stem Bolt	NC
					<input type="checkbox"/>	142CS3	K0073 Caster Pin Locks	\$195

**AXLE AND WHEELS**

**Axle Plate**

<input type="checkbox"/>	142AL2		Standard Adjustable	NC
<input type="checkbox"/>	142AL71	E0959	Amputee Adjustable	\$85

**Axle**

<input type="checkbox"/>	142AX1		Quick Release	STD
<input type="checkbox"/>	142AX11		Quad Quick Release	\$125

**Anti-Tip**

<input type="checkbox"/>	142AT1	E0971	Rear Anti-Tip	\$140
--------------------------	--------	-------	---------------	-------

**Rear Wheel Type**

<input type="checkbox"/>	142RW3		20" Mag	NC
<input type="checkbox"/>	142RW74		20" Lite Spoke	NC
<input type="checkbox"/>	142RW8		22" Mag	NC
<input type="checkbox"/>	142RW75		22" Lite Spoke	NC
<input type="checkbox"/>	142RW1		24" Mag	NC
<input type="checkbox"/>	142RW54		24" Spoke	\$100

**Rear Wheel Tire**

<input type="checkbox"/>	142RT1		Pneumatic <i>(Requires the selection of Rear Wheel Tire/Insert)</i>	NC
<input type="checkbox"/>	142RT29		Full Polyurethane <i>(avail w/Mag only)</i>	\$50
<input type="checkbox"/>	142RT4		Low Polyurethane <i>(avail w/24" Mag only)</i>	\$50

**Rear Wheel Tire Tube/Insert**

*Requires a selection of Pneumatic*

<input type="checkbox"/>	142RZ2	E2213	Pneumatic Airless Insert	\$80
<input type="checkbox"/>	142RT1		Pneumatic Tube	NC

**Hand Rim**

<input type="checkbox"/>	142HR1		Aluminum Anodized	NC
<input type="checkbox"/>	142HR4		Plastic Coated	\$95
<input type="checkbox"/>	142HR7		Omit	NC

**Wheel Lock**

<input type="checkbox"/>	142WL1		High Mount Push	NC
<input type="checkbox"/>	142WL2		High Mount Pull	NC

**Extension Handle**

<input type="checkbox"/>	142EH1	E0961	6" Pair	\$75
--------------------------	--------	-------	---------	------

**Grade Aid**

*Requires a selection of Rear Wheel Tire/Insert*

<input type="checkbox"/>	142GA1	E0974	Grade Aid	\$150
--------------------------	--------	-------	-----------	-------

**ACCESSORIES**

**Backpack**

<input type="checkbox"/>	142A1		Utility Pack - Black	\$65
<input type="checkbox"/>	142A2		Kids Black Backpack - Zippie	\$60

**Seat Pouch**

<input type="checkbox"/>	142A8		Black	\$50
--------------------------	-------	--	-------	------

**Miscellaneous**

<input type="checkbox"/>	142A49	E2207	Crutch Holder	\$195
<input type="checkbox"/>	142A83		Tool Kit	\$45
<input type="checkbox"/>	142A54	K0038	Extra Leg Strap	\$30

**Positioning Belt**

<input type="checkbox"/>	142A141	E0978	Belt 1.5" Auto Buckle w/Clamps	\$70
<input type="checkbox"/>	142A117	E0978	Belt 2" Aircraft Buckle w/Clamps	\$110
<input type="checkbox"/>	142A118	E0978	Belt 2" Aircraft Buckle w/Clamps	\$120
<input type="checkbox"/>	142A112	E0978	1" Dual Pull Side Squeeze Pad w /Clmp	\$85
<input type="checkbox"/>	142A814	E0978	1" Comfort Fit APB5116	\$120

**Spoke Guard**

<input type="checkbox"/>	142A24	K0065	Clear	\$115
<input type="checkbox"/>	142A23	K0065	Black	\$115

**Stroller Handles**

<input type="checkbox"/>	142SR2		Regular Stroller Handles	\$125
<input type="checkbox"/>	142SR3		Tall Stroller Handles	\$125

**Side Guard**

*Not available with height adj. flip back armrest*

<input type="checkbox"/>	142SG4	K0108	Plastic Kids	\$140
<input type="checkbox"/>	142SG3	K0108	Plastic Regular	\$140

**Upper Extremity Support Tray**

<input type="checkbox"/>	142A29	E0950	X-Small <i>(Suggested 10"-12")</i>	\$255
<input type="checkbox"/>	142A30	E0950	Small <i>(Suggested 13"-14")</i>	\$255
<input type="checkbox"/>	142A31	E0950	Medium <i>(Suggested 15"-17")</i>	\$255
<input type="checkbox"/>	142A31	E0950	Large <i>(Suggested 18"-20")</i>	\$255

**Lower Leg Length Matrix**

Hanger	Extension	<i>Seat Sling to Top of Footrest</i>	
		Min	Max
60° Swing In/Out	Short	12.5"	16.5"
	Medium	13.0"	18.0"
	Long	15.0"	19.5"
70° Swing In/Out Std & Hemi	Short	13.0"	17.0"
	Medium	13.5"	18.5"
	Long	15.5"	20.5"
60° Swing In/Out w/Platform Flip Up	Short	11.5"	14.5"
	Long	14.0"	18.5"
70° Swing In/Out w/Platform Flip Up	Short	12.5"	15.5"
	Long	15.5"	20.0"
*70° Front Mount	No Ext	6.0"	11.5"
	2" Ext	6.0"	14.0"
	4" Ext	6.0"	16.0"
*80° Front Mount	Short	6.0"	12.0"
	Medium	6.0"	14.5"
	Long	6.0"	16.5"
90° Front Mount	No Ext	5.5"	11.5"
ALR	Medium	15.0"	17.5"
	Long	17.5"	20.0"
ELR	Short	14.0"	18.0"
	Long	15.0"	20.0"

\*2" & 4" Extensions shipped separate, not mounted

**ANGLE ADJUSTABLE BACKREST CHART**

Back Tube Height	LOW	MEDIUM	TALL
Pivot Regular set at 3"	12" <b>13"</b> 14"	15" <b>16"</b> 17"	17" <b>18"</b> 19"
Pivot Regular set at 4"	13" <b>14"</b> 15"	16" <b>17"</b> 18"	18" <b>19"</b> 20"
Pivot Regular set at 5"	14" <b>15"</b> 16"	17" <b>18"</b> 19"	19" <b>20"</b> 21"
Pivot Tall set at 5"	14" <b>15"</b> 16"	17" <b>18"</b> 19"	19" <b>20"</b> 21"
Pivot Tall set at 6"	15" <b>16"</b> 17"	18" <b>19"</b> 20"	N/A
Pivot Tall set at 7"	16" <b>17"</b> 18"	19" <b>20"</b> 21"	N/A
Adjustable Pivot Point Preset Measurement from Seat Sling			
<b>If back height option is not chosen, middle option will be preset.</b>			



Sunrise Medical (US) LLC  
2842 Business Park Ave. · Fresno, CA 93727 · USA