

SERIAL NO: _____

AI

ALUMINIUM 7020



Customer Contact Details

Customer Name: _____

RGK Genius: _____

Address: _____

Phone: _____

Email: _____

Ship to Customer Address: ☐

Shipping Name: _____

Shipping Address: _____

Phone: _____

Email: _____

Additional Information

Notes: _____

Current Wheelchair: _____

Customer Weight: _____

Customer Height: _____

Customer Disability: _____

Other

VAT Exemption Declaration (UK Only):

I declare that I am chronically sick or disabled by reason of:

and that I am receiving from RGK Wheelchairs Ltd, the wheelchair specified herein and I claim relief from Value Added Tax under group 12 of schedule 8 to Value Added Tax Act 1994.

Signature: _____ Date: _____

RGK Contact Authorisation:

From time to time, RGK would like to send you exclusive information, news and offers on everything RGK, by email or post. We will not pass any of your details to any third party.

I would like to be contacted by RGK: ☐

Signature: _____ Date: _____

SERIAL NO: _____

Guide:

☒ Base product

☒ Required option

☐ Optional Extra



Critical Information



Information Needed

PRODUCT STYLE

☐ ELAXZ010 ELITE
\$5,534.00

☐ ELAXA062 Standard Frame

☐ ELAXA063 Standard Frame with 25mm seat rail
\$80.00

☐ ELAXA280 Heavy Duty Frame
\$632.00 (Spinergy R10 wheels also required)

☐ ELAXA145 Rugby League Upgrade
\$501.00

Max User Weight 125kg (275lbs)

Max User Weight ELAXA280 Heavy duty frame 158kg (348lbs)



PLEASE NOTE: ALL RGK CHAIRS ARE HAND BUILT AND HAVE A MANUFACTURING MEASUREMENT TOLERANCE OF +/- 5MM

CAD APPROVAL

☐ XW009 QUOTE ONLY (NO CAD FOR APPROVAL)
\$0.00

☐ XW012 QUOTE CAD drawing for approval
\$483.00 (2 customer amendments ONLY)

☐ XW011 ORDER (NO CAD FOR APPROVAL)
\$0.00

☐ XW010 ORDER CAD drawing for approval
\$299.00 (2 customer amendments ONLY)

CAD approved by: ☐ Customer ☐ RGK Genius

KEY MEASUREMENTS

☒ XY010 Seat width (Inside to inside of sideguards) (Max - 500mm)



mm in

☒ XY020 Seat length (Max - 500mm)



mm in

☐ XY030 Extra seat length (Provides additional leg support & hand grip)



mm in

☒ XY050 Front frame height (Min 360mm - Max 630mm)



mm in

☒ XY060 Rear frame height (Min 360mm - Max 630mm)



mm in

☒ XY070 Frame length (Min = XY020 + XY030 + 140mm)





mm in

☐ XA020 Ergonomic seat (Standard = 120mm)
\$96.00

Specify length



mm in




FRONT FRAME AND SEAT

Front Frame Options (Choose One)

☐ XY080 Seat taper design

i Specify taper width
 mm in



☐ XA130 V-seat design
\$255.00

i Specify taper width
 mm in


**NOTE: INSIDE TO INSIDE
WHEEL (XY470) REQUIRED**



Advanced Front Frame Options

☐ XY082 Seat taper position
\$95.00

i Specify measurement
 mm in



☐ XA131 V-seat and taper design
\$345.00


i Specify inside width a
front of sideguard (A)

mm in

Specify taper width (B)


i mm in

**NOTE: INSIDE TO INSIDE
WHEEL (XY470) REQUIRED**



☐ XA136 Parallel taper design
\$95.00
(Not with Seat Taper XY080)

i Specify taper width



☐ XA135 Asymmetrical taper design
\$403.00


i Specify taper width (A)
 mm in

i Specify offset (B)
(Inside of sideguard to
inside of taper width)

mm in

i Specify offset taper side

☐ Left ☐ Right



Seat Options

☐ XR040 Aluminum seat plate
\$160.00
(riveted to frame)


i Specify length of plate
 mm in
(Std XY020 - 20mm)

☐ XR050 Carbon fibre seat plate
\$407.00
(riveted to frame)

i Specify length of plate
 mm in
(Std XY020 - 20mm)

☐ XA260 Drop seat design
\$222.00
(only with either XR040 or 050)

i Specify drop
 mm in
(Std = 20mm)





FOOTREST

Fixed Footrest

- ☐ ELAXB080 Fixed position footrest

Height and Position

- ☐ Front footrest height
 mm in
 (min = 35mm) 
- ☐ Rear footrest height
 mm in
 (min = 35mm) 
- ☐ Footrest position
 (front of footplate to castor pot)
 mm in 

Footrest Size (Not ABS)

- ☐ Footrest width
 mm in
 (Std= 200mm) 
- ☐ Footrest length
 mm in
 (Std = 250mm) 



Footrest Plate

- ☐ XB180 Aluminium footplate
 (2.5mm)
- ☐ XB190 Carbon fibre footplate
\$265.00
- ☐ XB100 ABS moulded foot tray
\$141.00
 Specify length
 mm in Max = 380mm (Full)
 Std = 330mm
- ☐ ELAXB140 Spacer kit
\$22.00
 (Allows for 25mm increase in height)
- ☐ ELAXB145 Multi position footrest tubes
\$155.00 

Adjustable Footrest

- ☐ ELAXB130 Clamp on adjustable position footrest
\$243.00
 (100mm range)
- ☐ ELAXB085 Adjustable position footrest
\$148.00
 (50mm range)

Height and Position

- ☐ Front footrest height
 mm in
 (min = 40mm) 
- ☐ Rear footrest height
 mm in
 (min = 40mm) 
- ☐ Footrest position
 (front of footplate to castor pot)
 mm in 

Footrest Size (Not ABS)

- ☐ Footrest width
 mm in
 (Std= 200mm) 
- ☐ Footrest length
 mm in
 (Std = 250mm) 

Footrest Plate

- ☐ XB180 Aluminium footplate
 (5mm)
- ☐ XB100 ABS moulded foot tray
\$155.00
 Specify length
 mm in Max = 380mm (Full)
 Std = 330mm

Footrest Plate Options


- ☐ XB240 Aluminium raised edges
\$91.00
 (Cost per Edge)
☐ Front ☐ Back ☐ Left ☐ Right
 Specify height
 mm in (Max = 60mm)

- ☐ XB250 Carbon Fibre raised edge
\$350.00
 (1 side only)
☐ Front ☐ Back
 Specify height
 mm in (Max = 45mm)

Castor Pin Width

- ☒ XY330 Castor pin width
 mm in **MAX 320mm**

☐ Optimise measurements




CAMBER BAR, WING AND ANTI-TIP

Centre of Gravity

☐ ELAXY510 Fixed centre of gravity


Specify measurement
 mm in

(Note: dimension is taken from front of backrest to centre of axle)



☐ ELAXY515 Backrest COG adjustment
\$190.00
 (Not with welded sides)
 (+/-20mm)

Specify measurement
 mm in




Specify Camber Angle


☒ XY390 Camber angle

Wing Options

☐ ELAXA110 Round design basketball wing



☐ ELAXA125 Arrow design basketball wing
\$201.00



☐ ELAXA120 Straight sided design basketball wing
\$201.00



Anti-tip options

☐ ELAXJ040 Single fixed anti tip

☐ ELAXJ060 Double fixed anti tip
\$285.00
 (With high mount joining bar)

Specify castor pin width
 mm in



☐ Optimise measurements
 (Wide as possible - std if not indicated)


SIDEGUARDS

Sideguard Style


- | | | | |
|--|---|---|---|
| ☐ ELAXE230
\$430.00 | Welded aluminium sport
(Not with backrest COG Adj - ELAXY510) | ☐ XE030
\$213.00 | Screw on aluminium |
| ☐ ELAXE250
\$201.00 | Welded aluminium tubes only
(Not with backrest COG Adj - ELAXY510) | ☐ XE050
\$475.00 | Screw on carbon fibre
(Not recommended for basketball use) |
| ☐ XE310
\$202.00 | High design upgrade
(Only with option ELAXE230) | | |



Sideguard Design (Choose One)



☐ Straight design



☐ Round design



☐ Follow tyre

Sideguard Measurements


☐ Sideguard length
 mm in

☐ Sideguard height (back)

mm in

☐ Sideguard height (front)
(optional for max height players)

mm in



☐ Follow Tyre Design

☐ 25mm above tyre

☐ Individual above tyre
 mm in

☐ Level with tyre




BACKREST

Backrest Height (Choose One)


☐ XY230 Fixed backrest height

mm in



☐ XY245 Adjustable backrest height
(+/-25mm)


mm in



Ergonomic Options

☐ XY360 Backrest angle


Individual Deg
(Max -10/+10 Deg)



☐ XG100 Ergonomic backrest
\$96.00


Vertical height
 mm in

Rearward position
 mm in




☐ XG110 Inward taper
\$96.00
(not with ELAXY245)

Specify width or per side
 mm in



☐ XG120 Outward flare
\$96.00
(not with ELAXY245)

Specify width or per side
 mm in



CASTORS AND REAR WHEELS

Castor Options

- ☒ XC160 Sport 61mm castor
(Black)
- ☒ XC200 Sport micro fork

- ☐ XC181 Bones Reds bearings
(All castor bearings)
- ☐ XC182 Bones Swiss Ceramic bearings
(All castor bearings)
\$335.00

Rear Wheel Size

- ☐ XY420 24" Wheels ☐ XY430 25" Wheels ☐ XY440 26" Wheels ☐ XY450 700C Wheels


Rear Wheel Gap

- ☐ XY460 Rear wheel gap
- 15 20 25 30
- Individual mm in



OR

- ☐ XY470 Inside to inside of wheel
NOTE: REQUIRED FOR V SEAT AND TAPER
- mm in



Rear Wheel Options

- ☐ XK030 Equalizer 36/48 spoke
- ☐ XK080 Spinergy SLX 24 spoke
\$902.00

- ☐ XK330 Spinergy **SLX R10** heavy duty 24 spoke
\$1,589.00 (18 deg camber angle only)
- ☐ XK101 Spinergy **XSIX** 24 spoke + caps
\$951.00

Spinergy Colour Options:

- ☐ Black spokes ☐ White spokes ☐ Yellow spokes
- ☐ Red spokes ☐ Blue spokes ☐ Pink spokes
- ☐ Purple spokes ☐ Orange spokes ☐ Green spokes

- ☐ XK180 Individual spoke design
\$212.00 (Specify colour)

Inside:

Outside:

Spinergy Anodise Options

- ☐ XK190 Black hub (std option) ☐ XK200 Silver hub **\$111.00** ☐ XK201 Red hub/rims **\$329.00** ☐ XK205 Individual hub/rims **\$381.00**
- ☐ XK202 Purple hub/rims **\$329.00** ☐ XK203 Blue hub/rims **\$329.00** ☐ XK204 Gold hub/rims **\$329.00**
- Select colour:

Supply Without Wheels

- ☐ XK040 Supply without Wheels
-\$222.00

Wheel Type: Tyre Type:

Handrim Position / Options

- ☐ XY480 Narrow set handrim
- ☐ XL090 Remove excess mounting tabs
\$122.00 (recommended for small wheel gaps)
- ☐ XY485 Wide set handrim
- ☐ XL080 Super narrow set handrim
\$245.00 (only with XL040 and includes XL090)

Handrim Style

- ☐ XL010 Aluminium hard anodised **6 tab**
- ☐ XL070 Aluminium hard anodised **12 tab**
\$122.00 (24" not available)
- ☐ XL040 Titanium 6 tab
\$290.00 (Not with 700C)
- ☐ XL050 Stainless steel 6 tab
\$152.00

Tyre Options

- ☐ XM070 Kenda Kaliente ☐ Red ☐ Blue ☐ XM095 Kenda Wide
\$43.00 (Grey Only)
- ☐ XM150 IRC Exeracer Pro
\$100.00 (Red Only)
- ☐ XM080 Tufo MS3 tubeless
\$100.00 (Blue only)

Spare Wheel

- ☐ ELAXK500 Spare wheel Spinergy
\$1,018.00 (Built up + axle pin)
- ☐ XK520 Spare wheel upcharge - Titanium Rims
\$237.00
- ☐ XK522 Spare wheel upcharge - Remove Excess Tabs
\$63.00
- ☐ XK523 Spare wheel upcharge - Individual Hubs/Rims
\$147.00
- ☐ XK524 Spare wheel upcharge - Individual Spoke Colour
\$107.00

- ☐ XY500 Schraeder inner tube valve



- ☐ XY490 Presta inner tube valve
(Standard if no option selected)



FINISH OPTIONS

Frame & Component finish

☐ XN010 Powder coated frame☐ XN150 Powder coated components (Frame colour if not specified)Frame Comp. Gloss

- | | | |
|--------------------------|--------------------------|------------------------|
| ☐ | ☐ | Black gloss 103 |
| ☐ | ☐ | White gloss 100 |
| ☐ | ☐ | Grey gloss 132 |
| ☐ | ☐ | Cream gloss 135 |
| ☐ | ☐ | Red gloss 108 |
| ☐ | ☐ | Yellow gloss 120 |
| ☐ | ☐ | Forest green gloss 143 |
| ☐ | ☐ | Green gloss 107 |
| ☐ | ☐ | Dark blue gloss 151 |
| ☐ | ☐ | Royal blue gloss 154 |
| ☐ | ☐ | Blue gloss 174 |
| ☐ | ☐ | Duckegg blue gloss 141 |
| ☐ | ☐ | Pink gloss 188 |
| ☐ | ☐ | Orange gloss 180 |
| ☐ | ☐ | Purple gloss 157 |

Frame Comp. Matte

- | | | |
|--------------------------|--------------------------|---------------------|
| ☐ | ☐ | Black matt 105 |
| ☐ | ☐ | White matt 101 |
| ☐ | ☐ | Grey matt 131 |
| ☐ | ☐ | Light grey matt 133 |
| ☐ | ☐ | Navy blue matt 152 |

Frame Comp. Special

- | | | |
|--------------------------|--------------------------|-------------------------|
| ☐ | ☐ | Chrome Mirror 300 |
| ☐ | ☐ | Sparkle silver 170 |
| ☐ | ☐ | Metallic black matt 104 |
| ☐ | ☐ | Metallic grey matt 130 |
| ☐ | ☐ | Black sandpaper 201 |
| ☐ | ☐ | Grey sandpaper 202 |
| ☐ | ☐ | White sandpaper 200 |

☐ ELAXN030 Brushed frame
\$399.00 (All parts black unless XN160 chosen)☐ ELAXN050 Polished frame
\$1,054.00 (All parts black unless XN170 chosen)☐ XN200 Individual paint colour (State Colour/RAL code)
\$197.00☐ XN160 Brushed components
\$151.00 (not required with fixed footrest & backrest)☐ XN170 Polished components
\$222.00 (not required with fixed footrest & backrest)

Accent Colour

Note: Accent colour includes, forks, castor caps, decals.

- | | | | |
|--------------------------|---------------|--------------------------|---------------|
| ☐ | Black accent | ☐ | Silver accent |
| ☐ | Red accent | ☐ | Blue accent |
| ☐ | Purple accent | | |

☐ XN250 Individual accent
\$197.00Select Colour:

UPHOLSTERY

Seat Upholstery

☐ XP040 Seat sling strap adjustable☐ No seat sling (Seat Plate)

Backrest Upholstery

☒ XP055 Airtech sport backrest strap adjustable
(Black)☐ XP070 Backrest upgrade to leather corners
\$138.00 (Black)

Backrest Stitching Options

- | | | | | | | | |
|--------------------------|--------|--------------------------|--------|--------------------------|--------|--------------------------|-------------|
| ☐ | Silver | ☐ | Black | ☐ | Blue | ☐ | Individual: |
| ☐ | Red | ☐ | Orange | ☐ | Purple | <input type="text"/> | |

RGK Logo Colour

- | | | | | | | | |
|--------------------------|-------------|--------------------------|--------|--------------------------|--------|--------------------------|-------------|
| ☐ | White (Std) | ☐ | Green | ☐ | Blue | ☐ | Individual: |
| ☐ | Red | ☐ | Orange | ☐ | Purple | <input type="text"/> | |

Model Logo Colour

- | | | | | | | | |
|--------------------------|-------------|--------------------------|--------|--------------------------|--------|--------------------------|-------------|
| ☐ | White (Std) | ☐ | Green | ☐ | Blue | ☐ | Individual: |
| ☐ | Red | ☐ | Orange | ☐ | Purple | <input type="text"/> | |

☐ XP170 Individual backrest embroidery
\$156.00 ADD DETAILS ON PAGE 10

CUSHION

Cushion Foam Options

- ☐ **XP080**
\$96.00 Cushion foam
(select density and size below)
- ☐ Soft (CMS040) ☐ Medium (CMS050)
- ☐ 50mm ☐ 75mm ☐ 100mm ☐ Other mm in
- ☐ **XP340**
\$90.00 Hard Plastazote cushion upgrade
(Indicate max height)
- ☐ 580mm ☐ 630mm

Cushion Cover Options





- ☐ **XP085** Nylon cushion cover
(Black)
- ☐ **XP090** Air Lite cushion cover
(Black)
- ☐ **XP380**
\$102.00 Neoprene cushion cover
(Black)

STRAPPING

Velcro straps

- ☐ **XS050**
\$137.00 100mm Pro adjustable velcro strap
- ☐ **XS055**
\$137.00 50mm Pro adjustable velcro strap
- ☐ **XS057**
\$137.00 50mm Pro knee strap
(measurement from tube to tube over knees)

Velcro strap position

- ☐ **Strap 1**
-  
- Specify strap length mm in
- Specify distance from back post mm in
- ☐ **Strap 2**
-  
- Specify strap length mm in
- Specify distance from back post mm in

- ☐ **XS020**
\$93.00 Sport ratchet strap (Click strap)
(11" x 2" pad)
- ☐ **XS021**
\$94.00 Sport ratchet strap (Single)
- ☐ **XS025**
\$243.00 Pelvi Loc 2DA ratchet strap
(For thighs and knees)

 Select size:
☐ Medium (12") ☐ Large (14")

- ☐ **XS030**
\$260.00 Pelvi Loc 3DA ratchet strap
(For hips)



 Select size:
☐ S (8") ☐ M (12") ☐ L (14") ☐ XL (17")

- ☐ **XS035**
\$155.00 M2 ratchet strap
(For thighs and knees)

 ☐ 8" ☐ 10" ☐ 12"

- ☐ **XS036**
\$155.00 M2 ratchet strap (Single)
(For Feet)

 ☐ 8" ☐ 10" ☐ 12"

- ☐ **Strap 1**
- 
- Specify distance from back post mm in
- ☐ **Strap 2**
- 
- Specify distance from back post mm in

- ☐ **XS040**
\$97.00 Fixed mounting points
(one set per click strap) ☐ x1 ☐ x2 ☐ x3

- ☐ **XS045**
\$97.00 Clamped mounting points
(one set per click strap) ☐ x1 ☐ x2 ☐ x3

FRAME OPTIONS

- ☐ **XR060**
\$116.00 Calf strap mounting loops



- ☐ **XR065**
\$116.00 Strap location loops



ACCESSORIES

- ☐ **XR030**
\$114.00 Triple wheel bag
(Black)
- ☐ **XP130**
\$104.00 Frame protectors full length
- ☐ **XA150**
\$74.00 Plastic tube fitted to offensive wing

- ☐ **XS070** Foot strap (select if required)
- ☐ **XS060** Calf strap (select if required)
- ☒ **XP400** Back bar protector
- ☒ **XR090** Tool kit

ADDITIONAL NOTES

INDIVIDUAL BACKREST EMBROIDERY

PLEASE PROVIDE JPEG IMAGE AT TIME OF ORDER AND NOTE POSITION BELOW



BRACE BAR

SPEC NOTES:

☐

OFF SPEC REQUEST*

(please specify below the page number(s) and details to which the off spec request applies.

☐

**OFF SPEC REQUEST
APPROVED** (for internal
RGK office only)

X

RGK Internal Sign-off

Off spec requests require a feasibility study which may increase order price and lead time. Any dimension/ option/ frame request not listed on the order form will be treated as an off spec request.*

Estimated MSRP*

\$ 0.00