

AI

ALUMINIUM 7020



Customer Contact Details

Customer Name: _____

RGK Genius: _____

Address: _____

Phone: _____

Email: _____

Ship to Customer Address: ☐

Shipping Name: _____

Shipping Address: _____

Phone: _____

Email: _____

Additional Information

Notes: _____

Current Wheelchair: _____

Customer Weight: _____

Customer Height: _____

Customer Disability: _____

Other

VAT Exemption Declaration (UK Only):

I declare that I am chronically sick or disabled by reason of:

and that I am receiving from RGK Wheelchairs Ltd, the wheelchair specified herein and I claim relief from Value Added Tax under group 12 of schedule 8 to Value Added Tax Act 1994.

Signature: _____ Date: _____

RGK Contact Authorisation:

From time to time, RGK would like to send you exclusive information, news and offers on everything RGK, by email or post. We will not pass any of your details to any third party.

I would like to be contacted by RGK: ☐

Signature: _____ Date: _____

SERIAL NO: _____

Guide:

☒ Base product

☒ Required option

☐ Optional Extra



Critical Information



Information Needed

PRODUCT STYLE

☐ EXAXZ010 ELITE X
\$6,641.00

☒ XA600 Heat treated frame

☐ EXAXA062 Standard

☐ EXAXA063 Standard with 25mm seatrail
\$80.00

Max User Weight 125kg (275lbs)



PLEASE NOTE: ALL RGK CHAIRS ARE HAND BUILT AND HAVE A MANUFACTURING MEASUREMENT TOLERANCE OF +/- 5MM

CAD APPROVAL

CAD approved by: ☐ Customer ☐ RGK Genius

☐ XW009 QUOTE ONLY (NO CAD FOR APPROVAL)
\$0.00

☐ XW011 ORDER (NO CAD FOR APPROVAL)
\$0.00

☐ XW012 QUOTE CAD drawing for approval
\$483.00 (2 customer amendments ONLY)

☐ XW010 ORDER CAD drawing for approval
\$299.00 (2 customer amendments ONLY)

KEY MEASUREMENTS

☒ XY010 Seat width (Inside to inside of sideguards)



mm in

(Max - 500mm)

☒ XY020 Seat length



mm in

(Max - 500mm)

☐ XY030 Extra seat length (Provides additional leg support & hand grip)



mm in

☒ XY050 Front frame height



mm in

(Max - 630mm)
(Min - 360mm)

☒ XY060 Rear frame height



mm in

(Max - 630mm)
(Min - 360mm)

☒ XY070 Frame length



mm in

(Min = XY020 + XY030 + 140mm)


☐ XA020 Ergonomic seat
\$96.00

Specify length



mm in

(Standard = 120mm)



FRONT FRAME AND SEAT

Front Frame Options

☐ XY080 Seat taper design


i Specify taper width
 mm in



☐ XA130 V-seat design
\$255.00

i Specify taper width
 mm in


**NOTE: INSIDE TO INSIDE
 WHEEL (XY470) REQUIRED**



Advanced Front Frame Options

☐ XY082 Seat taper position
\$95.00

i Specify measurement
 mm in




☐ XA131 V-seat and taper design
\$345.00

i Specify inside width at
 front of sideguard (A)
 mm in


i Specify taper width (B)
 mm in

**NOTE: INSIDE TO INSIDE
 WHEEL (XY470) REQUIRED**



☐ XA136 Parallel taper design
\$95.00
 (Not with Seat Taper XY080)

i Specify taper width




☐ XA135 Asymmetrical taper design
\$403.00

i Specify taper width (A)
 mm in

i Specify offset (B)
 (Inside of sideguard to
 inside of taper width)
 mm in

Specify offset taper side
☐ Left ☐ Right




Seat Options

☐ XR050 Carbon fibre seat plate
\$407.00
 (riveted to frame)

i Specify length of plate
 mm in (Std XY020 - 20mm)

☐ XA260 Drop seat design
\$222.00
 (only with either XR040 or 050)

i Specify drop
 mm in (Std = 20mm)



FOOTREST

Footrest Type


☒ EXAXB080 Fixed position footrest

☐ EXAXB145 Multi position footrest tubes
\$161.00



Footrest Measurements

☐ Front footrest height
0 mm 0 in (min = 35mm)




☐ Footrest width
0 mm 0 in (Std = 200mm)




☐ Rear footrest height
0 mm 0 in (min = 35mm)



☐ Footrest length
0 mm 0 in (Std = 250mm)



☐ Footrest position
(front of footplate to castor pot)
0 mm 0 in



Footrest Plate

☐ XB190 Carbon fibre footplate
\$265.00

☐ XB100 ABS moulded foot tray
\$155.00

☐ EXAXB140 Spacer kit
\$22.00 (Allows for 25mm increase in height)

Specify length
0 mm 0 in Max = 380mm (Full)
Std = 330mm

Footrest Plate Options

☐ XB240 Aluminium raised edges
\$91.00 (Cost per Edge)

☐ Front ☐ Back ☐ Left ☐ Right

Specify height
0 mm 0 in (Max = 60mm)

☐ XB250 Carbon Fibre raised edge
\$350.00 (1 side only)

☐ Front ☐ Back


Specify height
0 mm 0 in (Max = 45mm)

Castor Pin Width

☒ XY330 Castor pin width

0 mm 0 in MAX 320mm

☐ Optimise measurements



CAMBER BAR, WING AND ANTI-TIP

Centre of Gravity and Camber Angle

☒ EXAXY510 Fixed centre of gravity

Specify measurement
0 mm 0 in

(Note: dimension is taken from front of backrest to centre of axle)



☒ XY390 Camber angle

18 0 20 0



Wing Options

☐ EXAXA110 Round design basketball wing




☐ EXAXA125 Arrow design basketball wing
\$201.00



☐ EXAXA120 Straight sided design basketball wing
\$201.00



☐ EXAXA129 Individual design basketball wing
\$201.00



Anti-tip options

☐ EXAXJ040 Single fixed anti tip

☐ EXAXJ060 Double fixed anti tip
\$285.00

SIDEGUARDS

Sideguard Style


☐ EXAXE230 Welded aluminium sport
\$430.00

☐ XE310 High design upgrade
\$202.00




☐ EXAXE250 Welded aluminium sport with Carbon Fibre
\$951.00

Sideguard Design (Choose One)



☐ Straight design




☐ Round design




☐ Follow tyre

Sideguard Measurements

☐ Sideguard length
0 mm 0 in



☐ Sideguard height (back)
0 mm 0 in



☐ Sideguard height (front)
(optional)
0 mm 0 in



☐ Follow Tyre Design

☐ 25mm above tyre



☐ Individual above tyre
0 mm 0 in



☐ Level with tyre

BACKREST


Backrest Height

☒ XY230 Fixed backrest height



Specify height

0 mm 0 in




Ergonomic Options

☒ XY360 Backrest angle

-5 -3 0 +3 +5

Individual
(Max -10/+10 Deg)



☐ XG100 Ergonomic backrest
\$96.00

Vertical height

0 mm in



Rearward position


0 mm in



☐ XG110 Inward taper
\$96.00
(inside of backrest posts)

Specify width / per side


0 mm in



☐ XG120 Outward flare
\$96.00
(inside of backrest posts)

Specify width / per side

0 mm in



CASTOR AND REAR WHEELS

Castor Options

- | | |
|---|---|
| ☐ XC160 Sport 61mm castor
(Black) | ☐ XC181 Bones Reds bearings
(All castor bearings) |
| ☐ XC200 Sport micro fork | ☐ XC182 Bones Swiss Ceramic bearings
(All castor bearings)
\$335.00 |

Rear Wheel Size

- | | | | |
|---|---|---|--|
| ☐ XY420 24" Wheels | ☐ XY430 25" Wheels | ☐ XY440 26" Wheels | ☐ XY450 700C Wheels |
|---|---|---|--|

Rear Wheel Gap

☐ XY460 Rear wheel gap  ☐ 15 ☐ 20 ☐ 25 ☐ 30 Individual 0 mm in	OR	☐ XY470 Inside to inside of wheel 0 mm in NOTE: REQUIRED FOR V SEAT AND TAPER
---	----	--

Rear Wheel Options

- | | |
|---|---|
| ☐ XK080 Spinergy SLX 24 spoke
\$902.00 | ☐ XK330 Spinergy SLX R10 heavy duty 24 spoke
(18 deg camber angle only)
\$1,589.00 |
| ☐ XK101 Spinergy XSLX 24 spoke + Caps
\$951.00 | ☐ XK521 Spinergy Ceramic bearing upgrade
\$289.00 |

Spinergy Colour Options:



- | | | |
|--|--|--|
| ☐ Black spokes | ☐ White spokes | ☐ Yellow spokes |
| ☐ Red spokes | ☐ Blue spokes | ☐ Pink spokes |
| ☐ Purple spokes | ☐ Orange spokes | ☐ Green spokes |

- ☐ XK180 Individual spoke design
(Specify colour)
\$212.00

Inside:

Outside:

Spinergy Anodise Options

- | | | | |
|---|---|---|---|
| ☐ XK190 Black hub
(std option) | ☐ XK200 Silver hub
\$111.00 | ☐ XK201 Red hub/rims
\$329.00 | ☐ XK205 Individual hub/rims
Select colour:
\$381.00 |
| ☐ XK202 Purple hub/rims
\$329.00 | ☐ XK203 Blue hub/rims
\$329.00 | ☐ XK204 Gold hub/rims
\$329.00 | <input type="text"/> |

Supply Without Wheels

- ☐ XK040 Supply without Wheels
-\$222.00
- 

Handrim Position / Options

- | | |
|---|--|
| ☐ XY480 Narrow set handrim | ☐ XY485 Wide set handrim |
| ☐ XL090 Remove excess mounting tabs
(recommended for small wheel gaps)
\$122.00 | ☐ XL080 Super narrow set handrim
(only with XL040 and includes XL090)
\$233.00 |

Handrim Style

- | | |
|--|---|
| ☐ XL010 Aluminium hard anodised 6 tab | ☐ XL040 Titanium 6 tab
(Not with 700C)
\$290.00 |
| ☐ XL070 Aluminium hard anodised 12 tab
(24" not available)
\$122.00 | ☐ XL050 Stainless steel 6 tab
\$152.00 |

Tyre Options

- | | | |
|--|--|--|
| ☐ XM070 Kenda Kaliente | ☐ Red ☐ Blue | ☐ XM095 Kenda Wide
(Grey Only)
\$43.00 |
| ☐ XM150 IRC Exeracer Pro
(Red Only)
\$100.00 | | |
| ☐ XM080 Tufo MS3 tubeless
(Blue only)
\$100.00 | | |

Spare Wheel

- | | |
|--|--|
| ☐ EXAXK500 Spare wheel Spinergy
(Built up + axle pin)
\$1,018.00 | ☐ XK523 Spare wheel upcharge - Individual Hubs/Rims
\$147.00 |
| ☐ XK520 Spare wheel upcharge - Titanium Rims
\$237.00 | ☐ XK524 Spare wheel upcharge - Individual Spoke Colour
\$107.00 |
| ☐ XK522 Spare wheel upcharge - Remove Excess Tabs
\$63.00 | |

- ☐ XY500 Schraeder inner tube valve



- ☐ XY490 Presta inner tube valve
(Standard if no option selected)



FINISH OPTIONS

Frame finish

☐ XN010 Powder coated frame

Frame Comp

Gloss

- ☐ ☐ Black gloss 103
- ☐ ☐ White gloss 100
- ☐ ☐ Grey gloss 132
- ☐ ☐ Cream gloss 135
- ☐ ☐ Red gloss 108
- ☐ ☐ Yellow gloss 120
- ☐ ☐ Forest green gloss 143
- ☐ ☐ Green gloss 107
- ☐ ☐ Dark blue gloss 151
- ☐ ☐ Royal blue gloss 154
- ☐ ☐ Blue gloss 174
- ☐ ☐ Duckegg blue gloss 141
- ☐ ☐ Pink gloss 188
- ☐ ☐ Orange gloss 180
- ☐ ☐ Purple gloss 157

Frame Comp

Matt

- ☐ ☐ Black matt 105
- ☐ ☐ White matt 101
- ☐ ☐ Grey matt 131
- ☐ ☐ Light grey matt 133
- ☐ ☐ Navy blue matt 152

Frame

Special

- ☐ ☐ Chrome Mirror 300
- ☐ ☐ Sparkle silver 170
- ☐ ☐ Metallic black matt 104
- ☐ ☐ Metallic grey matt 130
- ☐ ☐ Black sandpaper 201
- ☐ ☐ Grey sandpaper 202
- ☐ ☐ White sandpaper 200

☐ XN200 Individual paint colour (State Colour/RAL code)
\$197.00

☐ ☐

☐ EXAXN030 Brushed frame
\$399.00

☐ XN160 Brushed components
\$151.00 (not required with fixed footrest & backrest)

☐ EXAXN050 Polished frame
\$1,054.00

☐ XN170 Polished components
\$233.00 (not required with fixed footrest & backrest)

Accent Colour

Note: Accent colour includes, forks, castor caps, decals, ELITE logo and stitching. Cushion, toe strap, frame protectors and RGK logo are



- ☐ Black accent
- ☐ Red accent
- ☐ Purple accent

- ☐ Silver accent
- ☐ Blue accent
- ☐ Carbon accent

(caps only, black for others)

☐ XN250 Individual accent
\$197.00

UPHOLSTERY

Seat Upholstery

☐ EXAXP040 Seat sling strap adjustable

☐ No seat sling (Seat Plate)

Backrest Upholstery

☒ EXAXP055 Airtech sport backrest strap adjustable
(comes with Ti buckles)

☐ XP070 Backrest upgrade to leather corners
\$138.00 (Black)

Backrest Stitching Options

- ☐ Silver ☐ Black ☐ Blue ☐ Individual:
- ☐ Red ☐ Orange ☐ Purple

RGK Logo Colour

- ☐ White (Std) ☐ Green ☐ Blue ☐ Individual:
- ☐ Red ☐ Orange ☐ Purple

Model Logo Colour

- ☐ White (Std) ☐ Green ☐ Blue ☐ Individual:
- ☐ Red ☐ Orange ☐ Purple

Backrest Accessories

☐ XP170 Individual backrest embroidery
\$156.00 ADD DETAILS ON PAGE 9

CUSHION

Cushion foam options

- ☐ XP080 \$96.00 Cushion foam
(select density and size below)
- ☐ Soft (CMS040) ☐ Medium (CMS050)
- ☐ 50mm ☐ 75mm ☐ 100mm ☐ Other mm in
- ☐ XP340 \$90.00 Hard Plastazote cushion upgrade
(Indicate max height)
- ☐ 580mm ☐ 630mm

Cushion Cover Upgrades

- ☐ XP085 Nylon cushion cover (Black) ☐ XP090 \$74.00 Air Lite cushion cover (Black) ☐ XP380 \$102.00 Neoprene cushion cover (Black)


STRAPS

Velcro straps


- ☐ XS050 \$137.00 100mm Pro adjustable velcro strap
- ☐ XS055 \$137.00 50mm Pro adjustable velcro strap
- ☐ XS057 \$137.00 50mm Pro adjustable velcro **Knee** strap
(measurement from tube to tube over knees)

Velcro strap position

☐ Strap 1




Specify strap length

 mm in



Specify distance from back post

 mm in

☐ Strap 2



Specify strap length

 mm in


Specify distance from back post

 mm in

Ratchet Straps

- ☐ XS020 \$93.00 Sport ratchet strap (Click strap)
(11" x 2" pad)
- ☐ XS021 \$93.00 Sport ratchet strap (Single)
- ☐ XS025 \$243.00 Pelvi Loc 2DA ratchet strap
(For thighs and knees)

Select size:

☐ Medium (12") ☐ Large (14")

- ☐ XS030 \$260.00 Pelvi Loc 3DA ratchet strap
(For hips)

Select size:

☐ S (8") ☐ M (12") ☐ L (14") ☐ XL (17")

- ☐ XS035 \$155.00 M2 ratchet strap
(For thighs and knees)


☐ 8" ☐ 10" ☐ 12"

- ☐ XS036 \$155.00 M2 ratchet strap (Single)
(For Feet)

☐ 8" ☐ 10" ☐ 12"

Ratchet strap position


☐ Strap 1



Specify distance from back post

 mm in

☐ Strap 2



Specify distance from back post

 mm in

- ☐ \$97.00 Fixed mounting points
(one set per click strap)

☐ x1 ☐ x2 ☐ x3

FRAME OPTIONS

- ☐ XR060 \$116.00 Calf strap mounting loops



- ☐ XR065 \$116.00 Strap location loops



ACCESSORIES

- ☐ XR030 \$114.00 Triple wheel bag
(Black)
- ☐ XP130 \$104.00 Frame protectors full length
- ☐ XA150 \$74.00 Plastic tube fitted to offensive wing

- ☐ XS070 Foot strap (select if required)
- ☐ XS060 Calf strap (select if required)
- ☒ XP400 Back bar protector
- ☒ XR090 Tool kit

ADDITIONAL NOTES

INDIVIDUAL BACKREST EMBROIDERY

PLEASE PROVIDE JPEG IMAGE AT TIME OF ORDER AND NOTE POSITION BELOW



BRACE BAR

SPEC NOTES:

☐

OFF SPEC REQUEST*

(please specify below the page number(s) and details to which the off spec request applies.

☐

OFF SPEC REQUEST

APPROVED *(for internal RGK office only)*

X

RGK Internal Sign-off

Off spec requests require a feasibility study which may increase order price and lead time. Any dimension/ option/ frame request not listed on the order form will be treated as an off spec request.*

Estimated MSRP*