

All-terrain 4x4 power wheelchair

EIPW60



Customer				Ship To			
Account				Name			
Account #		Date		Address			
☐ Quote	☐ Order	PO #		City		State	
Buyer				ZIP		Phone	
Mark For				Email / Fax			

Please submit orders to: magicmobility@sunmed.com

Seating Specifications


HCPC	N/A	Typical WC Weight w/ Batteries	Starting at 320 lbs. Weights vary depending on options
Power Seating	50° of CG Tilt / Power Recline / 12" Seat Elevator	Seat Type	MPS and Rehab Style
Seat Size Range	Width: 13-23" Depth: 13-24"	STF Height Range	18.5", 19.5", 20.5"

Drive Base Specifications

Drive Wheel Position	4 Wheel Drive	Drive Control	R-Net PM 120A controller
User Weight Capacity	400 lbs. 340 lbs. with Elevate or Tilt	Drive Motors	4-pole 6mph
X8 Base Size	Width 28" Length 40.5" Turning Radius 56"	Drive Wheels	14" All Terrain Drive Wheels
		Tires	All Terrain Tires
		Charger	Dual Mode
		Battery Size	Group 24
		Color Accents	6 shroud colors and 6 color accents
		Ground Clearance	4"

Client Measurements

a) Top of head to seat	_____	in
b) Top of shoulder to seat	_____	in
c) Underarm to seat	_____	in
d) Elbow to seat	_____	in
e) Top of leg to seat	_____	in
f) Shoulder width	_____	in
g) Chest width	_____	in
i) Hip width	_____	in
j) Actual sitting depth	_____	in
m) Foot length	_____	in
n) Knee to foot	_____	in
o) Seat width	_____	in



Notes

✓ OPTION	HCPCS	DESCRIPTION	RETAIL		
Seating Options			\$		
☐ 600SF4		REHAB Seat and Back - Back Canes and adjustable flat seat pan	NC		
☐ 600SF5		MPS Seat and Back - Captain style seating	NC		
☐ 600SF6		MPS Back and REHAB Seat - Captain style backrest and adjustable flat seat pan	NC		
Base					
Extreme X8 Base		End Item Number: EIPW60	\$		
☐ 600MW79		Magic Mobility - Extreme X8 No Power Seating Options	\$ 20,037.00		
☐ 600MW80		Magic Mobility - Extreme X8 Single Power Seating Options	\$ 20,037.00		
☐ 600MW81		Magic Mobility - Extreme X8 Multi Power Seating Options	\$ 20,037.00		
☐ 600MW82		Magic Mobility - Extreme X8 No Power Seating Options - HD	\$ 20,037.00		
☐ 600MW83		Magic Mobility - Extreme X8 Single Power Seating Options - HD	\$ 20,037.00		
☐ 600MW84		Magic Mobility - Extreme X8 Multi Power Seating Options - HD	\$ 20,037.00		
☐ 600MW85		Magic Mobility - Extreme X8 No Power Seating Options - MPS (Captain)	\$ 20,037.00		
☐ 600MW86		Magic Mobility - Extreme X8 Single Power Seating Options - MPS (Captain)	\$ 20,037.00		
☐ 600MW87		Magic Mobility - Extreme X8 Multi Power Seating Options - MPS (Captain)	\$ 20,037.00		
☐ 600MW88		Magic Mobility - Extreme X8 No Power Seating Options - HD MPS (Captain)	\$ 20,037.00		
☐ 600MW89		Magic Mobility - Extreme X8 Single Power Seating Options - HD MPS (Captain)	\$ 20,037.00		
☐ 600MW90		Magic Mobility - Extreme X8 Multi Power Seating Options - HD MPS (Captain)	\$ 20,037.00		
Base Accessories					
☐ 6008A91		Jack for Drive Tire	\$ 213.00		
Tie Downs					
☐ 600TO52		Crash Tested X8 Tie Down Kit	\$ 830.00		
Additional Base Function					
☐ 600WL93		Directional Lock X8 Base - Used to reverse in straight line.	\$ 2,105.00		
Optional Drive Tire Options					
☐ 6008A101		Tire Sealant for 4 Pneumatic Drive Tires	\$ 80.00		
Shroud Color					
☐ 600C185	Midnight Black Shroud Color	NC	☐ 600C197	Clover Green Shroud Color	NC
☐ 600C99	Ocean Blue Shroud Color	NC	☐ 600C198	Tangerine Shroud Color	NC
☐ 600C196	Fire Engine Red Shroud Color	NC	☐ 600C199	Camouflage Shroud Color	\$ 638.00
Seat To Floor Height					
X8 Seat to Floor Heights start at 18.5" +/- 0.25"					
☐ 600HT44		Lowest Seat To Floor Height	NC		
☐ 600HT45		Seat To Floor Height +1"	NC		
☐ 600HT46		Seat To Floor Height +2"	NC		
Batteries					
☐ 600AG25	E2363	Group 24 Gel Battery (Pair) - Installed	\$ 1,138.00		
Battery Charger					
☐ 600BJ28		10A Battery Charger	NC		
☐ 600BJ38		12A Battery Charger	\$ 495.00		

Seating: Seat Functions			
✓ OPTION	HCPCS	DESCRIPTION	RETAIL
No Power Seat Options			\$
☐ 600ST86		No Power Seating System - 400 lb. weight capacity	\$ 4,345.00
☐ 600ST87		No Power Seating System with Power Functionality - 400 lb. weight capacity Required for lights, power recline, power legrest, steering lock, or future power seat functions	\$ 5,586.00
Power Tilt			\$
Utilized without Power Recline.			
☐ 600PA64	E1002	50° Power Tilt in Space - 340 lb. weight capacity	\$ 8,458.00
Power Seat Elevate			\$
Only available when Power Tilt or Power Tilt and Recline is selected.			
☐ 600PI16	E2298	12" Adjustable Seat Elevator - 340 lb. weight capacity	\$ 4,750.00
Power Recline			\$
Utilized without Power Tilt. Requires No Power Seating with Power Functionality.			
☐ 600PA65	E1003	150° Power Recline Only - 300 lb. weight capacity	\$ 7,225.00
Power Tilt and Recline			\$
☐ 600PA66	E1006	50° Tilt & 150° Power Recline - 300 lb. weight capacity	\$ 13,709.00

Seating: Seat Sizes			
MPS (Captain) Seat Cushion (width x depth)			\$
MPS Backrest can be paired with the REHAB flat pan. If flat pan is desired, omit the MPS cushion and select the desired REHAB flat pan.			
☐ 600SI24		16" x 18" Vinyl Seat	NC
☐ 600SI25		18" x 18" Vinyl Seat	NC
☐ 600SI26		18" x 20" Vinyl Seat	NC
☐ 600SI27		20" x 18" Vinyl Seat	NC
☐ 600SI28		20" x 20" Vinyl Seat	NC
☐ 600SI29		22" x 20" Vinyl Seat	\$ 486.00
REHAB Flat Pan Seat Sizes - Adjustable in width and depth for most sizes			\$

REHAB Seat Pan

Standard seat pan available 13-23" wide. Deep seat pan available 13"-23" wide.

☐ 600PN13		Standard Adjustable Flat Seat Pan - 13" to 19" deep	NC
☐ 600PN14		Deep Adjustable Flat Seat Pan - 20" to 22" deep	\$ 486.00
☐ 600PM15		Deep Adjustable Flat Seat Pan - 23" to 24" deep	\$ 858.00

REHAB Seat Width

☐ 13"	☐ 14"	☐ 15"	☐ 16"	☐ 17"	☐ 18"	☐ 19"	☐ 20"	☐ 21"	☐ 22"	☐ 23"
------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------

REHAB Seat Depth

☐ 13"	☐ 14"	☐ 15"	☐ 16"	☐ 17"	☐ 18"	☐ 19"	☐ 20"	☐ 21"	☐ 22"	☐ 23"	☐ 24"
------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------

REHAB Seating Cushion Options			\$
Cushion size will be determined by Rehab Seat size selected above			
☐ 600CU59	E2607	JAY® Ion® (Sizes 16" - 21" Wide)	\$ 460.00
☐ 600CU60	E2608	JAY® Ion® (Size 22" Wide)	\$ 578.00
☐ 600CU61	E2607	JAY® Union® (Sizes 16" - 21" Wide)	\$ 520.00
☐ 600CU62	E2608	JAY® Union® (Size 22" Wide)	\$ 635.00
☐ 600CU64	E2622	JAY® Fusion® (Sizes 16" - 21" Wide)	\$ 514.00
☐ 600CU65	E2623	JAY® Fusion® (Size 22" Wide)	\$ 613.00

Seating: Seat Sizes			
✓ OPTION	HCPCS	DESCRIPTION	RETAIL
Support Pads and Hardware (must pick pad and mount)			\$
Lateral Thigh Support Pads			
☐ 600KF22	E0953	Lateral Thigh Support Pad Short - Right (5x5)	\$ 147.00
☐ 600KF23	E0953	Lateral Thigh Support Pad Short - Left (5x5)	\$ 147.00
Lateral Thigh Support Mounts			
☐ 600KF26	E1028	Lateral Thigh Removable Mount - Right	\$ 292.00
☐ 600KF27	E1028	Lateral Thigh Removable Mount - Left	\$ 292.00
Seating: Seat Back Options			
MPS Backrest Options - Comes standard with Manual Recline			\$
Backrest Width			
☐ 600BN46		Standard Backrest Width - 16"-20" Wide	NC
☐ 600BN47		Wide Backrest Width - 22"-23" Wide	\$ 353.00
Backrest Height			
☐ 600B80		Standard Backrest Height - 22" Tall	NC
☐ 600B81		Tall Backrest Height - 26" Tall	\$ 586.00
Backrest Material			
☐ 600BU106		Cloth Backrest	NC
☐ 600BU107		Vinyl Backrest	\$ 319.00
MPS Back Accessories			\$
☐ 600PH33		Push Handle - Needed for O2 holder	\$ 167.00
☐ 6009A13		Accessory Mount - Needed for aftermarket headrest	\$ 127.00
REHAB Seating Backrest Options			\$
Backrest will have the width of the MPS seat or flat pan. Please choose your back height. Backrest will come with 2 point hardware with a 500 lb. weight capacity.			
JAY J3 Back Contour			
☐ 600JP01	E2313	J3® PA Back for 16"-20" Wide - 2.2" of lateral thoracic support	\$ 625.00
☐ 600JP02	E2614	J3® PA Back for 22" Wide - 2.2" of lateral thoracic support	\$ 928.00
☐ 600JP11	E2615	J3® PL Back for 16"-20" Wide - 3" of lateral thoracic support	\$ 729.00
☐ 600JP13	E2620	J3® PD Back for 16"-20" Wide (18" and 21" shell height only) - 6" of lateral thoracic support	\$ 875.00
JAY® J3® Back Shell Options and Height Ranges			
☐ 600JR14		16.5" tall shell (Upper Thoracic Short)	NC
☐ 600JR15		18" tall shell (Upper Thoracic Med)	NC
☐ 600JR16		19.5" tall shell (Upper Thoracic Tall)	NC
☐ 600JR17		21" tall shell (Shoulder Height)	NC
☐ 600JR27		22.5" tall shell (Shoulder Height)	NC
☐ 600JR28		24" tall shell (Shoulder Height)	NC
REHAB Flip Forward Backrest Option			\$
☐ 600D23		REHAB Flip Forward Backrest Option	\$ 311.00

Seating: Seat Back Options			
✓ OPTION	HCPCS	DESCRIPTION	RETAIL
JAY® Lateral Thoracic Swing-Away Supports (Requires Pad and Mount) - J3 back only			\$
JAY® Lateral Thoracic Support Pads			
☐ 600JW35	E0956	Lateral Thoracic (L x H) 4" x 5" Pad - Right	\$ 147.00
☐ 600JW38	E0956	Lateral Thoracic (L x H) 4" x 5" Pad - Left	\$ 147.00
☐ 600JW37	E0956	Lateral Thoracic (L x H) 6" x 4" Pad - Right	\$ 147.00
☐ 600JW40	E0956	Lateral Thoracic (L x H) 6" x 4" Pad - Left	\$ 147.00
JAY® Lateral Thoracic Swing-Away Mount			
☐ 600JW41	E1034	Lateral Thoracic Swing Away Bracket - Right	\$ 292.00
☐ 600JW42	E1034	Lateral Thoracic Swing Away Bracket - Left	\$ 292.00
JAY® Anterior Thoracic Support			\$
☐ 600ZR10	E0960	Bodypoint Monoflex Medium Chest Strap - Pad 3.5" x 18.5" (total length 49.5")	\$ 158.00
☐ 600ZR11	E0960	Bodypoint Monoflex Large Chest Strap - Pad 3.5" x 21.5" (total length 52.5")	\$ 158.00
☐ 600K10	E0960	JAY® Anterior Thoracic Support - Medium ATS8302	\$ 177.00
☐ 600K11	E0960	JAY® Anterior Thoracic Support - Large ATS8303	\$ 177.00
Pelvic Positioning Options - Select One			\$
Pelvic Positioning			
☐ 600A848		Positioning Belt, Auto Style Buckle, 2" x 60"	NC
☐ 600A849	E0978	Positioning Belt, Aircraft Buckle, Padded, 2" x 50"	\$ 88.00
☐ 600A850	K0108	Bodypoint - Evoflex Pelvic Stabilizer Belt	\$ 203.00
☐ 600A851	K0108	Bodypoint - Evoflex Pelvic Stabilizer Belt with Extension	\$ 224.00
Headrest Options - MPS Style			\$
☐ 600NU12	E0955	Cloth Headrest - MPS seat back only (12" x 4.5" x 6.25")	\$ 432.00
☐ 600NU13	E0955	Vinyl Headrest - MPS seat back only (12" x 4.5" x 6.25")	\$ 432.00
Headrest Options - REHAB Style			\$
Headrest Pad			
☐ 600HE39	E0955	Whitmyer® Plush Headrest Pad, 10"	\$ 239.00
☐ 600HE46	E0955	Whitmyer® Adjust-A-Plush Headrest Pad Medium (Adjustable from 7 to 16" Wide x 4" Tall)	\$ 367.00
Headrest Mount			
☐ 600MH33	E1033	Whitmyer® Linx Headrest Mount with Quick Release	\$ 192.00
☐ 600MH34	E1033	Whitmyer® Linx-2 Headrest Mount with Quick Release	\$ 234.00
Armrest Style Options			\$
☐ 600EG53		Flip-Up Removable Armrests (Post Style) - Attaches to base. Height adj 3" & angle adj 90°-155°	\$ 462.00
☐ 600EG54		Cantilever Armrests (Reclining Flip Back) - Attaches to back. Height adj 4" & angle adj 90°	\$ 928.00
☐ 600EG55		Cantilever 2 Point Quick Release Post - Attaches front of armrest to seat frame	\$ 1,232.00
Armrest Pad Options			\$
☐ 600AB81		10" Waterfall Armrest Pad - Right	NC
☐ 600AE81		10" Waterfall Armrest Pad - Left	NC
☐ 600AB82		14" Waterfall Armrest Pad - Right	NC
☐ 600AE82		14" Waterfall Armrest Pad - Left	NC

Legrest Options			
✓ OPTION	HCPCS	DESCRIPTION	RETAIL
Center Mount Fixed Legrest			\$
Center Mount Fixed Legrest			
☐ 600H210	K0040	Center Mount (14.5" - 17" Leg Length) with Flip-Up Footplate (11" x 10")	\$ 477.00
Center Mount Fixed Legrest - Accessories			
☐ 600DI59	E0995	Fixed Center Mount Calf Pads	\$ 604.00
☐ 600LS16	E0952	Foot Strap - 1.5" wide strap	\$ 93.00
Center Mount Power Elevating Legrest			\$
Center Mount Power Elevating Legrest (includes calf pads)			
☐ 600H215	E1012	Center Mount Power Elevating (11" - 17" Leg Length)- 1 Piece footplate (13.5" x 10")	\$ 4,384.00
☐ 600H216	E1012	Center Mount Power Elevating (11" - 17" Leg Length) - 2 Piece footplate (5" x 10.5" each)	\$ 4,250.00
Swing-Away Legrest			\$
Swing-Away Legrest			
☐ 600H211	K0040	Swing-Away Legrest (12" - 17" Leg Length) with flip-up footplates (6" x 8.5")	\$ 835.00
Swing-Away Power Elevating Legrest			\$
Swing-Away Power Elevating Legrest			
☐ 600H214	E1010	Power Elevating Legrest (12" - 17" Leg Length) with flip-up footplates (6" x 8.5")	\$ 4,783.00
S/A Power Legrest Positioning Options			
☐ 600DI58	K0042	2" Metal Heel Guards for 2 Piece Footplates - May not fit with all seat sizes	\$ 56.00
☐ 600HL15	E0951	Heel Loops - Attaches to the footplate. Length 11.5-13" and Width 1"	\$ 132.00
Accessories			\$
☐ 6008A93		Waterproof Wheelchair Cover	\$ 173.00
☐ 6008A94		Waterproof Joystick Cover	\$ 54.00
☐ 6008A95		Medical Necessity Bag - Velcro's to armrest (10.5" x 6" x 1.25")	\$ 40.00
☐ 600A221		Deluxe Backpack	\$ 113.00
☐ 600A615		Cup Holder - Self-leveling. Attaches to armrest opposite joystick	\$ 187.00
☐ 600A719	E2208	Articulating Oxygen Holder	\$ 1,230.00
☐ 600A612	E2207	Cane Holder - Includes push handles	\$ 791.00
☐ 6009A16		Luggage Rack for X8 - Attaches to base behind the chair (10.5" x 13")	\$ 393.00
☐ 6009A18		Lighting Kit for X8 - 2 headlights and 2 red marker lights with turn signals Not available with flip-down knee supports	\$ 1,395.00
☐ 6008A97		Charge Mate - 12 Volt/Car outlet and USB/cell phone charger with max draw of 6 amps	\$ 630.00
☐ 600A613		Cell Phone Holder - Adjustable to fit most cell phones	\$ 345.00
☐ 600A614		Fishing Rod Holder	\$ 345.00
Be Adaptive Products			
☐ 6008A98		Be Adaptive Limited Mobility Shooting Mount - Gun or crossbow	\$ 1,760.00
☐ 6008A99		Be Adaptive Trigger Mechanism	\$ 975.00
☐ 6008A100		Be Adaptive Limited Mobility Fishing Mount	\$ 1,760.00

Electronics									
✓ OPTION		HCPCS		DESCRIPTION				RETAIL	
Joystick									\$
☐ 600J121				Standard LED Joystick (w/o Screen Display)				NC	
☐ 600J116				Color Joystick (2.3" Color LCD Screen Display, CJSM1)				NC	
☐ 600J118				High Resolution Color Joystick,(CJSM2+BT) 3.5" Color LCD Display				\$ 1,230.00	
Joystick Mount									\$
☐ 600JM34				Fixed on Armrest - Depth adjustable 8.25"				NC	
☐ 600JM45		E1032		Swing-Away Joystick Mount - Includes latch to lock in the forward position				\$ 665.00	
☐ 600JM46		E1032		Midline Joystick Mount - 2 position mount (Rotation 0° - 180°)				\$ 1,608.00	
Joystick Mount Location									\$
☐ 600HM16				Right Hand Armrest				NC	
☐ 600HM17				Left Hand Armrest				NC	
Display Language Selection R-net									\$
☐ 600L1	English	NC		☐ 600L3	French	NC		☐ 600L7	Spanish NC
Joystick Handle Type									\$
Comes standard with PG Cone-Shaped Joystick Handle									
☐ 600JH13		E2323		Ball Handle				\$ 147.00	
☐ 600JH26		E2323		Bodypoint Grooved Mushroom Joystick Handle				\$ 191.00	
☐ 600JH24		E2323		Bodypoint – Rubber Dome Joystick Handle				\$ 191.00	
☐ 600JH25		E2323		Bodypoint – U-Shaped Joystick Handle (4 inch)				\$ 191.00	
Expandable Power Module									\$
Must select both items									
☐ 600EX05		E2377		Expandable 120A Controller - R-net				\$ 786.00	
☐ 600EX06		E2313		Harness for 120A Expandable Controller				\$ 710.00	
Power Seat Function									\$
☐ 600EO37		E2310		Thru-Drive Control for 1 Actuator Must be selected if ordering 1 actuator function				\$ 2,306.00	
☐ 600EO38		E2311		Thru-Drive Control for 2+ Actuators Must be selected if ordering 2 or more actuator functions				\$ 4,174.00	
Power Seat Control Option									\$
☐ 600QG11				Ctrl+5 Switch Box with Voice Activation				\$ 1,459.00	
☐ 600QG12				Ctrl+5 Switch Box with Voice Activation (Multi-power only)				NC	
Ctrl+5 Switch Type									\$
☐ 600EY22				Push Button				NC	
☐ 600EY23		K0108		Toggle Switches				\$ 350.00	
Power Seat Control Option - Ctrl+5 Mounting location									\$
☐ 600EY19	Right Armrest			NC		☐ 600EY20	Left Armrest NC		
Attendant Controls									\$
☐ 6000A05		E2331		Dual Control Module for Attendant - Includes accessory mount and push handle				\$ 1,487.00	
Attendant Control Mounting Location									\$
☐ 600DN5	Right Back Cane			NC		☐ 600DN6	Left Back Cane NC		
R-net Omni2 - 3.5" LCD Display w/ Built-In Infrared & Bluetooth									\$
☐ 600ON14				R-net Omni2 Only Input Control Module				\$ 2,072.00	
☐ 600ON15				R-net Omni2 Bundled Input Control Module Available w/ Switch-It Input Devices				NC	
R-net Omni2 - Mounting Location									\$
☐ 600DS23	Right Hand Armrest			NC		☐ 600DS24	Left Hand Armrest NC		
Head Control Packages									\$
☐ 600EU25		E2328		Switch-It DualPro® Proportional Head Array - Facial Pads				\$ 6,668.00	
☐ 600EU26		E2328		Switch-It DualPro® Proportional Head Array - Spot Pads				\$ 6,668.00	