



Customer				Ship To			
Account				Name			
Account #		Date		Address			
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #		City		State	
Buyer				ZIP	Phone		
Mark For				Email / Fax			

Please submit orders to: magicmobility@sunmed.com

Seating Specifications


HPCP	N/A	Typical WC Weight w/ Batteries	Starting at 320 lbs. Weights vary depending on options
Power Seating	50° of CG Tilt / Power Recline / 12" Seat Elevator	Seat Type	MPS and Rehab Style
Seat Size Range	Width: 13-23" Depth: 13-24"	STF Height Range	18.5", 19.5", 20.5"

Drive Base Specifications

Drive Wheel Position	4 Wheel Drive	Drive Control	R-Net PM 120A controller
User Weight Capacity	400 lbs. 340 lbs. with Elevate or Tilt	Drive Motors	4-pole 6mph
X8 Base Size	Width 28" Length 40.5" Turning Radius 56"	Drive Wheels	14" All Terrain Drive Wheels
		Tires	All Terrain Tires
		Charger	Dual Mode
		Battery Size	Group 24
		Color Accents	6 shroud colors and 6 color accents
		Ground Clearance	4"

Client Measurements

- a) Top of head to seat in
- b) Top of shoulder to seat in
- c) Underarm to seat in
- d) Elbow to seat in
- e) Top of leg to seat in
- f) Shoulder width in
- g) Chest width in
- i) Hip width in
- j) Actual sitting depth in
- m) Foot length in
- n) Knee to foot in
- o) Seat width in



Notes

✓ OPTION	HCPCS	DESCRIPTION	RETAIL		
Seating Options			\$		
<input type="checkbox"/> 600SF4		REHAB Seat and Back - Back Canes and adjustable flat seat pan	NC		
<input type="checkbox"/> 600SF5		MPS Seat and Back - Captain style seating	NC		
<input type="checkbox"/> 600SF6		MPS Back and REHAB Seat - Captain style backrest and adjustable flat seat pan	NC		
Base					
Extreme X8 Base		End Item Number: EIPW60	\$		
<input type="checkbox"/> 600MW79		Magic Mobility - Extreme X8 No Power Seating Options	\$ 19,644.00		
<input type="checkbox"/> 600MW80		Magic Mobility - Extreme X8 Single Power Seating Options	\$ 19,644.00		
<input type="checkbox"/> 600MW81		Magic Mobility - Extreme X8 Multi Power Seating Options	\$ 19,644.00		
<input type="checkbox"/> 600MW82		Magic Mobility - Extreme X8 No Power Seating Options - HD	\$ 19,644.00		
<input type="checkbox"/> 600MW83		Magic Mobility - Extreme X8 Single Power Seating Options - HD	\$ 19,644.00		
<input type="checkbox"/> 600MW84		Magic Mobility - Extreme X8 Multi Power Seating Options - HD	\$ 19,644.00		
<input type="checkbox"/> 600MW85		Magic Mobility - Extreme X8 No Power Seating Options - MPS (Captain)	\$ 19,644.00		
<input type="checkbox"/> 600MW86		Magic Mobility - Extreme X8 Single Power Seating Options - MPS (Captain)	\$ 19,644.00		
<input type="checkbox"/> 600MW87		Magic Mobility - Extreme X8 Multi Power Seating Options - MPS (Captain)	\$ 19,644.00		
<input type="checkbox"/> 600MW88		Magic Mobility - Extreme X8 No Power Seating Options - HD MPS (Captain)	\$ 19,644.00		
<input type="checkbox"/> 600MW89		Magic Mobility - Extreme X8 Single Power Seating Options - HD MPS (Captain)	\$ 19,644.00		
<input type="checkbox"/> 600MW90		Magic Mobility - Extreme X8 Multi Power Seating Options - HD MPS (Captain)	\$ 19,644.00		
Base Accessories			\$		
<input type="checkbox"/> 6008A91		Jack for Drive Tire	\$ 203.00		
Tie Downs			\$		
<input type="checkbox"/> 600TO52		Crash Tested X8 Tie Down Kit	\$ 790.00		
Additional Base Function			\$		
<input type="checkbox"/> 600WL93		Directional Lock X8 Base - Used to reverse in straight line.	\$ 2,005.00		
Optional Drive Tire Options			\$		
<input type="checkbox"/> 6008A101		Tire Sealant for 4 Pneumatic Drive Tires	\$ 76.00		
Shroud Color			\$		
<input type="checkbox"/> 600C185	Midnight Black Shroud Color	NC	<input type="checkbox"/> 600C197	Clover Green Shroud Color	NC
<input type="checkbox"/> 600C99	Ocean Blue Shroud Color	NC	<input type="checkbox"/> 600C198	Tangerine Shroud Color	NC
<input type="checkbox"/> 600C196	Fire Engine Red Shroud Color	NC	<input type="checkbox"/> 600C199	Camouflage Shroud Color	\$ 608.00
Seat To Floor Height			\$		
X8 Seat to Floor Heights start at 18.5" +/- 0.25"					
<input type="checkbox"/> 600HT44		Lowest Seat To Floor Height	NC		
<input type="checkbox"/> 600HT45		Seat To Floor Height +1"	NC		
<input type="checkbox"/> 600HT46		Seat To Floor Height +2"	NC		
Batteries			\$		
<input type="checkbox"/> 600AG25	E2363	Group 24 Gel Battery (Pair) - Installed	\$ 1,116.00		
Battery Charger			\$		
<input type="checkbox"/> 600BJ28		10A Battery Charger	NC		

Seating: Seat Functions			
✓ OPTION	HCPCS	DESCRIPTION	RETAIL
No Power Seat Options			\$
<input type="checkbox"/> 600ST86		No Power Seating System - 400 lb. weight capacity	\$ 4,138.00
<input type="checkbox"/> 600ST87		No Power Seating System with Power Functionality - 400 lb. weight capacity Required for lights, power recline, power legrest, steering lock, or future power seat functions	\$ 5,320.00
Power Tilt			\$
Utilized without Power Recline.			
<input type="checkbox"/> 600PA64	E1002	50° Power Tilt in Space - 340 lb. weight capacity	\$ 8,292.00
Power Seat Elevate			\$
Only available when Power Tilt or Power Tilt and Recline is selected.			
<input type="checkbox"/> 600PI16	E2300	12" Adjustable Seat Elevator - 340 lb. weight capacity	\$ 4,750.00
Power Recline			\$
Utilized without Power Tilt. Requires No Power Seating with Power Functionality.			
<input type="checkbox"/> 600PA65	E1003	150° Power Recline Only - 300 lb. weight capacity	\$ 6,750.00
Power Tilt and Recline			\$
<input type="checkbox"/> 600PA66	E1006	50° Tilt & 150° Power Recline - 300 lb. weight capacity	\$ 12,750.00

Seating: Seat Sizes			
MPS (Captain) Seat Cushion (width x depth)			\$
MPS Backrest can be paired with the REHAB flat pan. If flat pan is desired, omit the MPS cushion and select the desired REHAB flat pan.			
<input type="checkbox"/> 600SI24		16" x 18" Vinyl Seat	NC
<input type="checkbox"/> 600SI25		18" x 18" Vinyl Seat	NC
<input type="checkbox"/> 600SI26		18" x 20" Vinyl Seat	NC
<input type="checkbox"/> 600SI27		20" x 18" Vinyl Seat	NC
<input type="checkbox"/> 600SI28		20" x 20" Vinyl Seat	NC
<input type="checkbox"/> 600SI29		22" x 20" Vinyl Seat	\$ 463.00
REHAB Flat Pan Seat Sizes - Adjustable in width and depth for most sizes			\$

REHAB Seat Pan

Standard seat pan available 13-23" wide. Deep seat pan available 13"-23" wide.

<input type="checkbox"/> 600PN13		Standard Adjustable Flat Seat Pan - 13" to 19" deep	NC
<input type="checkbox"/> 600PN14		Deep Adjustable Flat Seat Pan - 20" to 22" deep	\$ 463.00
<input type="checkbox"/> 600PM15		Deep Adjustable Flat Seat Pan - 23" to 24" deep	\$ 817.00

REHAB Seat Width

- 13"
 14"
 15"
 16"
 17"
 18"
 19"
 20"
 21"
 22"
 23"

REHAB Seat Depth

- 13"
 14"
 15"
 16"
 17"
 18"
 19"
 20"
 21"
 22"
 23"
 24"

REHAB Seating Cushion Options			
Cushion size will be determined by Rehab Seat size selected above			
<input type="checkbox"/> 600CU59	E2607	JAY® Ion® (Sizes 16" - 21" Wide)	\$ 451.00
<input type="checkbox"/> 600CU60	E2608	JAY® Ion® (Size 22" Wide)	\$ 567.00
<input type="checkbox"/> 600CU61	E2607	JAY® Union® (Sizes 16" - 21" Wide)	\$ 510.00
<input type="checkbox"/> 600CU62	E2608	JAY® Union® (Size 22" Wide)	\$ 623.00
<input type="checkbox"/> 600CU64	E2622	JAY® Fusion® (Sizes 16" - 21" Wide)	\$ 504.00
<input type="checkbox"/> 600CU65	E2623	JAY® Fusion® (Size 22" Wide)	\$ 601.00

✓ OPTION	HCPCS	DESCRIPTION	RETAIL
Lateral Thigh Support Pads and Hardware (must pick pad and mount)			\$
Lateral Thigh Support Pads			
<input type="checkbox"/> 600KF22	E0953	Lateral Thigh Support Pad Short - Right (5"x5")	\$ 144.00
<input type="checkbox"/> 600KF23	E0953	Lateral Thigh Support Pad Short - Left (5"x5")	\$ 144.00
Lateral Thigh Support Mounts			
<input type="checkbox"/> 600KF26	E1028	Lateral Thigh Removable Mount - Right	\$ 286.00
<input type="checkbox"/> 600KF27	E1028	Lateral Thigh Removable Mount - Left	\$ 286.00
Seating: Seat Back Options			
MPS Backrest Options - Comes standard with Manual Recline			\$
Backrest Width			
<input type="checkbox"/> 600BN46		Standard Backrest Width - 16-20" Wide	NC
<input type="checkbox"/> 600BN47		Wide Backrest Width - 22-23" Wide	\$ 336.00
Backrest Height			
<input type="checkbox"/> 600B80		Standard Backrest Height - 22" Tall	NC
<input type="checkbox"/> 600B81		Tall Backrest Height - 26" Tall	\$ 558.00
Backrest Material			
<input type="checkbox"/> 600BU106		Cloth Backrest	NC
<input type="checkbox"/> 600BU107		Vinyl Backrest	\$ 304.00
MPS Back Accessories			\$
<input type="checkbox"/> 600PH33		Push Handle - Needed for O2 holder	\$ 159.00
<input type="checkbox"/> 6009A13		Accessory Mount - Needed for aftermarket headrest	\$ 121.00
REHAB Seating Backrest Options			\$
Backrest will have the width of the MPS seat or flat pan. Please choose your back height. Backrest will come with 2 point hardware with a 500 lb. weight capacity.			
JAY J3 Back Contour			
<input type="checkbox"/> 600JP01	E2313	J3® PA Back for 16"-20" Wide - 2.2" of lateral thoracic support	\$ 613.00
<input type="checkbox"/> 600JP02	E2614	J3® PA Back for 22" Wide - 2.2" of lateral thoracic support	\$ 910.00
<input type="checkbox"/> 600JP11	E2615	J3® PL Back for 16"-20" Wide - 3" of lateral thoracic support	\$ 715.00
<input type="checkbox"/> 600JP13	E2620	J3® PD Back for 16"-20" Wide (18" and 21" shell height only) - 6" of lateral thoracic support	\$ 858.00
JAY® J3® Back Shell Options and Height Ranges			
<input type="checkbox"/> 600JR14		16.5" tall shell (Upper Thoracic Short)	NC
<input type="checkbox"/> 600JR15		18" tall shell (Upper Thoracic Med)	NC
<input type="checkbox"/> 600JR16		19.5" tall shell (Upper Thoracic Tall)	NC
<input type="checkbox"/> 600JR17		21" tall shell (Shoulder Height)	NC
<input type="checkbox"/> 600JR27		22.5" tall shell (Shoulder Height)	NC
<input type="checkbox"/> 600JR28		24" tall shell (Shoulder Height)	NC
REHAB Flip Forward Backrest Option			\$
<input type="checkbox"/> 600D23		REHAB Flip Forward Backrest Option	\$ 296.00
JAY® Lateral Thoracic Swing-Away Supports (Requires Pad and Mount) - J3 back only			\$
JAY® Lateral Thoracic Support Pads			
<input type="checkbox"/> 600JW35	E0956	Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Right	\$ 144.00
<input type="checkbox"/> 600JW38	E0956	Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Left	\$ 144.00
<input type="checkbox"/> 600JW37	E0956	Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Right	\$ 144.00
<input type="checkbox"/> 600JW40	E0956	Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Left	\$ 144.00
JAY® Lateral Thoracic Swing-Away Mount			
<input type="checkbox"/> 600JW41	E1028	Lateral Thoracic Swing Away Bracket - Right	\$ 286.00
<input type="checkbox"/> 600JW42	E1028	Lateral Thoracic Swing Away Bracket - Left	\$ 286.00

JAY® Anterior Thoracic Support				\$
<input type="checkbox"/>	600ZR10	E0960	Bodypoint Monoflex Medium Chest Strap - Pad 3.5" x 18.5" (total length 49.5")	\$ 155.00
<input type="checkbox"/>	600ZR11	E0960	Bodypoint Monoflex Large Chest Strap - Pad 3.5" x 21.5" (total length 52.5")	\$ 155.00
<input type="checkbox"/>	600K10	E0960	JAY® Anterior Thoracic Support - Medium ATS8302	\$ 174.00
<input type="checkbox"/>	600K11	E0960	JAY® Anterior Thoracic Support - Large ATS8303	\$ 174.00
Pelvic Positioning Options - Select One				\$
Pelvic Positioning				
<input type="checkbox"/>	600A848		Positioning Belt, Auto Style Buckle, 50 mm 1600 mm / 2"	NC
<input type="checkbox"/>	600A849	E0978	Positioning Belt, Aircraft Buckle, Padded 50 mm x 1300 mm / 2"	\$ 86.00
<input type="checkbox"/>	600A850	K0108	Bodypoint - Evoflex Pelvic Stabilizer Belt	\$ 199.00
<input type="checkbox"/>	600A851	K0108	Bodypoint - Evoflex Pelvic Stabilizer Belt with Extension	\$ 220.00
Headrest Options - MPS Style				\$
<input type="checkbox"/>	600NU12	E0955	Cloth Headrest - MPS seat back only (12" x 4.5" x 6.25")	\$ 411.00
<input type="checkbox"/>	600NU13	E0955	Vinyl Headrest - MPS seat back only (12" x 4.5" x 6.25")	\$ 411.00
Headrest Options - REHAB Style				\$
Headrest Pad				
<input type="checkbox"/>	600HE39	E0955	Whitmyer® Plush Headrest Pad 254 mm / 10"	\$ 234.00
<input type="checkbox"/>	600HE46	E0955	Whitmyer® Adjust-A-Plush Headrest Pad Medium (Adjustable from 7 to 16" Wide x 4" Tall)	\$ 360.00
Headrest Mount				
<input type="checkbox"/>	600MH33	E1028	Whitmyer® Linx Headrest Mount with Quick Release	\$ 188.00
<input type="checkbox"/>	600MH34	E1028	Whitmyer® Linx-2 Headrest Mount with Quick Release	\$ 229.00
Armrest Style Options				\$
<input type="checkbox"/>	600EG53		Flip-Up Removable Armrests (Post Style) - Attaches to base. Height adj 3" & angle adj 90°-155°	\$ 440.00
<input type="checkbox"/>	600EG54		Cantilever Armrests (Reclining Flip Back) - Attaches to back. Height adj 4" & angle adj 90°	\$ 884.00
<input type="checkbox"/>	600EG55		Cantilever 2 Point Quick Release Post - Attaches front of armrest to seat frame	\$ 1,173.00
Armrest Pad Options				\$
<input type="checkbox"/>	600AB80		4" x 16" Tapered Pads - Right	\$ 54.00
<input type="checkbox"/>	600AE80		4" x 16" Tapered Pads - Left	\$ 54.00
<input type="checkbox"/>	600AB81		10" Waterfall Armrest Pad - Right	NC
<input type="checkbox"/>	600AE81		10" Waterfall Armrest Pad - Left	NC
<input type="checkbox"/>	600AB82		14" Waterfall Armrest Pad - Right	NC
<input type="checkbox"/>	600AE82		14" Waterfall Armrest Pad - Left	NC
Legrest Options				
<input checked="" type="checkbox"/>	OPTION	HCPCS	DESCRIPTION	RETAIL
Center Mount Fixed Legrest				\$
Center Mount Fixed Legrest				
<input type="checkbox"/>	600H210	K0040	Center Mount Flip-Up Footplate	\$ 454.00
Center Mount Fixed Legrest - Calfpads				
<input type="checkbox"/>	600DI59	E0995	Fixed Center Mount Calf Pads	\$ 575.00
Center Mount Power Elevating Legrest				\$
Center Mount Power Elevating Legrest				
<input type="checkbox"/>	600H215	E1012	X8 - Center Mount Power Elevating - 1 Piece footplate - Includes calf pads (11" x 11.5")	\$ 4,175.00
<input type="checkbox"/>	600H216	E1012	X8 - Center Mount Power Elevating - 2 Piece footplate - Includes calf pads (5" x 10.5")	\$ 4,048.00

Swing-Away Legrest				\$							
Swing-Away Legrest											
<input type="checkbox"/>	600H211	K0040	Side Mount Swing-Away Legrest - 2 piece flip-up footplate (8.5" x 6"). 13" height adjustable	\$ 795.00							
S/A Non-Power Positioning Options											
<input type="checkbox"/>	600LS15	E0951	Heel and Side Support - Attaches to the footplate stem and footplate	\$ 657.00							
<input type="checkbox"/>	600LS16	E0952	Foot Strap - 1.5" wide strap	\$ 89.00							
Swing-Away Power Elevating Legrest				\$							
Swing-Away Power Elevating Legrest											
<input type="checkbox"/>	600H214	E1010	Side Mount Power Elevating Legrest - Swing away legrest with flip-up 2" height adj footplate	\$ 4,555.00							
S/A Power Legrest Positioning Options											
<input type="checkbox"/>	600HL15	E0951	Heel Loops - Attaches to the footplate. Length 11.5-13" and Width 1"	\$ 126.00							
Accessories				\$							
<input type="checkbox"/>	6008A93		Waterproof Wheelchair Cover	\$ 165.00							
<input type="checkbox"/>	6008A94		Waterproof Joystick Cover	\$ 51.00							
<input type="checkbox"/>	6008A95		Medical Necessity Bag - Velcro's to armrest (10.5" x 6" x 1.25")	\$ 38.00							
<input type="checkbox"/>	600A221		Deluxe Backpack	\$ 108.00							
<input type="checkbox"/>	600A615		Cup Holder - Self-leveling. Attaches to armrest opposite joystick	\$ 178.00							
<input type="checkbox"/>	600A719	E2208	Articulating Oxygen Holder	\$ 1,171.00							
<input type="checkbox"/>	600A612	E2207	Cane Holder - Includes push handles	\$ 753.00							
<input type="checkbox"/>	6009A16		Luggage Rack for X8 - Attaches to base behind the chair (10.5" x 13")	\$ 374.00							
<input type="checkbox"/>	6009A18		Lighting Kit for X8 - 2 headlights and 2 red marker lights with turn signals Not available with flip-down knee supports	\$ 1,329.00							
<input type="checkbox"/>	6008A97		Charge Mate - 12 Volt/Car outlet and USB/cell phone charger with max draw of 6 amps	\$ 600.00							
<input type="checkbox"/>	600A613		Cell Phone Holder - Adjustable to fit most cell phones	\$ 329.00							
<input type="checkbox"/>	600A614		Fishing Rod Holder	\$ 329.00							
Be Adaptive Products											
<input type="checkbox"/>	6008A98		Be Adaptive Limited Mobility Shooting Mount - Gun or crossbow	\$ 1,676.00							
<input type="checkbox"/>	6008A99		Be Adaptive Trigger Mechanism	\$ 929.00							
<input type="checkbox"/>	6008A100		Be Adaptive Limited Mobility Fishing Mount	\$ 1,676.00							
Electronics											
<input checked="" type="checkbox"/>	OPTION	HCPCS	DESCRIPTION	RETAIL							
Joystick				\$							
<input type="checkbox"/>	600J121		Standard LED Joystick (w/o Screen Display)	NC							
<input type="checkbox"/>	600J116		Color Joystick (2.3" Color LCD Screen Display, CJSM1)	NC							
<input type="checkbox"/>	600J118		High Resolution Color Joystick,(CJSM2+BT) 3.5" Color LCD Display	\$ 1,206.00							
Joystick Mount				\$							
<input type="checkbox"/>	600JM34		Fixed on Armrest - Depth adjustable 8.25"	NC							
<input type="checkbox"/>	600JM45	E1028	Swing-Away Joystick Mount - Includes latch to lock in the forward position	\$ 633.00							
<input type="checkbox"/>	600JM46	E1028	Midline Joystick Mount - 2 position mount (Rotation 0° - 180°)	\$ 1,531.00							
Joystick Mount Location				\$							
<input type="checkbox"/>	600HM16		Right Hand Armrest	NC							
<input type="checkbox"/>	600HM17		Left Hand Armrest	NC							
Display Language Selection R-net				\$							
<input type="checkbox"/>	600L1	English	NC	<input type="checkbox"/>	600L3	French	NC	<input type="checkbox"/>	600L7	Spanish	NC
Joystick Handle Type				\$							
Comes standard with PG Cone-Shaped Joystick Handle											
<input type="checkbox"/>	600JH13	E2323	Ball Handle	\$ 144.00							
<input type="checkbox"/>	600JH26	E2323	Bodypoint Grooved Mushroom Joystick Handle	\$ 187.00							
<input type="checkbox"/>	600JH24	E2323	Bodypoint – Rubber Dome Joystick Handle	\$ 187.00							
<input type="checkbox"/>	600JH25	E2323	Bodypoint – U-Shaped Joystick Handle (4 inch)	\$ 187.00							

Expandable Power Module					\$	
Must select both items						
<input type="checkbox"/>	600EX05	E2377	Expandable 120A Controller - R-net		\$ 749.00	
<input type="checkbox"/>	600EX06	E2313	Harness for 120A Expandable Controller		\$ 676.00	
Power Seat Function					\$	
<input type="checkbox"/>	600EO37	E2310	THRU DRIVE CNTRL 1 ACTUATOR- Must be selected if ordering 1 actuator function		\$ 2,196.00	
<input type="checkbox"/>	600EO38	E2311	THRU DRIVE CNTRL 2+ ACTUATORS- Must be selected if ordering 2 or more actuator functions		\$ 3,975.00	
Power Seat Control Option					\$	
<input type="checkbox"/>	600QG11		Ctrl+5 Switch Box with Free Assignable Functions		\$ 1,430.00	
<input type="checkbox"/>	600QG12		Ctrl+5 Switch Box with Free Assignable Functions (Multi-power only)		NC	
Ctrl+5 Switch Type					\$	
<input type="checkbox"/>	600EY22		Push Button		NC	
<input type="checkbox"/>	600EY23	K0108	Toggle Switches		\$ 343.00	
Power Seat Control Option - Mounting location					\$	
<input type="checkbox"/>	600EY19	CTRL+5 MNT ARMREST RIGHT	NC	<input type="checkbox"/>	600EY20 CTRL+5 MNT ARMREST LEFT	NC
Attendant Controls					\$	
<input type="checkbox"/>	6000A05	E2331	Dual Control Module for Attendant - Includes accessory mount and push handle		\$ 1,487.00	
Attendant Control Mounting Location					\$	
<input type="checkbox"/>	600DN5	Right Back Cane	NC	<input type="checkbox"/>	600DN6 Left Back Cane	NC
R-net Omni-2 - 3.5" LCD Display w/ Built-In Infrared & Bluetooth					\$	
<input type="checkbox"/>	6000N14		R-net Omni-2 Only Input Control Module		\$ 2,031.00	
<input type="checkbox"/>	6000N15		R-net Omni2 Bundled Input Control Module Available w/ Switch-It Input Devices		NC	
R-net Omni-2 -Mounting Location					\$	
<input type="checkbox"/>	600DS23	Right Hand Armrest	NC	<input type="checkbox"/>	600DS24 Left Hand Armrest	NC
Head Control Packages					\$	
<input type="checkbox"/>	600EU25	E2328	Switch-It DualPro® Proportional Head Array - Facial Pads		\$ 6,537.00	
<input type="checkbox"/>	600EU26	E2328	Switch-It DualPro® Proportional Head Array - Spot Pads		\$ 6,537.00	