



Quickie M6™ February 2023

EIF1

Mark For:

Date: _____
 Dealer Acct #: _____
 Dealer: _____
 Dealer Contact: _____
 Dealer Address: _____
 Dealer City: _____ ST: _____ ZIP: _____
 Dealer Phone: () _____ Fax: () _____
 Confirmation Email: _____
 Confirm Via: Fax Email

Submitting for: Quote Order

PO#: _____
ADDITIONAL SHIPPING INFORMATION
 Ship To: _____
 Attention: _____
 Address: _____
 Address: _____
 Ship To City: _____ ST: _____ ZIP: _____
 Ship To Phone: () _____ Fax: () _____

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

Option #	HCPCS	Definition	Price	Option #	HCPCS	Definition	Price
----------	-------	------------	-------	----------	-------	------------	-------

1. FRAME

Base Model

<input type="checkbox"/>	115WM21	K0007	Quickie M6	\$3,414.00
--------------------------	---------	-------	------------	------------

2. SIZE / SEATING

Frame Width

<input type="checkbox"/>	115W13	E2201	22"	\$559.00
<input type="checkbox"/>	115W14	E2202	24"	\$559.00
<input type="checkbox"/>	115W15	E2202	26"	\$614.00
<input type="checkbox"/>	115W16	K0108	28"	\$614.00
<input type="checkbox"/>	115W17	K0108	30"	\$614.00

Frame Length

<input type="checkbox"/>	115FL9		Regular (18"-19" depth)	NC
<input type="checkbox"/>	115FL18		Long (18"-22" depth)	NC

Seat Depth

<input type="checkbox"/>	115SD11		18"	NC
<input type="checkbox"/>	115SD12		19"	NC
<input type="checkbox"/>	115SD13	E2203	20" (Long frame only)	\$614.00
<input type="checkbox"/>	115SD14	E2203	21" (Long frame only)	\$614.00
<input type="checkbox"/>	115SD33	E2204	22" (Long frame only)	\$1,006.00

Seat Height

<input type="checkbox"/>	115HT65		17" Includes low adjustable axle plate and shorter fork. Only available with 5" caster	\$335.00
<input type="checkbox"/>	115HT11		18.5" Only available w/ 5" casters	NC
<input type="checkbox"/>	115HT14		20"	NC

Cushion

<input type="checkbox"/>	115CU31	K0669	2" Cushion	\$168.00
--------------------------	---------	-------	------------	----------

JAY Cushions (Optional)

To order any JAY cushions on-chair, please see the JAY price list.
 Cushion name or part number must be filled out with desired width and depth

<input type="checkbox"/>	115JC20	JAY Bundled Cushion		Price List
		Cushion Name or Part #		
		Width: _____	Depth: _____	

3. CASTERS

Castors

<input type="checkbox"/>	115CW6		5" Soft Roll	\$96.00
<input type="checkbox"/>	115CW36	E2214 E2215	8" x 2" Pneumatic	\$162.00
<input type="checkbox"/>	115CW37	E2214 K0108	8" x 2" Pneumatic w/Airless Insert	\$252.00

Seat Height Matrix

Caster	Front/Rear
5"	17" (Different Axle & Fork)
5"	18.5"
8"	20"

4. REAR WHEELS

Rear Wheel Type

<input type="checkbox"/>	115RW2		24" Spoke Wheel	NC
<input type="checkbox"/>	115RW29		24" Spinergy Wheel	\$1,185.00

Adjustable Axle Plate

<input type="checkbox"/>	115AL2		Standard (for 18.5", 20" Seat Ht.)	NC
<input type="checkbox"/>	115AL28		Low (only for 17" Seat Ht.)	NC

Rear Wheel Tire

<input type="checkbox"/>	115RT1	E2211 E2212	Pneumatic Tire	NC
<input type="checkbox"/>	115RT2	E2211 E2213	Pneumatic Airless Insert Tire	\$96.00
<input type="checkbox"/>	115RT25	E2218	Primo Orion Poly	\$121.00

Axle

<input type="checkbox"/>	115AX1		Quick Release	NC
--------------------------	--------	--	---------------	----

Rear Tire & Wheel Matrix

Tire	Rear Wheel	
	Spoke wheel	Spinergy wheel
Pneumatic	24	24
Pneumatic w/ Airless Insert	24	24
Primo Orion Poly	24	24

5. WHEEL LOCKS				
Wheel Lock			Extension Handle	
<input type="checkbox"/>	115WL1	High Mount Push Wheel Lock	NC	
<input type="checkbox"/>	115WL12	Do Not Mount Wheel Lock	NC	
<input type="checkbox"/>	115WL10	Omit Wheel Lock	NC	
<input type="checkbox"/>	115EH1	E0961 6" Pair		\$84.00

6. HANGERS AND FOOTREST				
Front Hangers			Extension Tubes	
<input type="checkbox"/>	115H1	60° Swing Away	NC	
<input type="checkbox"/>	115H146	70° Heavy-Duty Ft. Mount Lift Off <i>Only available with 5" caster and ang adj footplates. NA w/ extension tubes</i>		\$106.00
<input type="checkbox"/>	115H132	80° Heavy-Duty Ft. Mount Lift Off <i>Only available with 5" caster and ang adj footplates. NA w/ extension tubes</i>		\$106.00
<input type="checkbox"/>	115H41	E0990 ELR - Heavy-Duty Reinforced <i>Only available with aluminum footplates</i>		\$665.00
<input type="checkbox"/>	115E1	Short (13"-15.5")	NC	
<input type="checkbox"/>	115E2	Medium (15.5"-18")	NC	
<input type="checkbox"/>	115E5	Long (18"-20.5") NA w/ 5" caster	NC	

Footplates				
<i>Heel loop or leg strap not included, See Footplate Options Section</i>				
<input type="checkbox"/>	115F1	Composite	NC	
<input type="checkbox"/>	115F25	K0040 Adult Angle Adj		\$212.00
<input type="checkbox"/>	115F5	Aluminum	NC	
Footplate Options				
<input type="checkbox"/>	115HL1	E0951 Heel Loop		\$45.00

Lower Leg Length Range Matrix				
Hanger	Ext Tube	Footplate		
		Composite	Alum Angle Adj	Aluminum
60°	Short	11.50" - 13.25"	11.75" - 14.00"	11.25" - 13.50"
	Medium	13.75" - 16.00"	14.00" - 16.25"	13.50" - 14.00"
	Long	16.00" - 16.50"	16.00" - 17.75"	NA
70°	No Ext Tube	NA	9.75" - 13.50"	NA
	80°	No Ext Tube	NA	10.75" - 14.75"
ELR	Short	NA	NA	13.00" - 17.00"
	Medium	NA	NA	14.25" - 17.00"

All values are ± 0.25" and based on flat seat heights and do not account for the minimum 2" ground clearance or footplate thickness. Full range not available with all seat heights. ELR values based on 70° hanger angle

7. BACKREST				
Backrest			Back Height Preset	
<input type="checkbox"/>	115BT34	8° w/ Push Handles <i>Adjustable (16-3/4" to 20-1/8")</i>	NC	
<input type="checkbox"/>	115B91	16-3/4"		NC
<input type="checkbox"/>	115B96	17-7/8"		NC
<input type="checkbox"/>	115B39	19"		NC
<input type="checkbox"/>	115B98	20-1/8"		NC

8. ARMREST				
Armrest				
<input type="checkbox"/>	115AA33	E0973 Dual Post Height Adj - Desk 10"		\$274.00
<input type="checkbox"/>	115AA34	E0973 Dual Post Height Adj - Full 14"		\$274.00

9. KEY OPTIONS & ACCESSORIES				
Rear Anti-Tip			Positioning Belt	
<input type="checkbox"/>	115AT1	E0971 Rear Anti-Tip		\$168.00
<input type="checkbox"/>	115AT3	Omit Rear Anti-Tip	NC	
<input type="checkbox"/>	115A115	E0978 Belt 2" Aircraft Buckle Long		\$123.00
<input type="checkbox"/>	115A116	E0978 Belt 2" Aircraft Buckle Padded Long		\$133.00
<input type="checkbox"/>	115A142	E0978 Belt 1.5" Auto Buckle		\$78.00
Miscellaneous				
<input type="checkbox"/>	115A83	Tool Kit		\$61.00

10. COLOR				
Standard Frame Colors			Matte Frame Colors	
<input type="checkbox"/>	115C140	Aztec Gold	NC	
<input type="checkbox"/>	115C2	Black (Not avail as base w/Zebra)	NC	
<input type="checkbox"/>	115C141	Black Cherry	NC	
<input type="checkbox"/>	115C10	Black Opal	NC	
<input type="checkbox"/>	115C38	Blue Opal	NC	
<input type="checkbox"/>	115C36	Candy Blue	NC	
<input type="checkbox"/>	115C14	Candy Purple	NC	
<input type="checkbox"/>	115C35	Candy Red	NC	
<input type="checkbox"/>	115C142	Electric Blue	NC	
<input type="checkbox"/>	115C143	Evergreen	NC	
<input type="checkbox"/>	115C126	Glow	NC	
<input type="checkbox"/>	115C144	Green Apple	NC	
<input type="checkbox"/>	115C147	Hot Sparkle Pink	NC	
<input type="checkbox"/>	115C148	Mauve Pink	NC	
<input type="checkbox"/>	115C103	Rootbeer	NC	
<input type="checkbox"/>	115C114	Sparkle Silver	NC	
<input type="checkbox"/>	115C145	Sunrise Orange	NC	
<input type="checkbox"/>	115C113	Titanium Color Paint	NC	
<input type="checkbox"/>	115C8	Yellow	NC	
<input type="checkbox"/>	115C112	Desert Camouflage		\$441.00
<input type="checkbox"/>	115C43	Mossy Oak Camouflage		\$441.00
<input type="checkbox"/>	115C42	Stars & Stripes		\$441.00
<input type="checkbox"/>	115C108	Zebra (Choose above base color) Std with White unless above base paint is chosen NA w/ Blck, Matte Black, Blck Opal or Titanium)		\$441.00



Sunrise Medical (US) LLC
2842 Business Park Ave. · Fresno, CA 93727 · USA