



Pal Classroom Seat

August 2020

Mark For:

Date: _____

Dealer Acct #: _____

Dealer: _____

Dealer Contact: _____

Dealer Address: _____

Dealer City: _____ ST: _____ ZIP: _____

Dealer Phone: _____ Fax: _____

Confirmation Email: _____

Confirm Via: Fax EmailSubmitting for: Quote Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____

Attention: _____

Address: _____

Address: _____

Ship To City: _____ ST: _____ ZIP: _____

Ship To Phone: _____ Fax: _____

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

Pal Classroom Seat

Size 1 Weight Capacity: 66 lbs | Size 2 Weight Capacity: 77 lbs | Size 3 Weight Capacity: 88 lbs | Size 4 Weight Capacity: 110 lbs

Pal Classroom Seat

Selection Required. Includes multi-functional casters.

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	135-1600-03	Size 1 PAL Classroom Seat - Blue		\$971
<input type="checkbox"/>	135-1600-01	Size 1 PAL Classroom Seat - Green		\$971
<input type="checkbox"/>	135-2600-03	Size 2 PAL Classroom Seat - Blue		\$1,066
<input type="checkbox"/>	135-2600-01	Size 2 PAL Classroom Seat - Green		\$1,066
<input type="checkbox"/>	135-3600-03	Size 3 PAL Classroom Seat - Blue		\$1,166
<input type="checkbox"/>	135-3600-01	Size 3 PAL Classroom Seat - Green		\$1,166
<input type="checkbox"/>	135-4600-03	Size 4 PAL Classroom Seat - Blue		\$1,358
<input type="checkbox"/>	135-4600-01	Size 4 PAL Classroom Seat - Green		\$1,358

Side Pad

Must match seat size.

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	135-1635-03	Size 1 PAL Side Pads 1" Thick - Blue (Pair)		\$117
<input type="checkbox"/>	135-1635-01	Size 1 PAL Side Pads 1" Thick - Green (Pair)		\$117
<input type="checkbox"/>	135-2635-03	Size 2 PAL Side Pads 1" Thick - Blue (Pair)		\$117
<input type="checkbox"/>	135-2635-01	Size 2 PAL Side Pads 1" Thick - Green (Pair)		\$117
<input type="checkbox"/>	135-3635-03	Size 3 and 4 PAL Side Pads 1" Thick - Blue (Pair)		\$117
<input type="checkbox"/>	135-3635-01	Size 3 and 4 PAL Side Pads 1" Thick - Green (Pair)		\$117
<input type="checkbox"/>	135-1654-03	Size 1 PAL Side Pads 2" Thick - Blue (Pair)		\$117
<input type="checkbox"/>	135-1654-01	Size 1 PAL Side Pads 2" Thick - Green (Pair)		\$117
<input type="checkbox"/>	135-2654-03	Size 2 PAL Side Pads 2" Thick - Blue (Pair)		\$117
<input type="checkbox"/>	135-2654-01	Size 2 PAL Side Pads 2" Thick - Green (Pair)		\$117
<input type="checkbox"/>	135-3654-03	Size 3 and 4 PAL Side Pads 2" Thick - Blue (Pair)		\$117
<input type="checkbox"/>	135-3654-01	Size 3 and 4 PAL Side Pads 2" Thick - Green (Pair)		\$117

Ramped Seat Cushion*Must match seat size.*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	135-1651-03	Size 1 PAL Ramped Seat Cushion - Blue		\$79
<input type="checkbox"/>	135-1651-01	Size 1 PAL Ramped Seat Cushion - Green		\$79
<input type="checkbox"/>	135-2651-03	Size 2 PAL Ramped Seat Cushion - Blue		\$79
<input type="checkbox"/>	135-2651-01	Size 2 PAL Ramped Seat Cushion - Green		\$79
<input type="checkbox"/>	135-3651-03	Size 3 PAL Ramped Seat Cushion - Blue		\$79
<input type="checkbox"/>	135-3651-01	Size 3 PAL Ramped Seat Cushion - Green		\$79
<input type="checkbox"/>	135-4651-03	Size 4 PAL Ramped Seat Cushion - Blue		\$79
<input type="checkbox"/>	135-4651-01	Size 4 PAL Ramped Seat Cushion - Green		\$79

Tray and Tray Pad*Must match seat size.*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	135-1620	Size 1 PAL Tray		\$213
<input type="checkbox"/>	135-2620	Size 2 PAL Tray		\$213
<input type="checkbox"/>	135-3620	Size 3 PAL Tray		\$213
<input type="checkbox"/>	135-4620	Size 4 PAL Tray		\$213
<input type="checkbox"/>	135-1656	Size 1 PAL Tray Pad		\$98
<input type="checkbox"/>	135-2656	Size 2 PAL Tray Pad		\$98
<input type="checkbox"/>	135-3656	Size 3 PAL Tray Pad		\$98
<input type="checkbox"/>	135-4656	Size 4 PAL Tray Pad		\$98

Pommel*Must match seat Size.*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	135-1602-03	Size 1 PAL Pommel - Blue		\$154
<input type="checkbox"/>	135-1602-01	Size 1 PAL Pommel - Green		\$154
<input type="checkbox"/>	135-3602-03	Size 2 PAL Pommel - Blue		\$154
<input type="checkbox"/>	135-3602-01	Size 2 PAL Pommel - Green		\$154

Footrest*Must match seat Size.*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	135-1603	Size 1 PAL Footrest		\$117
<input type="checkbox"/>	135-2603	Size 2 PAL Footrest		\$117
<input type="checkbox"/>	135-3603	Size 3 PAL Footrest		\$117
<input type="checkbox"/>	135-4603	Size 4 PAL Footrest		\$117

Rocker / Stabilizer*Not available with PAL ski with Casters*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	135-1604	Size 1 PAL Rocker/Stabilizer (Pair)		\$117
<input type="checkbox"/>	135-2604	Size 2 PAL Rocker/Stabilizer (Pair)		\$117
<input type="checkbox"/>	135-3604	Size 3 PAL Rocker/Stabilizer (Pair)		\$117
<input type="checkbox"/>	135-4604	Size 4 PAL Rocker/Stabilizer (Pair)		\$117

Caster*Not available with Rocker/Stabilizer*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	135-1610	Size 1 PAL Ski with Casters (Pair)		\$231
<input type="checkbox"/>	135-2610	Size 2 PAL Ski with Casters (Pair)		\$231
<input type="checkbox"/>	135-3610	Size 3 PAL Ski with Casters (Pair)		\$231
<input type="checkbox"/>	135-4610	Size 4 PAL Ski with Casters (Pair)		\$231

Hi Seat Frame*Available only with Size 1 PAL Classroom Seat*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	135-1630	Size 1 PAL Hi Seat Frame		\$290

Assembly

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	ASSY-LKY	Product Assembly Fee		\$120

Technical Specifications - Imperial

Specification	Size 1	Size 2	Size 3	Size 4
Age	1 - 3 years	2 - 5 years	4 - 8 years	7 - 12 years
Child Height	27.6" - 37.4"	35.4" - 41.3"	39.4" - 51.2"	47" - 63"
Max User Weight	66 lbs	77 lbs	88 lbs	110 lbs
Seat to Footrest	5.5" - 9.4"	6.5" - 11.2"	9.1" - 13.8"	11.2" - 16"
Seat Depth	7.1" - 9.4"	8.7" - 11.4"	9.8" - 13"	12" - 16"
Backrest Height	10.4"	12.4"	15"	17.7"
Backrest Angle	95° - 105°	95° - 105°	95° - 105°	95° - 105°
Seat Width	5.9" - 9.8"	6.7" - 10.6"	9.1" - 13"	11.1" - 15.2"
Seat Weight	9.6 lbs	13.2 lbs	17.9 lbs	22 lbs
Tray Size	20"W x 18"D	21.3"W x 19.7"D	23.6"W x 21.7"D	28.7"W x 23.6"D
Hi Seat Frame	14.6"W x 18.7"D x 23.2"H	N/A	N/A	N/A

Technical Specifications - Metric

Specification	Size 1	Size 2	Size 3	Size 4
Age	1 - 3 years	2 - 5 years	4 - 8 years	7 - 12 years
Child Height	700 - 950 mm	900 - 1050 mm	1000 - 1300 mm	1200 - 1600 mm
Max User Weight	30 kg	35 kgs	40 kgs	50 kgs
Seat to Footrest	40 - 240 mm	165 - 285 mm	230 - 350 mm	280 - 400 mm
Seat Depth	80 - 240 mm	220 - 290 mm	250 - 330 mm	300 - 400 mm
Backrest Height	265 mm	315mm	380 mm	450 mm
Backrest Angle	95° - 105°	95° - 105°	95° - 105°	95° - 105°
Seat Width	50 - 250 mm	170 - 270 mm	230 - 330 mm	280 - 380 mm
Seat Weight	4.35kg	6 kgs	8.15 kgs	10 kgs
Tray Size	508mm W x 457mm D	541mm W x 500mm D	599mm W x 551mm D	729mmW x 599mmD
Hi Seat Frame	370mmW x 475mmD x 589mmH	N/A	N/A	N/A



Sunrise Medical (US) LLC2842 Business Park Ave. • Fresno, CA 93727 • USA