

Magic Mobility Extreme X8

March 2022



Order Form

Account:	
Date:	
PO#:	
Buyer Name:	
Marked For:	

Ship To:	
Name:	
Address:	
City, State, Zip:	
Phone:	


Please submit order to: magicmobility@sunmed.com

Sunrise Medical (US) LLC
2842 Business Park
Fresno, CA 93727
1-800-333-4000



Client Measurements

a)	Top of head to seat	Inch
b)	Top of shoulder to seat	Inch
c)	Underarm to seat	Inch
d)	Elbow to seat	Inch
e)	Top of leg to seat	Inch
f)	Shoulder width	Inch
g)	Chest width	Inch
i)	Hip width	Inch
j)	Actual sitting depth	Inch
m)	Foot length	Inch
n)	Knee to foot	Inch
o)	Seat Width	Inch



Notes:

Seating Specifications

HCPC	N/A	Seat Type	MPS and Rehab Style
Typical WC Weight w/ Batteries	Starting at 320 lbs. Weights vary depending on options	Seat Size Range	Width 13"-23" Depth 13"-24"
		STF Height Range	18.5", 19.5", 20.5"
Power Seating	50 degrees of CG Tilt / Power Recline / 12" seat elevator		

Drive Base Specifications

Drive Wheel Position	4 wheel Drive	Drive Wheels	14" All Terrain Drive Wheels
User Weight Capacity	400 lbs. 340 lbs. with Elevate and/or Tilt	Tires	All Terrain Tires
Base Size	Width 28", Length 40.5", Turn Radius 56"	Charger	Dual Mode Charger
Drive Control	R-Net PM 120A controller	Battery Size	Group 24
Drive Motors	4-pole 6mph	Color Accents	6 shroud colors and 6 color accents
		Ground Clearance	4"

OPTION	HCPC	Description	Retail
Seating Options			
<input type="checkbox"/>	MM-RHBSB	REHAB Seat & Back (Back canes and Adjustable Flat Seat Pan)	N/C
<input type="checkbox"/>	MM-MPSSB	MPS Seat & Back (Captain Style Seating)	N/C
<input type="checkbox"/>	MM-MPSBRHBS	MPS Back & REHAB Seat (Captain Style Backrest and Adjustable Flat Seat Pan)	N/C

Base

OPTION	HCPC	Description	Retail
Extreme X8 Base End Item Number: E1MM			
<input type="checkbox"/>	MM-X8NP	Magic Mobility - Extreme X8 No Power Seating Options	\$16,617.00
<input type="checkbox"/>	MM-X8SP	Magic Mobility - Extreme X8 Single Power Seating Options	\$16,617.00
<input type="checkbox"/>	MM-X8MP	Magic Mobility - Extreme X8 Multi Power Seating Options	\$16,617.00
<input type="checkbox"/>	MM-X8NPHD	Magic Mobility - Extreme X8 No Power Seating Options - HD	\$16,617.00
<input type="checkbox"/>	MM-X8SPHD	Magic Mobility - Extreme X8 Single Power Seating Options - HD	\$16,617.00
<input type="checkbox"/>	MM-X8MPHD	Magic Mobility - Extreme X8 Multi Power Seating Options - HD	\$16,617.00
<input type="checkbox"/>	MM-X8NPC	Magic Mobility - Extreme X8 No Power Seating Options - MPS (Captain)	\$16,617.00
<input type="checkbox"/>	MM-X8SPC	Magic Mobility - Extreme X8 Single Power Seating Options - MPS (Captain)	\$16,617.00
<input type="checkbox"/>	MM-X8MPC	Magic Mobility - Extreme X8 Multi Power Seating Options - MPS (Captain)	\$16,617.00
<input type="checkbox"/>	MM-X8NPCHD	Magic Mobility - Extreme X8 No Power Seating Options - HD MPS (Captain)	\$16,617.00
<input type="checkbox"/>	MM-X8SPCHD	Magic Mobility - Extreme X8 Single Power Seating Options - HD MPS (Captain)	\$16,617.00
<input type="checkbox"/>	MM-X8MPCHD	Magic Mobility - Extreme X8 Multi Power Seating Options - HD MPS (Captain)	\$16,617.00

Shroud Color			
<input type="checkbox"/>	MM-BKSHRD	Midnight Black Shroud Color	N/C
<input type="checkbox"/>	MM-BLSHRD	Ocean Blue Shroud Color	N/C
<input type="checkbox"/>	MM-RDSHRD	Fire Engine Red Shroud Color	N/C
<input type="checkbox"/>	MM-GNSHRD	Clover Green Shroud Color	N/C
<input type="checkbox"/>	MM-TGSHRD	Tangerine Shroud Color	N/C
<input type="checkbox"/>	MM-CMSHRD	Camouflage Shroud Color	\$514.00

Batteries			
<input type="checkbox"/>	MM-BAT24G	E2363 Group 24 Gel Battery installed - Pair	\$888.00
<input type="checkbox"/>	MM-BATCHG	Battery Charger	N/C

Drive Wheels			
<input type="checkbox"/>	MM-ATDW	All Terrain Drive Wheels and Tires	N/C

Tie Downs			
<input type="checkbox"/>	MM-X8TRNST	Crash Tested X8 Tie Down Kit	\$669.00

Additional Base Function			
<input type="checkbox"/>	MM-02944	Directional Lock X8 Base - used to reverse in a straight line	\$1,696.00
<i>Note: Must be ordered with SFCK or MFCK</i>			

Optional Drive Tire Options			
<input type="checkbox"/>	MM-TS4	Tire Sealant for 4 Pneumatic Drive Tires	\$64.00

Base Accessories			
<input type="checkbox"/>	MM-01595	Jack for Drive Tire	\$171.00

Seating System

OPTION	HCPC	Description	Retail
Seating System			

No Power Seat Options			
<input type="checkbox"/>	MM-NP	No Power Seating System (400 lb. Weight Capacity)	\$3,500.00
<input type="checkbox"/>	MM-NPFO	No Power Seating System with Power Functionality (400 lb. Weight Capacity)	\$4,500.00
<i>Note: required for lights, power recline, power legrest, steering lock, or future power seat functions.</i>			

Power Tilt - Utilized without Power Recline			
<input type="checkbox"/>	MM-TLT	E1002 50° Power Tilt in Space - 340 lb. Weight Capacity	\$6,474.00

Power Seat Elevate - Only available when Power Tilt or Power Tilt & Recline is selected			
<input type="checkbox"/>	MM-ELVT	E2300 12° Adjustable Seat Elevator - 340 lb. Weight Capacity	\$3,750.00

Power Recline - Utilized without Power Tilt - Requires No Power Seating with Power Functionality			
<input type="checkbox"/>	MM-PRMPS	E1003 150° Power Recline for MPS Seating - 400 lb. Weight Capacity	\$4,606.00
<input type="checkbox"/>	MM-PRRHB	E1003 150° Power Recline w/ Shear Reduction for REHAB Seating - 400 lb. Weight Capacity	\$4,606.00

Power Tilt and Recline			
<input type="checkbox"/>	MM-TLTPRMPS	E1006 50° Tilt & 150° Power Recline - MPS Seating - 340 lb. Weight Capacity	\$10,224.00
<input type="checkbox"/>	MM-TLTPRRHB	E1007 50° Tilt & 150° PR w/ Shear Reduction - REHAB Seating - 340 lb. Weight Capacity	\$11,765.00

Power Seat Function			
<input type="checkbox"/>	MM-SFCK	E2310 Single Seat Function Control Kit - Thru Drive	\$1,498.00
<i>Note: MUST be selected if ordering 1 actuator functions</i>			
<input type="checkbox"/>	MM-MFCK	E2311 Multiple Seat Function Control Kit - Thru Drive	\$3,050.00
<i>Note: MUST be selected if ordering 2 or more actuator functions</i>			

Seat To Floor Height

X8 Seat to Floor Heights start at 18.5" +/- 0.25"

<input type="checkbox"/> MM-STFH	Lowest Seat To Floor Height	N/C
<input type="checkbox"/> MM-STFH1	Seat To Floor Height +1"	N/C
<input type="checkbox"/> MM-STFH2	Seat To Floor Height +2"	N/C

Power Seat Options

<input type="checkbox"/> MM-CRTL5JS	Ctrl+5 - 5 Button Box for controlling power seat functions - Includes USB port - Joystick Side	\$803.00
<input type="checkbox"/> MM-CRTL5OS	Ctrl+5 - 5 Button Box for controlling power seat functions - Includes USB port - Opposite Side	\$803.00

Seat Sizes (width x depth)

MPS (Captain) Seat Cushion (width x depth)

Note: MPS Backrest can be paired with the REHAB flat pan - If Flat pan is desired omit the MPS cushion and select the desired REHAB flat pan

<input type="checkbox"/> MM-01673	16" x 16" Vinyl Seat	N/C
<input type="checkbox"/> MM-01674	16" x 18" Vinyl Seat	N/C
<input type="checkbox"/> MM-01675	18" x 16" Vinyl Seat	N/C
<input type="checkbox"/> MM-01676	18" x 18" Vinyl Seat	N/C
<input type="checkbox"/> MM-01678	18" x 20" Vinyl Seat	N/C
<input type="checkbox"/> MM-01679	20" x 18" Vinyl Seat	N/C
<input type="checkbox"/> MM-01680	20" x 20" Vinyl Seat	N/C
<input type="checkbox"/> MM-01681	22" x 18" Vinyl Seat	\$391.00
<input type="checkbox"/> MM-01682	22" x 20" Vinyl Seat	\$391.00

REHAB Seat Sizes - Adjustable in width & depth for most sizes

REHAB Seat Pan

<input type="checkbox"/> MM-02792STD	Standard Adjustable Flat Seat Pan - 13" to 19" deep	N/C
<i>Note: Standard seat pan available 13"-23" wide</i>		
<input type="checkbox"/> MM-02792DP	Deep Adjustable Flat Seat Pan - 20" to 21" deep	\$391.00
<i>Note: Deep seat pan available 13"-23" wide</i>		

Deep Seat Pan Extension Kits

<input type="checkbox"/> MM-05690	23" Deep Seat Pan Extension Kit - Requires selection of MM-02792DP	\$690.00
<input type="checkbox"/> MM-05741	24" Deep Seat Pan Extension Kit - Requires selection of MM-02792DP	\$690.00

REHAB Seat Width

- 13" 14" 15" 16" 17" 18" 19" 20" 21" 22" 23"

REHAB Seat Depth

- 13" 14" 15" 16" 17" 18" 19" 20" 21" 22" 23" 24"

Seatback Options - select one

MPS Backrest Options - Comes standard with Manual Recline

Breathable Cloth MPS Backrest

<input type="checkbox"/> MM-01360	Cloth Backrest - 16"-20" Wide & 22" Tall	N/C
<input type="checkbox"/> MM-01520	Cloth Tall Backrest - 16"-20" Wide & 26" Tall	\$471.00
<input type="checkbox"/> MM-01361	Cloth Wide Backrest - 22"-23" Wide & 22" Tall	\$284.00
<input type="checkbox"/> MM-04598	Cloth Wide & Tall Backrest - 22"-23" Wide & 26" Tall	\$562.00

Wipeable Vinyl MPS Backrest

<input type="checkbox"/> MM-01360VI	Vinyl Backrest - 16"-20" Wide & 22" Tall	\$257.00
<input type="checkbox"/> MM-01520VI	Vinyl Tall Backrest - 16"-20" Wide & 26" Tall	\$674.00
<input type="checkbox"/> MM-01361VI	Vinyl Wide Backrest - 22"-23" Wide & 22" Tall	\$342.00
<input type="checkbox"/> MM-04598VI	Vinyl Tall & Wide Backrest - 22"-23" Wide & 26" Tall	\$797.00

MPS Back Accessories

<input type="checkbox"/> MM-04784	Push Handle - Needed for O2 holder, Backpack, Cane holder	\$134.00
<input type="checkbox"/> MM-02533	Accessory Mount - Needed for aftermarket headrest and/or attendant control	\$102.00

REHAB Seating Backrest Options

Backrest will have the width of the MPS seat or flat pan. Please choose your back height

Backrest will come with 2 point hardware with a 500 lb. capacity

<input type="checkbox"/> J3PA	E2613	JAY J3 PA - Posterior - 2.2" of lateral thoracic support - goes to 20" wide	\$592.00
Select Back Height	<input type="checkbox"/> 16.5" <input type="checkbox"/> 18" <input type="checkbox"/> 19.5" <input type="checkbox"/> 21" <input type="checkbox"/> 22.5" <input type="checkbox"/> 24"		
<input type="checkbox"/> J3PAW	E2614	JAY J3 PAW - Posterior - 2.2" of lateral thoracic support - goes to 22" wide	\$810.00
Select Back Height	<input type="checkbox"/> 18" <input type="checkbox"/> 21"		
<input type="checkbox"/> J3PL	E2615	JAY J3 PL - Posterior Lateral - 3" of lateral thoracic support - goes to 20" wide	\$683.00
Select Back Height	<input type="checkbox"/> 16.5" <input type="checkbox"/> 18" <input type="checkbox"/> 19.5" <input type="checkbox"/> 21" <input type="checkbox"/> 22.5" <input type="checkbox"/> 24"		
<input type="checkbox"/> J3PD	E2620	JAY J3 PD - Posterior Deep - 6" of lateral thoracic support - goes to 20" wide	\$811.00
Select Back Height	<input type="checkbox"/> 16.5" <input type="checkbox"/> 18" <input type="checkbox"/> 19.5" <input type="checkbox"/> 21" <input type="checkbox"/> 22.5" <input type="checkbox"/> 24"		

MPS Backrest Options For use with REHAB Back Canes

<input type="checkbox"/> MM-03996	Cloth MPS Removable Backrest (Works with REHAB back canes)	\$1,669.00
<input type="checkbox"/> MM-03996VI	Vinyl MPS Removable Backrest (Works with REHAB back canes)	\$1,669.00

REHAB Flip Forward Backrest Option			
<input type="checkbox"/>	MM-FLPFWD	REHAB Flip Forward Backrest Option	\$250.00
REHAB Seating Cushion Options			
<input type="checkbox"/>	JION	E2607 JAY Ion	\$405.00
<input type="checkbox"/>	JU	E2607 JAY Union	\$454.00
<input type="checkbox"/>	JFUSION	E2622 JAY Fusion	\$450.00
Headrest Options - MPS Style			
<input type="checkbox"/>	MM-02770	E0955 Cloth Headrest - 12" x 4.5" x 6.25" - MPS seat back only	\$348.00
<input type="checkbox"/>	MM-02770VI	E0955 Vinyl Headrest - 12" x 4.5" x 6.25" - MPS seat back only	\$348.00
Headrest Options - REHAB Style - Select both			
<input type="checkbox"/>	J3PLSH10	E0955 WHITMYER PLUSH Headrest Pad - 10" x 5" - requires J3HEDX & 2533	\$164.00
<input type="checkbox"/>	J3HEDX	E1028 WHITMYER AXYS Headrest w/ Removable Hardware - depth adj .5" behind backrest and 7" forward	\$224.00
Armrest Style Options			
<input type="checkbox"/>	MM-ARMR	Removable Armrests - Post style (Attaches to the Seat Pan - Height Adjustable 2")	N/C
<input type="checkbox"/>	MM-ARMFR	Flip-Up Removable Armrests - Post style (Attaches to the base-height adj 3" & angle adj 90° - 155°)	\$372.00
<input type="checkbox"/>	MM-ARMC	Cantilever Armrests - Reclining Flip Back (Attaches to the back - Height Adjustable 4" & Angle adj 90°)	\$748.00
<input type="checkbox"/>	MM-ARMC2	Cantilever 2 Point Quick Release Post (Attaches front of armrest to seat frame)	\$992.00
Armrest Pad Options			
<input type="checkbox"/>	MM-ARMPAD01	3" x 12" Desk Length Pads	N/C
<input type="checkbox"/>	MM-ARMPAD02	4" x 16" Tapered Pads	\$91.00
<input type="checkbox"/>	MM-ARMPAD03	3" x 12" Rounded Leatherette Pad	\$246.00
<input type="checkbox"/>	MM-ARMPAD04	3" x 14" Rounded Leatherette Pad	\$246.00
Lateral Support Options - REHAB Seat Only			
<input type="checkbox"/>	J3SAMNTR	E1028 Single Swing away Lateral Hardware - Right	\$269.00
<input type="checkbox"/>	J3SAMNTL	E1028 Single Swing away Lateral Hardware - Left	\$269.00
<input type="checkbox"/>	J3LAT6X4R	E0956 Single 6x4 Lateral Pad - Right	\$128.00
<input type="checkbox"/>	J3LAT6X4L	E0956 Single 6x4 Lateral Pad - Left	\$128.00
Positioning Options			
Positioning Hardware - Compatible with REHAB Seat pan only			
<input type="checkbox"/>	MM-FDHIP	E1028 Flip Down Hip Guide Hardware - fixed to the seat pan near the back	\$578.00
<input type="checkbox"/>	MM-RMHIP	E1028 Removable Hip Guide Hardware - fixed to the seat pan near the back	\$470.00
<input type="checkbox"/>	MM-RMKNE	E1028 Removable Knee Adductors Hardware - fixed to the seat pan near the front	\$470.00
<input type="checkbox"/>	MM-FDKNE	E1028 Flip Down Knee Adductors - fixed to the seat pan near the front	\$578.00
Positioning Pad			
<input type="checkbox"/>	MM-4X6PAD	E0956 4" x 6" Positioning Pad	\$354.00
<input type="checkbox"/>	MM-4X10PAD	E0956 4" x 10" Positioning Pad	\$354.00
Positioning Pad Accessories			
<input type="checkbox"/>	MM-4X6GEL	E0956 4" x 6" Gel Cover	\$257.00
Chest Harness Options			
<input type="checkbox"/>	MM-SH120M-B2	E0960 Bodypoint Monoflex Medium Chest Strap - pad 3.5" x 18.5" total length 49.5"	\$144.00
<input type="checkbox"/>	MM-SH120L-B2	E0960 Bodypoint Monoflex Large Chest Strap - pad 3.5" x 21.5" total length 52.5"	\$144.00
<input type="checkbox"/>	MM-50-CH110	E0960 Healthwares Padded Butterfly Positioning Chest Strap - 7" at Chest	\$166.00
<input type="checkbox"/>	MM-SH331L-B1	E0960 Bodypoint Stayflex Butterfly Harness - 3.5" at Chest	\$225.00
Pelvic Positioning Options select one			
Pelvic Positioning			
<input type="checkbox"/>	MM-33140	E0978 Therafin Push Button Seat Restraint - 2" x 64"	N/C
<input type="checkbox"/>	MM-33148	E0978 Therafin Extra Large Push Button Seat Restraint - 2" x 80"	\$54.00
<input type="checkbox"/>	MM-ABSB	E0978 Therafin Airline Buckle Seat Restraint - 2" x 80"	\$139.00
<input type="checkbox"/>	MM-EB205-L62	K0108 Bodypoint Evoflex Pelvic Stabilizer - sitting hip width 14" - 18" requires Hw312-KIT	\$102.00
<input type="checkbox"/>	MM-WC19SB6K	E0978 Crash Tested Seat Belt Kit - 600mm	\$893.00
<input type="checkbox"/>	MM-WC19SB7K	E0978 Crash Tested Seat Belt Kit - 700mm	\$893.00
Pelvic Positioning Options			
<input type="checkbox"/>	MM-HW312-KIT	Hardware for Bodypoint Pelvic Stabilizer	\$75.00
<input type="checkbox"/>	MM-EBE100	Extension for Bodypoint Evoflex Pelvic Stabilizer	\$45.00

Legrest Options				
Fixed Legrest				
Fixed Center Mount Legrest				
<input type="checkbox"/>	MM-040531	K0040	Center Mount with Flip-Up Footplate for X8 height adj - 1 piece Footplate	\$385.00
<input type="checkbox"/>	MM-040532	K0040	Center Mount with Flip-Up Footplate for X8 height adj - 2 piece Footplate	\$385.00
Swing Away Legrest				
<input type="checkbox"/>	MM-03528	K0040	Side Mount Swing Away Leg Rest - w/ 2 piece flip-up footplate 8.5" x 6", 13" height adj	\$672.00
Fixed Center Mount Legrest Calf Pad Options				
<input type="checkbox"/>	MM-05547SM	E0995	Fixed Center Mount Calf Pads - w/ breathable fabric - Small (13"-15")	\$487.00
<input type="checkbox"/>	MM-05547ME	E0995	Fixed Center Mount Calf Pads - w/ breathable fabric - Medium (16"-19")	\$487.00
<input type="checkbox"/>	MM-05547LG	E0995	Fixed Center Mount Calf Pads - w/ breathable fabric - Large (20"-23")	\$487.00
Fixed Legrest Positioning Strap Options				
<input type="checkbox"/>	MM-SP102M	K0039	Bodypoint Medium Calf Strap - side mount only length 36" & width 3.5"	\$102.00
<input type="checkbox"/>	MM-SP102L	K0039	Bodypoint Large Calf Strap - side mount only length 40" & width 3.5"	\$118.00
<input type="checkbox"/>	MM-04597	E0951	Heel & Side Support - Attaches to the footplate stem and footplate	\$556.00
<input type="checkbox"/>	MM-50-SR230	E0952	Foot Strap - 1.5" wide strap	\$75.00
Power Elevating Legrest				
Center Mount Power Elevating Legrest				
<input type="checkbox"/>	MM-05765	E1012	Center Mount Power Elevating Legrest - 1 pc footplate 11" x 11.5" - Includes Calf Pads	\$3,531.00
<input type="checkbox"/>	MM-06639	E1012	Center Mount Power Elevating Legrest - 2 pc footplate 5" x 10.5" - Includes Calf Pads	\$3,424.00
Swing Away Power Elevating Legrest				
<input type="checkbox"/>	MM-03526	E1010	Side Mount Power Elevating Legrest* - Swing away legrest w/ Flip Up 2" height adj footplate	\$3,852.00
*Must be ordered with upgraded electronics MM-SFCK or MM-MFCK depending on number of options selected				
Power Legrest Positioning Options				
<input type="checkbox"/>	MM-05324	K0042	2" Metal Heel Guards for 2 Piece Footplates	\$717.00
<i>Note: May not fit with all seat sizes - requires MM-06639 or MM-03526</i>				
<input type="checkbox"/>	MM-827A	E0952	Foot Straps for 2 Piece Footplates - 11"-14" (pair)	\$75.00
<input type="checkbox"/>	MM-HL	E0951	Heel Loops - Attaches to the footplate - length 11.5" - 13" & width 1"	\$107.00
Electronics				
OPTION	HCPC	Description	Retail	
Joystick				
<input type="checkbox"/>	MM-PG-JSMLEDL	LED R-net Joystick - No LCD display or mono jacks	N/C	
<input type="checkbox"/>	MM-PG-JSMD	Color LCD Joystick - with mono jack ports for external switches - CJSM1	\$337.00	
<input type="checkbox"/>	MM-PG-JSM2IRBT	Hi-Resolution Color LCD Joystick with Bluetooth and Infra-Red - CJSM2	\$503.00	
Expandable Power Module - Must Select both items				
<input type="checkbox"/>	MM-PM120	E2377	R-net 120 amp Expandable Controller	\$696.00
<input type="checkbox"/>	MM-PG-SL	E2313	R-net Harness for Expandable Controller	\$487.00
Joystick Handle Type - Comes standard with PG Cone-Shaped Joystick Handle				
<input type="checkbox"/>	MM-PC101B	E2323	Bodypoint U-Shaped Joystick Handle - 3"	\$161.00
<input type="checkbox"/>	MM-PC102B	E2323	Bodypoint U-Shaped Joystick Handle - 4"	\$161.00
<input type="checkbox"/>	MM-PC110B	E2323	Bodypoint Grooved Mushroom Joystick Handle	\$150.00
Joystick Mount				
<input type="checkbox"/>	MM-03069R		Fixed Joystick Mount - Depth Adj 8.25" - Right	N/C
<input type="checkbox"/>	MM-03069L		Fixed Joystick Mount - Depth Adj 8.25" - Left	N/C
<input type="checkbox"/>	MM-04607	E1028	Swing Away Joystick Mount - Right (Includes latch to lock in the forward position)	\$535.00
<input type="checkbox"/>	MM-04608	E1028	Swing Away Joystick Mount - Left (Includes latch to lock in the forward position)	\$535.00
<input type="checkbox"/>	MM-JMMIDR	E1028	Midline S/A Joystick Mount - Rotation 0° - 180° - Right	\$1,295.00
<input type="checkbox"/>	MM-JMMIDL	E1028	Midline S/A Joystick Mount - Rotation 0° - 180° - Left	\$1,295.00
Joystick Accessories				
<input type="checkbox"/>	MM-03330AM	E2331	Attendant Control Joystick for MPS Seating with mounting kit	\$1,792.00
<input type="checkbox"/>	MM-03329AM	E2331	Attendant Control Joystick for REHAB Seating with mounting kit	\$1,792.00
<input type="checkbox"/>	139920		Tash Buddy Button - Diameter 2.5" & Height .7" Outer edge force 3.4 oz. Center force 7 oz.	\$166.00
<input type="checkbox"/>	MM-ESBLK		Egg Switch - single function push button that can scroll through modes or profiles - requires CPS9	\$139.00
<input type="checkbox"/>	MM-CPS9-L-HW		Swing Away Egg Switch Hardware - Stealth	\$428.00
Omni 2				
<input type="checkbox"/>	MM-05436R		OMNI 2 Display - for alternative drive systems - Right	\$2,579.00
<input type="checkbox"/>	MM-05436L		OMNI 2 Display - for alternative drive systems - Left	\$2,579.00
Alternative Drive Controls				
<input type="checkbox"/>	MM-HDR-4M-S-C	E2327	Switch-It Head Array - three switch head array	\$5,671.00
<input type="checkbox"/>	MM-03335	E2373	Chin Control for MPS Seating System w/ mount	\$3,585.00
<input type="checkbox"/>	MM-03540	E2373	Chin Control for REHAB Seating System w/ mount	\$3,585.00

Accessories

OPTION	HCPC	Description	Retail
Accessories			
<input type="checkbox"/> MM-01605		Waterproof Wheelchair Cover	\$139.00
<input type="checkbox"/> MM-01606		Waterproof Joystick Cover	\$43.00
<input type="checkbox"/> MM-01603		Medical Necessity Bag - Velcro's to armrest - 10.5" x 6" x 1.25"	\$32.00
<input type="checkbox"/> MM-001B		Deluxe Backpack - 13.5" x 15" x 4"	\$91.00
<i>Note: Requires Push handle (MM-04784) on MPS seat to attach backrest to chair</i>			
<input type="checkbox"/> MM-PP0052		Cup Holder - Attaches to the armrest opposite joystick - self leveling	\$150.00
<input type="checkbox"/> MM-06807	E2208	Articulating Oxygen Holder for MPS Seating System	\$990.00
<input type="checkbox"/> MM-06808	E2208	Articulating Oxygen Holder for Rehab Seating System	\$990.00
<input type="checkbox"/> MM-05275	E2207	Cane Holder for MPS Seating System	\$637.00
<input type="checkbox"/> MM-05276	E2207	Cane Holder for REHAB Seating System	\$637.00
<input type="checkbox"/> MM-03198		Luggage Rack for X8 - Attaches to base behind the chair 10.5" x 13"	\$316.00
<input type="checkbox"/> MM-02857		Lighting Kit for X8 - 2 headlights, 2 red marker lights w/ turn signals-NA with Flip Down Knee Supports	\$1,124.00
<input type="checkbox"/> MM-05420		Lift brackets for X8	\$920.00
<input type="checkbox"/> MM-04705		Charge Mate - 12 Volt/Car outlet and USB port/cell phone charger, with a max draw of 6 amps	\$508.00
<input type="checkbox"/> MM-PP0133		Cell Phone Holder - Adj. to fit most cell phones	\$278.00
<input type="checkbox"/> MM-RAM117BU		Fishing Rod Holder	\$278.00
<input type="checkbox"/> MM-DDT		Hand-held Programmer	\$1,530.00
Be Adaptive Products			
<input type="checkbox"/> MM-LM100		Be Adaptive Limited Mobility Shooting Mount - gun or crossbow	\$1,418.00
<input type="checkbox"/> MM-TM100		Be Adaptive Trigger Mechanism	\$786.00
<input type="checkbox"/> MM-FM200		Be Adaptive Limited Mobility Fishing Mount	\$1,418.00