

LECKEY



Mygo Stander

Infinite positioning to help kids stand
with confidence

Supporting Children

Stand out from the crowd in the true 3-in-1 stander

Mygo Stander provides the best standing therapy solution on the market – it is the only 3-in-1 stander which delivers truly excellent positioning whether a child is in prone, supine or upright. Mygo Stander doesn't compromise on quality, ensuring children get the support and comfort they need to enjoy standing and the benefits this brings.

OUTSTANDING ADJUSTABILITY TO ACCOMMODATE HIP AND KNEE CONTRACTURES

- Mygo Stander's footplate and knee supports uniquely move together, taking the strain off the hamstrings and other contracted muscles, so that children can stand in a fully supported, natural and comfortable position.
- This linked movement opens up standing therapy to children who were previously unsuitable because of their contractures.



DESIGNED TO IMPROVE UPPER LIMB FUNCTION, INCREASING ACTIVITY AND INDEPENDENCE

- The ergonomically designed chest pad replicates the shape of the rib cage allowing the shoulder girdle to freely protract and retract in prone standing. The supine shoulder support also stabilizes the shoulder girdle.
- Combined with Mygo Stander's lower tray height and a more upright head and chest position, the child's upper limb function is significantly improved. The child has space to move, supporting freer arm movement and an increased ability to engage in everyday activities such as feeding themselves, playing, self-care tasks and communicating with others.



MODULAR STANDER WHICH CAN BE TAILORED TO MEET SPECIFIC NEEDS

- Mygo Stander can be easily configured to evolve with the child as they grow and their postural support needs and therapy goals change.
- An extensive range of options and accessories enable an individualized fit, ensuring the best possible posture is achieved e.g., pelvic de-rotation belts, chest/hip harnesses, knee straps and a variety of laterals.



Technical Information

	Size 1	Size 2
Age (Approx.)	4-10	8-14
Max User Weight	50kgs / 110lbs	60kgs / 132lbs
User Height	Min 1080mm / 42.5" Max 1450mm / 57.1"	Min 1350mm / 53.1" Max 1700mm / 66.9"
Chest Support Height (footplate to armpit)	Min 760mm / 29.9" Max 1080mm / 42.5"	Min 980mm / 38.6" Max 1300mm / 51.2"
Chest Width (distance between laterals)	Min 160mm / 6.3" Max 280mm / 11"	Min 200mm / 7.9" Max 320mm / 12.6"
Hip Support Height (footplate to mid buttocks)	Min 460mm / 18.1" Max 730mm / 28.7"	Min 620mm / 24.4" Max 875mm / 34.4"
Hip Width (distance between hip guides)	Min 220mm / 8.7" Max 330mm / 13"	Min 280mm / 11" Max 390mm / 15.4"
Footplate Angle Plantarflexion \ Dorsiflexion	10° \ 10°	10° \ 10°
Maximum Knee Contracture Accommodated	25°	25°
Tray Height Range (from top of hip pad)	0mm – 150mm / 0"-6"	0mm – 150mm / 0"-6"
Tray Angle Adjustment Prone / Supine	45° \ 45°	45° \ 45°
Product Angle Range Prone / Supine	0°-80° / 0°-80°	0°-80° / 0°-80°
Stander Footprint Length Width	978mm / 38.5" 630mm / 24.8"	1130mm / 44.5" 650mm / 25.6"

To arrange a product demonstration/assessment or for more information, please contact us.