



Quickie GP Swing Away

August 2026

EIGPSA

Mark For:

Date: _____
Dealer Acct #: _____
Dealer: _____
Dealer Contact: _____
Dealer Address: _____
Dealer City: _____ ST: _____ ZIP: _____
Dealer Phone: () _____ Fax: () _____
Confirmation Email: _____
Confirm Via: ☐ Fax ☐ Email

Submitting for:

☐ Quote

☐ Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____
Attention: _____
Address: _____
Address: _____
Ship To City: _____ ST: _____ ZIP: _____
Ship To Phone: () _____ Fax: () _____

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
---	----------	-------	------------	-------	---	----------	-------	------------	-------

1. FRAME

For sizes 14" or less please use HCPCS E1235

☐	025WM26	K0005	Quickie GP Swing-Away	\$2,467.00
--------------------------	---------	-------	-----------------------	------------

2. SIZE / SEATING

Seat Width

Seat Width 1/2" Narrower

☐	025W3	12"	N/A w/ adj. upholstery	NC
☐	025W4	13"	N/A w/ adj. upholstery	NC
☐	025W5	14"		NC
☐	025W6	15"		NC
☐	025W7	16"		NC
☐	025W8	17"		NC
☐	025W9	18"		NC
☐	025W41	19"		\$147.00
☐	025W11	E2201	20"	\$583.00

Seat Depth

☐	025SD5		12"		NC
☐	025SD6		13"		NC
☐	025SD7		14"		NC
☐	025SD8		15"		NC
☐	025SD9		16"		NC
☐	025SD10		17"		NC
☐	025SD11		18"		NC
☐	025SD74		19"		\$147.00
☐	025SD13	E2203	20"		\$641.00
☐	025SD14	E2203	21"		\$641.00
☐	025SD33	E2204	22"		\$1,051.00

Cushion

☐	025CU4	Omit Seat Sling	-\$15.00
--------------------------	--------	-----------------	----------

JAY Bundled Cushions - Optional

To order any JAY cushions on-chair, please see JAY Price List. Cushion name or part number must be filled out with desired width & depth.

☐	025JC20	JAY Bundled Cushion Cushion Name or Part #: _____ Width: _____ Depth: _____	Price List
--------------------------	---------	---	------------

JAY Personalized Seating- Optional

To order any JAY personalized seating on-chair, please see desired order form

☐	025QC02	Jay SureFit	See OF
☐	025QC03	Jay Fit	See OF

3. CASTERS

Caster Wheels

☐	025CW1		3" Micro	NC
☐	025CW3		4" Low Profile Polyurethane	NC
☐	025CW31		4" Micro Lighted	\$82.00
☐	025CW5		5" Low Profile Polyurethane	NC
☐	025CW7	E2219	5" x 1.5" Semi-Pneumatic <i>only avail. for use in bottom fork holes</i>	\$100.00
☐	025CW6	E2219	5" Soft Roll Semi-Pneumatic	\$105.00
☐	025CW33	E2214 E2215	6" Pneumatic	\$88.00
☐	025CW34		6" Polyurethane	NC
☐	025CW35	E2219	6" x 1.5" Semi-Pneumatic <i>only avail. for use in bottom fork holes</i>	\$100.00
☐	025CW21	E2214 E2215	8" Pneumatic	\$100.00
☐	025CW32		8" Polyurethane	NC

Caster Wheels - Continued

☐	025CW36	E2214 E2215	8" x 2" Pneumatic	\$169.00
☐	025CW37	E2214 K0108	8" x 2" Pneumatic w/Airless Insert	\$263.00

Caster Fork

☐	025CS5	3/4" Longer Fork Stem Bolt	NC
☐	025CS9	1-1/2" Longer Fork Stem Bolt	NC

✓	Option #	HPCPS	Definition	Price	✓	Option #	HPCPS	Definition	Price
4. REAR WHEELS									

Xtender Power Assist				
<i>RSH 15"-19" with 22" wheels, RSH 16"-20" with 24" wheels. Spare batteries must match set up. Anti-Tips included. Dual Drive mode available w/B4ME</i>				
☐	025XT6	E0986	22" Xtender Power Assist JWX-2 NICKEL	\$10,404.00
☐	025XT7	E0986	22" Xtender Power Assist JWX-2 LITHIUM	\$11,390.00
☐	025XT4	E0986	24" Xtender Power Assist JWX-2 NICKEL	\$10,404.00
☐	025XT5	E0986	24" Xtender Power Assist JWX-2 LITHIUM	\$11,390.00

Please select one of the following handrim options

☐	025XH5		Plastic Coated	NC
☐	025XH6		Stainless Steel	NC

Xtender - Other Options				
☐	025XA10		Spare NiMH Power Pack - NICKEL	\$2,128.00
☐	025XA11		Spare Li-ion Power Pack - LITHIUM	\$2,594.00

Rear Wheel Size				
☐	025RW8		22" Mag	NC
☐	025RW9		22" Spoke	\$234.00
☐	025RW1		24" Mag	NC
☐	025RW2		24" Spoke	NC
☐	025RW11		26" Spoke <i>N/A with Xtender</i>	\$234.00

Rear Wheel Tire				
☐	025RT1		Pneumatic	NC
☐	025RT2	E2211 E2213	Pneumatic Airless Insert N/A with 26" wheels & Xtender	\$100.00
☐	025RT3		Full Poly N/A with 26" wheels	NC
☐	025RT6	E2211 E2212	High-Pressure Clincher <i>Only avail. On 24" or 26" spoke rims</i>	\$110.00

Hand Rim				
<i>N/A with Xtender Power Assist. Plastic coated not available on all sizes</i>				
☐	025HR1		Aluminum Anodized	NC
☐	025HR4		Plastic Coated	\$129.00
☐	025HR7		Omit Hand Rim	NC

Projection Hand Rim				
<i>N/A with plastic coated hand rims.</i>				
☐	025HR8		8 Oblique	\$379.00

Axle				
☐	025AX4		Quick Release Steel	NC
☐	025AX11		Quad Quick Release	\$147.00

Axle Plate				
☐	025AL2		Standard Adjustable	NC
☐	025AL6	E0959	Amputee Adjustable <i>N/A with padded swing-away armrests or stroller handles</i>	\$205.00
☐	025AL13		Double Slotted Adjustable	\$70.00

5. WHEEL LOCKS

Wheel Lock				
☐	025WL1		High Mount Push	NC
☐	025WL2		High Mount Pull	NC
☐	025WL5		High Mount Scissor <i>N/A w/ grade aids or extension handles</i>	\$147.00
☐	025WL12		Do Not Mount	NC
☐	025WL10		Omit	NC

Extension Handle				
☐	025EH1	E0961	6" Pair	\$88.00

Grade Aid				
☐	025GA1	E0974	Grade Aid <i>N/A w/ poly, high-pressure clincher tires</i>	\$175.00

6. HANGER AND FOOTREST

Hangers				
☐	025H145		60° Swing In/Out Hanger	NC
☐	025H139		70° Swing In/Out Hanger	NC
☐	025H158		90° Swing In/Out Front Mount <i>Only available Ang adj adult & kids footplates</i>	NC
☐	025H159	K0053	Articulating Hanger <i>N/A w/ Platform Flip Up and Alum Angle Adj Front Mount footplates.</i>	\$431.00
☐	025H160	E0990	Elevating Hanger - Ext Mount <i>N/A w/ Platform Flip Up and Alum Angle Adj Front Mount footplates.</i>	\$367.00
☐	025H161	E0990	Elevating Hanger - Front Mount <i>Only available w/ front mount alum angle adj footplates and 14"-22" width</i>	\$367.00

Footplates				
<i>Heel loop or leg strap not included, See Footplate Options Section</i>				
☐	025F1		Composite Footplate 14"-20" <i>N/A w/ 90° hangers and ELRs</i>	NC
☐	025F3	K0040	Adult Angle Adjustable	\$221.00
☐	025F4	K0040	Kids Angle Adjustable <i>Only available on 12"-16" width</i>	\$221.00
☐	025F48	K0040	Platform Flip Up <i>Available w/ 60° and 70° hangers</i>	\$221.00
☐	025F5		Aluminum <i>Only available w/ ELRs</i>	NC
☐	025F7		Omit Footplate & Ext Tube	-\$15.00

6. HANGER AND FOOTREST (CONT)

Seat Sling to Top of Footrest

Hanger	Extension	Seat Sling to Top of Footrest in inches		*Minimum Allowable FSH
		Min	Max	
60° Swing In/Out	Short	11.5	15	15
	Medium	11.5	17	15
	Long	14	18.5	17
70° Swing In/Out	Short	13	18	15
	Medium	13	20	15.5
	Long	16	20.5	17.5
60° Swing In/Out w/Platform Flip Up	Short	11.5	14.5	14
	Long	14	18.5	16
70° Swing In/Out w/Platform Flip Up	Short	13	15.5	15
	Long	15.5	20	17.5
90° Front Mount	No Ext	5.5	11.5	14
ALR	Short	16.5	20	15
	Medium	17.5	20.5	16
	Long	19	20.5	17
ELR Ext. Mount	Short	15.5	19.5	15
	Medium	15.5	20.5	16
	Long	17	20.5	17
ELR Front Mount	No Ext	5.5	9.5	13

Footrest Extension Tubes

See Chart on the bottom of the page for lengths

☐	025E1	Short	N/A w/ ELR or ALR	NC
☐	025E2	Medium		NC
☐	025E5	Long		NC
☐	025E16	High Mount	Only avail. with adult and kids angle adjustable footplates	\$88.00
☐	025E12	Short ELR		NC
☐	025E15	Long ELR		NC

Footrest Options

☐	025HL1	E0951	Heel Loop	\$47.00
☐	025HL5	K0038	Leg Strap	\$47.00
☐	025HL4		Omit Heel Loop	NC
☐	025IG2		Impact Guard Neoprene (Only avail. On 60° or 70° hangers)	\$105.00

7. BACKREST

Back Type

☐	025BT1	Folding Backrest (Adjustable)		
☐	025BA1	11" - 14" Back Height		\$75.00
☐	025BA2	12" - 16" Back Height		\$75.00
☐	025BA3	14" - 18" Back Height		\$75.00
☐	025BA4	16" - 20" Back Height		\$75.00
☐	025BT3	Non-Folding Backrest (Fixed)		
		N/A with push handles or adj. upholstery		
☐	025BB1	8" Back Height		NC
☐	025BB2	9" Back Height		NC
☐	025BB3	10" Back Height		NC
☐	025BB4	11" Back Height		NC
☐	025BB5	12" Back Height		NC
☐	025BB6	13" Back Height		NC
☐	025BB7	14" Back Height		NC
☐	025BB8	15" Back Height		NC
☐	025BB9	16" Back Height		NC

Back Upholstery

☐	025BU1		Standard	NC
☐	025BU8	E2611	Adjustable	\$205.00
☐	025BU16		Reinforced (N/A w/ adj. upholstery)	\$36.00
☐	025BU14		Omit Back Upholstery	-\$15.00

Backrest Options

☐	025BT33		Folding w/ Integral Push Handles	NC
☐	025A85		Bolt-On Push Handle (Folding only) N/A with Xtender Power Assist	\$134.00
☐	025SR2		Regular Stroller Handle N/A with amputee axle plate	\$158.00
☐	025SR3		Tall Stroller Handle N/A with amputee axle plate	\$158.00

JAY Back (Optional)

☐	025JB01		J3 Back (See J3 Order Form)	Price List
--------------------------	---------	--	-----------------------------	------------

8. ARMREST / SIDEGUARDS

Armrest - Padded Swing Away

All Padded Swing-Away N/A with amputee axle plate or Xtender Power Assist

☐	025AA19	E0973	Padded S-A Medium	\$263.00
☐	025AA16	E0973	Padded S-A Adult	\$263.00

Single Post

☐	025AA9	E0973	Single Post Ht. Adjust. Low Must be removed to fold back down	\$333.00
☐	025AM1	E0973	Desk Length - Classic (10")	NC
☐	025AM2	E0973	Full Length - Classic (14")	NC
☐	025AA8	E0973	Single Post Ht. Adjust. Std Must be removed to fold back down	\$333.00
☐	025AM1	E0973	Desk Length - Classic (10")	NC
☐	025AM2	E0973	Full Length - Classic (14")	NC

Side Guard

N/A w/ height adjustable arms

☐	025SG3	K0108	Plastic Regular	\$169.00
☐	025SG4	K0108	Plastic Kids	\$169.00

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
9. KEY OPTIONS AND ACCESSORIES									
Rear Anti-Tip					Spoke Guards				
☐	025AT1	E0971	Rear Anti-Tip	\$175.00	<i>N/A w/ Xtender Power Assist</i>				
Backpack					☐	025A23	K0065	Black Spoke Guards	\$134.00
☐	025A1		Quickie Backpack Black	\$88.00	☐	025A24	K0065	Clear Spoke Guards	\$134.00
☐	025A2		Quickie Kids Backpack Black	\$70.00	Crutch Holder				
Seat Pouch					☐	025A52	E2207	Right Side Mount	\$227.00
☐	025A8		Quickie Seat Pouch Black	\$58.00	☐	025A53	E2207	Left Side Mount	\$227.00
Positioning Belts					Extra Leg Strap				
☐	025A142	E0978	Belt 1.5" Auto Buckle	\$82.00	☐	025A55	K0038	Extra Leg Strap Black	\$36.00
☐	025A12	E0978	Belt 2" Aircraft Buckle	\$129.00	Miscellaneous				
☐	025A13	E0978	Belt 2" Aircraft Buckle Padded	\$139.00	☐	025A83		Tool Kit	\$63.00
					☐	025A999		Quickie Water Bottle & Holder	\$63.00
					Tray Table				
					☐	025A31	E0950	Tray Table	\$298.00

10. COLOR

NEW!

Standard Frame Colors				
☐	025C214		Black (Gloss)	NC
☐	025C215		Black (Matte)	NC
☐	025C211		Celestial Sapphire	NC
☐	025C212		Desert Sage	NC
☐	025C213		Desert Tan	NC
☐	025C216		Metallic Shadow	NC
☐	025C217		Midnight Blue	NC
☐	025C218		Mystic Purple	NC
☐	025C219		Petal Pink	NC
☐	025C220		Pink Power	NC
☐	025C221		Red Pepper	NC
☐	025C222		Royal Blue	NC
☐	025C223		Salute Green	NC
☐	025C224		Scarlet Noir	NC
☐	025C225		Silver Ore	NC
☐	025C226		Solar Orange	NC
☐	025C227		Starlight White	NC
☐	025C228		Storm Gray	NC
☐	025C229		Turbo Green	NC
☐	025C230		Velocity Yellow	NC

Personalized Frame Color				
☐	025C112		Desert Camouflage	\$460.00
☐	025C43		Mossy Oak Camouflage	\$460.00
☐	025C42		Stars & Stripes	\$460.00
☐	025C108		Zebra (Choose above base color) Std with White unless above base paint is chosen N/A w/ Black, Matte Black)	\$460.00

Built For Me

To inquire about any unique personalized Quickie products, contact:
Phone: (800) 333 - 4000, Option #3, Fax: (800) 333 9011

built4me

Commonly Requested Options Include:

☒ Heavy Duty (350 lb)

☒ One Arm Drive

☒ Unilateral Wheel Locks

☒ Seat Pan

☒ Lifestyle Options (Paint, Anodizing, etc...)



Sunrise Medical (US) LLC
12002 Volunteer Blvd. Mount Juliet, TN 37122

Customer Service: 800-333-4000 Fax: 800-300-7502 Email: orders@sunmed.com www.sunrisemedical.com

© 8/2026 MK-100023 REV O