



JAY GO Cushion Price List



Effective January 2020

E2603 - Skin Protection Wheelchair Cushion

E2604 - Skin Protection Wheelchair Cushion 22" and over in width

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.


How to Order:

Part #	JAY GO Cushions	Price	HCPC
<input checked="" type="checkbox"/> JGO1818	JAY GO Cushion 18" x 18"	\$200	E2603

Each JAY GO part number is simply based upon the size you are ordering. You may also reference the size matrix on page 1 of this price list.



Part # Example



JAY GO Cushion Cover

The JAY GO cover is designed using dartex™ fabric which is known to be both stretchable and water resistant. This cover's sizing and sewing pattern is designed to reduce shear and peak pressure buildup. The material itself provides stretch which is an uncommon property in most incontinent products.



Figure A



Please visit www.dartexcoatings.com for additional information on the use of dartex coatings

JAY GO Cushion Foam

The JAY GO foam base is designed to contain an accommodating seating surface in a thicker, softer dual layered foam base. The optimal combination of two different foam densities provides skin protection.



Cross section view

		JAY GO Size Matrix											
		WIDTH											
Depth		14	15	16	17	18	19	20	21	22	23	24	
	14												
	15												
	16												
	17												
	18												
	19												
	20												
	21												
	22												
	23												
24													

NO Modified sizes are available

Size is available Size is not available

JAY GO Cushion Sealed Zipper

The JAY GO zipper and cover provide protection to the foam base against odor because the location is low and away from fluid. The water resistant material, combined with the Aquaguard™ zipper and anti wicking thread provides a thorough protection system for a cushion base.



Figure B

Optimized Aquaguard Zipper location

The JAY GO Cushion is an excellent compliment to the **JAY GO Back**.

Please see the **JAY GO Back** order form (part # 119284) for additional information.



Each JAY GO cushion includes an outer cover, dual layered foam base, plastic bag with label and owner's manual

Ordering Information PHONE (800) 333-4000, Fax (800) 300-7502, WEB: www.sunrisemedical.com

Specifications Subject to change without notice. Copyright © 2020 Sunrise Medical (8.11) Jay GO Cushion 119288 Rev B Page 1 of 2

JAY GO Cushion Price List

Effective January 2020



E2603 - Skin Protection Wheelchair Cushion

E2604 - Skin Protection Wheelchair Cushion 22" and over in width

E2603 and E2604 is based upon the code verification by PDAC or our interpretation of the codes. Proper use of HCPC codes does not insure coverage or payment. For coverage information, please verify the appropriate payor's coverage

Part #	Description	Price	HCPC CODE	Part #	Description	Price	HCPC CODE
<input type="checkbox"/> JGO1414	JAY GO Cushion 14" x 14"	\$200	E2603	<input type="checkbox"/> JGO2121	JAY GO Cushion 21" x 21"	\$200	E2603
<input type="checkbox"/> JGO1415	JAY GO Cushion 14" x 15"	\$200	E2603	<input type="checkbox"/> JGO2122	JAY GO Cushion 21" x 22"	\$200	E2603
<input type="checkbox"/> JGO1416	JAY GO Cushion 14" x 16"	\$200	E2603	<input type="checkbox"/> JGO2123	JAY GO Cushion 21" x 23"	\$200	E2603
<input type="checkbox"/> JGO1417	JAY GO Cushion 14" x 17"	\$200	E2603	<input type="checkbox"/> JGO2124	JAY GO Cushion 21" x 24"	\$200	E2603
<input type="checkbox"/> JGO1418	JAY GO Cushion 14" x 18"	\$200	E2603	<input type="checkbox"/> JGO2217	JAY GO Cushion 22" x 17"	\$300	E2604
<input type="checkbox"/> JGO1419	JAY GO Cushion 14" x 19"	\$200	E2603	<input type="checkbox"/> JGO2218	JAY GO Cushion 22" x 18"	\$300	E2604
<input type="checkbox"/> JGO1420	JAY GO Cushion 14" x 20"	\$200	E2603	<input type="checkbox"/> JGO2219	JAY GO Cushion 22" x 19"	\$300	E2604
<input type="checkbox"/> JGO1514	JAY GO Cushion 15" x 14"	\$200	E2603	<input type="checkbox"/> JGO2220	JAY GO Cushion 22" x 20"	\$300	E2604
<input type="checkbox"/> JGO1515	JAY GO Cushion 15" x 15"	\$200	E2603	<input type="checkbox"/> JGO2221	JAY GO Cushion 22" x 21"	\$300	E2604
<input type="checkbox"/> JGO1516	JAY GO Cushion 15" x 16"	\$200	E2603	<input type="checkbox"/> JGO2222	JAY GO Cushion 22" x 22"	\$300	E2604
<input type="checkbox"/> JGO1517	JAY GO Cushion 15" x 17"	\$200	E2603	<input type="checkbox"/> JGO2223	JAY GO Cushion 22" x 23"	\$300	E2604
<input type="checkbox"/> JGO1518	JAY GO Cushion 15" x 18"	\$200	E2603	<input type="checkbox"/> JGO2224	JAY GO Cushion 22" x 24"	\$300	E2604
<input type="checkbox"/> JGO1519	JAY GO Cushion 15" x 19"	\$200	E2603	<input type="checkbox"/> JGO2319	JAY GO Cushion 23" x 19"	\$300	E2604
<input type="checkbox"/> JGO1520	JAY GO Cushion 15" x 20"	\$200	E2603	<input type="checkbox"/> JGO2320	JAY GO Cushion 23" x 20"	\$300	E2604
<input type="checkbox"/> JGO1521	JAY GO Cushion 15" x 21"	\$200	E2603	<input type="checkbox"/> JGO2321	JAY GO Cushion 23" x 21"	\$300	E2604
<input type="checkbox"/> JGO1614	JAY GO Cushion 16" x 14"	\$200	E2603	<input type="checkbox"/> JGO2322	JAY GO Cushion 23" x 22"	\$300	E2604
<input type="checkbox"/> JGO1615	JAY GO Cushion 16" x 15"	\$200	E2603	<input type="checkbox"/> JGO2323	JAY GO Cushion 23" x 23"	\$300	E2604
<input type="checkbox"/> JGO1616	JAY GO Cushion 16" x 16"	\$200	E2603	<input type="checkbox"/> JGO2324	JAY GO Cushion 23" x 24"	\$300	E2604
<input type="checkbox"/> JGO1617	JAY GO Cushion 16" x 17"	\$200	E2603	<input type="checkbox"/> JGO2419	JAY GO Cushion 24" x 19"	\$300	E2604
<input type="checkbox"/> JGO1618	JAY GO Cushion 16" x 18"	\$200	E2603	<input type="checkbox"/> JGO2420	JAY GO Cushion 24" x 20"	\$300	E2604
<input type="checkbox"/> JGO1619	JAY GO Cushion 16" x 19"	\$200	E2603	<input type="checkbox"/> JGO2421	JAY GO Cushion 24" x 21"	\$300	E2604
<input type="checkbox"/> JGO1620	JAY GO Cushion 16" x 20"	\$200	E2603	<input type="checkbox"/> JGO2422	JAY GO Cushion 24" x 22"	\$300	E2604
<input type="checkbox"/> JGO1621	JAY GO Cushion 16" x 21"	\$200	E2603	<input type="checkbox"/> JGO2423	JAY GO Cushion 24" x 23"	\$300	E2604
<input type="checkbox"/> JGO1622	JAY GO Cushion 16" x 22"	\$200	E2603	<input type="checkbox"/> JGO2424	JAY GO Cushion 24" x 24"	\$300	E2604
<input type="checkbox"/> JGO1714	JAY GO Cushion 17" x 14"	\$200	E2603				
<input type="checkbox"/> JGO1715	JAY GO Cushion 17" x 15"	\$200	E2603				
<input type="checkbox"/> JGO1716	JAY GO Cushion 17" x 16"	\$200	E2603				
<input type="checkbox"/> JGO1717	JAY GO Cushion 17" x 17"	\$200	E2603				
<input type="checkbox"/> JGO1718	JAY GO Cushion 17" x 18"	\$200	E2603				
<input type="checkbox"/> JGO1719	JAY GO Cushion 17" x 19"	\$200	E2603				
<input type="checkbox"/> JGO1720	JAY GO Cushion 17" x 20"	\$200	E2603				
<input type="checkbox"/> JGO1721	JAY GO Cushion 17" x 21"	\$200	E2603				
<input type="checkbox"/> JGO1722	JAY GO Cushion 17" x 22"	\$200	E2603				
<input type="checkbox"/> JGO1723	JAY GO Cushion 17" x 23"	\$200	E2603				
<input type="checkbox"/> JGO1814	JAY GO Cushion 18" x 14"	\$200	E2603				
<input type="checkbox"/> JGO1815	JAY GO Cushion 18" x 15"	\$200	E2603				
<input type="checkbox"/> JGO1816	JAY GO Cushion 18" x 16"	\$200	E2603				
<input type="checkbox"/> JGO1817	JAY GO Cushion 18" x 17"	\$200	E2603				
<input type="checkbox"/> JGO1818	JAY GO Cushion 18" x 18"	\$200	E2603				
<input type="checkbox"/> JGO1819	JAY GO Cushion 18" x 19"	\$200	E2603				
<input type="checkbox"/> JGO1820	JAY GO Cushion 18" x 20"	\$200	E2603				
<input type="checkbox"/> JGO1821	JAY GO Cushion 18" x 21"	\$200	E2603				
<input type="checkbox"/> JGO1822	JAY GO Cushion 18" x 22"	\$200	E2603				
<input type="checkbox"/> JGO1823	JAY GO Cushion 18" x 23"	\$200	E2603				
<input type="checkbox"/> JGO1824	JAY GO Cushion 18" x 24"	\$200	E2603				
<input type="checkbox"/> JGO1915	JAY GO Cushion 19" x 15"	\$200	E2603				
<input type="checkbox"/> JGO1916	JAY GO Cushion 19" x 16"	\$200	E2603				
<input type="checkbox"/> JGO1917	JAY GO Cushion 19" x 17"	\$200	E2603				
<input type="checkbox"/> JGO1918	JAY GO Cushion 19" x 18"	\$200	E2603				
<input type="checkbox"/> JGO1919	JAY GO Cushion 19" x 19"	\$200	E2603				
<input type="checkbox"/> JGO1920	JAY GO Cushion 19" x 20"	\$200	E2603				
<input type="checkbox"/> JGO1921	JAY GO Cushion 19" x 21"	\$200	E2603				
<input type="checkbox"/> JGO1922	JAY GO Cushion 19" x 22"	\$200	E2603				
<input type="checkbox"/> JGO1923	JAY GO Cushion 19" x 23"	\$200	E2603				
<input type="checkbox"/> JGO1924	JAY GO Cushion 19" x 24"	\$200	E2603				
<input type="checkbox"/> JGO2015	JAY GO Cushion 20" x 15"	\$200	E2603				
<input type="checkbox"/> JGO2016	JAY GO Cushion 20" x 16"	\$200	E2603				
<input type="checkbox"/> JGO2017	JAY GO Cushion 20" x 17"	\$200	E2603				
<input type="checkbox"/> JGO2018	JAY GO Cushion 20" x 18"	\$200	E2603				
<input type="checkbox"/> JGO2019	JAY GO Cushion 20" x 19"	\$200	E2603				
<input type="checkbox"/> JGO2020	JAY GO Cushion 20" x 20"	\$200	E2603				
<input type="checkbox"/> JGO2021	JAY GO Cushion 20" x 21"	\$200	E2603				
<input type="checkbox"/> JGO2022	JAY GO Cushion 20" x 22"	\$200	E2603				
<input type="checkbox"/> JGO2023	JAY GO Cushion 20" x 23"	\$200	E2603				
<input type="checkbox"/> JGO2024	JAY GO Cushion 20" x 24"	\$200	E2603				
<input type="checkbox"/> JGO2117	JAY GO Cushion 21" x 17"	\$200	E2603				
<input type="checkbox"/> JGO2118	JAY GO Cushion 21" x 18"	\$200	E2603				
<input type="checkbox"/> JGO2119	JAY GO Cushion 21" x 19"	\$200	E2603				
<input type="checkbox"/> JGO2120	JAY GO Cushion 21" x 20"	\$200	E2603				

JAY GO MODIFICATIONS

MOD	Description	Price	MOD CODE
<input type="checkbox"/>	Custom Emboidery Must specify location and font? Location _____ Font _____ Includes up to 32 characters in one single spot Write specific text in "NOTES" box below	\$25	EMBROIDERY
<input type="checkbox"/>	Notch for Seat Rails 1.25" from bottom of cushion removed 1.25" from side of cushion removed 	NC	NR
<input type="checkbox"/>	Notch for Posts (Back Canes) 1" from side of cushion removed Must specify 1,2 or 3" from back of cushion  Notch depth _____ No modification to inner or outer cover	NC	NP
<input type="checkbox"/>	Anterior or Posterior Slope Must specify posterior or anterior and height Foam wedge glued underneath cushion, covers sewn to match Anterior or Posterior _____ Height increase _____	\$25	WGP or WGA
<input type="checkbox"/>	Leg Length Cut	NC	LLCR/LLCL
<input type="checkbox"/>	Contracture Cut (1" x 45°)	NC	CCT
<input type="checkbox"/>	Leg Trough Cut (1" std)	NC	LTC
<input type="checkbox"/>	Custom Drop Seat	N/A	CDS
<input type="checkbox"/>	Medial Support Cut	NC	MSC
<input type="checkbox"/>	Color Cover	\$25	CC
<input type="checkbox"/>	Scrotal Cut	N/A	SC

Notes:

Please visit SunParts Online for spare parts information



Ordering Information PHONE (800) 333-4000, Fax (800) 300-7502, WEB: www.sunrisemedical.com

Specifications Subject to change without notice. Copyright © 2020 Sunrise Medical (8.11) Jay GO Cushion 119288 Rev B Page 2 of 2