



Kid Kart Xpress
January 2022
 EIXP

Date: _____ Dealer Acct #: _____ Dealer: _____ Dealer Contact: _____ Dealer Address: _____ Dealer City: _____ ST: _____ ZIP: _____ Dealer Phone: () _____ Fax: () _____ Confirm Via: <input type="checkbox"/> Fax <input type="checkbox"/> Email Mark For: _____	SUBMITTING FOR: <input type="checkbox"/> Quote <input type="checkbox"/> Order PO#: _____ Ship To: _____ Attention: _____ Address: _____ Address: _____ Ship To City: _____ ST: _____ ZIP: _____ Ship To Phone: () _____ Fax: () _____
--	---

The HCPCS CODES herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative HCPCS CODES, including misc. codes to ensure access for their enrollees. The use of HCPCS CODES does not ensure coverage or payment.

1. FRAME

Option #	HCPCS	Definition	Price
<input type="checkbox"/> 103BS1	E1232	Kid Kart Xpress with Tilt <i>(incl. standard planar seat/back, top shoulder measurements <11")</i>	\$3,155
<input type="checkbox"/> 103BS5	E1234	Kid Kart Xpress with Tilt no seating	\$2,982
Transit System			
<input type="checkbox"/> 103KJ3	K0108	Transit Option <i>(only available at time of order; not retrofittable)</i>	\$322

Frame Options			
<input type="checkbox"/> 103PG2		Inner Back Shell - Tall <i>(seat to top of shoulder measurements greater than 11")</i>	NC
<input type="checkbox"/> 103KJ4	E1014	Recline Feature	\$215
<input type="checkbox"/> 103A76	E1029	Ventilator and Battery Tray <i>(requires 14" and 8" wheel package)</i>	\$473

2. SIZE/SEATING
 Choose Standard OR Adjustable

STANDARD			
<input type="checkbox"/> 103CU28		Standard Cushion w/o Lateral Supports	NC
<input type="checkbox"/> 103CU11		Standard Cushion	NC
<input type="checkbox"/> 103CU12	E2294	Wedge Cushion	\$242
<input type="checkbox"/> 103CU13	E2294	Ledge Cushion - Medium Density Contoured Foam - 6" to 12.5" Hip Width - 3" Lateral Hip to Thigh Supports <i>Not avail. w/o lateral support option</i> - Slight Medial Thigh Contour - Fixed Ischial Ledge - 2" Seat Depth Growth Kit - Removable Neoprene Cover	\$242

PERSONALIZED ADJUSTABLE			
<input type="checkbox"/> 103ST3	E2294	Adjustable Contour Cushion - Medium Density Contoured Foam - Adjustable Pre-Ischial Ledge - Adjustable 1" High Foam Wedges for Lateral Hip and Thigh Support - Medial Knee and Thigh Support - 3" Seat Depth Growth Kit - Removable Neoprene Cover	\$343

Seat Depth Options			
<input type="checkbox"/> 103SD26		Short 5.5" through 9"	NC
<input type="checkbox"/> 103SD28		Long 7.5" through 11"	NC
Client Actual Seat Depth:		_____ inches	
Set Seat Depth at:		_____ inches	

Seat Depth Options			
<input type="checkbox"/> 103SD27		Short 5.5" through 10"	NC
<input type="checkbox"/> 103SD29		Long 7.5" through 11"	NC
Foam Options			
<input type="checkbox"/> 103FM1		Conform Soft Foam	NC
<input type="checkbox"/> 103FM2		Conform Medium/Soft Foam	NC
<input type="checkbox"/> 103FM3		Conform Medium Foam	NC
Client Actual Seat Depth:		_____ inches	
Set Seat Depth at:		_____ inches	

3. WHEELS AND CASTERS

<input type="checkbox"/> 103WQ05		10" Semi-Pneu Rear w/ 7" Casters	NC
<input type="checkbox"/> 103WQ06		12" Pneumatic Rear w/ 7" Casters	NC
<input type="checkbox"/> 103WQ07	E2211 E2213	12" Pneu Rear w/ Airless Inserts with 7" Casters	\$81
<input type="checkbox"/> 103WQ08	E2211 E2212	14" Pneu Rear with 8" Casters	\$145

<input type="checkbox"/> 103WQ09	E2211 E2213	14" Pneu Rear w/ Airless Inserts with 8" Casters	\$225
<input type="checkbox"/> 103WQ10		14" Pneu Rear w/ 8" Casters w/ Vent option	NC
<input type="checkbox"/> 103WQ11	E2211 E2213	14" Pneu Rear w/ Airless Inserts with 8" Casters w/ Vent Option	\$81

4. BACKREST
 Choose 1" OR 2.5"

" T " Back Cushion 1" Thick			
<input type="checkbox"/> 103BT35		Short 8" to 10"	NC
<input type="checkbox"/> 103BT28		Medium 10" to 12"	NC
<input type="checkbox"/> 103BT29		Tall 12" to 14"	NC
<input type="checkbox"/> 103BT56		Extra Tall 14" to 16"	NC

" T " Back Cushion 2.5" Thick			
<i>2.5" "T" Back used for short cushion seat depths less than 7.5" and long cushion seat depths less than 9.5"</i>			
<input type="checkbox"/> 103BT36	K0108	Short 8" to 10"	\$97
<input type="checkbox"/> 103BT30	K0108	Medium 10" to 12"	\$97
<input type="checkbox"/> 103BT31	K0108	Tall 12" to 14"	\$97
<input type="checkbox"/> 103BT57	K0108	Extra Tall 14" to 16"	\$97

5. SUPPORTS

**Priced as each, selecting these options indicates Qty 2. Hardware Bracket must be chosen with any support option as pair.*

Lateral Hip Support

Not available w/standard cushion.

<input type="checkbox"/>	103KC4/6	E0956	SM (4" x 4", 3/8" thick)*	\$107
<input type="checkbox"/>	103KC3/5	E0956	LG (5-3/8" x 4-3/4", 1/2" thick)*	\$107
<input type="checkbox"/>	103HE68	E1028	Hardware Hip Support Brackets	\$54

Lateral Thigh Support

Not available w/standard cushion.

<input type="checkbox"/>	103KF4/6	E0953	SM (4" x 4", 3/8" thick)*	\$107
<input type="checkbox"/>	103KF3/5	E0953	LG (5-3/8" x 4-3/4", 1/2" thick)*	\$107
<input type="checkbox"/>	103HE69	E1028	Hardware Lateral Thigh Brackets	\$54

Medial Support

<input type="checkbox"/>	103KH1	E0957	SM (3"H x 2-1/2"W x 4.5"L)	\$188
<input type="checkbox"/>	103KH2	E0957	LG (3-3/4"H x 4"W x 6"L)	\$188
<input type="checkbox"/>	103HE70	E1028	Hardware Medial Support Bracket	\$54

Lateral Trunk Support

Fixed Laterals (Adjustable 1"-10" widths)

<input type="checkbox"/>	103KE24		Fixed Angle Adjustable - Right	NC
<input type="checkbox"/>	103KE25		Fixed Angle Adjustable - Left	NC
<input type="checkbox"/>	103KD15	E0956	Large Straight Pad 3" x 6" - Right	\$107
<input type="checkbox"/>	103KD19	E0956	Large Straight Pad 3" x 6" - Left	\$107
<input type="checkbox"/>	103KD13	E0956	Small Straight Pad 2" x 6" - Right	\$107
<input type="checkbox"/>	103KD17	E0956	Small Straight Pad 2" x 6" - Left	\$107

Lateral Trunk Support - Cont.

Swing Away Laterals (Adjustable 1"-8" widths)

<input type="checkbox"/>	103KE26		Swing Away Bracket Adjustable-Right	\$188
<input type="checkbox"/>	103KE27		Swing Away Bracket Adjustable-Left	\$188
<input type="checkbox"/>	103KD15	E0956	Large Straight Pad 3" x 6" - Right	\$107
<input type="checkbox"/>	103KD19	E0956	Large Straight Pad 3" x 6" - Left	\$107
<input type="checkbox"/>	103KD16	E0956	Large Curved Pad 3" x 6" - Right	\$107
<input type="checkbox"/>	103KD20	E0956	Large Curved Pad 3" x 6" - Left	\$107
<input type="checkbox"/>	103KD13	E0956	Small Straight Pad 2" x 6" - Right	\$107
<input type="checkbox"/>	103KD17	E0956	Small Straight Pad 2" x 6" - Left	\$107
<input type="checkbox"/>	103KD14	E0956	Small Curved Pad 2" x 6" - Right	\$107
<input type="checkbox"/>	103KD18	E0956	Small Curved Pad 2" x 6" - Left	\$107

Anterior Trunk Support

Dimensions are seat to shoulder height.

<input type="checkbox"/>	103K11	E0960	XX Small Butterfly Pad (< 10")	\$118
<input type="checkbox"/>	103K12	E0960	X Small Butterfly Pad (10" - 12")	\$118
<input type="checkbox"/>	103K13	E0960	Small Butterfly Pad (12" - 14")	\$118
<input type="checkbox"/>	103K14	E0960	Medium Butterfly Pad (14" - 16")	\$118
<input type="checkbox"/>	103K15	E0960	Large Butterfly Pad (>16")	\$118

Trunk Support Options

<input type="checkbox"/>	103SS1	E0960	Shoulder Strap Pads	\$27
--------------------------	--------	-------	---------------------	------

6. KEY OPTIONS

Head Support Options

When choosing a headrest, Hardware Bracket must be chosen as well

<input type="checkbox"/>	103HE25	E0955	Whitmyer 8" Plush	\$220
<input type="checkbox"/>	103HE26	E0955	Whitmyer 10" Plush	\$220
<input type="checkbox"/>	103HE27	E0955	Whitmyer SM Contour Cradle	\$220
<input type="checkbox"/>	103HE28	E0955	Whitmyer Med Contour Cradle	\$220
<input type="checkbox"/>	103HE29	E0955	Whitmyer Sm Adjust-a-Plush	\$289
<input type="checkbox"/>	103HQ03	K0108	Whitmyer DFS Squared for Single pad, must order Whitmyer headrest	\$128
<input type="checkbox"/>	103HE31	E1028	Hardware Bracket Only (can be chosen w/o headrest)	\$113

Shoe Holders

<input type="checkbox"/>	103A68	E0952 E0951	6" Long Shoe Holder - Pair	\$128
<input type="checkbox"/>	103A69	E0952 E0951	8" Long Shoe Holder - Pair	\$128

Foot Straps

<input type="checkbox"/>	103A79	E0952	16" Foot Straps with Covers - Pair	\$38
<input type="checkbox"/>	103A80	E0952	18" Foot straps with Covers - Pair	\$38

Upper Extremity Support Tray

<input type="checkbox"/>	103A41	E0950	ABS (Light Grey)	\$135
<input type="checkbox"/>	103A39	E0950	Polycarbonate (Clear)	\$204

Pelvic Positioning Belt

Must use when ordering Anterior Trunk Supports.

<input type="checkbox"/>	103A14	E0978	1" Padded	\$48
<input type="checkbox"/>	103A15	E0978	1-1/2" Padded	\$48
<input type="checkbox"/>	103A263	E0978	1" Padded Auto Style	\$59
<input type="checkbox"/>	103A256	E0978	1-1/2" Padded Auto Style	\$59
<input type="checkbox"/>	103PP1		Set at 45°	NC
<input type="checkbox"/>	103PP2		Set at 90° (90° for Adjustable Contour Cushion Only)	NC

Floor Sitter Feeder Base

<input type="checkbox"/>	103A81	E1399	Folding Fixed Height Feeder Base	\$725
<input type="checkbox"/>	103A78	E1399	Folding Booster Base (Ht. Adj.)	\$671

Life Support Accessories

<input type="checkbox"/>	103A462		Rain Cape	\$188
<input type="checkbox"/>	103A74	K0105	IV Support Pole	\$156
<input type="checkbox"/>	103A711	E2208	O2 Tank Holder	\$156
<input type="checkbox"/>	103A75	K0105 E2208	IV / O2 Combination	\$236

7. COLORS

Package includes: Silver Frame and Black Upholstery

Color Package

Canopy, shell and catch all

<input type="checkbox"/>	103CL9	Purple	NC
<input type="checkbox"/>	103CL10	Soft Pink	NC
<input type="checkbox"/>	103CL11	Green	NC
<input type="checkbox"/>	103CL12	Navy Blue	NC
<input type="checkbox"/>	103CL13	Orange	NC
<input type="checkbox"/>	103CL14	Butterflies	NC
<input type="checkbox"/>	103CL15	Red Race Flags	NC

Canopy Pattern

<input type="checkbox"/>	103UF01	Cross Over	NC
<input type="checkbox"/>	103UF02	Center Stripe	NC



Sunrise Medical (US) LLC

2842 Business Park Ave. · Fresno, CA 93727 · USA